# ESR CONTROL CUBE

**REF** 10435

**REF** 10436

DIESSE Diagnostica  
Senese S.p.A.  
Via delle Rose, 10  
53035 Monteriggioni (Siena)  
Italy

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## 11. BIBLIOGRAFIA

1. Clinical and Laboratories Standard Institute, H26-A2, Validation, Verification and Quality assurance of automated Haematology analyzer. Approved standard – second edition.

## INSTRUCTIONS FOR USE

### ESR CONTROL CUBE

**REF** 10435      4 x 9 ml      (2 x Level I + 2 x Level II)

**REF** 10436      2 x 9 ml      (1 x Level I + 1 x Level II)

### For the hematological control of the determination of the Erythrocyte Sedimentation Rate (ESR)

#### For in vitro diagnostic use only

#### 1. INTENDED USE

Hematological control of the determination of the Erythrocyte Sedimentation Rate (ESR).

#### 2. INTRODUCTION AND PRINCIPLE OF THE METHOD

It is good laboratory practice to use control materials to check the precision of the instruments and methods used.

The ESR Control Cube is a material to be used as a sample in the control of precision of the VES-Matic Cube models and of the Cube Next models (MINI-CUBE e Cube 30 touch).

#### 3. REAGENTS

Contents: Suspension of stabilized human red blood cells in a preservative medium.

Preparation: ready for use

Stability: Store the product at 2-10°C.

- The unopened product is stable up to the expiry date and is stable for 95 days at 2-30°C after opening.  
Close the vial carefully after use and store in an upright position.
- Once prepared in the tube, the product is stable for 7 days at 18-30°C.

#### 4. ADDITIONAL MATERIALS

Bar-code labels reporting lot number, expiration date and range of the expected results for both levels (VES-Matic Cube 200/80; VES-Matic Cube 30, MINI-CUBE/Cube 30 touch).

#### 5. PRECAUTIONS

##### FOR IN VITRO DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and found negative with FDA-approved methods, for the presence of HbsAg and for the presence of anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee of the absence of infective agents, all materials of human origin must be considered potentially infectious.

All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

**Waste disposal:** material of biological origin once used must be treated as infectious residuals and eliminated according to law.

#### 6. WARNINGS

- Check that the batch number on the vial corresponds to the lot number reported on the Certificate of Analysis (CoA) enclosed with the present Instructions for Use.
- Avoid exposure to direct sunlight. After use, clean the screw part of the cap and the vial, and close the cap tightly; store in an upright position.
- DO NOT FREEZE AND DO NOT EXPOSE TO TEMPERATURES OVER 30°C.**
- Good function of the product is guaranteed only if it is stored and used as described in the present Instructions for Use.
- It may occur that in some laboratories the results obtained do not fall within the declared range. In these cases it is advisable to perform several tests following the laboratory's own method and thus establish a personal range of acceptability. The laboratory can thus consider this range valid until there is a change of lot, or a change in the execution method for the test.**
- The supernatant in the unmixed vial may appear clear red; this is normal and does not indicate that the product is unsuitable for use.
- A loss of colour of the product may be a sign of degradation. **Do not use the product if a loss of quality is suspected.**

#### 7. PROCEDURE

##### IMPORTANT: Mix the vial contents carefully

- Remove the vial from the refrigerator and take to room temperature (20-30 minutes).
- To mix:
  - Prior to first time use carefully resuspend the product contained in the unopened vial using vortex (up to 60 seconds).
  - Hold the vial in vertical position and roll it between the palms of the hands for 15-20 seconds.

- Keep mixing the vial, turning it completely and rapidly upside down, for 20 times at least.

- d. Check the bottom of the vial to verify that the product is completely resuspended and repeat steps 2b and 2c if necessary.  
Allow bubbles to disperse and mix again by inversion 8-10 times immediately before the sampling.
  - e. Divide the sample in aliquots immediately after mixing.
3. Transfer 2.5 – 3.0 mL of suspension into an EDTA (Top Lavender) tube.
  4. Always apply to the tube the bar-code label relative to the type of ESR Control Cube transferred into it: Normal/Abnormal, VES-Matic Cube 200/80; VES-Matic Cube 30, MINI-CUBE/Cube 30 touch.  
**Carefully check to make sure to use the correct labels for the instrument.**
  5. Insert the tube in the instrument like any other sample (see the operating manual of the instrument).
  6. Compare the result found with the value reported on the CoA.

## 8. EXPECTED RESULTS

The ranges of acceptability for the different instruments are reported on the enclosed CoA.

If the results of the control are not within the range reported on the CoA but are reproducible, the relative working conditions may be the cause of the values not falling within the given range. In this case it is advisable to perform several runs following the laboratory's usual method and thus determine the relative acceptable range.

If the results of the control are out of the range reported in the CoA and are not reproducible, it is advisable to repeat the control using a new batch.

It must be remembered that the ESR Control Cube is for use to check the precision of the system (see the paragraph on WARNINGS) and not as a reference standard.

## 9. PERFORMANCE

The minimum and maximum values reported in the table on the CoA are derived from replicate analysis on automated methods.

## 10. LIMITATIONS

The ESR Control Cube is not intended for use as a reference standard




The ESR Control Cube cannot be used to control any hematological parameters other than the ESR.

## 11. REFERENCES

1. Clinical and Laboratories Standard Institute, H26-A2, Validation, Verification and Quality assurance of automated Haematology analyzer. Approved standard – second edition.

**DIESSE Diagnostica Senese S.p.A.**  
Via delle Rose 10  
53035 Monteriggioni (Siena)  
Italy

	EN ES IT PT	Date of manufacture Fecha de fabricación Data di fabbricazione Data de fabrico	GK FR RO DE	Ημερομηνία Παραγωγής Date de fabrication Data fabricatiei Herstellungsdatum
	EN ES IT PT	Use By Fecha de caducidad Utilizzare entro Prazo de validade	GK FR RO DE	Ημερομηνία λήξης Utiliser jusque A se folosi pana la Verwendbar bis
	EN ES IT PT	Caution, consult accompanying documents Atención, ver instrucciones de uso Attenzione, vedere le istruzioni per l'uso Atenção, consulte a documentação incluída	GK FR RO DE	Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα Attention voir notice d'instructions Atentie, consultați documentele insotitoare Achtung, Begleitdokumente beachten
	EN ES IT PT	Manufacturer Fabricante Fabbricante Fabricante	GK FR RO DE	Κατασκευαστής Fabricant Productator Hersteller
	EN ES IT PT	Temperature limitation Límite de temperatura Limiti di temperatura Limites de temperatura	GK FR RO DE	Περιορισμοί θερμοκρασίας Limites de température Limita da temperatura Temperaturbegrenzung
	EN ES IT PT	Consult Instructions for Use Consulte las instrucciones de uso Consultare le istruzioni per l'uso Consulte as instruções de utilização	GK FR RO DE	Συμβουλευτείτε τις οδηγίες χρήσης Consulter les instructions d'utilisation Pentru utilizare consultați instructiunile Gebrauchsanweisung beachten
	EN ES IT PT	Biological risks Riesgo biológico Rischio biologico Risco biológico	GK FR RO DE	Βιολογικοί κίνδυνοι Risques biologiques Risk biologic Biogefährdung
	EN ES IT PT	Catalogue number Número de catálogo Numero di catalogo Referência de catálogo	GK FR RO DE	Αριθμός καταλόγου Référence du catalogue Numar de catalog Bestellnummer
	EN ES IT PT	In Vitro Diagnostic Medical Device Producto sanitario para diagnóstico in vitro Dispositivo medico-diagnostico in vitro Dispositivo médico para diagnóstico in vitro	GK FR RO DE	In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν Dispositif médical de diagnostic in vitro Dizpositiv medical pentru diagnosticare <i>in vitro</i> In-Vitro-Diagnostikum
	EN ES IT PT	Batch code Código de lote Codice del lotto Código do lote	GK FR RO DE	Αριθμός Παρτίδας Code du lot Lot Chargenbezeichnung
	EN ES IT PT	Level I Nivel I Livello I Nível I	GK FR RO DE	Επίπεδο I Niveau I Nivel I Stufe I
	EN ES IT PT	Level II Nivel II Livello II Nível II	GK FR RO DE	Επίπεδο II Niveau II Nivel II Stufe II