SEDICHECK CUBE 30

REF 10439

DIESSE Diagnostica
Senese S.p.A.
Via delle Rose, 10
53035 Monteriggioni (Siena)
Italy

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SEDICHECK CUBE 30

REF 10439 4 x 3 ml (2 x Level I + 2 x Level II)

For the functionality control of the VES-Matic Cube 30, MINI-CUBE and Cube 30 touch instruments

For in vitro diagnostic use only

1. INTENDED USE

Functionality control of the VES-Matic Cube 30, MINI-CUBE and Cube 30 touch instruments.

2. INTRODUCTION AND PRINCIPLE OF THE METHOD

It is good laboratory practice to use materials of control to check the precision of the methods used.

The product SediCheck Cube 30 is a material of control to be used to monitor the precision of the Ves-Matic Cube 30, MINI-CUBE and Cube 30 touch instruments.

3. WARNINGS

1. Check that the batch number on the tube corresponds to the lot number reported on the Certificate of Analysis (CoA) enclosed with the present Instructions for Use.
2. **DO NOT FREEZE AND DO NOT EXPOSE TO TEMPERATURES OVER 30°C.**
3. Good function of the product is guaranteed only if it is stored and used as described in the present Instructions for Use.
4. To resuspend the product turn the tube upside down repeatedly.
5. **VES-Matic Cube 30: USE ONLY POSITION 1 AND 2 FOR THE CONTROL.**

4. REAGENTS

Contents:

Level I - stabilized suspension of photo-absorbing particles.

Preservative: glycerol > 40%.

Level II - stabilized suspension of photo-absorbing particles.

Preservative: ethanol 20%.

Preparation: ready for use

Stability: Store at 2 - 30°C.

The unopened product is stable up to the expiry date.

5. PROCEDURE

1. **During first use turn the tube upside down repeatedly until complete resuspension. Check that the bottom and the walls of the tube do not contain residuals of the product that are not resuspended. Keep shaking for 2 minutes at least. DO NOT USE THE VORTEX.**

2. Insert the tube in the instrument like any other sample (see the operating manual of the instrument), **but only using first and second position (VES-Matic Cube 30).**
3. Compare the result found with the value reported on the CoA. If the result obtained is within the range reported on the CoA and on the barcode the instrument will print the data without any alarm. Out of range values will be reported with an asterisk.
4. Do not use the same tube for more than 14 consecutive days.

6. EXPECTED RESULTS

If the results of the control are not within the range reported on the CoA but are reproducible, the relative working conditions may be the cause of the values not falling within the given range. In this case it is advisable to perform several runs following the laboratory's usual method and thus determine the relative acceptable range.

If the results of the control are out of the range reported in the CoA and are not reproducible, it is advisable to repeat the control using a new batch.

It must be remembered that SediCheck Cube 30 is used to check the precision of the system.

7. LIMITATIONS

SediCheck Cube 30 is not intended for use as a reference standard.



SediCheck Cube 30 cannot be used to control any hematological parameters other than the ESR.

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DIESE Diagnostica Senese S.p.A.