CHORUS

**CYTOMEGALOVIRUS
IgG AVIDITY**

REF 81092

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INSTRUCTIONS FOR USE

CHORUS CYTOMEGALOVIRUS IgG AVIDITY

**For the qualitative determination of the avidity of
IgG-class antibodies to Cytomegalovirus**

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the qualitative determination of the avidity of IgG-class antibodies to Cytomegalovirus in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Cytomegalovirus is a herpes virus transmitted by close human contact. No symptoms of infection are apparent in the majority of cases. However, the virus is very dangerous and may be fatal in immunodepressed patients. Serum-negative female patients who become infected during pregnancy may transmit the disease to the fetus. In 95% of cases this occurs without symptoms, but some neonates may present jaundice, hepatosplenomegaly and retarded psycho-motorial development. For this reason it is of great importance to determine the immunitary state of the patient and check for serum conversion.

Sometimes it is not sufficient to test for IgM antibodies, as the presence of this class may be due to the persistence of IgM or asymptomatic reinfection by Cytomegalovirus, without risk for the fetus.

For this reason it is useful to assay the avidity of the IgG antibodies, in as far as it has been demonstrated that avidity is low in the acute phase of infection and increases over a period of time. The presence of low avidity antibodies is therefore an indication of a recent or current infection. The avidity of IgG antibodies in the serum can be assayed with this ELISA kit.

3. PRINCIPLE OF THE METHOD

The Chorus Cytomegalovirus IgG Avidity device is ready to use for the detection of the avidity of IgG antibodies against Cytomegalovirus, in the Chorus/ Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay).

A dissociating agent causes the detachment of human anti-Cytomegalovirus IgG which are adhered to the antigen bound on the solid phase. This detachment is inversely proportional to the avidity; the percentual ratio is calculated compared to a non treated well, and the avidity value of the IgG is obtained.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.

The results are expressed in percentage of avidity.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.


4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.
5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 12 tests.

DD DEVICES 2 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: Avidity Buffer

Contents: Urea in saline buffer solution containing phenol (0.05%) and Bronidox (0.02%).

Position 6-5: Microplate wells coated with Cytomegalovirus

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing phenol 0.05%, Bronidox 0.02% and indicator to reveal the presence of the serum.

Position 2: CONJUGATE

Contents: anti-human IgG monoclonal antibodies labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

In which the operator must place the undiluted serum

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CONTROL HIGH High Avidity Control Serum 1 x 0.3 ml

Contents: Diluted human serum, containing high avidity antibodies and preservative. Liquid ready for use.

SAMPLE DILUENT Sample Diluent 1 x 1.5 ml

Contents: IgG negative serum to be used as sample diluent for those which give a result over range.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 – 83608
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
HIGH AVIDITY	
CONTROL SERUM	8 weeks at 2/8°C
SAMPLE DILUENT	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.

- Dispense 50 µl of undiluted test serum in well no. 1 of each device.
- Place the devices in the Chorus/Chorus TRIO instrument. Perform the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the test must be repeated.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in percentage of avidity.

A percentage higher than 40% indicates that the IgG have high avidity. Values below 30% indicate low avidity. Intermediate values indicate medium avidity IgG.

When the sign NC appears on the instrument the IU/ml of the sample are below 1.5 IU/ml.

When the sign DIL appears on the instrument the IU/ml are over 10 IU/ml; in this case dilute the sample (e.g. 1:5) in Sample Diluent and repeat the test.

Repeat the test in case of a medium avidity result.

It is advisable to repeat the test after 2-3 weeks, in order to evaluate a possible variation of percentage of avidity.

11. LIMITATIONS

The test is applicable to samples containing at least 2 IU/ml of IgG anti-Cytomegalovirus and less than 10 IU/ml.

The presence of high avidity IgG (particularly in transplanted patients), even if it does not absolutely exclude the possibility of a recent infection, is strongly indicative of an infection contracted 5 months earlier at least.

On the contrary, the presence of low avidity IgG is strongly indicative of an infection contracted during the last 4 months.

In the case of newborn babies, a low avidity is not an indication of a current acute infection.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. CROSS-REACTIONS

23 samples, positive to Epstein-Barr, Herpes simplex, Mumps and Rubella were tested.

No significant cross-reactions were found.

13. METHOD COMPARISON

In an experimentation 40 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		Low	High	Total
	Diesse	8	1	9
	High	0	31	31
	Total	8	32	40

Percent Positive Agreement (~Diagnostic Sensitivity):

100% CI_{95%}: 67.5.-99.6.

Percent Negative Agreement: (~Diagnostic Specificity):

96.9.% CI_{95%}: 84.2.-99.3

The agreement between the two methods is excellent with a Cohen's Kappa of 0.94.

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