CHORUS

Borrelia IgG CSF

DIESSE

REF 81048

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INSTRUCTIONS FOR USE

CHORUS Borrelia IgG CSF

For the qualitative determination of anti-Borrelia burgdorferi sensu latu IgG antibodies in cerebrospinal fluid

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Third generation immunoenzymatic method for the qualitative determination of IgG class antibodies against *Borrelia burgdorferi sensu latu* in human cerebrospinal fluid, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Lyme disease (Borreliosis) is a multi-systemic infectious disease caused by *Borrelia burgdorferi* bacteria belonging to the Spirochaete family.

The infection is transmitted to humans by the bite of ticks of the Ixodes genus and it can be divided in three clinical phases: early localized infection, early disseminated infection, late disseminated infection. During the first phase about two-thirds of patients develop a skin erythema together with flu-like symptoms. During the second phase, which occurs after weeks or months, the infection spreads through the bloodstream and the lymphatic system causing muscle-skeletal and neurological disorders. The late phase occurs after months or years and can provoke acrodermatitis chronica atrophicans (ACA), chronic neuroborreliosis, polyarthrititis.

IgG class antibodies appear in the cerebrospinal fluid starting from the second phase of the disease and persist also during the third one. Their detection is necessary for the diagnosis of early and late neuroborreliosis.

3. PRINCIPLE OF THE METHOD

The Chorus *Borrelia* IgG CSF device is a third generation immunoenzymatic method ready to use for the detection of IgG antibodies against *Borrelia burgdorferi sensu latu* in cerebrospinal fluid with high diagnostic sensitivity and specificity, in the Chorus/Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted sample.

After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgG antibodies conjugated to horseradish peroxidase. The unbound conjugate is eliminated, and the peroxidase substrate added. The blue colour which develops is proportional to the concentration of specific antibodies present in the sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments. The results are expressed in Index (OD sample/OD cut-off).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.


5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use the device after the expiry date.
11. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 12 tests.

DD DEVICES 2 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with selected parts of *Borrelia burgdorferi* sensu lato specific antigens VlsE, p83, internal flagellin (p41i), p39, OspC and p17.

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8).

Position 3: SAMPLE DILUENT

Contents: Proteic buffer solution at pH 7.1.

Position 2: CONJUGATE

Contents: anti-human IgG antibodies labelled with peroxidase.

Position 1: EMPTY WELL

In which the operator must place the undiluted sample

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.300 ml

Contents: Liquid, ready for use, containing anti-*Borrelia burgdorferi* IgG antibodies and preservative.

CONTROL + POSITIVE CONTROL 1 x 0.550 ml

Contents: Liquid, ready for use, containing anti-*Borrelia burgdorferi* IgG antibodies and preservative.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.

- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the positive control (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of cerebrospinal fluid handled with all precautions dictated by good laboratory practice.

The fresh sample may be stored for 7 days at 2/8°C, or frozen for longer periods at -20°C.

Do not keep the samples in auto-defrosting freezers.

Defrosted samples must be shaken carefully before use.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense in well no. 1 of each device:

SAMPLE	130 µl/device
CALIBRATOR	130 µl/device
POSITIVE CONTROL	130 µl/device

4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the positive control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in
IO-09/378-C IFU 81048 – Ed. 25.01.2016

Index (OD sample/OD cut-off).

The test on the examined sample can be interpreted as follows:

POSITIVE: when the result is > 1.2 .

NEGATIVE: when the result is < 0.8 .

DOUBTFUL/EQUIVOCAL: for all values between 0.8 and 1.2.

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new sample.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. METHOD COMPARISON

In an experimentation 55 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

Diesse		Reference		
		+	-	Total
	+	12	0	12
	-	0	43	43
	Total	12	43	55

Percent Positive Agreement (~Diagnostic Sensitivity):

100% CI_{95%}: 75.7-99.7.

Percent Negative Agreement: (~Diagnostic Specificity):

100% CI_{95%}: 91.8.-99.9.

The agreement between the two methods is excellent with a Cohen's Kappa of 1.0.

13. PRECISION AND REPEATABILITY

Sample	Within run	
	Mean (Index)	CV%
1	0.2	-
2	0.2	25.0*
3	3.5	4.6
4	2.4	4.2

Sample	Between run	
	Mean (Index)	CV%
1	0.2	25.0*
2	0.2	25.0*
3	3.1	7.7
4	2.1	11.0

Sample	Between lots		Between Instruments	
	Mean (Index)	CV%	Mean (Index)	CV%
1	0.2	20.0*	0.2	40.0*
2	0.2	30.0*	0.2	40.0*
3	3.6	8.9	3.3	10.9
4	2.4	6.3	2.2	11.8

* Artifact caused by the known fault of Variation Coefficient which becomes extremely sensitive to even very small changes in the mean when the mean value is near zero.

14. REFERENCES

1. L.A.Magnarelli, J.F.Anderson, R.C.Johnson. Cross-Reactivity in Serological Tests for Lyme disease and Other Spirochetal Infections. The Journal of Infectious disease. 1987 ; 156(1) :183-8.
2. R.A.Kalish, G.McHugh, J.Granquist, B.Shea, R.Ruthazer, A.C.Steere. Persistence of Immunoglobulin M or G Antibody Responses to Borrelia burgdorferi 10-20 Years after Active Lyme Disease. Clinical Infectious Disease. 2001 ; 33 :780-5
3. Dandache P, Nadelman RB, Erythema migrans? in Infect Dis Clin North Am, vol. 22, n° 2, 2008, pp. 235-60, DOI:10.1016/j.idc.2007.12.012, PMID 18452799.
4. B.-W.Pfister, L.Zöller, V.Brade, H.Eiffert, U.B. Göbel, G.Stanek in cooperation with H.-W.Pfister. MIQ 12/2000 Lyme Borreliosis ; Microbiological diagnosis
5. Andrew Szczepanski and Jorge L. Benach. Lyme Borreliosis: Host Responses to Borrelia burgdorferi. Microbiological Reviews, Mar. 1991, p. 21-34
6. Austin Vaz, Lisa Glickestein, Jodie A. Field, Gail Mchugh, Vijay K. Sikand, Nitin Damle, and Allen C. Steere. Cellular and Humoral Immune Responses to Borrelia burgdorferi antigens in Patients with Culture-Positive Early Lyme Disease. Infection and Immunity, Dec. 2001, p. 7437-7444

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Protected by the European Patents No.0894143 and No.1012181 licensed to TestLine Clinical Diagnostics (Brno, CZ).

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- L.A.Magnarelli, J.F.Anderson, R.C.Johnson. Cross-Reactivity in Serological Tests for Lyme disease and Other Spirochetal Infections. The Journal of Infectious disease. 1987 ; 156(1) :183-8.
- R.A.Kalish, G.McHugh, J.Granquist, B.Shea, R.Ruthazer, A.C.Steere. Persistence of Immunoglobulin M or G Antibody Responses to Borrelia burgdorferi 10-20 Years after Active Lyme Disease. Clinical Infectious Disease. 2001 ; 33 :780-5
- Dandache P, Nadelman RB, Erythema migrans? in Infect Dis Clin North Am, vol. 22, n° 2, 2008, pp. 235-60, DOI:10.1016/j.idc.2007.12.012, PMID 18452799.
- B.-W.Pfister, L.Zöller, V.Brade, H.Eiffert, U.B. Göbel, G.Stanek in cooperation with H.-W.Pfister. MIQ 12/2000 Lyme Borreliosis ; Microbiological diagnosis
- Andrew Szczepanski and Jorge L. Benach. Lyme Borreliosis: Host Responses to Borrelia burgdorferi. Microbiological Reviews, Mar. 1991, p. 21-34
- Austin Vaz, Lisa Glickestein, Jodie A. Field, Gail Mchugh, Vijay K. Sikand, Nitin Damle, and Allen C. Steere. Cellular and Humoral Immune Responses to Borrelia burgdorferi antigens in Patients with Culture-Positive Early Lyme Disease. Infection and Immunity, Dec. 2001, p. 7437-7444

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28. LITERATUR

1. L.A.Magnarelli, J.F.Anderson, R.C.Johnson. Cross-Reactivity in Serological Tests for Lyme disease and Other Spirochetal Infections. The Journal of Infectious disease. 1987 ; 156(1) :183-8.
2. R.A.Kalish, G.McHugh, J.Granquist, B.Shea, R.Ruthazer, A.C.Steere. Persistence of Immunoglobulin M or G Antibody Responses to Borrelia burgdorferi 10-20 Years

after Active Lyme Disease. Clinical Infectious Disease. 2001 ; 33 :780-5

3. Dandache P, Nadelman RB, Erythema migrans? in Infect Dis Clin North Am, vol. 22, n° 2, 2008, pp. 235–60, DOI:10.1016/j.idc.2007.12.012, PMID 18452799.
4. B.-W.Pfister, L.Zöller, V.Brade, H.Eiffert, U.B. Göbel, G.Stanek in cooperation with H.-W.Pfister. MIQ 12/2000 Lyme Borreliosis ; Microbiological diagnosis
5. Andrew Szczepanski and Jorge L. Benach. Lyme Borreliosis: Host Responses to Borrelia burgdorferi. Microbiological Reviews, Mar. 1991, p. 21-34
6. Austin Vaz, Lisa Glickestein, Jodie A. Field, Gail Mchugh, Vijay K. Sikand, Nitin Damle, and Allen C. Steere. Cellular and Humoral Immune Responses to Borrelia burgdorferi antigens in Patients with Culture-Positive Early Lyme Disease. Infection and Immunity, Dec. 2001, p. 7437–7444

DIESSE Diagnostica Senese S.p.A.
Strada dei Laghi 39
53035 Monteriggioni (SIENA)
Italien

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1. L.A.Magnarelli, J.F.Anderson, R.C.Johnson. Cross-Reactivity in Serological Tests for Lyme disease and Other Spirochetal Infections. The Journal of Infectious disease. 1987 ; 156(1) :183-8.
2. R.A.Kalish, G.McHugh, J.Granquist, B.Shea, R.Ruthazer, A.C.Steere. Persistence of Immunoglobulin M or G Antibody Responses to Borrelia burgdorferi 10-20 Years after Active Lyme Disease. Clinical Infectious Disease. 2001 ; 33 :780-5

3. Dandache P, Nadelman RB, Erythema migrans? in Infect Dis Clin North Am, vol. 22, n° 2, 2008, pp. 235–60, DOI:10.1016/j.idc.2007.12.012, PMID 18452799.
4. B.-W.Pfister, L.Zöller, V.Brade, H.Eiffert, U.B. Göbel, G.Stanek in cooperation with H.-W.Pfister. MIQ 12/2000 Lyme Borreliosis ; Microbiological diagnosis
5. Andrew Szczepanski and Jorge L. Benach. Lyme Borreliosis: Host Responses to Borrelia burgdorferi. Microbiological Reviews, Mar. 1991, p. 21-34
6. Austin Vaz, Lisa Glickestein, Jodie A. Field, Gail Mchugh, Vijay K. Sikand, Nitin Damle, and Allen C. Steere. Cellular and Humoral Immune Responses to Borrelia burgdorferi antigens in Patients with Culture-Positive Early Lyme Disease. Infection and Immunity, Dec. 2001, p. 7437–7444

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14. BIBLIOGRAFIA

- L.A.Magnarelli, J.F.Anderson, R.C.Johnson. Cross-Reactivity in Serological Tests for Lyme disease and Other Spirochetal Infections. The Journal of Infectious disease. 1987; 156(1):183-8.
- R.A.Kalish, G.McHugh, J.Granquist, B.Shea, R.Ruthazer, A.C.Steere. Persistence of Immunoglobulin M or G Antibody Responses to Borrelia burgdorferi 10-20 Years after Active Lyme Disease. Clinical Infectious Disease. 2001; 33:780-5
- Dandache P, Nadelman RB, Erythema migrans? in Infect Dis Clin North Am, vol. 22, nº 2, 2008, pp. 235-60, DOI:10.1016/j.idc.2007.12.012, PMID 18452799.
- B.-W.Pfister, L.Zöller, V.Brade, H.Eiffert, U.B. Göbel, G.Stanek in cooperation with H.-W.Pfister. MIQ 12/2000 Lyme Borreliosis; Microbiological diagnosis
- Andrew Szczepanski and Jorge L. Benach. Lyme Borreliosis: Host Responses to Borrelia burgdorferi. Microbiological Reviews, Mar. 1991, p. 21-34
- Austin Vaz, Lisa Glickestein, Jodie A. Field, Gail Mchugh, Vijay K. Sikand, Nitin Damle, and Allen C. Steere. Cellular and Humoral Immune Responses to Borrelia burgdorferi antigens in Patients with Culture-Positive Early Lyme Disease. Infection and Immunity, Dec. 2001, p. 7437-7444

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