**CHORUS**

**U1-70 RNP**

**REF 86088**

**REF 86088/12**

DIESSE Diagnostica Senese  
S.p.A.  
Via delle Rose, 10  
53035 Monteriggioni (Siena)  
Italy

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## **INSTRUCTIONS FOR USE**

### **CHORUS U1-70 RNP**

#### **For the semiquantitative determination of IgG antibodies anti-U1-70**

#### **For *In Vitro* Diagnostic Use Only**

##### **1. INTENDED USE**

Immunoenzymatic method for the semiquantitative determination of IgG class antibodies against U1-70 in human serum, using a disposable device applied on the chorus and Chorus TRIO instruments.

##### **2. INTRODUCTION**

The U1-snRNP complex is a small nuclear ribonucleoprotein particle composed of uridine rich (thus U) small nuclear RNA and a set of proteins: the 70 kDa U1-specific ribonucleoprotein plus proteins A and C (all formerly summarized as RNPs) and the Sm (Smith antigen which comprises eight proteins: B/B', D1, D2, D3, E, F and G. Because of its protein components, the Sm and RNPs, the complex is often named RNP/Sm complex. The U1-sn RNP complex is a part of the splicosomal complex, that facilitates the processing of pre-mRNA to mature mRNA in the nucleus.

Antibodies against the 70 kDa protein of the U1-snRNP complex belong to the heterogenous group of anti-nuclear antibodies (ANA), which are associated with various autoimmune diseases. They are directed against various proteins of the nucleus. Indirect immunofluorescence test (IFT) on eucaryotic cells like HeLa has been the established method for the detection of ANAs. Single antibody specificities are distinguished by fluorescence patterns but more specific testing by ELISAs employing the target antigens are available too for a simple and reliable differentiation of ANAs.

Antibodies against the 70 kDa U1-specific protein are found in 95% of mixed connective tissue diseases (MCTD) but do also occur in systemic lupus erythematosus (SLE) with a prevalence of 40%. Isolated occurrence of anti-70 kDa antibodies is typical for the Sharp syndrome, whilst antibodies against Sm are highly specific for systemic lupus erythematosus (SLE) and thus are included in diagnostic and classification criteria for SLE. Anti-SM antibodies are found in 20-30% of patients with SLE.

##### **3. PRINCIPLE OF THE METHOD**

The Chorus U1-70 device is ready to use for the detection of IgG antibodies against U1-70, in the Chorus/Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay). The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human serum. After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human immunoglobulins antibodies conjugated to horse radish peroxidase.

The unbound conjugate is eliminated and the peroxidase substrate is added. The colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test in the Chorus/Chorus TRIO instruments.

The results are expressed in Arbitrary Units (AU/ml).

##### **4. WARNINGS AND PRECAUTIONS**

###### **FOR *IN VITRO* DIAGNOSTIC USE ONLY**

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.

**Waste disposal: serum samples, calibrators and strips, once used, must be treated as infectious residuals and eliminated according to law.**

###### **Health and Safety Information**

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

**Analytical Precautions**

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged. Do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use and the Instrument Operating Manual must be carefully followed.
5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

## 5. KIT COMPOSITION AND REAGENT PREPARATION


The kit is sufficient for 36 tests (REF 86088).

The kit is sufficient for 12 tests (REF 86088/12).

**DD DEVICES**

6 packages each containing 6 devices (REF 86088).

2 packages each containing 6 devices (REF 86088/12).

Description:

**Position 8:** Space for application of bar code label

**Position 7:** Empty

**Position 6:** MICROPLATE WELL

Coated with human recombinant 70 kDa U1-snRNP

**Position 5:** Uncoated MICROPLATE WELL

**Position 4: TMB SUBSTRATE**

Contents: Tetramethylbenzidine 0.26 mg/mL and H<sub>2</sub>O<sub>2</sub> 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

**Position 3: SAMPLE DILUENT**

Contents: Proteic solution in Tris, saline, with sodium azide

**Position 2: CONJUGATE**

Contents: anti-human IgG monoclonal antibodies labeled with horse radish peroxidase.

**Position 1: EMPTY WELL**

in which undiluted serum must be added

**Use: equilibrate a package at room temperature**, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

**CALIBRATOR** CALIBRATOR

Contents: Diluted human serum. Liquid, ready for use.

**CONTROL +** POSITIVE CONTROL

Contents: Diluted human serum. Liquid, ready for use.

**MATERIALS REQUIRED BUT NOT PROVIDED:**

- WASHING BUFFER AUTOIMMUNITY (REF 86004)
- CLEANING SOLUTION 2000 (REF 83609)
- SANITIZING SOLUTION (REF 83604 - 83608)
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT (REF 83607)
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

**6. STORAGE AND STABILITY OF REAGENTS**

**Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).**

**The expiry date is printed on each component and on the kit label.**

**Reagents have a limited stability after opening:**

|                  |                  |
|------------------|------------------|
| DEVICES          | 8 weeks at 2/8°C |
| CALIBRATOR       | 8 weeks at 2/8°C |
| POSITIVE CONTROL | 8 weeks at 2/8°C |

**7. SPECIMEN COLLECTION AND STORAGE**

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

## 8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the instrument Operating Manual.

## 9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554  
Fax: 0039 0577 366605  
email: scientificsupport@diesse.it

## 10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in Arbitrary Units (AU/ml) calculated on the basis of a lot-dependent graph stored in the instrument.

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is > 18 AU/ml

NEGATIVE: when the result is < 12 AU/ml

DOUBTFUL/EQUIVOCAL: for all values between 12 and 18 AU/ml

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

## 11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

## 12. CALIBRATION RANGE

Calibration range 3.0-300.0 AU/ml.

For samples > 300.0 AU/ml retest the diluted sample in the Negative Control/Sample Diluent (PF83607-not supplied with the kit).

## 13. CROSS-REACTIONS

3 samples, positive to tTG, gliadin, cardiolipin, Tg and TPO were tested.

No significant cross-reactions were found.

## 14. METHOD COMPARISON

In an experimentation 57 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

|        |      | Reference |    |      |
|--------|------|-----------|----|------|
|        |      | +         | -  | Tot. |
|        |      |           |    |      |
| Diesse | +    | 12        | 1  | 13   |
|        | -    | 0         | 44 | 44   |
|        | Tot. | 12        | 45 | 57   |

Percent Positive Agreement (~Diagnostic Sensitivity):

100% CI<sub>95%</sub>: 75.7-99.7

Percent Negative Agreement (~Diagnostic Specificity):

97.8% CI<sub>95%</sub>: 88.3-99.6

The agreement between the two methods is excellent with a Cohen's Kappa of 0.94.

## 15. PRECISION AND REPEATABILITY

| Sample | Within run   |     | Between run  |     |
|--------|--------------|-----|--------------|-----|
|        | Mean (AU/ml) | CV% | Mean (AU/ml) | CV% |
| 1      | 61.3         | 5.6 | 59.3         | 4.8 |
| 2      | 32.2         | 8.0 | 31.3         | 4.3 |
| 3      | 15.1         | 7.2 | 15.5         | 3.2 |

| Sample | Between lots and Between instruments |     |
|--------|--------------------------------------|-----|
|        | Mean (AU/ml)                         | CV% |
| 1      | 59.9                                 | 1.2 |
| 2      | 30.5                                 | 2.4 |
| 3      | 15.0                                 | 2.4 |

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**DIESSE Diagnostica Senese S.p.A.**