# CHORUS

## SPECIFIC IgE

**DIESSE**

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Disponibili su richiesta.

### 14. BIBLIOGRAFIA

1. Geha RS Human IgE. J. Allergy Clin. Immunol. 1994; 74:109-120.
2. Gordon RR et al. Immunoglobulin E and eczema- asthma syndrome in early childhood. Lancet 1982; 72-74
3. Norman PS. The clinical significance of IgE. Hospital Practice 1975 Aug; 10(8): 41-49
4. Johansson SGO, Lancet 1967; 2:951
5. Zefferstron and Johansson SGO, Allergy 198; 36:537
6. Wathrich B. Clin. Allergy 1978; 8:214
7. Barbee RA et al. Distribution of IgE in a community population sample: correlation with age, sex and allergen skin test reactivity. J. Allergy Clin. Immunology 1981; 68: 106-111
8. Kjellman IM, Anderson SGO, Roth A. Serum IgE levels in healthy children. Clin. Allergy 1976; 65:51-59
9. Halpern G.M. Markers of human allergic disease. J. Clin. Immunoassay 1983 Jun; 6 (2): 131-139
10. Hamilton RG, Adkinson NF. Quantitation of allergen-specific IgE in serum using the radioallergosorbent test. J. Clin. Immunoassay 1983 Jun; 6 (12): 147-154
11. Homberger HA, Yunginger JW. Laboratory testing in the diagnosis and management of allergic diseases. Clin. Lab. Annual 1983; 2:351-388
12. Bukley RH. Immunopharmacology of Allergic Disease, 1975; 253: 474
13. Capmon A, Dessaint JP, Capron M at al. Nature 1975; 253:474
14. Jahonsson SGO, Bennick H, Berg T. Clin. Immunol. 1972; p1

## INSTRUCTIONS FOR USE

### CHORUS SPECIFIC IgE

**For the semiquantitative determination of  
allergen-specific IgE antibodies**

**For *In Vitro* Diagnostic Use Only**

#### 1. INTENDED USE

Immunoenzymatic method for the semiquantitative determination of allergen-specific IgE class antibodies in human serum, using a disposable device applied on the Chorus TRIO instruments.

#### 2. INTRODUCTION

The allergen-specific IgE antibodies circulate in the blood of allergic patients like subjects suffering from hay fever, atopic eczema or asthma.

The research of allergen-specific IgE antibodies, in combination with other clinical data, is useful to ensure a proper identification of the substances responsible for allergies.

#### 3. PRINCIPLE OF THE METHOD

The Chorus devices are ready to use for the detection of allergen-specific IgE antibodies, in the Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay). Anti-human IgE monoclonal antibody is bound to the solid phase. The specific immunoglobulins are bound to the antibody through incubation with diluted human serum. After washing to eliminate the proteins which have not reacted, incubation is performed with the biotinylated conjugate. The biotinylated conjugate is eliminated and the Streptavidin-Peroxidase conjugate is added. The unbound conjugate is eliminated and the peroxidase substrate is added.

The enzymatic reaction is stopped by the adding of the Stop Solution, which makes the solution turn to yellow. The colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test in the Chorus TRIO instruments.

The results are expressed in Units (KU/L).

#### 4. WARNINGS AND PRECAUTIONS

##### FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory

practice should be followed when handling reagents and samples.

Calibrator and positive control contain horse serum, obtained from healthy animals. Since no method can offer complete assurance about the absence of infectious agents, it is recommended to handle with care the components containing serum.

**Waste disposal:** serum samples, calibrators and strips, once used, must be treated as infectious residuals and eliminated according to law.

#### Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

#### Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged. Do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus TRIO instrument; the Instructions for Use and the Instrument Operating Manual must be carefully followed.

**The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website (<http://www.diesse.it/en/Support/Download/strumento:39/>)**

5. Check that the Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.


8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

## 5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 12 tests.

**DD** DEVICES 2 packages each containing 6 devices

Description:

**Position 8:** Space for application of bar code label

**Position 7:** ALLERGEN-BIOTIN

Contents: allergen extract conjugated with biotin in TRIS buffer pH 8.1 containing preservatives.

**Position 6:** MICROPLATE WELL

Coated with anti-human IgE monoclonal antibodies

**Position 5:** Uncoated MICROPLATE WELL

**Position 4:** TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H<sub>2</sub>O<sub>2</sub> 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

**Position 3:** STOP SOLUTION

Contents: Sulfuric acid solution 0.3 M

**Position 2:** CONJUGATE

Contents: Streptavidin-Peroxidase in a red buffer solution (pH 5.5) containing preservative.

**Position 1:** EMPTY WELL

In which undiluted serum must be added

**Use:** equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

**CALIBRATOR** CALIBRATOR 1 x 0.250 ml

Contents: Horse serum containing human IgE antibodies and preservative. Liquid, ready for use.

**CONTROL +** POSITIVE CONTROL 1 x 0.450 ml

Contents: Horse serum containing human IgE antibodies and preservative. Liquid, ready for use.

## MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Chorus TRIO Instrument

- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

## 6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

|                  |                  |
|------------------|------------------|
| DEVICES          | 8 weeks at 2/8°C |
| CALIBRATOR       | 8 weeks at 2/8°C |
| POSITIVE CONTROL | 8 weeks at 2/8°C |

## 7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 2 days at 2/8°C, or frozen for longer periods at -20°C.

Freeze-thawing cycles have to be avoided.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

## 8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 75-100 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the instrument Operating Manual.

## 9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554  
 Fax: 0039 0577 366605  
 email: scientificsupport@diesse.it

## 10. INTERPRETATION OF THE RESULTS

The Chorus TRIO instrument expresses the result in Units (KU/L), calculated on the basis of a lot-dependent graph stored in the instrument and in allergen-specific IgE classes.

The test on the examined serum can be interpreted as follows:

| KU/L        | Allergen-specific IgE class | Interpretation           |
|-------------|-----------------------------|--------------------------|
| < 0.40      | 0                           | Absent or not detectable |
| 0.40 – 0.69 | 1                           | Low                      |
| 0.70 – 3.5  | 2                           | Moderate                 |
| 3.6 – 17.9  | 3                           | High                     |
| 18.0 – 49.9 | 4                           | Very high                |
| 50.0 - > 50 | 5                           | Very high                |

## 11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

Because there is no international standardization, both for the specific IgE units and for the reference allergens, the specific IgE concentrations obtained with various methods can be quite different. Therefore, it is advisable to take care when comparing or seeking to correlate different methods.

The specific IgE affinity constants may differ from allergen to allergen, as a consequence be aware that identical results for different allergens may not necessarily involve an equivalent clinical situation.

## 12. CALIBRATION RANGE

Calibration range 0.40 – 50.0 KU/L.

For samples > 50.0 KU/L retest the diluted sample in the Negative Control/Sample Diluent (PF83607-not supplied with the kit).

## 13. PERFORMANCE CHARACTERISTICS

Available on request.

## 14. REFERENCES

1. Geha RS Human IgE. J. Allergy Clin. Immunol. 1994; 74:109-120.
2. Gordon RR et al. Immunoglobulin E and eczema-asthma syndrome in early childhood. Lancet 1982; 72-74
3. Norman PS. The clinical significance of IgE. Hospital Practice 1975 Aug; 10(8): 41-49
4. Johansson SGO, Lancet 1967; 2:951
5. Zefferstron and Johansson SGO, Allergy 198; 36:537
6. Wathrich B. Clin. Allergy 1978; 8:214
7. Barbee RA et al. Distribution of IgE in a community population sample: correlation with age, sex and allergen skin test reactivity. J. Allergy Clin. Immunology 1981; 68: 106-111
8. Kjellman IM, Anderson SGO, Roth A. Serum IgE levels in healthy children. Clin. Allergy 1976; 65:51-59
9. Halpern G.M. Markers of human allergic disease. J. Clin. Immunoassay 1983 Jun; 6 (2): 131-139

10. Hamilton RG, Adkinson NF. Quantitation of allergen-specific IgE in serum using the radioallergosorbent test. J. Clin. Immunoassay 1983 Jun; 6 (12): 147-154
11. Homberger HA, Yunginger JW. Laboratory testing in the diagnosis and management of allergic diseases. Clin. Lab. Annual 1983; 2:351-388
12. Bukley RH. Immunopharmacology of Allergic Disease, 1975; 253: 474
13. Capmon A, Dessaint JP, Capron M et al. Nature 1975; 253:474
14. Jahonsson SGO, Bennick H, Berg T. Clin. Immunol. 1972; p1

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1. Geha RS Human IgE. J. Allergy Clin. Immunol. 1994; 74:109-120.
2. Gordon RR et al. Immunoglobulin E and eczema- asthma syndrome in early childhood. Lancet 1982; 72-74
3. Norman PS. The clinical significance of IgE. Hospital Practice 1975 Aug; 10(8): 41-49
4. Johansson SGO, Lancet 1967; 2:951
5. Zefferstron and Johansson SGO, Allergy 198; 36:537
6. Wathrich B. Clin. Allergy 1978; 8:214
7. Barbee RA et al. Distribution of IgE in a community population sample: correlation with age, sex and allergen skin test reactivity. J. Allergy Clin. Immunology 1981; 68: 106-111
8. Kjellman IM, Anderson SGO, Roth A. Serum IgE levels in healthy children. Clin. Allergy 1976; 65:51-59
9. Halpern G.M. Markers of human allergic disease. J. Clin. Immunoassay 1983 Jun; 6 (2): 131-139
10. Hamilton RG, Adkinson NF. Quantitation of allergen-specific IgE in serum using the radioallergosorbent test. J. Clin. Immunoassay 1983 Jun; 6 (12): 147-154
11. Homberger HA, Yunginger JW. Laboratory testing in the diagnosis and management of allergic diseases. Clin. Lab. Annual 1983; 2:351-388
12. Bukley RH. Immunopharmacology of Allergic Disease, 1975; 253: 474
13. Capmon A, Dessaint JP, Capron M et al. Nature 1975; 253:474
14. Jahonsson SGO, Bennick H, Berg T. Clin. Immunol. 1972; p1

|                                                                                     |                |                                                                                                                               |                |                                                                                                                                           |
|-------------------------------------------------------------------------------------|----------------|-------------------------------------------------------------------------------------------------------------------------------|----------------|-------------------------------------------------------------------------------------------------------------------------------------------|
|    | EN<br>ES<br>IT | Date of manufacture<br>Fecha de fabricación<br>Data di fabbricazione                                                          | FR<br>EL<br>PT | Date de fabrication<br>Ημερομηνία Παραγωγής<br>Data de fabrico                                                                            |
|    | EN<br>ES<br>IT | Use By<br>Fecha de caducidad<br>Utilizzare entro                                                                              | FR<br>EL<br>PT | Utiliser jusque<br>Ημερομηνία λήξης<br>Prazo de validade                                                                                  |
|    | EN<br>ES<br>IT | Do not reuse<br>No reutilizar<br>Non riutilizzare                                                                             | FR<br>EL<br>PT | Ne pas réutiliser<br>Μην κάνετε επαναληπτική χρήση<br>Não reutilizar                                                                      |
|    | EN<br>ES<br>IT | Caution, consult accompanying documents<br>Atención, ver instrucciones de uso<br>Attenzione, vedere le istruzioni per l'uso   | FR<br>EL<br>PT | Attention voir notice d'instructions<br>Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα<br>Atenção, consulte a documentação incluída       |
|    | EN<br>ES<br>IT | Manufacturer<br>Fabricante<br>Fabbicante                                                                                      | FR<br>EL<br>PT | Fabricant<br>Κατασκευαστής<br>Fabricante                                                                                                  |
|    | EN<br>ES<br>IT | Contains sufficient for <n> tests<br>Contenido suficiente para <n> ensayos<br>Contenuto sufficiente per "n" saggi             | FR<br>EL<br>PT | Contenu suffisant pour "n" tests<br>Περιεχόμενο επαρκές για «n» εξετάσεις<br>Conteúdo suficiente para "n" ensaios                         |
|    | EN<br>ES<br>IT | Temperature limitation<br>Límite de temperatura<br>Limiti di temperatura                                                      | FR<br>EL<br>PT | Limites de température<br>Περιορισμοί θερμοκρασίας<br>Limites de temperatura                                                              |
|    | EN<br>ES<br>IT | Consult Instructions for Use<br>Consulte las instrucciones de uso<br>Consultare le istruzioni per l'uso                       | FR<br>EL<br>PT | Consulter les instructions d'utilisation<br>Συμβουλευτείτε τις οδηγίες χρήσης<br>Consulte as instruções de utilização                     |
|    | EN<br>ES<br>IT | Biological risks<br>Riesgo biológico<br>Rischio biologico                                                                     | FR<br>EL<br>PT | Risques biologiques<br>Βιολογικοί κίνδυνοι<br>Risco biológico                                                                             |
|   | EN<br>ES<br>IT | Catalogue number<br>Número de catálogo<br>Numero di catalogo                                                                  | FR<br>EL<br>PT | Référence du catalogue<br>Αριθμός καταλόγου<br>Referência de catálogo                                                                     |
|  | EN<br>ES<br>IT | In Vitro Diagnostic Medical Device<br>Producto sanitario para diagnóstico in vitro<br>Dispositivo medico-diagnostico in vitro | FR<br>EL<br>PT | Dispositif médical de diagnostic in vitro<br>In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν<br>Dispositivo médico para diagnóstico in vitro |
|  | EN<br>ES<br>IT | Batch code<br>Código de lote<br>Codice del lotto                                                                              | FR<br>EL<br>PT | Code du lot<br>Αριθμός Παρτίδας<br>Código do lote                                                                                         |
|  | EN<br>ES<br>IT | CE marking of conformity<br>Marcado CE de conformidad<br>Marcatura CE di conformità                                           | FR<br>EL<br>PT | Marquage de conformité CE<br>Σημανση συμμορφωσης CE<br>Marcação CE de conformidade                                                        |