# CHORUS

## Diphtheria IgG

**REF 81266**

**REF 81266/12**

DIESSE Diagnostica Senese  
S.p.A.  
Via delle Rose, 10  
53035 Monteriggioni (Siena)  
Italy

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## 16. BIBLIOGRAPHIE

- Bell, F. et al. Combined diphtheria, tetanus, pertussis and Haemophilus influenzae type b vaccines for primary immunisation; Archives of Disease in Childhood; 75: 298-303 (1996).
- Zasada, A. et al. Seroprevalence of diphtheria toxoid IgG antibodies in children, adolescents and adults in Poland; BMC Infectious Diseases; 13: 551(2013.)

3. Frank, J.W. et al. Diphtheria overimmunization in children. Can. Med. Assoc. J., 141: 1241 (1989).
4. Lau, R.C. Detection of diphtheria toxin antibodies in human sera by ELISA. J. Hyg., 96: 415 (1986).
5. Melville-Smith, M. et al. Estimation of diphtheria antitoxin in human sera by ELISA. J. Med. Microbiol., 25: 279 (1988).
6. World Health Organization. The immunological basis for immunization series; module 2: Diphtheria – Update 2009. WHO Library Cataloguing (2009).

**DIESSE Diagnostica Senese S.p.A.**  
Via delle Rose 10  
53035 Monteriggioni (Siena)  
Italie

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1. Bell, F. et al. Combined diphtheria, tetanus, pertussis and Haemophilus influenzae type b vaccines for primary immunisation; Archives of Disease in Childhood; 75: 298-303 (1996).
2. Zasada, A. et al. Seroprevalence of diphtheria toxoid IgG antibodies in children, adolescents and adults in Poland; BMC Infectious Diseases; 13: 551(2013.)
3. Frank, J.W. et al. Diphtheria overimmunization in children. Can. Med. Assoc. J., 141: 1241 (1989).
4. Lau, R.C. Detection of diphtheria toxin antibodies in human sera by ELISA. J. Hyg., 96: 415 (1986).
5. Melville-Smith, M. et al. Estimation of diphtheria antitoxin in human sera by ELISA. J. Med. Microbiol., 25: 279 (1988).
6. World Health Organization. The immunological basis for immunization series; module 2: Diphtheria – Update 2009. WHO Library Cataloguing (2009).

**DIESSE Diagnostica Senese S.p.A.**  
Via delle Rose 10  
53035 Monteriggioni (Siena)  
Italy

## CHORUS Diphtheria IgG

**For the quantitative determination of anti-Diphtheria toxin IgG antibodies**

**For *In Vitro* Diagnostic Use Only**

### 1. INTENDED USE

Immunoenzymatic method for the quantitative determination of IgG class antibodies against Diphtheria toxin in human serum and plasma, using a disposable device applied on the Chorus and Chorus TRIO instruments.

### 2. INTRODUCTION

Diphtheria is a bacterial infectious disease which appears predominantly during the childhood. The disease leads particularly to an inflammation of the pharynx, larynx and nasal mucosa. Additionally, bacterial toxins cause via long-distance effect damages of the heart, circulation and CNS. Only the toxigenic strains are pathogenic.

The etiologic agent is the *Corynebacterium diphtheriae*.

Depending on the stage of disease, the three types 'slight, middle and serious' can be distinguished.

The appearance of Diphtheria shows a seasonal prevalence with the greatest incidence in winter. Especially non-vaccinated children will be infected.

In most cases children will be vaccinated (DTP = Diphtheria-Tetanus-Pertussis) after the third month of life. The state of immunity can be monitored by determining the IgG concentration.

### 3. PRINCIPLE OF THE METHOD

The Chorus Diphtheria IgG device is ready to use for the detection of IgG antibodies against Diphtheria toxin, in the Chorus/ Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human sample.

After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgG antibodies conjugated to horseradish peroxidase. The unbound conjugate is eliminated, and the peroxidase substrate added. The blue colour which develops is proportional to the concentration of specific antibodies present in the sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.

The results are expressed in milli-International Units (mIU/ml) calculated in reference to "Who Reference Preparation 04/496".

## 4. WARNINGS AND PRECAUTIONS

### FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

**Waste disposal:** samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

#### Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

#### Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.

**The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the**

table published on Diesse website (<http://www.diesse.it/en/Support/Download/strumento:39/>)

5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

#### 5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 81266).


The kit is sufficient for 12 tests (REF 81266/12).

#### DD DEVICES

6 packages each containing 6 devices (REF 81266).

2 packages each containing 6 devices (REF 81266/12).

#### Description:

**Position 8:** Space for application of bar code label

**Position 7:** Empty

**Position 6:** MICROPLATE WELL

Coated with Diphtheria Toxoid antigen

**Position 5:** MICROPLATE WELL

Uncoated

**Position 4:** TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H<sub>2</sub>O<sub>2</sub> 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

**Position 3:** SAMPLE DILUENT

Contents: Proteic solution containing sodium azide 0.09%

**Position 2:** CONJUGATE

Contents: anti-human IgG monoclonal antibodies labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

**Position 1:** EMPTY WELL

In which the operator must place the undiluted sample

**Use:** equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

#### CALIBRATOR CALIBRATOR 1 x 0.175 ml

Contents: Diluted human serum containing anti-Diphtheria toxin IgG antibodies and preservative. Liquid, ready for use.

#### CONTROL + POSITIVE CONTROL 1 x 0.425 ml

Contents: Diluted human serum containing anti-Diphtheria toxin IgG antibodies and preservative. Liquid, ready for use.

#### MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER (REF 83606)
- CLEANING SOLUTION 2000 (REF 83609)
- SANITIZING SOLUTION (REF 83604 - 83608)
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

#### 6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the positive control (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

#### Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

#### 7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum and plasma collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh sample may be stored for 7 days at 2/8°C, or frozen for longer periods at -20°C.

Freeze-thawing cycles have to be avoided.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

#### 8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 50 µl of undiluted test sample in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

#### 9. TEST VALIDATION

Use the positive control to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the positive

control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the positive control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554  
 Fax: 0039 0577 366605  
 email: scientificsupport@diesse.it

#### 10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in milli-International Units (mIU/ml) (Who Reference Preparation 04/496), calculated on the basis of a batch-dependent curve stored in the instrument.

The test on the examined sample can be interpreted as follows:

<10 mIU/mL	Basic immunisation recommended
10-100 mIU/mL	Booster vaccination recommended
>100 mIU/mL	Good immunity

**Attention: do not consider the qualitative report N; D/E; P on the instrument.**

#### 11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

#### 12. CALIBRATION RANGE

Calibration range 10-2000 mIU/ml

#### 13. ANALYTICAL SPECIFICITY

3 samples were spiked with the following potentially interfering factors and then tested:

Bilirubin (1.9 mg/dl – 30 mg/dl)  
 Triglycerides (31 mg/dl – 500 mg/dl)  
 Hemoglobin (0.5 mg/ml – 8 mg/ml)

The presence in the sample of the interfering substances described above does not affect the test result.

#### 14. METHOD COMPARISON

In an experimentation 119 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		Immune	Not immune	Total
Diesse	Immune	95	0	95
	Not immune	0	24	24
	Total	95	24	119

Percent Positive Agreement (~Diagnostic Sensitivity):

100% CI<sub>95%</sub>: 96.1-99.9

Percent Negative Agreement: (~Diagnostic Specificity):

100% CI<sub>95%</sub>: 86.2-99.8

The agreement between the two methods is excellent with a Cohen's Kappa of 1.0.

#### 15. PRECISION AND REPEATABILITY

Sample	Within run	
	Mean (mIU/ml)	CV%
1	83.5	6.5
2	277.3	4.7
3	325.2	9.1
4	396.6	12.9
5	1668.7	12.1

Sample	Between run	
	Mean (mIU/ml)	CV%
1	76.8	11.9
2	301.6	6.9
3	311.9	7.7
4	324.7	12.6
5	1625.8	12.2

Sample	Between lots		Between Instruments	
	Mean (mIU/ml)	CV%	Mean (mIU/ml)	CV%
1	82.6	12.0	76.8	13.6
2	283.7	10.1	313.4	8.6
3	303.1	4.2	293.7	9.4
4	385.8	9.7	315.6	11.0
5	1554.9	9.0	1478.7	13.2




#### 16. REFERENCES

- Bell, F. et al. Combined diphtheria, tetanus, pertussis and Haemophilus influenzae type b vaccines for primary immunisation; Archives of Disease in Childhood; 75: 298-303 (1996).
- Zasada, A. et al. Seroprevalence of diphtheria toxoid IgG antibodies in children, adolescents and adults in Poland; BMC Infectious Diseases; 13: 551(2013.)
- Frank, J.W. et al. Diphtheria overimmunization in children. Can. Med. Assoc. J., 141: 1241 (1989).
- Lau, R.C. Detection of diphtheria toxin antibodies in human sera by ELISA. J. Hyg., 96: 415 (1986).
- Melville-Smith, M. et al. Estimation of diphtheria antitoxin in human sera by ELISA. J. Med. Microbiol., 25: 279 (1988).
- World Health Organization. The immunological basis for immunization series; module 2: Diphtheria – Update 2009. WHO Library Cataloguing (2009).

**DIESE Diagnostica Senese S.p.A.**  
 Via delle Rose 10  
 53035 Monteriggioni (Siena)  
 Italy

	EN Date of manufacture ES Fecha de fabricación IT Data di fabbricazione	FR Date de fabrication GR Ημερομηνία Παραγωγής PT Data de fabrico
	EN Use By ES Fecha de caducidad IT Utilizzare entro	FR Utiliser jusque GR Ημερομηνία λήξης PT Prazo de validade
	EN Do not reuse ES No reutilizar IT Non riutilizzare	FR Ne pas réutiliser GR Μην κάνετε επαναληπτική χρήση PT Não reutilizar
	EN Caution, consult accompanying documents ES Atención, ver instrucciones de uso IT Attenzione, vedere le istruzioni per l'uso	FR Attention voir notice d'instructions GR Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα PT Atenção, consulte a documentação incluída
	EN Manufacturer ES Fabricante IT Fabricante	FR Fabricant GR Κατασκευαστής PT Fabricante
	EN Contains sufficient for <n> tests ES Contenido suficiente para <n> ensayos IT Contenuto sufficiente per "n" saggi	FR Contenu suffisant pour "n" tests GR Περιεχόμενο επαρκές για «n» εξετάσεις PT Conteúdo suficiente para "n" ensaios
	EN Temperature limitation ES Límite de temperatura IT Limiti di temperatura	FR Limites de température GR Περιορισμοί θερμοκρασίας PT Limites de temperatura
	EN Consult Instructions for Use ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso	FR Consulter les instructions d'utilisation GR Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização
	EN Biological risks ES Riesgo biológico IT Rischio biologico	FR Risques biologiques GR Βιολογικοί κίνδυνοι PT Risco biológico
	EN Catalogue number ES Número de catálogo IT Numero di catalogo	FR Référence du catalogue GR Αριθμός καταλόγου PT Referência de catálogo
	EN In Vitro Diagnostic Medical Device ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro	FR Dispositif médical de diagnostic in vitro GR In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro
	EN Batch code ES Código de lote IT Codice del lotto	FR Code du lot GR Αριθμός Παρτίδας PT Código do lote