CHORUS**DIESSE**
DIESSE**Poliovirus IgG**DIESSE Diagnostica Senese
S.p.A.Via delle Rose, 10
53035 Monteriggioni (Siena)
Italy**REF** 81264**REF** 81264/12

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16. BIBLIOGRAFIA

- Cohen JI, Chapter 175: Enteroviruses and Reoviruses, in Kasper DL, Braunwald E, Fauci AS, et al. (eds.) (a cura di), Harrison's Principles of Internal Medicine, 16th, McGraw-Hill Professional, 2004, p. 1144, ISBN 0-07-140235-7.
- Atkinson W, Hamborsky J, McIntyre L, Wolfe S (eds.), Poliomyelitis (PDF), in Epidemiology and Prevention of Vaccine-Preventable Diseases (The Pink Book), 11th, Washington DC, Public Health Foundation, 2009, pp. 231-44.
- Carlsson B, Zaman S, Mellander L, Jalil F, Hanson LA. Secretary and serum immunoglobulin class-specific antibodies to poliovirus after vaccination.
- WHO (1990): Global eradication of poliomyelitis by the year 2000. World Health Assembly Resolution WHA 41.28 Handbook of resolutions and decisions of the World Health Assembly and the Executive Board, vol III, 2nd edit, (1985-1989), Geneva

5. Brown B, Oberste MS, Maher K, Pallansch MA. Complete Genomic Sequencing Shows that Polioviruses and Members of Human Enterovirus Species C Are Closely Related in the Noncapsid Coding Region. Journal of virology, 2003, 8973-8984.

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

INSTRUCTIONS FOR USE

CHORUS Poliovirus IgG

For the qualitative determination of anti-Poliovirus IgG antibodies

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the qualitative determination of IgG class antibodies against Poliovirus in human serum and plasma, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Poliomyelitis is a severe infection of the Central Nervous System, which primarily affects the spinal cord neurons and can lead to paralysis and death.

The infection is caused by three types of Poliovirus (1, 2 and 3), belonging to the Enterovirus genus.

The infection can occur by fecal-oral route, by ingesting contaminated water or food, or saliva and droplets emitted by sick people during coughing and sneezing.

In more than 90% of the infections, the patient does not suffer any symptoms; in the remaining cases, the symptoms are fever, nausea, vomiting, diarrhoea and sore throat; very rarely the paralysis is seen.

There exists no treatment of the disease, only symptomatic therapy, which can partially reduce the effects of the infection, and prevention measures through vaccination.

In many Asiatic countries poliomyelitis is still endemic, but WHO has started ambitious projects to eradicate the disease.

The diagnosis of poliomyelitis can be performed by isolation and detection of the virus or by determination of antibodies in the blood using neutralization tests and, more recently, ELISA technique.

The detection, in the tested sample, of IgG class antibodies against the three types of Poliovirus can be due to a previous infection or to immunization by vaccination.

3. PRINCIPLE OF THE METHOD

The Chorus Poliovirus IgG device is ready to use for the detection of IgG antibodies against Poliovirus, in the Chorus/Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human sample.

After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgG antibodies conjugated to horseradish peroxidase. The unbound conjugate is eliminated, and the peroxidase substrate added. The blue colour which develops is

proportional to the concentration of specific antibodies present in the sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.

The results are expressed in Index (OD sample/OD cut-off).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully

followed and the Instrument Operating Manual must be consulted.

5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 81264).


The kit is sufficient for 12 tests (REF 81264/12).

DD DEVICES

6 packages each containing 6 devices (REF 81264).

2 packages each containing 6 devices (REF 81264/12).

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with Poliovirus antigen (type 1, type 2 and type 3)

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: saline proteic solution with Proclin (0.1%)

Position 2: CONJUGATE

Contents: anti-human IgG antibodies labelled with horseradish peroxidase, in proteic solution containing preservative.

Position 1: EMPTY WELL

In which the operator must place the undiluted sample

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.175 ml

Contents: Diluted human serum containing IgG antibodies anti-Poliovirus and preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.425 ml

Diluted human serum containing IgG antibodies anti-Poliovirus and preservative. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER (REF 83606)
- CLEANING SOLUTION 2000 (REF 83609)
- SANITIZING SOLUTION (REF 83604 - 83608)
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the positive control (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum or plasma collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh sample may be stored for 7 days at 2/8°C, or frozen for longer periods at -20°C.

Freeze-thawing cycles have to be avoided.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 50 µl of undiluted test sample in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the results obtained. It should be used as reported in the Instrument

Operating Manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the positive control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in Index (OD sample/OD cut-off).

The test on the examined sample can be interpreted as follows:

POSITIVE: when the result is > 1.2

NEGATIVE: when the result is < 0.8

DOUBTFUL/EQUIVOCAL: for all values between 0.8 and 1.2

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new sample after 2-4 weeks and test the new sample with the previous one.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. ANALYTICAL SPECIFICITY

3 samples were spiked with the following potentially interfering factors and then tested:

Bilirubin (1.9 mg/dl – 30 mg/dl)
 Triglycerides (31 mg/dl – 500 mg/dl)
 Hemoglobin (0.5 mg/ml – 4 mg/ml)

The presence in the sample of the interfering substances described above does not affect the test result.

13. CROSS-REACTIONS

It is not possible to exclude cross-reactions of antibodies directed against enterovirus different from Poliovirus.

14. METHOD COMPARISON

In an experimentation 81 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	49	0	49
	-	4	28	32
	Total	53	28	81

Percent Positive Agreement (~Diagnostic Sensitivity):

92.5% CI_{95%}: 82.1-97.0

Percent Negative Agreement: (~Diagnostic Specificity):

100% CI_{95%}: 87.9-99.8

The agreement between the two methods is excellent with a Cohen's Kappa of 0.89.

15. PRECISION AND REPEATABILITY

Sample	Within run	
	Mean (Index)	CV%
1	0.2	40.0*
2	0.1	50.0*
3	2.5	9.2
4	2.7	5.2

Sample	Between run	
	Mean (Index)	CV%
1	0.2	20.0*
2	0.1	50.0*
3	2.2	7.3
4	2.4	6.7

Sample	Between lots		Between Instruments	
	Mean (Index)	CV%	Mean (Index)	CV%
1	0.2	20.0*	0.2	25.0*
2	0.1	40.0*	0.2	30.0*
3	2.5	8.4	2.4	12.9
4	2.5	13.2	2.5	12.0

* Artifact caused by the known fault of Variation Coefficient which becomes extremely sensitive to even very small changes in the mean when the mean value is near zero.

16. REFERENCES

- Cohen JI, Chapter 175: Enteroviruses and Reoviruses, in Kasper DL, Braunwald E, Fauci AS, et al. (eds.) (a cura di), Harrison's Principles of Internal Medicine, 16th, McGraw-Hill Professional, 2004, p. 1144, ISBN 0-07-140235-7.
- Atkinson W, Hamborsky J, McIntyre L, Wolfe S (eds.), Poliomyelitis (PDF), in Epidemiology and Prevention of Vaccine-Preventable Diseases (The Pink Book), 11th, Washington DC, Public Health Foundation, 2009, pp. 231-44.
- Carlsson B, Zaman S, Mellander L, Jalil F, Hanson LA. Secretary and serum immunoglobulin class-specific antibodies to poliovirus after vaccination.
- WHO (1990): Global eradication of poliomyelitis by the year 2000. World Health Assembly Resolution WHA 41.28, Handbook of resolutions and decisions of the World Health Assembly and the Executive Board, vol III, 2nd edit, (1985-1989), Geneva
- Brown B, Oberste MS, Maher K, Pallansch MA. Complete Genomic Sequencing Shows that Polioviruses and Members of Human Enterovirus Species C Are Closely Related in the Noncapsid Coding Region. Journal of virology, 2003, 8973-8984.

DIESSE Diagnostica Senese S.p.A.
 Via delle Rose 10
 53035 Monteriggioni (Siena)
 Italy

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4. WHO (1990): Global eradication of poliomyelitis by the year 2000. World Health Assembly Resolution WHA 41.28, Handbook of resolutions and decisions of the World Health Assembly and the Executive Board, vol III, 2nd edit, (1985–1989), Geneva
5. Brown B, Oberste MS, Maher K, Pallansch MA. Complete Genomic Sequencing Shows that Polioviruses and Members of Human Enterovirus Species C Are Closely Related in the Noncapsid Coding Region. Journal of virology, 2003, 8973-8984.

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

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Percent Positive Agreement (~Sensibilidade Diagnóstica):
 92.5% CI_{95%}: 82.1-97.0

Percent Negative Agreement: (~Especificidade Diagnóstica):
 100% CI_{95%}: 87.9-99.8

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16. BIBLIOGRAFIA

1. Cohen JI, Chapter 175: Enteroviruses and Reoviruses, in Kasper DL, Braunwald E, Fauci AS, et al. (eds.) (a cura di), Harrison's Principles of Internal Medicine, 16th, McGraw-Hill Professional, 2004, p. 1144, ISBN 0-07-140235-7.
2. Atkinson W, Hamborsky J, McIntyre L, Wolfe S (eds.), Poliomyelitis (PDF), in Epidemiology and Prevention of Vaccine-Preventable Diseases (The Pink Book), 11th, Washington DC, Public Health Foundation, 2009, pp. 231-44.
3. Carlsson B, Zaman S, Mellander L, Jalil F, Hanson LA.



Secretary and serum immunoglobulin class-specific antibodies to poliovirus after vaccination.

4. WHO (1990): Global eradication of poliomyelitis by the year 2000. World Health Assembly Resolution WHA 41.28 Handbook of resolutions and decisions of the World Health Assembly and the Executive Board, vol III, 2nd edit, (1985–1989), Geneva
5. Brown B, Oberste MS, Maher K, Pallansch MA. Complete Genomic Sequencing Shows that Polioviruses and Members of Human Enterovirus Species C Are Closely Related in the Noncapsid Coding Region. Journal of virology, 2003, 8973-8984.

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

	EN Date of manufacture ES Fecha de fabricación IT Data di fabbricazione	FR Date de fabrication GR Ημερομηνία Παραγωγής PT Data de fabrico
	EN Use By ES Fecha de caducidad IT Utilizzare entro	FR Utiliser jusque GR Ημερομηνία λήξης PT Prazo de validade
	EN Do not reuse ES No reutilizar IT Non riutilizzare	FR Ne pas réutiliser GR Μην κάνετε επαναληπτική χρήση PT Não reutilizar
	EN Caution, consult accompanying documents ES Atención, ver instrucciones de uso IT Attenzione, vedere le istruzioni per l'uso	FR Attention voir notice d'instructions GR Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα PT Atenção, consulte a documentação incluída
	EN Manufacturer ES Fabricante IT Fabricante	FR Fabricant GR Κατασκευαστής PT Fabricante
	EN Contains sufficient for <n> tests ES Contenido suficiente para <n> ensayos IT Contenuto sufficiente per "n" saggi	FR Contenu suffisant pour "n" tests GR Περιεχόμενο επαρκές για «n» εξετάσεις PT Conteúdo suficiente para "n" ensaios
	EN Temperature limitation ES Límite de temperatura IT Limiti di temperatura	FR Limites de température GR Περιορισμοί θερμοκρασίας PT Limites de temperatura
	EN Consult Instructions for Use ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso	FR Consulter les instructions d'utilisation GR Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização
	EN Biological risks ES Riesgo biológico IT Rischio biologico	FR Risques biologiques GR Βιολογικοί κίνδυνοι PT Risco biológico
	EN Catalogue number ES Número de catálogo IT Numero di catalogo	FR Référence du catalogue GR Αριθμός καταλόγου PT Referência de catálogo
	EN In Vitro Diagnostic Medical Device ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro	FR Dispositif médical de diagnostic in vitro GR In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro
	EN Batch code ES Código de lote IT Codice del lotto	FR Code du lot GR Αριθμός Παρτίδας PT Código do lote