CHORUS

Tick-Borne Encephalitis Virus IgG

REF 81276

**DIESSE**

DIESSE Diagnostica Senese
S.p.A.
Via delle Rose, 10
53035 Monteriggioni (Siena)
Italy

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15. BIBLIOGRAFIA

- K. L. Mansfield, N. Johnson, L. P. Phipps, J. R. Stephenson, A. R. Fooks and T. Solomon. Tick-borne encephalitis virus – a review of an emerging zoonosis Journal of General Virology (2009), 90, 1781–1794
- Anu Ja'a'skela'inen, Xiuqi Han, Matthias Niedrig, Antti Vaheri, and Olli Vapalahti. Diagnosis of Tick-Borne Encephalitis by a -Capture Immunoglobulin M-Enzyme Immunoassay Based on Secreted Recombinant Antigen Produced in Insect Cells. JOURNAL OF CLINICAL MICROBIOLOGY, Sept. 2003, p. 4336–4342
- Mansfield K.L., Horton D.L., Johnson N., Li Li, Barrett A.D.T., Smith D.J., Galbraith S.E., Solomon T and Fooks A.R. Flavivirus-induced antibody cross-reactivity. Journal of General Virology (2011), 92, 2821-2829.

DIESSE Diagnostica Senese S.p.A.
 Via delle Rose 10
 53035 Monteriggioni (Siena)
 Italy

INSTRUCTIONS FOR USE

CHORUS Tick-Borne Encephalitis Virus IgG

For the semiquantitative determination of anti-Tick-Borne Encephalitis Virus IgG antibodies

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the semiquantitative determination of IgG class antibodies against Tick-Borne Encephalitis Virus in human serum and citrate plasma, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Tick-borne encephalitis (TBE) is a viral infectious disease, which affects Central Nervous System, caused by an arbovirus from the family *Flaviviridae*. Reservoir hosts of the virus are forest animals and the transmission, to humans and animals, occurs through tick bites or, exceptionally, by consumption of infected heat untreated milk. Up to 70% of infections are clinically inapparent or with mild symptoms; in the remaining 30% of cases, after an incubation period of 3-14 days, non-specific flu-like symptoms occur. After some days of remission, the second phase of the disease develops. This phase is characterized by aggravated headaches, ophtalmoplegia, vomiting, malaise, meningeal symptoms, cranial nerve paralysis and paresis of the limbs. An acute stage of tick-borne encephalitis lasts 1 to 3 weeks. More severe course of the disease, often with lasting ill effects, can be observed in senior patients.

The diagnosis is based on anamnesis, clinical picture and results of laboratory tests such as the detection of specific IgM and IgG antibodies. The IgG antibodies protect the organism against a new infection and they can be detected over a long period (several years) after the past infection or vaccination.

3. PRINCIPLE OF THE METHOD

The Chorus Tick-Borne Encephalitis Virus IgG device is ready to use for the detection of IgG antibodies against TBE-Virus, in the Chorus/Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human sample.

After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgG antibodies conjugated to horseradish peroxidase. The unbound conjugate is eliminated, and the peroxidase substrate added. The blue colour which develops is proportional to the concentration of specific antibodies present in the sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments. The results are expressed in Arbitrary Units (AU/ml).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.

The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website (<http://www.diesse.it/en/Support/Download/strumento:39/>)


5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices.

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with Tick-borne encephalitis virus purified and inactivated native antigen

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic buffer solution.

Position 2: CONJUGATE

Contents: anti-human IgG antibodies labelled with horseradish peroxidase.

Position 1: EMPTY WELL

In which the operator must place the undiluted sample.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.450 ml

Contents: anti-TBE Virus IgG antibodies and preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.900 ml

Contents: anti-TBE Virus IgG antibodies and preservative. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the positive control (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum or citrate plasma collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh sample may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense the following volumes in well no. 1 of each device:

SAMPLE	50 µl/device
CALIBRATOR	130 µl/device
POSITIVE CONTROL	130 µl/device

4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the positive control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in Arbitrary Units (AU/ml) calculated on the basis of a lot-dependent graph stored in the instrument.

The test on the examined sample can be interpreted as follows:

POSITIVE: when the result is > 22

NEGATIVE: when the result is < 18

DOUBTFUL/EQUIVOCAL: for all values between 18 and 22

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new sample.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. CROSS-REACTIONS

It is known the possible cross-reaction with anti-Flavivirus antibodies.

13. METHOD COMPARISON

In an experimentation 94 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	34	2	36
	-	2	56	58
	Total	36	58	94

Percent Positive Agreement (~Diagnostic Sensitivity):

94.4% CI_{95%}: 81.8 - 98.4

Percent Negative Agreement: (~Diagnostic Specificity):

96.6% CI_{95%}: 88.2 - 99.0

The agreement between the two methods is excellent with a Cohen's Kappa of 0.92.

14. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision	
	Mean (AU/ml)	CV%	Mean (AU/ml)	CV%
1	43.5	4.9	44.7	3.2

2	22.6	2.5	22.6	2.5
3	13.3	5.0	14.9	5.2

Sample	Precision between batches		Precision between instruments	
	Mean (AU/ml)	CV%	Mean (AU/ml)	CV%
1	46.1	4.6	46.3	4.3
2	25.6	8.8	22.2	5.1
3	15.1	4.8	15.6	6.2

15. REFERENCES

1. K. L. Mansfield, N. Johnson, L. P. Phipps, J. R. Stephenson, A. R. Fooks¹ and T. Solomon. Tick-borne encephalitis virus – a review of an emerging zoonosis *Journal of General Virology* (2009), 90, 1781–1794
2. Anu Ja"askela"inen, Xiuqi Han, Matthias Niedrig, Antti Vaheri, and Olli Vapalahti. Diagnosis of Tick-Borne Encephalitis by a -Capture Immunoglobulin M-Enzyme Immunoassay Based on Secreted Recombinant Antigen Produced in Insect Cells. *JOURNAL OF CLINICAL MICROBIOLOGY*, Sept. 2003, p. 4336–4342
3. Mansfield K.L., Horton D.L., Johnson N.,¹ Li Li, Barrett A.D.T., Smith D.J., Galbraith S.E., Solomon T and Fooks A.R. Flavivirus-induced antibody cross-reactivity. *Journal of General Virology* (2011), 92, 2821-2829.

DIESE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

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CHORUS Tick-Borne Encephalitis Virus IgG

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23. LITERATUR

1. K. L. Mansfield, N. Johnson, L. P. Phipps, J. R. Stephenson, A. R. Fooks¹ and T. Solomon. Tick-borne encephalitis virus – a review of an emerging zoonosis Journal of General Virology (2009), 90, 1781–1794.
2. Anu Ja"askela"inen, Xiuqi Han, Matthias Niedrig, Antti Vaheri, and Olli Vapalahti. Diagnosis of Tick-Borne Encephalitis by a -Capture Immunoglobulin M-Enzyme Immunoassay Based on Secreted Recombinant Antigen Produced in Insect Cells. JOURNAL OF CLINICAL MICROBIOLOGY, Sept. 2003, p. 4336–4342.

3. Mansfield K.L., Horton D.L., Johnson N.,¹ Li Li, Barrett A.D.T., Smith D.J., Galbraith S.E., Solomon T and Fooks A.R. Flavivirus-induced antibody cross-reactivity. Journal of General Virology (2011), 92, 2821-2829.

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italien

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1. K. L. Mansfield, N. Johnson, L. P. Phipps, J. R. Stephenson, A. R. Fooks¹ and T. Solomon. Tick-borne encephalitis virus – a review of an emerging zoonosis *Journal of General Virology* (2009), 90, 1781–1794.
2. Anu Jaakkola¹, Xiuqi Han, Matthias Niedrig, Antti Vaheri, and Olli Vapalahti. Diagnosis of Tick-Borne Encephalitis by a -Capture Immunoglobulin M-Enzyme Immunoassay Based on Secreted Recombinant Antigen Produced in Insect Cells. *JOURNAL OF CLINICAL MICROBIOLOGY*, Sept. 2003, p. 4336–4342.
3. Mansfield K.L., Horton D.L., Johnson N., Li Li, Barrett A.D.T., Smith D.J., Galbraith S.E., Solomon T and Fooks A.R. Flavivirus-induced antibody cross-reactivity. *Journal of General Virology* (2011), 92, 2821-2829.

DIESSE Diagnostica Senese S.p.A.
 Via delle Rose 10
 53035 Monteriggioni (Siena)
 Italy

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15. BIBLIOGRAFÍA

1. K. L. Mansfield, N. Johnson, L. P. Phipps, J. R. Stephenson, A. R. Fooks¹ and T. Solomon. Tick-borne encephalitis virus – a review of an emerging zoonosis Journal of General Virology (2009), 90, 1781–1794
2. Anu Ja"ä"skela"inen, Xiuqi Han, Matthias Niedrig, Antti Vaheri, and Olli Vapalahti. Diagnosis of Tick-Borne Encephalitis by a -Capture Immunoglobulin M-Enzyme Immunoassay Based on Secreted Recombinant Antigen Produced in Insect Cells. JOURNAL OF CLINICAL MICROBIOLOGY, Sept. 2003, p. 4336–4342
3. Mansfield K.L., Horton D.L., Johnson N., Li Li, Barrett A.D.T., Smith D.J., Galbraith S.E., Solomon T and Fooks A.R. Flavivirus-induced antibody cross-reactivity. Journal of General Virology (2011), 92, 2821-2829.

DIESSE Diagnostica Senese S.p.A.
 Via delle Rose 10
 53035 Monteriggioni (Siena)
 Italy

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15. BIBLIOGRAPHIE

- K. L. Mansfield, N. Johnson, L. P. Phipps, J. R. Stephenson, A. R. Fooks¹ and T. Solomon. Tick-borne encephalitis virus – a review of an emerging zoonosis *Journal of General Virology* (2009), 90, 1781–1794.
- Anu Ja"ä"skela"inen, Xiuqi Han, Matthias Niedrig, Antti Vaheri, and Olli Vapalahti. Diagnosis of Tick-Borne Encephalitis by a -Capture Immunoglobulin M-Enzyme Immunoassay Based on Secreted Recombinant Antigen Produced in Insect Cells. *JOURNAL OF CLINICAL MICROBIOLOGY*, Sept. 2003, p. 4336–4342.
- Mansfield K.L., Horton D.L., Johnson N., Li Li, Barrett A.D.T., Smith D.J., Galbraith S.E., Solomon T and Fooks A.R. Flavivirus-induced antibody cross-reactivity. *Journal of General Virology* (2011), 92, 2821–2829.

DIESSE Diagnostica Senese S.p.A.
 Via delle Rose 10
 53035 Monteriggioni (Siena)
 Italie

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15. BIBLIOGRAFIA

1. K. L. Mansfield, N. Johnson, L. P. Phipps, J. R. Stephenson, A. R. Fooks¹ and T. Solomon. Tick-borne encephalitis virus – a review of an emerging zoonosis Journal of General Virology (2009), 90, 1781–1794
2. Anu Ja"ä"skela"inen, Xiuqi Han, Matthias Niedrig, Antti Vaheri, and Olli Vapalahti. Diagnosis of Tick-Borne Encephalitis by a -Capture Immunoglobulin M-Enzyme Immunoassay Based on Secreted Recombinant Antigen Produced in Insect Cells. JOURNAL OF CLINICAL MICROBIOLOGY, Sept. 2003, p. 4336–4342
3. Mansfield K.L., Horton D.L., Johnson N., Li Li, Barrett A.D.T., Smith D.J., Galbraith S.E., Solomon T and Fooks A.R. Flavivirus-induced antibody cross-reactivity. Journal of General Virology (2011), 92, 2821–2829.

DIESSE Diagnostica Senese S.p.A.
 Via delle Rose 10
 53035 Monteriggioni (Siena)
 Italy

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CHORUS**Tick-Borne Encephalitis Virus IgG**

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15. BIBLIOGRAFIE

- K. L. Mansfield, N. Johnson, L. P. Phipps, J. R. Stephenson, A. R. Fooks¹ and T. Solomon. Tick-borne encephalitis virus – a review of an emerging zoonosis Journal of General Virology (2009), 90, 1781–1794.
- Anu Ja"ä"skela"inen, Xiuqi Han, Matthias Niedrig, Antti Vaheri, and Olli Vapalahti. Diagnosis of Tick-Borne Encephalitis by a -Capture Immunoglobulin M-Enzyme Immunoassay Based on Secreted Recombinant Antigen Produced in Insect Cells. JOURNAL OF CLINICAL MICROBIOLOGY, Sept. 2003, p. 4336–4342.
- Mansfield K.L., Horton D.L., Johnson N.,¹ Li Li, Barrett A.D.T., Smith D.J., Galbraith S.E., Solomon T and Fooks A.R. Flavivirus-induced antibody cross-reactivity. Journal of General Virology (2011), 92, 2821–2829.

DIESSE Diagnostica Senese S.p.A.
 Via delle Rose 10
 53035 Monteriggioni (Siena)
 Italy

	EN Use By ES Fecha de caducidad IT Utilizzare entro	FR Utiliser jusque GR Ημερομηνία λήξης PT Prazo de validade
	EN Do not reuse ES No reutilizar IT Non riutilizzare	FR Ne pas réutiliser GR Μην κάνετε επαναληπτική χρήση PT Não reutilizar
	EN Caution, consult accompanying documents ES Atención, ver instrucciones de uso IT Attenzione, vedere le istruzioni per l'uso	FR Attention voir notice d'instructions GR Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα PT Atenção, consulte a documentação incluída
	EN Manufacturer ES Fabricante IT Fabbicante	FR Fabricant GR Κατασκευαστής PT Fabricante
	EN Contains sufficient for <n> tests ES Contenido suficiente para <n> ensayos IT Contenuto sufficiente per "n" saggi	FR Contenu suffisant pour "n" tests GR Περιεχόμενο επαρκές για «n» εξετάσεις PT Conteúdo suficiente para "n" ensaios
	EN Temperature limitation ES Límite de temperatura IT Limiti di temperatura	FR Limites de température GR Περιορισμοί θερμοκρασίας PT Limites de temperatura
	EN Consult Instructions for Use ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso	FR Consulter les instructions d'utilisation GR Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização
	EN Biological risks ES Riesgo biológico IT Rischio biologico	FR Risques biologiques GR Βιολογικοί κίνδυνοι PT Risco biológico
	EN Catalogue number ES Número de catálogo IT Numero di catalogo	FR Référence du catalogue GR Αριθμός καταλόγου PT Referência de catálogo
	EN In Vitro Diagnostic Medical Device ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro	FR Dispositif médical de diagnostic in vitro GR In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro
	EN Batch code ES Código de lote IT Codice del lotto	FR Code du lot GR Αριθμός Παρτίδας PT Código do lote