CHORUS

Clostridium difficile GDH

DIESSE

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5. Carman RJ, Wickham KN, Chen L, Lawrence AM, Boone JH, Wilkins TD, Kerkering TM, Lyerly DM. Glutamate dehydrogenase is highly conserved among *Clostridium difficile* ribotypes. J Clin Microbiol. 2012; 50(4): 1425-6.
6. Lyerly D.M., Ball D.W., Toth J and Wilkins T.D. Characterization of Cross-Reactive Proteins Detected by Culturette Brand Rapid Latex Test for *Clostridium difficile*. Journal of Clinical Microbiology. (Mar. 1988) p. 397-400.

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INSTRUCTIONS FOR USE

CHORUS Clostridium difficile GDH

For the qualitative determination of Clostridium difficile GDH in stools.

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the qualitative determination of Clostridium difficile GDH in stools, using a disposable device applied on the Chorus TRIO instruments.

2. INTRODUCTION

The gram-positive anaerobic Clostridium difficile is the leading causative agent of antibiotic-associated diarrhea and pseudomembranous colitis. This pathogen is capable of causing diseases that could be severe or fatal if not diagnosed and treated on time. The exposure to antibiotics is the major risk factor for C. difficile infection. The infection can develop if the normal gastrointestinal flora is disrupted by an antibiotic therapy and a person acquires spores that produce C. difficile, typically via fecal-oral route. C. difficile's key virulence factors are toxin A and toxin B. These toxins show high sequence and functional homology.

All strains of C. difficile produce high levels of GDH. Therefore, C. difficile's GDH enzyme is considered a very good antigen marker for the detection of this organism.

The Chorus Clostridium difficile GDH test is a highly specific and sensitive diagnostic kit for the determination of GDH in stool specimens. The presence of toxins in the GDH positive samples is confirmed by using the Chorus Clostridium difficile A/B Toxins (REF 81170).

3. PRINCIPLE OF THE METHOD

The Chorus Clostridium difficile GDH device is ready to use for the detection of C. difficile GDH present in human feces, in the Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). The specific monoclonal anti-C. difficile GDH antibodies are bound to the solid phase. The specific antigen is bound to the monoclonal antibody and to the anti- C. difficile GDH monoclonal antibody conjugated to horseradish peroxidase through incubation with the samples and conjugate.

After washing to eliminate the proteins which have not reacted, the peroxidase substrate is added.

The enzymatic reaction is stopped by the adding of the Stop Solution, which makes the solution turn to yellow.

The colour which develops is proportional to the concentration of the C. difficile GDH present in the sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus TRIO instruments.

The results are expressed in Index (OD sample/OD cut-off).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.

Waste disposal: samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus TRIO instrument, the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.
The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website (<http://www.diesse.it/en/Support/Download/strumento:39/>)
5. Check that the Chorus TRIO instrument is set up correctly (see Operating Manual).


6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use the device after the expiry date.
11. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 12 tests.

DD DEVICES 2 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with monoclonal antibodies specific for *C. difficile*

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: STOP SOLUTION

Contents: Sulfuric acid solution 0.3 M

Position 2: CONJUGATE

Contents: ready-to-use solution containing anti-*C. difficile* GDH monoclonal antibodies labelled with horseradish peroxidase and preservative.

Position 1: EMPTY WELL

In which the operator must place the sample.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.300 ml

Contents: Solution containing GDH recombinant enzyme and preservative.

Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.450 ml

Contents: Solution containing GDH recombinant enzyme and preservative.

Liquid, ready for use.

STOOL DILUENT/NEGATIVE CONTROL

C. DIFFICILE STOOL DILUENT – NEGATIVE CONTROL

1 x 6 ml

Contents: ready-to-use buffer solution and preservative.

Use: Stool Diluent (see chapter 8 Assay Procedure).

The Diluent can also be used as negative control solution.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-1000 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the positive control (see chapter 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The type of sample is composed of feces handled with all precautions dictated by good laboratory practice.

Keep the samples at 2/8°C and test within 48 hours after collection. If testing cannot be performed within 48 hours, store the samples at -20°C or lower.

Minimize specimen freezing thawing which may result in loss of activity.

8. ASSAY PROCEDURE

1. Set up one dilution tube for each specimen to be tested. 1.5 ml tubes are recommended for this purpose. Add 400 µl of *C. difficile* Stool Diluent/Negative Control to each tube.
2. **Formed samples:** use a wooden applicator stick or a disposable teaspoon to transfer the fecal specimen to the tube. Transfer approximately 0.1 to 0.15 g of specimen (about the size of a small pea) into 400 µl of Stool Diluent/Negative Control, removing all the material from the collector.
Liquid samples: transfer 150 µl of specimen to the tube. Make sure the liquid specimens are evenly suspended.
3. **Thoroughly mix (vortex) the fecal specimen for 15 seconds.**
4. Let the tubes stand in vertical position for 10 minutes.
5. Centrifuge the tubes at 1000 g for 30 seconds.
6. Store the diluted samples at 2/8°C until the test is performed.
7. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
8. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
9. Dispense 130 µl of sample in well no. 1 of each device; at each change of batch, use a device for the calibrator.

10. Place the devices in the Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the positive control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus TRIO instrument expresses the result in Index (OD sample/OD cut-off).

The test on the examined sample can be interpreted as follows:

POSITIVE: when the result is > 0.5
NEGATIVE: when the result is < 0.5
DOUBTFUL/EQUIVOCAL: for value 0.5

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new sample.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

ATTENTION: Chorus *Clostridium difficile* GDH test (REF 81168) can only be performed in a dedicated run or together with Chorus *Clostridium difficile* A/B Toxins (REF 81170) and Chorus *Helicobacter pylori* Ag (REF 81063). Otherwise none of the programmed tests will be performed.

12. CROSS-REACTIONS

In the determination of GDH, it is known the enzyme interference between toxigenic and non-toxigenic organisms as well as the anti-GDH antibodies interference of other *Clostridium* species.

13. METHOD COMPARISON

In an experimentation 48 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	23	0	23
	-	0	25	25
	Total	23	25	48

Percent Positive Agreement (~Diagnostic Sensitivity):
100% CI_{95%}: 85.7-99.8

Percent Negative Agreement (~Diagnostic Specificity):
100% CI_{95%}: 86.7-99.8

The agreement between the two methods is excellent with a Cohen's Kappa of 1.0.

14. PRECISION AND REPEATABILITY

Sample	Within-run Precision	
	Mean (Index)	CV%
1	3.1	3.2
2	2.3	3.5
3	0.3	13.3

Sample	Between-run precision	
	Mean (Index)	CV%
1	3.1	2.6
2	2.0	3.0
3	0.3	-

Sample	Precision between batches		Precision between instruments	
	Mean (Index)	CV%	Mean (Index)	CV%
1	3.0	3.0	3.4	4.1
2	2.0	4.0	2.5	7.6
3	0.3	-	0.3	10.0

15. REFERENCES

- Owens RC Jr, Donskey CJ, Gaynes RP, Loo VG, Muto CA. Antimicrobial-associated risk factors for *Clostridium difficile* infection. Clin. Infect Dis. 2008; 46 Suppl. 1:S19-31
- Savidge TC, Pan WH, Newman P, O'Brien M, Anton PM, Pothoulakis C. *Clostridium difficile* toxin B is an inflammatory enterotoxin in human intestine. Gastroenterology 2003; 125(2): 413-20
- Pituch H, van den Braak N, van Leeuwen W, van Belkum A, Martirosian G, Obuch-Woszczatyński P, Łuczak M, Meisel-Mikolajczyk F. Clonal dissemination of a toxin-A-negative/toxin-B-positive *Clostridium difficile* strain from patients with antibiotic-associated diarrhea in Poland. Clin Microbiol Infect. 2001; 7(8): 442-6
- Lyerly DM, Barroso LA, Wilkins TD. Identification of the latex test-reactive protein of *Clostridium difficile* as glutamate dehydrogenase. J Clin Microbiol. 1991; 29(11): 2639-42.
- Carman RJ, Wickham KN, Chen L, Lawrence AM, Boone JH, Wilkins TD, Kerkerling TM, Lyerly DM. Glutamate dehydrogenase is highly conserved among *Clostridium difficile* ribotypes. J Clin Microbiol. 2012; 50(4): 1425-6.
- Lyerly D.M., Ball D.W., Toth J and Wilkins T.D. Characterization of Cross-Reactive Proteins Detected by Culturette Brand Rapid Latex Test for *Clostridium difficile*. Journal of Clinical Microbiology. (Mar. 1988) p. 397-400.

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1. Owens RC Jr, Donskey CJ, Gaynes RP, Loo VG, Muto CA. Antimicrobial-associated risk factors for *Clostridium difficile* infection. Clin. Infect Dis. 2008; 46 Suppl. 1:S19-31
2. Savidge TC, Pan WH, Newman P, O'Brien M, Anton PM, Pothoulakis C. *Clostridium difficile* toxin B is an inflammatory enterotoxin in human intestine. Gastroenterology 2003; 125(2): 413-20
3. Pituch H, van den Braak N, van Leeuwen W, van Belkum A, Martirosian G, Obuch-Woszczatyński P, Łuczak M, Meisel-Mikolajczyk F. Clonal dissemination of a toxin-A-negative/toxin-B-positive *Clostridium difficile* strain from patients with antibiotic-associated diarrhea in Poland. Clin Microbiol Infect. 2001; 7(8): 442-6
4. Lysterly DM, Barroso LA, Wilkins TD. Identification of the latex test-reactive protein of *Clostridium difficile* as glutamate dehydrogenase. J Clin Microbiol. 1991; 29(11): 2639-42.
5. Carman RJ, Wickham KN, Chen L, Lawrence AM, Boone JH, Wilkins TD, Kerkerling TM, Lysterly DM. Glutamate dehydrogenase is highly conserved among *Clostridium difficile* ribotypes. J Clin Microbiol. 2012; 50(4): 1425-6.
6. Lysterly D.M., Ball D.W., Toth J and Wilkins T.D. Characterization of Cross-Reactive Proteins Detected by Culturette Brand Rapid Latex Test for *Clostridium difficile*. Journal of Clinical Microbiology. (Mar. 1988) p. 397-400.

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30. LITERATUR

1. Owens RC Jr, Donskey CJ, Gaynes RP, Loo VG, Muto CA. Antimicrobial-associated risk factors for *Clostridium difficile* infection. Clin. Infect Dis. 2008; 46 Suppl. 1:S19-31
2. Savidge TC, Pan WH, Newman P, O'brien M, Anton PM, Pothoulakis C. *Clostridium difficile* toxin B is an inflammatory enterotoxin in human intestine. Gastroenterology 2003; 125(2): 413-20
3. Pituch H, van den Braak N, van Leeuwen W, van Belkum A, Martirosian G, Obuch-Woszczatyński P, Łuczak M, Meisel-Mikołajczyk F. Clonal dissemination of a toxin-A-negative/toxin-B-positive *Clostridium difficile* strain from patients with antibiotic-associated diarrhea in Poland. Clin Microbiol Infect. 2001; 7(8): 442-6
4. Lyerly DM, Barroso LA, Wilkins TD. Identification of the latex test-reactive protein of *Clostridium difficile* as glutamate dehydrogenase. J Clin Microbiol. 1991; 29(11): 2639-42.
5. Carman RJ, Wickham KN, Chen L, Lawrence AM, Boone JH, Wilkins TD, Kerkering TM, Lyerly DM. Glutamate dehydrogenase is highly conserved among *Clostridium difficile* ribotypes. J Clin Microbiol. 2012; 50(4): 1425-6.
6. Lyerly D.M., Ball D.W., Toth J and Wilkins T.D. Characterization of Cross-Reactive Proteins Detected by Culturette Brand Rapid Latex Test for *Clostridium difficile*. Journal of Clinical Microbiology. (Mar. 1988) p. 397-400.

DIESSE Diagnostica Senese S.p.A.
Strada dei Laghi 39
53035 Monteriggioni (SIENA)
Italien

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- Owens RC Jr, Donskey CJ, Gaynes RP, Loo VG, Muto CA. Antimicrobial-associated risk factors for *Clostridium difficile* infection. Clin. Infect Dis. 2008; 46 Suppl. 1:S19-31
- Savidge TC, Pan WH, Newman P, O'brien M, Anton PM, Pothoulakis C. *Clostridium difficile* toxin B is an inflammatory enterotoxin in human intestine. Gastroenterology 2003; 125(2): 413-20
- Pituch H, van den Braak N, van Leeuwen W, van Belkum A, Martirosian G, Obuch-Woszczatyński P, Łuczak M, Meisel-Mikołajczyk F. Clonal dissemination of a toxin-A-negative/toxin-B-positive *Clostridium difficile* strain from patients with antibiotic-associated diarrhea in Poland. Clin Microbiol Infect. 2001; 7(8): 442-6
- Lyerly DM, Barroso LA, Wilkins TD. Identification of the latex test-reactive protein of *Clostridium difficile* as glutamate dehydrogenase. J Clin Microbiol. 1991; 29(11): 2639-42.
- Carman RJ, Wickham KN, Chen L, Lawrence AM, Boone JH, Wilkins TD, Kerkering TM, Lyerly DM. Glutamate dehydrogenase is highly conserved among *Clostridium difficile* ribotypes. J Clin Microbiol. 2012; 50(4): 1425-6.
- Lyerly D.M., Ball D.W., Toth J and Wilkins T.D. Characterization of Cross-Reactive Proteins Detected by Culturette Brand Rapid Latex Test for *Clostridium difficile*. Journal of Clinical Microbiology. (Mar. 1988) p. 397-400.

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15. BIBLIOGRAFÍA

1. Owens RC Jr, Donskey CJ, Gaynes RP, Loo VG, Muto CA. Antimicrobial-associated risk factors for *Clostridium difficile* infection. Clin. Infect Dis. 2008; 46 Suppl. 1:S19-31
2. Savidge TC, Pan WH, Newman P, O'brien M, Anton PM, Pothoulakis C. *Clostridium difficile* toxin B is an inflammatory enterotoxin in human intestine. Gastroenterology 2003; 125(2): 413-20
3. Pituch H, van den Braak N, van Leeuwen W, van Belkum A, Martirosian G, Obuch-Woszczatyński P, Łuczak M, Meisel-Mikołajczyk F. Clonal dissemination of a toxin-A-negative/toxin-B-positive *Clostridium difficile* strain from patients with antibiotic-associated diarrhea in Poland. Clin Microbiol Infect. 2001; 7(8): 442-6
4. Lysterly DM, Barroso LA, Wilkins TD. Identification of the latex test-reactive protein of *Clostridium difficile* as glutamate dehydrogenase. J Clin Microbiol. 1991; 29(11): 2639-42.
5. Carman RJ, Wickham KN, Chen L, Lawrence AM, Boone JH, Wilkins TD, Kerkering TM, Lysterly DM. Glutamate dehydrogenase is highly conserved among *Clostridium difficile* ribotypes. J Clin Microbiol. 2012; 50(4): 1425-6.
6. Lysterly D.M., Ball D.W., Toth J and Wilkins T.D. Characterization of Cross-Reactive Proteins Detected by Culturette Brand Rapid Latex Test for *Clostridium difficile*. Journal of Clinical Microbiology. (Mar. 1988) p. 397-400.

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15. BIBLIOGRAPHIE

1. Owens RC Jr, Donskey CJ, Gaynes RP, Loo VG, Muto CA. Antimicrobial-associated risk factors for *Clostridium*

- difficile* infection. Clin. Infect Dis. 2008; 46 Suppl. 1:S19-31
2. Savidge TC, Pan WH, Newman P, O'brien M, Anton PM, Pothoulakis C. *Clostridium difficile* toxin B is an inflammatory enterotoxin in human intestine. Gastroenterology 2003; 125(2): 413-20
 3. Pituch H, van den Braak N, van Leeuwen W, van Belkum A, Martirosian G, Obuch-Woszczatyński P, Łuczak M, Meisel-Mikołajczyk F. Clonal dissemination of a toxin-A-negative/toxin-B-positive *Clostridium difficile* strain from patients with antibiotic-associated diarrhea in Poland. Clin Microbiol Infect. 2001; 7(8): 442-6
 4. Lysterly DM, Barroso LA, Wilkins TD. Identification of the latex test-reactive protein of *Clostridium difficile* as glutamate dehydrogenase. J Clin Microbiol. 1991; 29(11): 2639-42.
 5. Carman RJ, Wickham KN, Chen L, Lawrence AM, Boone JH, Wilkins TD, Kerkerling TM, Lysterly DM. Glutamate dehydrogenase is highly conserved among *Clostridium difficile* ribotypes. J Clin Microbiol. 2012; 50(4): 1425-6.
 6. Lysterly D.M., Ball D.W., Toth J and Wilkins T.D. Characterization of Cross-Reactive Proteins Detected by Culturette Brand Rapid Latex Test for *Clostridium difficile*. Journal of Clinical Microbiology. (Mar. 1988) p. 397-400.

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- glutamate dehydrogenase. J Clin Microbiol. 1991; 29(11): 2639-42.
5. Carman RJ, Wickham KN, Chen L, Lawrence AM, Boone JH, Wilkins TD, Kerkering TM, Lysterly DM. Glutamate dehydrogenase is highly conserved among *Clostridium difficile* ribotypes. J Clin Microbiol. 2012; 50(4): 1425-6.
 6. Lysterly D.M., Ball D.W., Toth J and Wilkins T.D. Characterization of Cross-Reactive Proteins Detected by Culturette Brand Rapid Latex Test for *Clostridium difficile*. Journal of Clinical Microbiology. (Mar. 1988) p. 397-400.

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