ENZY-WELL

CHIKUNGUNYA IgM

REF 91134 (96 tests)

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16. **BIOBLOGRAFIA**

1. Abrahm & Shridharan "Chikungunya Virus Infection a Resurgent Scourge" Indian J. med. Res. 126, December 2007, pp 502-504
2. M. Chhabra, V. Mittal, D. Bhattacharya, UVS Rana, S Lal "Chikungunya Fever: a re-emerging viral infection" Indian J. of Medical Microbiology, 2008- 26(1): pp 5-12

INSTRUCTIONS FOR USE

ENZY-WELL CHIKUNGUNYA IgM

REF 91134

(English)

1. INTENDED USE

IMMUNOENZYMATIC METHOD FOR THE QUALITATIVE DETERMINATION OF IgM-CLASS ANTIBODIES TO CHIKUNGUNYA VIRUS IN HUMAN SERUM, AS AN AID IN THE DIAGNOSIS OF CHIKUNGUNYA VIRUS INFECTION.

2. SUMMARY AND EXPLANATION OF THE TEST

Chikungunya Virus, belonging to Togaviridae family, genus Alphavirus, causes a disease called “Chikungunya Fever” spread by bite of *Aedes albopictus* (tiger mosquito) e *Aedes aegypti* mosquito.

Chikungunya Virus often produces a mild illness that may be confused with Dengue Virus, indeed symptoms include fever, headache, arthralgia, myalgia and rash.

Although serious complications are uncommon arthralgia can be debilitating and may persist for months, or even years, after infection.

The disease is endemic to Africa, India and South-East Asia, while its detection outside of these countries is normally related to trips to those areas.

The potential risks of transmission through donations of blood or organs should also be considered. For this reason it is of great importance to conduct the serological tests for the diagnosis and for blood and organs donations.

3. PRINCIPLE OF THE TEST

The test is based on the ELISA technique (Enzyme linked Immunosorbent Assay).

The antigen, composed of inactivated chikungunya virus, is bound to the solid phase (8-well strips). The specific immunoglobulins are bound to the antigen through incubation with dilute human serum.

After washings to eliminate the immunoglobulins which have not reacted, incubation is performed with the conjugate, composed of human IgM monoclonal antibodies conjugated to horse radish peroxidase (HRP).

The unbound conjugate is eliminated, and the peroxidase substrate added.

The blue colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

When the enzymatic reaction is interrupted by the addition of a sulphuric acid solution, the yellow colour which develops can be easily read using a microplate reader.

4. KIT CONTENTS AND REAGENT PREPARATION

Bring to room temperature before use.

MT PLATE MICROPLATE 12x8 wells coated with inactivated Chikungunya virus.

Use: open the package on the opposite side from the code (CH followed by the lot number) which is useful for identification purposes, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

CONJ CONJUGATE 1x15 mL

Contents: human anti-IgM monoclonal antibodies labelled with HRP, in phosphate buffer solution containing phenol 0.05% and Bronidox 0.02%. Ready for use without further dilution.

CONTROL - NEGATIVE CONTROL (PF93933) 1x1.6 mL

Contents: Human serum not containing anti-Chikungunya antibodies, diluted in Phosphate buffer 0.01 mol/L, with BSA 1% and sodium azide 0.09%. Liquid, ready for use without further dilution.

CONTROL + POSITIVE CONTROL (PF 91833) 1x1.6 mL

Contents: Human serum, containing anti-Chikungunya IgM antibodies, diluted in Phosphate buffer 0.01 mol/L, with BSA 1% and sodium azide 0.09%. Liquid, ready for use.

CUT-OFF CONTROL (PF 91834) 1x2 mL

Contents: Human serum, containing anti-Chikungunya IgM antibodies, diluted in Phosphate buffer 0.01 mol/L, with BSA 1% and sodium azide 0.09%. Liquid, ready for use.

WASH BUF 10x WASH BUFFER 10X (PF93603) 1x100 mL INTERCHANGEABLE BETWEEN LOTS

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5%.

Preparation: dilute the required volume 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SORBENT CHM SORBENT (PF30173) 1x6 mL

Contents: Proteic solution in phosphate buffer containing goat IgG with sodium azide 0.09% and methyl orange as dye.

DILUENT CHK SAMPLE DILUENT (PF30266) 1x60 mL INTERCHANGEABLE BETWEEN LOTS

To be used to dilute samples.

Contents: Proteic solution in phosphate buffer with sodium azide 0.09% containing methyl orange as dye.

SUBS TMB SUBSTRATE (PF93619) 1x12 mL. Ready for use INTERCHANGEABLE BETWEEN LOTS

Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilised in citrate buffer 0.05 mol/L (pH 3.8).

H₂SO₄ 0.3 M STOP SOLUTION (PF93602) 1x16 mL INTERCHANGEABLE BETWEEN LOTS

H₂SO₄ 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

MATERIALS REQUIRED BUT NOT PROVIDED

- Incubator at 37°C
- Microplate reader (wave length 450 or 450/620 nm, with linearity up to OD ≥ 2000)
- Microplate washer (preferable) able to dispense volumes in the range 225-375 µL
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µL solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
MICROPLATE	6 weeks at 2/8°C, polythene bag
CONTROLS	6 weeks at 2/8°C
CONJUGATE	6 weeks at 2/8°C
SUBSTRATE	up to the expiry date at 2/8°C, 1 week at 15-30°C in the dark
SAMPLE DILUENT	up to the expiry date at 2/8°C
WASH BUFFER	ready for use 2 weeks at 2/8°C, 5 days at 15/30°C.
STOP SOLUTION	up to the expiry date at 2/8°C

6. PRECAUTIONS**FOR IN VITRO DIAGNOSTIC USE ONLY**

Caution: This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Avoid the use of self-defrosting freezers for the storage of samples.
6. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
7. Do not allow wells to become dry during the assay procedure.
8. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
9. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
10. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO₂ incubators must not be used.
11. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
12. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.
13. The use of the kit on automatic instruments must be validated by the user.
14. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:
 - installation and particular requisites
 - operating principles, instructions, precautions and risks
 - manufacturer's specifications and instrument performance
 - servicing and maintenance.

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted

samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric or contaminated samples should be avoided.

8. **TEST PROCEDURE**

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Wash Buffer 10x (100 mL + 900 mL H₂O).
- Dilute samples 1:51 distributing 10 µL of serum into 0.5 mL of Diluent CHK.
- Dispense undiluted controls (100 µL per well) in duplicate at least.
- Dispense 50 µL of Sorbent CHM in each well where the examined sera will be placed and add 50 µL of each diluted sample.
- Leave one well for the blank (100 µL of substrate)

Wells are covered with protective film and incubated for 45 minutes at 37°C. After washing four times for 30 seconds (300 µL), add 100 µL of conjugate to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film. The plate is washed again 4 times, as described above. Finally, the substrate is distributed, 100 µL/well. After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µL of Stop Solution.

The adsorbance (O.D.) is read at 450 nm or 450/620 nm within 30 min. Take the reading again at 405 nm if there are OD higher than 2.000.

9. SCHEME OF TEST PROCEDURE

Manual Technique

- | | |
|--------|---|
| STEP 1 | Place 100 µL of undiluted controls in the strip wells. |
| | - |
| STEP 2 | Place 50 µL of Sorbent CHM to the wells for the samples |
| | - |
| STEP 3 | Add 50 µL of diluted serum 1:51 |
| | - |
| | Incubate for 45 min. at 37°C |
| | - |
| | Wash 4 times (300 µL) |
| | - |
| STEP 4 | Add 100 µL of conjugate to each well |
| | - |
| | Incubate for 45 min. at 37°C |
| | - |
| | Wash 4 times (300 µL) |
| | - |
| STEP 5 | Add 100 µL of Substrate to each well |
| | - |
| | Incubate for 15 min. at R.T. |
| | - |
| STEP 6 | Add 100 µL of Stop Solution |
| | - |
| | Read absorbance at 450 nm within 30 min |

10. **TEST VALIDATION**

Subtract the value of the blank (< 0.150) from all the other readings.

Negative control: the ratio between the negative control OD and cut-off OD must be < 0.6.

Positive control: the positive control must have an OD at least 1.5 times higher than the cut-off.

Cut-off control: the OD of the cut-off must be > 0.200.

11. **INTERPRETATION OF THE RESULTS**

Calculate the ratio O.D. sample/OD Cut-off (Index).

The sample is considered:

Positive: if the Index is > 1.1

Negative: if the Index is < 0.9

Doubtful: if the Index is between 0.9 and 1.1

12. LIMITATIONS OF THE PROCEDURE

All positive test results require careful interpretation since false positive reactions may occur.

Absence of anti-Chikungunya specific IgM does not exclude the possibility of infection.

The test results should always be used in conjunction with information available from the clinical evaluation and other diagnostic procedures.

Plasma or heat-inactivated serum should not be used since it is not known how they could affect the test.

Samples containing antibodies known to give cross-reactions with Chikungunya Virus have not been studied for the potential effects.

13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

In an experimentation 126 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	30	3	33
	-	1	92	93
	Total	31	95	126

Diagnostic Sensitivity: 96.8% (CI_{95%}: 83.8-99.3)

Diagnostic Specificity: 96.8% (CI_{95%}: 91.1-98.9)

14. PRECISION

Sample	Within-run Precision	
	Mean (Index)	CV%
1	1.0	7.0
2	3.9	5.4

Sample	Between-run Precision	
	Mean (Index)	CV%
1	0.4	12.5
2	0.7	12.9
3	1.5	14.0
4	5.2	11.2

15. TROUBLE SHOOTING GUIDE

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert paragraph 4 for correct code).
		Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow). Repeat test
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.
	Inadequate washing	Ensure that wash apparatus works well
Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate Color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.

	Inadequate volume of substrate added to the plate	Check pipette function.
--	---	-------------------------

16. **REFERENCES**

1. Abrahm & Shridharan “Chikungunya Virus Infection a Resurgent Scourge “Indian J. med. Res. 126,December 2007, pp 502-504
2. M. Chhabra, V. Mittal, D. Bhattacharya, UVS Rana, S Lal “Chikungunya Fever: a re-emerging viral infection” Indian J. of Medical Microbiology, 2008- 26(1): pp 5-12

	EN ES IT	Date of manufacture Fecha de fabricación Data di fabbricazione	FR GR PT	Date de fabrication Ημερομηνία Παραγωγής Data de fabrico
	EN ES IT	Use By Fecha de caducidad Utilizzare entro	FR GR PT	Utiliser jusque Ημερομηνία λήξης Prazo de validade
	EN ES IT	Caution, consult accompanying documents Atención, ver instrucciones de uso Attenzione, vedere le istruzioni per l'uso	FR GR PT	Attention voir notice d'instructions Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα Atenção, consulte a documentação incluída
	EN ES IT	Manufacturer Fabricante Fabbicante	FR GR PT	Fabricant Κατασκευαστής Fabricante
	EN ES IT	Contains sufficient for <n> tests Contenido suficiente para <n> ensayos Contenuto sufficiente per "n" saggi	FR GR PT	Contenu suffisant pour "n" tests Περιεχόμενο επαρκές για «n» εξετάσεις Conteúdo suficiente para "n" ensaios
	EN ES IT	Temperature limitation Límite de temperatura Limiti di temperatura	FR GR PT	Limites de température Περιορισμοί θερμοκρασίας Limites de temperatura
	EN ES IT	Consult Instructions for Use Consulte las instrucciones de uso Consultare le istruzioni per l'uso	FR GR PT	Consulter les instructions d'utilisation Συμβουλευτείτε τις οδηγίες χρήσης Consulte as instruções de utilização
	EN ES IT	Biological risks Riesgo biológico Rischio biologico	FR GR PT	Risques biologiques Βιολογικοί κίνδυνοι Risco biológico
	EN ES IT	Catalogue number Número de catálogo Numero di catalogo	FR GR PT	Référence du catalogue Αριθμός καταλόγου Referência de catálogo
	EN ES IT	In Vitro Diagnostic Medical Device Producto sanitario para diagnóstico in vitro Dispositivo medico-diagnostico in vitro	FR GR PT	Dispositif médical de diagnostic in vitro In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν Dispositivo médico para diagnóstico in vitro
	EN ES IT	Batch code Código de lote Codice del lotto	FR GR PT	Code du lot Αριθμός Παρτίδας Código do lote

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Diesse Diagnostica Senese SpA