CHORUS

DIESSE

ADENOVIRUS IgA

REF
81198
REF
81198/12

DIESSE Diagnostica Senese
S.p.A.

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CHORUS Adenovirus IgA

1. INTENDED USE

Immunoenzymatic method for the qualitative determination of IgA-class antibodies to *Adenovirus* in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

CHORUS Adenovirus IgA For the Qualitative Determination of IgA Antibodies to Adenovirus

For *In Vitro* Diagnostic Use Only

2. INTRODUCTION

Adenoviruses are nonenveloped icosahedral particles of 60-90 nm in diameter with fiber proteins on each spike that facilitate the attack against the target cells: they represent the largest nonenveloped viruses. Because of their large size they are able to be transported inside the target cells through pores (envelope fusion is not necessary). Up to now, there are 51 described serotypes in humans, divided in 6 species. Most infections with adenovirus result in infections of the respiratory tract that range from conjunctivitis to gastroenteritis. Adenovirus infections are common and frequent and are mainly contracted during childhood. The course of these illnesses is often latent, for this reason the virus can still be detected in the tonsils two years after the contraction of the infection. The virus is transmitted by saliva and stools and is passed from one person to another through mouth, nose and eye.

Most of the infections are asymptomatic. About 5% of colds contracted by children are caused by Adenovirus. The virus can spread among people in crowded environments causing epidemics. Adenovirus infections can be detected through direct verification (pathogenic antigen growth or biochemical identification) or through serological analysis that are represented by complement fixation tests (CFT) and immunoassay tests (ELISA). Viral preparations mainly containing genospecific epitopes are generally used. Currently the verification of the presence of subtypes in case of respiratory infections does not result as necessary from a clinical point of view.

Among the normal population a high seropositivity can be verified and is caused by infections from different subtypes contracted during childhood. In case of an acute adenoviral

infection the normal serological procedures record significant increases of the titer on two sera collected at a distance of 8-10 days from each other. With the development of the ELISA test procedures the evaluation of single sera is acquiring importance.

Even though the IgM determination has a diagnostic value in cases of primary infections, the IgA determination, together with the IgG determination, represents the main immunitary response for patients with Adenovirus infections.

3. PRINCIPLE OF THE ASSAY

The test is based on the ELISA principle (Enzyme Linked ImmunoSorbent Assay) which uses the reaction between the antibodies present in the tested sample and the immobilized antigen bound to solid phases of the polystyrene.

The immunoglobulins bind to the antigen through incubation with diluted human serum. After washings to eliminate the proteins which have not reacted, incubation with the conjugate (anti-human IgA monoclonal antibodies conjugated with horse radish peroxidase) is performed. The conjugate which has not been bound is eliminated and the peroxidase substrate added.

The blue color which develops is proportional to the concentration of specific antibodies present in the test serum.

The disposable devices contain all the reagents to perform the test when applied on the Chorus instruments. The result is expressed as an INDEX (ratio between the OD value of the sample and that of the Cut-off).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the CHORUS instrument.
4. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The conjugate and controls contains phenol
 - b) The substrate is acid

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.

5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
5. Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument.
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed samples or samples presenting microbial contamination may all constitute sources of error.
11. Before inserting the devices in the instrument, check that the reaction well does not contain foreign bodies.
12. Pipette the test serum (50 µL) in well 1 of the device (see figure).
13. Do not use the device after the expiry date.

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 81198).


The kit is sufficient for 12 tests (REF 81198/12).

DD DEVICES

6 packages each containing 6 devices (REF 81198).

2 packages each containing 6 devices (REF 81198/12).

Description of device:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with Adenovirus antigen

Position 5: MICROPLATE WELL

Uncoated

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/l citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing phenol 0.05%, Bronidox 0.02% and an indicator to reveal the presence of the serum.

Position 2: CONJUGATE

Contents: monoclonal antibodies anti-human IgA labeled with peroxidase in phosphate buffer containing phenol 0.05% and Bronidox 0.02%

Position 1: EMPTY WELL

In which the operator must place the undiluted serum

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.175 mL

Content: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin. Liquid, ready for use.

CONTROL + POSITIVE CONTROL

1 x 0.425 mL

Content: Diluted human serum, at known antibody concentration, containing Proclin and Gentamycin. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER (REF 83606)
- CLEANING SOLUTION 2000 (REF 83609)
- SANITIZING SOLUTION (REF 83604 - 83608)
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT (REF 83607)
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples cannot be used. The test cannot be applied to plasma.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions, points 1 and 8.
3. Dispense 50 µL of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected. If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus instrument expresses the result as an INDEX (lot-dependent - ratio between the OD value of the sample and that of the Cut-off) stored in the instrument.

The test serum can be interpreted as follows:

POSITIVE: when the ratio is > 1.1
NEGATIVE: when the ratio is < 0.9
DOUBTFUL/EQUIVOCAL: for all values between 0.9 and 1.1

In the case of a doubtful/equivocal result, repeat the test. If the test remains doubtful/equivocal, collect a new sample after 1-2 weeks.

11. LIMITATIONS

The test cannot be used as only method for a clinical diagnosis. Negative results may not exclude an eventual infection. Samples of serum taken during the acute stage of the infection could result negative with this method because the seroconversion can take weeks (3-4) until its reveal. In case of positive results, it is recommendable to analyze a second sample obtained 8-14 days later, in order to check a possible increase of the IgA-class antibodies. However, the test results should be interpreted with caution and used in conjunction with information available from the clinical evaluation and other diagnostic procedures.

12. ANALYTICAL SPECIFICITY

The following 11 samples containing potentially interfering substances were tested:

Triglycerides (n=2)
C-Reactive Protein (n=7)
Rheumatoid factor (n=2)

The presence in the tested serum of the interfering substances listed above did not alter the results of the test.

13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

81 samples were tested in an experimentation with DIESSE kit as well as with another commercial kit; the discordant samples were confirmed using an immunofluorescence test.

Below are the schematized data of the trial:

		Reference		
		+	-	Total
Diesse	+	16	3	19
	-	0	62	62
	Total	16	65	81

Diagnostic Sensitivity: 100.0 % CI_{95%}: 80.6 – 100.0
Diagnostic Specificity: 95.4 % CI_{95%}: 87.3 - 98.4

14. PRECISION

WITHIN-RUN PRECISION

	Replicates	Mean	SD	CV%
Lot 616-S	15	2.4	0.24	10.0
Lot 617-S	15	3.1	0.26	8.4
Lot 618-S	15	2.7	0.20	7.4

BETWEEN-RUN PRECISION

	Sample	Mean	SD	CV%
Lot 616-S	Standard 1	0.3	0.00	N/A
	Standard 2	0.9	0.10	11.1
	Standard 3	3.1	0.17	5.5
Lot 617-S	Standard 1	0.3	0.00	N/A
	Standard 2	1.0	0.15	15.0
	Standard 3	3.4	0.12	3.5
Lot 618-S	Standard 1	0.3	0.00	N/A
	Standard 2	1.0	0.17	17.0
	Standard 3	3.3	0.50	15.2

N/A – the CV% cannot be calculated because the mean corresponds to 0 (division by 0)

15. REFERENCES

1. Wu and Nemerow (2004). Virus yoga: the role of flexibility in virus host cell recognition. Trends Microbiol 12: 162–168. DOI:10.1016/j.tim.2004.02.005.
2. Meier and Greber (2004). *Adenovirus endocytosis*. J Gene Med 6: S152–S163. DOI:10.1002/jgm.553.
3. Hierholzer, J. C. und Tannock, G. A. Adenoviruses. in Rose, N. R., Friedman, H. und Fahey, J. L. (Eds.) Manual of Clinical Laboratory Immunology. American Society for Microbiology, 3rd ed. (1986) 527-531.
4. Pachucki C.T. Seminars in Respiratory Infections, 7, 46-53, 1992.
5. Harmon M.W. "Laboratory Diagnosis of viral infections" Ed. Lennette E.H., Smith T.S., 587-601, (1999).
6. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
7. R. Zieglermaier et al.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980)

16. GLOSSARY OF LABELING SYMBOLS

	Date of manufacture
	Use By
	Do not reuse
	Caution, consult accompanying documents
	Manufacturer
	Contains sufficient for <n> tests
	Temperature limitation
	Consult instructions for use
	Biological risks
	Catalogue number
	<i>In vitro</i> diagnostic medical device
	Batch code

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DIESSE Diagnostica Senese S.p.A.