CHORUS

HSV 2 SCREEN

REF 81019 (36 tests)

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15. BIBLIOGRAFIA

1. W.E. Rawls et al. Measurement of antibodies to herpesvirus types 1 and 2 in human sera. J. Immunol. 104: 599 (1970).
2. J.J. Gibson et al. A cross-sectional study of herpes simplex virus types 1 and 2 in college students: occurrence and determinants of infection. J. Inf. Dis. 162: 306 (1990).
3. R.M. Coleman et al. Determination of Herpes Simplex Virus antibodies by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 18: 287 (1983).

INSTRUCTIONS FOR USE

CHORUS HSV 2 SCREEN

REF 81019

(English)

1. INTENDED USE

QUALITATIVE DETERMINATION OF ANTI-HERPES SIMPLEX VIRUS (TYPE 2) IMMUNOGLOBULINS IN HUMAN SERUM, USING A DISPOSABLE DEVICE APPLIED ON THE CHORUS INSTRUMENT

2. INTRODUCTION

The Herpes simplex virus (HSV) is a member of the Herpesviridae family, of which two types are known: type 1 (HSV-1) and type 2 (HSV-2) which present slight antigenic differences. HSV-1 causes chiefly oral-facial lesions, while HSV-2 is mainly responsible for genital lesions, but this distinction is not binding, both types occasionally causing infection in either anatomical site. HSV may also cause a form of ocular cheratitis, and lesions of the central nervous system.

HSV can affect practically the whole population. The primary infection is often in a subclinical form and is rarely diagnosed. After a latency period of variable duration, reactivation may occur and viral replication may or may not give rise to clinical lesions. Infection contracted during birth is of particular interest, this being an important cause of morbidity and mortality. It is therefore important to determine the immunitary state of women during pregnancy in order to detect serum conversion.

3. PRINCIPLE OF THE METHOD


The test is based on the principle of analysis by competition. The anti-HSV-2 antibodies present in the sample compete with the marker (peroxydase/anti-HSV 2 monoclonal antibody conjugate) to occupy the few binding sites available on the antigen coated on the solid phase. The higher the antibody concentration in the sample, the lower the number of conjugated antibodies which are bound. The unbound components are eliminated by washing, and the bound enzymatic activity is determined colorimetrically by transformation of a chromogen substrate (blue-coloured reaction). The intensity of the colour which develops is inversely proportional to the antibody concentration in the sample. The disposable devices contain all the reagents to perform the test when applied on the Chorus instrument.

4. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Use: **equilibrate a package at room temperature**, open the package and remove the devices required; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

Description:

Position 8: Space for application of bar code label

Position 7: empty

Position 6: MICROPLATE WELL coated with Herpes Simplex Virus type 2

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: empty

Position 2: CONJUGATE

Contents: anti-HSV 2 monoclonal antibodies labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL in which the operator must place the undiluted serum

CALIBRATOR CALIBRATOR 1 x 0.3 mL

Contents: Newborn Calf Serum in phosphate buffer 0.01 mol/L with phenol 0.05% and Bronidox 0.02%, liquid, ready for use.

CONTROL + POSITIVE CONTROL

Contents: Freeze-dried human serum containing anti-HSV type 2 antibodies.

Use: reconstitute with the volume of deionised water reported on the label, **wait at least 5 minutes without shaking; then mix by inversion avoiding froth formation.**

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER [REF] 83606
- CLEANING SOLUTION 2000 [REF] 83609
- SANITIZING SOLUTION [REF] 83604 – 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT [REF] 83607
- Chorus / Chorus Trio Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see point 9, Test validation)

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening and/or preparation:

DEVICES	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C

6. PRECAUTIONS

ONLY FOR DIAGNOSTIC USE IN VITRO.

Attention:

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens. Wash hands thoroughly after placing the devices in the CHORUS instrument.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The conjugate contains phenol
 - b) The substrate is acid
 If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
4. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to

clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
5. Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument
7. Defective barcodes can be inserted manually in the instrument.
8. Do not expose the devices to strong light or to hypochlorite vapours during storage and use.
9. The use of strongly hemolyzed samples, or samples presenting microbial contamination may all constitute sources of error.
10. Before inserting the devices in the instrument, check that the reaction well does not contain foreign bodies.
11. Resuspend the freeze-dried control serum homogeneously according to the volumes reported on the labels.
12. Pipette the test serum (50 µl) in well 1 of the device (see figure).
13. Do not use the device after the expiry date.

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers.

Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples should be avoided.

The test cannot be applied to plasma.

8. PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in the Analytical Warnings points 1 and 8.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus instrument expresses the result as an INDEX (ratio between the OD value of the Cut-Off and that of the test sample).

The test serum can be interpreted as follows:

POSITIVE: when the ratio is > 1.2

NEGATIVE: when the ratio is <0.8

DOUBTFUL: when the ratio is between 0.8 and 1.2

In the case of a doubtful result, repeat the test. If the test remains doubtful, collect a new sample.

11. LIMITATIONS OF THE PROCEDURE

The test is not able to discriminate between the presence of IgG and that of IgM.

The test cannot be used on its own for clinical diagnostic purposes. A negative result does not exclude the possibility of the disease.

As in any immunoenzymatic method, results may be found which are not repeatedly reactive.

12. ANALYTICAL SPECIFICITY

The presence of antibodies against viruses such as CMV, Epstein Barr., Rubella and Mumps does not influence the result of the test

13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

56 serum samples were tested, in comparison with another commercial method to determine the presence of specific Ig. 33 negative and 15 positive samples were found, with agreement between the different methods. The results which were not in agreement gave 2 false negative and 6 false positive results.

The results are summarized in the following table.

Table 1

		REFERENCE	
		+	-
DIESSE	+	15	6
	-	2	33

The kit shows a sensitivity of 88% (CI 95%: 66–97) and a specificity of 85% (CI 95%: 70–93).

14. PRECISION**Table 2. Within-run precision**

	No. of replicates	Mean Index	St. Dev.	CV%
Lot no. 033	9	3.3	0.18	5
Lot no. 034	6	2.5	0.16	6
Lot no. 035	6	3.0	0.21	7

Table 3. Between-run precision and Precision between batches

Sample	Mean Index			Mean	St. Dev.	CV%
	Lot 033	Lot 034	Lot 035			
H2S 1	0.5	0.4	0.4	0.4	0.06	15
H2S 2	1.3	1.2	1.1	1.3	0.15	12
H2S 3	4.0	3.9	3.6	3.5	0.26	7

15. REFERENCES

1. W.E. Rawls et al. Measurement of antibodies to herpesvirus types 1 and 2 in human sera. J. Immunol. 104: 599 (1970).
2. J.J. Gibson et al. A cross-sectional study of herpes simplex virus types 1 and 2 in college students: occurrence and determinants of infection. J. Inf. Dis. 162: 306 (1990).
3. R.M. Coleman et al. Determination of Herpes Simplex Virus antibodies by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 18: 287 (1983).

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