CHORUS

INFLUENZA A

IgG

DIESSE
REF 81190
REF 81190/12

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CHORUS INFLUENZA A IgG

1. INTENDED USE

Immunoenzymatic method for the qualitative determination of IgG-class antibodies to *Influenza A Virus* in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

CHORUS INFLUENZA A IgG For the Qualitative Determination of IgG Antibodies to *Influenza A*

For In Vitro Diagnostic Use Only

2. INTRODUCTION

Influenza viruses are categorized in type A, B, and C on the basis of the antigenic differences of their nucleoprotein (NP) and matrix (M1) proteins. Only type A and B are perceived to be clinically relevant in humans.

Influenza virus infections are the causal agents of recurrent epidemics of acute respiratory diseases in humans. In particular, influenza is highly contagious, can spread easily and is responsible for a considerable morbidity and mortality each year. Elderly and compromised individuals are especially at risk of developing severe illness and complications.

During epidemics influenza viruses can cause 10.000-20.000 deaths in the elderly and in patients with chronic cardiovascular and pulmonary diseases.

The clinical symptoms of influenza are very similar to those associated with other respiratory viruses, often circulating in the community at the same time. The immune response to influenza virus infection is influenced by the previous exposure history of the host to the influenza antigens. Serum antibodies appear in the 2nd week after the onset of the illness, reach peak titers by 4 weeks, and persist for months to years before gradually declining.

Serological diagnosis of acute influenza infection is based upon the detection a fourfold or greater increase in the IgG titer between the acute-phase and the convalescent-phase sera. Thus, it requires collections of paired serum samples at a distance of about 2 weeks between the first and the second collection. Although a serological diagnosis through the detection of IgG is retrospective and therefore of limited usefulness in many clinical situations, it is considered essential for surveying influenza outbreaks, for assessing the efficacy of vaccines and specific antiviral therapies.

3. PRINCIPLE OF THE ASSAY

The test is based on the ELISA principle (Enzyme Linked ImmunoSorbent Assay).

Influenza A antigen is bound to the solid phase. Specific immunoglobulins bind to the antigen through incubation with diluted human serum. After washings to eliminate the proteins which have not reacted, incubation with the conjugate (anti-human IgG monoclonal antibodies conjugated with horse radish peroxidase) is performed. The conjugate which has not been bound is eliminated and the peroxidase substrate added.

The blue colour which develops is proportional to the concentration of specific antibodies present in the test serum.

The disposable devices contain all the reagents to perform the test when applied on the Chorus instruments.

The results are expressed in Index – ratio between OD value of the test sample and that of the Cut-Off.

4. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the CHORUS instrument.
4. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The conjugate contains phenol
 - b) The substrate is acid

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minutes exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills

unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
5. Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument and do not use devices with defective labels.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument.
9. Do not expose the devices to strong light or to hypochlorite vapours during storage and use.
10. The use of strongly hemolyzed samples or samples presenting microbial contamination may all constitute sources of error.
11. Before inserting the devices in the instrument, check that the reaction well does not contain foreign bodies.
12. Pipette the test serum (50 µl) in well 1 of the device (see figure).
13. Do not use the device after the expiry date.

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 81190).


The kit is sufficient for 12 tests (REF 81190/12).

DD DEVICES

6 packages each containing 6 devices (REF 81190).

2 packages each containing 6 devices (REF 81190/12).

Description of device:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with *Influenza A* antigen

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing phenol 0.05%, Bronidox 0.02% and an indicator to reveal the presence of the serum.

Position 2: CONJUGATE

Contents: monoclonal antibodies anti-human IgG labelled with peroxidase in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

In which the operator must place the undiluted serum

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.175 mL

Contents: Diluted human serum containing IgG antibodies anti-Influenza A and preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.425 mL

Contents: Diluted human serum containing IgG antibodies anti-Influenza A and preservative. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER (REF 83606)
- CLEANING SOLUTION 2000 (REF 83609)
- SANITIZING SOLUTION (REF 83604 - 83608)
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT (REF 83607)
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples cannot be used. The test cannot be applied to plasma.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions, points 1 and 8.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected. If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554

Fax: 0039 0577 366605

email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus instrument expresses the result as an INDEX (ratio between the OD value of the sample and that of the Cut-off).

The test serum can be interpreted as follows:

POSITIVE: when the ratio is > 1.1

NEGATIVE: when the ratio is < 0.9

DOUBTFUL/EQUIVOCAL: for all values between 0.9 and 1.1

In the case of a doubtful/equivocal result, repeat the test. If the test remains doubtful/equivocal, collect a new sample after 1-2 weeks.

11. LIMITATIONS

The test cannot be used as only method for a clinical diagnosis. Negative results may not exclude an eventual infection.

Samples of serum taken during the acute stage of the infection could result negative with this method because the seroconversion can take weeks (3-4) until its reveal. In this

case or in case of IgG positive results, it is recommended an IgM and IgA test and the analysis of a second sample obtained 8-14 days later, in order to check a possible increase of the IgG-class antibodies. However, the test results should be interpreted in conjunction with information available from the clinical evaluation and other diagnostic procedures.

IgG positive results in new born babies must be carefully interpreted since IgG-class antibodies can be transmitted from mother to fetus.

Serological Interpretation		
IgG	IgA/IgM	Diagnostic significance
+	-	Previous infection
-	+	Recent or current infection; a new collection after 1-2 weeks is recommended in order to highlight the IgG seroconversion
+	+	Possible recent infection The positivity of specific IgM and IgA can persist for a significant period

12. ANALYTICAL SPECIFICITY

Two serum samples (Serum A – at cut-off, Serum B – negative) were spiked with the following interfering factors:

Rheumatoid Factor

Bilirubin (5 mg/dl – 30mg/dl)

Triglycerides (1 mg/ml – 30 mg/ml)

Hemoglobin (1.2 mg/ml – 10 mg/ml)

The presence in the tested serum of the interfering substances listed above did not alter the results of the test.

13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

In an experimentation, 107 samples were analyzed with the Diesse kit as well as with another commercial method. Below are the schematized results of the trial:

		Reference		
		+	-	Total
Diesse	+	106	0	106
	-	0	1	1
	Total	106	1	107

Diagnostic Sensitivity: 100.0 % CI_{95%}: 96.5 – 100.0

Diagnostic Specificity: 100.0 % CI_{95%}: 20.1 – 100.0

14. PRECISION

WITHIN-RUN PRECISION

Sample	No. of replicates	Mean	SD	CV%
IAg-C 2	15	0.9	0.07	7.8

BETWEEN-RUNS PRECISION

Sample	Mean Index			Mean	SD	CV%
	Day 1	Day 2	Day 3			
IAG-C 1	0.1	0.1	0.1	0.1	0.00	-
IAG-C 2	0.9	1.0	0.8	0.9	0.10	11.1
IAG-C 3	2.2	2.3	2.1	2.2	0.10	4.5

15. REFERENCES

1. Voeten J.T.M et al. Journal of Clinical Microbiology, 36, 3527-3531, 1998.
2. Pachucki C.T.. Seminars in Respiratory Infections, 7, 46-53, 1992.
3. Harmon M.W. "Laboratory Diagnosis of viral infections" Ed. Lennette E.H. , Smith T.S., 587-601 , (1999).
4. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
5. R. Zieglermaier et al.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980)

16. GLOSSARY OF LABELING SYMBOLS

	Date of manufacture
	Use By
	Do not reuse
	Caution, consult accompanying documents
	Manufacturer
	Contains sufficient for <n> tests
	Temperature limitation
	Consult instructions for use
	Biological risks
	Catalogue number
	<i>In vitro</i> diagnostic medical device
	Batch code

Manufactured by
DIESSE Diagnostica Senese S.p.A.
 Strada dei Laghi 39
 53035 Monteriggioni (SIENA)
 Italy

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