ENZY-WELL

Epstein-Barr VCA IgM

REF 91056

DIESSE

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18. BIBLIOGRAFIA

1. A. Andersson, V. Vetter et al. Avidities of IgG directed against viral capsid antigen or early antigen: useful markers for significant Epstein-Barr Virus serology. J. Med. Virology 43: 238 (1994).
2. J. Middeldorp and P. Herbrink: Epstein Barr Virus specific marker molecules for early diagnosis of infectious mononucleosis. J. Virol. Methods 21: 133 (1988).
3. C. Valent Sumaya: Serological testing for Epstein Barr Virus - developments in interpretation. J. Inf. Dis. 151: 984 (1985).
4. J. Luka, R.C. Chase and G. Pearson: A sensitive enzyme-linked immunosorbent assay (ELISA) against the major EBV-associated antigens. I. Correlation between ELISA and immunofluorescence titers using purified antigens. J. Immunol. Methods 67: 145 (1984).
5. Van Grunsven WM, van Heerde EC, de Haard HJ, Spaan WJ, Middeldorp JM: Gene mapping and expression of two immunodominant Epstein-Barr virus capsid proteins. J Virol. 1993 Jul;67(7):3908-16.

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INSTRUCTIONS FOR USE

ENZY-WELL Epstein-Barr VCA IgM

**For the semiquantitative determination of IgM
antibodies anti- Epstein-Barr viral capsid antigen**

For In Vitro Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the semiquantitative determination of IgM class antibodies to Epstein-Barr viral capsid antigen in human serum.

2. INTRODUCTION

Epstein Barr Virus (EBV) is a herpesvirus which causes infectious mononucleosis (IM). It is also associated with Burkitt's lymphoma, nasopharyngeal carcinoma and lymphatic proliferative syndromes in immunodepressed patients. The virus is widespread throughout the world and 80-90% of the population is serum-positive.

The laboratory diagnosis of IM is traditionally performed by detecting heterophile antibodies which develop in the serum during the course of the infection, and which agglutinate horse erythrocytes. However, these antibodies may not always be present in patients affected by IM, particularly if below 14 years of age; furthermore, they may also persist for over a year after the infection. The determination of heterophile antibodies alone may therefore lead to an erroneous diagnosis. It is therefore important to determine the presence of antibodies towards the viral antigens. In particular, the detection of antibodies directed to the "Viral Capsid Antigen" (VCA) and the nuclear antigen (EBNA) is particularly useful.

During the course of IM, the IgM- and IgG-class antibodies to VCA appear early, while the IgG to EBNA develop later during the infection. The presence of IgM against VCA in the absence of IgG against EBNA therefore indicates that there is a current infection, while the presence of IgG against both VCA and EBNA is indicative of a prior infection.

3. PRINCIPLE OF THE TEST

The test is based on the ELISA technique (Enzyme Linked ImmunoSorbent Assay). The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with dilute human serum. After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgM monoclonal antibodies conjugated to horseradish peroxidase. The unbound conjugate is eliminated, and the peroxidase substrate added. The blue color which develops is proportional to the concentration of specific antibodies present in the serum sample. When the enzymatic reaction is interrupted by the addition of a sulphuric acid solution, a yellow coloring develops. The color, proportional to the quantity of specific antibodies present in the serum, can be easily read using a microplate reader.

4. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens and performing the assay.
3. Wash hands thoroughly when finished.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Non-disposable apparatus should be sterilized after use, the preferred method is to autoclave for 1 h at 121°C. All the disposable apparatus should be eliminated according to law.
6. 2 M hydrochloric acid, used for washing glassware, is corrosive and should be handled with appropriate care. In case of contact with skin or eyes, wash thoroughly with water.
7. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be sufficient to ensure effective decontamination.
8. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

1. The OD of Cut-Off, controls and samples can be slightly different among different plates. For such reason if during the same run strips from different plates are used, even if the lot is the same, it is necessary to repeat the determination of the Cut-Off.
2. Allow all reagents to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **Working at a correct temperature is important for the incubation of the strips. Check that the thermostat does not go below 35°C or over 39°C.**

3. Open the package containing the strips after 30 minutes at room temperature at least.
4. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
5. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
6. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or deionized water.
7. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
8. Do not allow wells to become dry during the test.
9. Do not keep the samples in auto-defrosting freezers.
10. Care should be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the substrate and conjugate.
11. Care should be taken to avoid touching the rim of the well with conjugate.
12. Do not "blow-out" from microplates.
13. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack to support the plates, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate Operating Manual for further details. CO₂ incubators must not be used.
14. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
15. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
16. Care should be taken to avoid contaminating the microplate wells with disposable gloves powder.
17. The use of the kit with automated equipment has to be validated by the user.
18. For each instrument used, read the Operating Manual, in particular to obtain additional information on the following points:
 - installation and particular requisites
 - operating principles, instructions, precautions and risks
 - manufacturer's specifications and instrument performances
 - technical servicing and maintenance

5. KIT COMPOSITION AND REAGENT PREPARATION

Reagents are sufficient for 96 determinations.

- Bring to room temperature before use.

MT PLATE MICROPLATE 12x8 (PF93054-C)

Content: 1 microplate (96 wells) coated with synthetic viral capsid antigen.

Use: open the package at the opposite end from the code (EM followed by the lot number) which is useful for identification purposes, remove the support and strips required. Place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

CONJ CONJUGATE 1x16 mL (PF93571)

Content: anti-human IgM monoclonal antibodies, labeled with peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%. Liquid, ready for use.

CAL CALIBRATOR 3x1.6 mL (PF91871 – PF30077 – PF92072)

Content: Diluted human serum, containing a known concentration of IgM antibodies, in phosphate buffer 0.01 mol/L with BSA 1% and Sodium Azide 0.09%.

Liquid, ready for use.

The titer in arbitrary units, reported on the label, is 20, 80, 160 AU/ml. Calibrator 1 containing 20 AU/ml corresponds to the cut-off.

Colour: the colour is proportional to the relative antibody titer.

CONTROL IgG - IgG NEGATIVE CONTROL 1x1.6 mL (PF93910)

INTERCHANGEABLE BETWEEN LOTS

Content: Human serum, not containing anti-IgM antibodies anti-Epstein-Barr viral capsid antigen, diluted in phosphate buffer 0.01 mol/L with BSA 1% and Sodium Azide 0.09%. Liquid, ready for use.

SAMP DIL 2 DILUENT 2 1x100 mL (PF93611)

INTERCHANGEABLE BETWEEN LOTS

Content: Proteic solution in phosphate buffer with Sodium Azide 0.09% and methyl orange as dye.

Use: To be used to dilute samples. Liquid, ready for use.

SAMP SORBENT 10 SAMPLE SORBENT 10 1x7 mL (PF30253)

INTERCHANGEABLE BETWEEN LOTS

Content: Proteic solution in phosphate buffer with Sodium Azide 0.09%, methyl orange as dye and anti-human IgG antibodies.

Use: To be used to dilute samples (in the microwell). Liquid, ready for use.

WASH BUF 10x WASH BUFFER 10X 1x100 mL (PF93603)

INTERCHANGEABLE BETWEEN LOTS

Content: Phosphate buffered saline (PBS), concentrated 10 times; contains Brij 0.5%.

Preparation: Dilute the required volume 1:10 with distilled or deionized water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SUBS TMB SUBSTRATE 1x12 mL (PF93619)

INTERCHANGEABLE BETWEEN LOTS

Content: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in citrate buffer 0.05 mol/L (pH 3.8). Liquid, ready for use.

H₂SO₄ 0.3 M STOP SOLUTION 1x16 mL (PF93602)

INTERCHANGEABLE BETWEEN LOTS

Content: H₂SO₄ 0.3 mol/L. Liquid, ready for use.

ADHESIVE FILMS (2).
POLYETHYLENE BAG (1).

MATERIALS REQUIRED BUT NOT PROVIDED

- Incubator at 37-40°C
- Microplate reader (wave length 450 or 450/620 nm, with linearity up to OD ≥ 2.000)
- Microplate washer (optional) able to dispense volumes in the range 225-375 µl
- Distilled or deionized water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl of solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2-8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation:

MICROPLATE	8 weeks at 2-8°C in polyethylene bag
CONJUGATE	8 weeks at 2-8°C
CALIBRATORS	8 weeks at 2-8°C
NEGATIVE CONTROL	8 weeks at 2-8°C
DILUENT 2	up to the expiry date at 2-8°C
SAMPLE SORBENT 10	up to the expiry date at 2-8°C
SUBSTRATE	up to the expiry date at 2-8°C; 1 week at 15-30°C; store in the dark
WASHING BUFFER	2 weeks at 2-8°C; 5 days at 15-30°C
STOP SOLUTION	up to the expiry date at 2-8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

PREPARATION

Bring all the reagents to room temperature (18-30°C) before use.

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Wash Buffer 10x (100 mL + 900 mL H₂O).
- Dilute samples 1:26 distributing 40 µL of serum into 1 mL of diluent.

PROCEDURE

1. Distribution of the samples:
Dispense 50 µL of Sample Sorbent 10 per well and 50 µL of diluted sample (duplicate testing is recommended). Dispense 100 µL of negative control and 100 µL of each calibrator per well.
The minimum requisite is 1 negative control and 1 control for each calibration level.
Leave one well for the blank (100 µL of substrate).
2. Incubation:
Incubate the plate covered with an adhesive protective film for 45 minutes at 37°C.
3. Washing:
Remove the adhesive film, aspirate the contents of all the wells and rinse 4 times, filling each well with 300 µL of washing solution. Wait 30 seconds between each wash.
4. Distribution of the conjugate:
Dispense 100 µL of conjugate in each well.
5. Conjugate Incubation:
Incubate the plate covered with an adhesive protective film for 45 minutes at 37°C.
6. Washing:
Remove the adhesive film, aspirate the contents of all the wells and rinse 4 times, filling each well with 300 µL of washing solution. Wait 30 seconds between each wash.
7. Distribution of the substrate:
Dispense 100 µL of substrate in each well.
8. Substrate incubation:
Incubate the plate for 15 minutes at room temperature.
9. Interruption of the reaction:
Dispense 100 µL of stop solution, in the same order as that followed for point 4.
10. Reading:
Read the O.D. at 450 nm or 450/620 nm within 30 minutes. Repeat the reading at 405 nm in case of O.D. > 2.000.

9. SCHEME OF TEST PROCEDURE FOR ENZY- WELL EPSTEIN-BARR VCA IgM

- STEP 1 Place 50 µL of Sample Sorbent 10 and 50 µL of diluted sample. Place 100 µL of negative control and of each calibrator in the wells of the strips. Mix well.
-
- Incubate for 45 min. at 37°C
-

- Wash 4 times (300 µL)
- STEP 2 Add 100 µL of conjugate to each well
- Incubate for 45 min. at 37°C
- Wash 4 times (300 µL)
- STEP 3 Add 100 µL of Substrate to each well
- Incubate for 15 min. at R.T.
- STEP 4 Add 100 µL of Stop Solution
- Read the O.D. at 450 nm or 450/620 nm within 30 min

10. TEST VALIDATION

Subtract the value of the blank (< 0.150) from all the other readings.

Negative control: the ratio between the negative control OD and cut-off OD (Calibrator 1) must be ≤ 0.6 .

Positive control (Calibrator 3): the OD must be ≥ 1.0

Cut-off control (Calibrator 1): the OD must be ≥ 0.2 .

Repeat the test, if one of the results of the control sera is not within the acceptability range. If the problem persists contact the Scientific Support.

Tel: 0039 0577 319554

Fax: 0039 0577 366605

email: scientificsupport@diesse.it

11. INTERPRETATION OF THE RESULTS

QUALITATIVE RESULTS

Calculate the ratio between the OD value of the sample and that of the Cut-off (Index).

The sample is considered:

POSITIVE: when the Index is > 1.2

DOUBTFUL: for all the values between 0.8 and 1.2

NEGATIVE: when the Index is < 0.8

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

SEMIQUANTITATIVE RESULTS

Report on a graph the OD of the calibrators, after subtracting the OD of the blank, for each calibrator titer (AU/ml).

The corresponding titer of the test sample can be obtained by extrapolation.

A standard curve must be performed for each run.

The sample is considered:

POSITIVE: when the result is > 24

DOUBTFUL: for all the values between 16 and 24

NEGATIVE: when the result is < 16

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample about 2 weeks later.

12. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

An accurate interpretation of an EBV infection must be based on the results of VCA IgG, VCA IgM, EBNA IgG, EA IgG ed EA IgM results.

Performance characteristics have not been evaluated in patients affected by nasopharyngeal carcinoma, Burkitt's lymphoma or other lymphadenopathies associated with EBV, apart from EBV-related mononucleosis.

13. ANALITICAL SPECIFICITY

5 samples (2 Negative, 1 Cut-Off and 2 Positive) were spiked with the following potentially interfering factors and then tested:

Rheumatoid factor (44-220 IU/ml)

Bilirubin (4.5-45 mg/dl)

Triglycerides (10-250 mg/dl)

Hemoglobin (5-30 mg/ml)

The presence in the serum sample of the interfering substances described above does not affect the test result.

14. CROSS-REACTIONS

59 samples, positive to Toxoplasma, Varicella, Cytomegalovirus, Herpes simplex, Rubella and ANA were tested.

No significant cross-reactions were found.

15. METHOD COMPARISON

In an experimentation 485 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	69	3	72
	-	1	412	413
	Total	70	415	485

Percent Positive Agreement (~Diagnostic Sensitivity):

98.6% CI_{95%}: 92.3 – 99.7

Percent Negative Agreement (~Diagnostic Specificity):

99.3% CI_{95%}: 97.9.-99.7

16. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision		Precision between batches	
	Mean (Index)	CV%	Mean (Index)	CV%	Mean (Index)	CV%
1	0.4	12.5	0.4	12.5	0.4	-
2	0.5	14.0	0.5	14.0	0.4	15.0
3	0.3	-	0.3	0.0	0.3	-
4	0.9	7.8	0.9	12.2	0.9	6.7
5	0.9	3.3	0.9	14.4	0.9	11.1
6	1.3	6.9	1.2	7.5	1.2	-
7	1.9	7.4	1.9	8.4	1.8	3.3
8	2.3	6.1	2.7	9.3	2.5	6.0
9	3.3	3.9	3.5	8.3	3.4	4.4
10	2.3	5.2	2.5	10.8	2.3	2.6

17. TROUBLESHOOTING

PROBLEM	POSSIBLE ERROR	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure. Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert).
		Check for moisture in unused plate. (Silica gel desiccant must be pale yellow). Repeat test.
Invalid run (all positive)	Contamination of substrate	Collect new aliquot of substrate.
	Inadequate washing	Ensure that plate washer works well
Poor precision	Inadequate aspiration of wells	Ensure that plate washer works well
	Pipetting error	Check pipette function
	Addition of reagents too slow	Avoid drying of the plate after washing step. Add reagents immediately

	Presence of air-bubbles	Avoid air-bubbles formation during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
	Inadequate volume of substrate added	Check pipette function

18. REFERENCES

1. A. Andersson, V. Vetter et al. Avidities of IgG directed against viral capsid antigen or early antigen: useful markers for significant Epstein-Barr Virus serology. J. Med. Virology 43: 238 (1994).
2. J. Middeldorp and P. Herbrink: Epstein Barr Virus specific marker molecules for early diagnosis of infectious mononucleosis. J. Virol. Methods 21: 133 (1988).
3. C. Valent Sumaya: Serological testing for Epstein Barr Virus - developments in interpretation. J. Inf. Dis. 151: 984 (1985).
4. J. Luka, R.C. Chase and G. Pearson: A sensitive enzyme-linked immunosorbent assay (ELISA) against the major EBV-associated antigens. I. Correlation between ELISA and immunofluorescence titers using purified antigens. J. Immunol. Methods 67: 145 (1984).
5. Van Grunsven WM, van Heerde EC, de Haard HJ, Spaan WJ, Middeldorp JM : Gene mapping and expression of two immunodominant Epstein-Barr virus capsid proteins. J Virol. 1993 Jul;67(7):3908-16.

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36. BIBLIOGRAPHIE

- A. Andersson, V. Vetter et al. Avidities of IgG directed against viral capsid antigen or early antigen: useful markers for significant Epstein-Barr Virus serology. J. Med. Virology 43: 238 (1994).
- J. Middeldorp and P. Herbrink: Epstein Barr Virus specific marker molecules for early diagnosis of infectious mononucleosis. J. Virol. Methods 21: 133 (1988).
- C. Valent Sumaya: Serological testing for Epstein Barr Virus - developments in interpretation. J. Inf. Dis. 151: 984 (1985).
- J. Luka, R.C. Chase and G. Pearson: A sensitive enzyme-linked immunosorbent assay (ELISA) against the major EBV-associated antigens. I. Correlation

between ELISA and immunofluorescence titers using purified antigens. J. Immunol. Methods 67: 145 (1984).

10. Van Grunsven WM, van Heerde EC, de Haard HJ, Spaan WJ, Middeldorp JM: Gene mapping and expression of two immunodominant Epstein-Barr virus capsid proteins. J Virol. 1993 Jul;67(7):3908-16.

DIESSE Diagnostica Senese
S.p.A.
Strada dei Laghi 39
53035 Monteriggioni (SIENA)
Italien

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1. A. Andersson, V. Vetter et al. Avidities of IgG directed against viral capsid antigen or early antigen: useful markers for significant Epstein-Barr Virus serology. J. Med. Virology 43: 238 (1994).
2. J. Middeldorp and P. Herbrink: Epstein Barr Virus specific marker molecules for early diagnosis of infectious mononucleosis. J. Virol. Methods 21: 133 (1988).
3. C. Valent Sumaya: Serological testing for Epstein Barr Virus - developments in interpretation. J. Inf. Dis. 151: 984 (1985).
4. J. Luka, R.C. Chase and G. Pearson: A sensitive enzyme-linked immunosorbent assay (ELISA) against the major EBV-associated antigens. I. Correlation between ELISA and immunofluorescence titers using purified antigens. J. Immunol. Methods 67: 145 (1984).
5. [Van Grunsven WM](#), [van Heerde EC](#), [de Haard HJ](#), [Spaan WJ](#), [Middeldorp JM](#): Gene mapping and expression of two immunodominant Epstein-Barr virus capsid proteins. [J Virol](#). 1993 Jul;67(7):3908-16.

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1. A. Andersson, V. Vetter et al. Avidities of IgG directed against viral capsid antigen or early antigen: useful markers for significant Epstein-Barr Virus serology. J. Med. Virology 43: 238 (1994).
2. J. Middeldorp and P. Herbrink: Epstein Barr Virus specific marker molecules for early diagnosis of infectious mononucleosis. J. Virol. Methods 21: 133 (1988).
3. C. Valent Sumaya: Serological testing for Epstein Barr Virus - developments in interpretation. J. Inf. Dis. 151: 984 (1985).
4. J. Luka, R.C. Chase and G. Pearson: A sensitive enzyme-linked immunosorbent assay (ELISA) against the major EBV-associated antigens. I. Correlation

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