CHORUS

TG High Sensitivity

**DIESSE****REF 86104**

DIESSE Diagnostica Senese
S.p.A.
Via delle Rose, 10
53035 Monteriggioni (Siena)
Italy

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17. BIBLIOGRAFIA

- Indrasena B.S.H. Use of thyroglobulin as a tumor marker. World J Biol Chem 2017 February 26; 8(1) ; 8(1): 81-85.
- Spencer C. A. Clinical Utility of Thyroglobulin Antibody (TgAb) Measurement for Patients with Differentiated Thyroid Cancers (DTC). J Clin Endocrinol Metab, December 2011, 69(12):3615-3627.
- Giovannella L., M.Clark P., Chiovato L., Duntas L., Elisei R., Feldt-Rasmussen U., Leenhardt L., Luster M., Schalin-Jantti C., Schott M., Seregini E., Rimmele H., Smith J., and A Verburg F. DIAGNOSIS OF ENDOCRINE DISEASE. Thyroglobulin measurement using highly sensitive assay in patients with differentiated thyroid cancer: a clinical position paper. European Journal of Endocrinology (2014) 171, R33-R46.
- Spencer C., LoPresti J., and Fatemi S. How sensitive (second-generation) thyroglobulin measurement is changing paradigms for monitoring patients with differentiated thyroid cancer, in the absence or presence of thyroglobulin autoantibodies. Curr Opin Endocrinol Diabetes Obes 2014, 21:394-404.
- Gunter E. Analyte stability and freeze-thaw information Specimen solution LCC.

DIESE Diagnostica Senese
S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

INSTRUCTIONS FOR USE

CHORUS TG High Sensitivity

**For the quantitative determination of human
thyroglobulin (Tg)**

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the quantitative determination of human thyroglobulin (Tg) in human serum, using a disposable device applied on the Chorus TRIO instruments.

2. INTRODUCTION

Thyroglobulin (Tg) is a glycoprotein of 660 kDa localized within the colloid of the thyroid follicle where it is synthesized under the influence of thyrotropin. It plays an essential role in the storage of iodine and acts as precursor for the synthesis of thyroid hormones T3 and T4. Elevated thyroglobulin serum concentrations have been reported in various thyroid diseases, such as hyperthyroidism, non-toxic goitre, thyroiditis and differentiated thyroid carcinoma.

The main clinical application for the Tg determination, however, is represented by the post-operative monitoring of differentiated thyroid carcinoma. The assay of Tg is used for the early detection or exclusion of metastases, tumor relapses and the follow-up of radioiodine treatments. Serum Tg is not detectable in patients who underwent total thyroidectomy and are free of metastases and tumour. These patients in true complete remission will not display Tg levels, even by endogenous TSH stimulation.

Consequently detectable Tg values in this group of patients are an important indication for still existing or newly developed neoplasia, particularly if these detectable Tg values are increasing under a TSH-suppressive thyroid hormone treatment (Tg profiles). For this reason a high sensitivity assay test is very useful.

In case of patients that need an extended period monitoring, it is recommended to keep using the same analytical test, inserting in the same run the samples previously tested as well.

The presence of anti-Tg antibodies in the serum of the patient can influence the Tg assay.

3. PRINCIPLE OF THE METHOD

The Chorus Tg High Sensitivity device is ready to use for the detection of human thyroglobulin (Tg), in the Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay). During two different incubations, the Tg in the sample is bound to two different anti-Tg monoclonal antibodies: one present in the solid phase and the other one conjugated with Biotin.

After washing to eliminate the components which have not reacted, incubation is performed with the conjugate, composed of Streptavidin conjugated to horse radish peroxidase.

The unbound conjugate is eliminated and the peroxidase substrate is added. The enzymatic reaction is subsequently stopped by adding the Stop Solution that turns the solution into yellow. The colour which develops is proportional to the concentration of hormone present in the serum sample.

The disposable devices contain all the reagents to perform the test in the Chorus TRIO instruments.

The results are expressed in ng/ml calculated in reference to BCR-457 reference material.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.

Waste disposal: serum samples, calibrators and strips, once used, must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**

2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged. Do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.

The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website (<http://www.diesse.it/en/Support/Download/strumento:39/>)


5. Check that the Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: STREPTAVIDIN-POD

Contents: Streptavidin conjugated to peroxidase, in phosphate buffer solution containing phenol 0.05% and Bronidox 0.02%.

Position 6: MICROPLATE WELL

Coated with anti-Tg monoclonal antibodies

Position 5: empty

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine and H₂O₂ stabilized in buffer

Position 3: STOP SOLUTION

Contents: sulphuric acid solution 0.3 M

Position 2: CONJUGATE

Contents: anti-Tg monoclonal antibodies labeled with Biotin, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

in which undiluted serum must be added

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.500 mL

Contents: Buffer containing Tg and preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL

1 x 1 mL

Contents: Buffer containing Tg and preservative. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

Fresh serum may be stored for 1 day at 2/8°C.

Heat-inactivation can rise to erroneous results.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.

- Dispense 130 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
- Place the devices in the Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the instrument Operating Manual.

9. TEST VALIDATION

Use the positive control serum to check the validity of the results obtained. It should be used as reported in the instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus TRIO instrument expresses the result in ng/ml (BCR-457 reference material), calculated on the basis of a lot-dependent graph stored in the instrument.

The test on the examined serum can be interpreted as follows:

Euthyroid	>3
Patients who underwent total thyroidectomy with suppressed TSH	<0.29 (functional sensitivity)

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. CALIBRATION RANGE

Calibration range 0.10-5.1 ng/ml

For samples > 5.1 ng/ml retest the diluted sample in the Negative Control/Sample Diluent (PF83607-not supplied with the kit).

13. REFERENCE RANGE

Among the normal population the expected values are between 3 and 60 ng/ml. Considering the type of test the expected values in the normal population are >3 ng/ml.

14. ANALYTICAL SPECIFICITY

A study was conducted using the critical reagents of the kit. 5 samples (3 Negative and 2 Positive) were spiked with the following potentially interfering factors and then tested:

Rheumatoid factor (1.25 IU/ml - 20 IU/ml)
Bilirubin (5 mg/dl - 100 mg/dl)
Triglycerides (50 mg/dl - 3000 mg/dl)
Hemoglobin (0.6 mg/dl - 10 mg/dl)

The presence in the serum sample of the interfering substances described above does not affect the test result.

15. METHOD COMPARISON

In an experimentation 104 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table :

Correlation	r	95%CI
Pearson	0.93	0.90-0.95
Spearman	0.87	0.82-0.91

The correlation between the two methods is very high (Pearson) and high (Spearman).

16. PRECISION AND REPEATABILITY

Sample	Within run		Between run	
	Mean (ng/ml)	CV%	Mean (ng/ml)	CV%
1	3.6	2.2	4.9	11.0
2	1.6	5.0	2.7	7.4
3	2.7	4.1	1.5	6.7
4	0.24	12.5	0.85	12.9

Sample	Between lots		Between Instruments	
	Mean (ng/ml)	CV%	Mean (ng/ml)	CV%
1	4.0	3.8	4.0	1.5
2	2.1	8.1	2.1	4.8
3	1.1	9.1	1.1	5.5
4	0.58	6.9	0.57	7.0

17. REFERENCES

- Indrasena B.S.H. Use of thyroglobulin as a tumor marker. World J Biol Chem 2017 February 26; 8(1); 8(1): 81-85.
- Spencer C. A. Clinical Utility of Thyroglobulin Antibody (TgAb) Measurement for Patients with Differentiated Thyroid Cancers (DTC). J Clin Endocrinol Metab, December 2011, 69(12):3615-3627.
- Giovannella L., M.Clark P., Chiovato L., Duntas L., Elisei R., Feldt-Rasmussen U., Leenhardt L., Luster M., Schalin-Jantti C., Schott M., Seregini E., Rimmele H., Smith J., and A Verburg F. DIAGNOSIS OF ENDOCRINE DISEASE. Thyroglobulin measurement using highly sensitive assay in patients with differentiated thyroid cancer: a clinical position paper. European Journal of Endocrinology (2014) 171, R33-R46.
- Spencer C., LoPresti J., and Fatemi S. How sensitive (second-generation) thyroglobulin measurement is changing paradigms for monitoring patients with differentiated thyroid cancer, in the absence or presence of thyroglobulin autoantibodies. Curr Opin Endocrinol Diabetes Obes 2014, 21:394-404.
- Gunter E. Analyte stability and freeze-thaw information Specimen solution LCC.

DIESE Diagnostica Senese
S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

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1. Indrasena B.S.H. Use of thyroglobulin as a tumor marker. World J Biol Chem 2017 February 26; 8(1) ; 8(1): 81-85.
2. Spencer C. A. Clinical Utility of Thyroglobulin Antibody (TgAb) Measurement for Patients with Differentiated Thyroid Cancers (DTC). J Clin Endocrinol Metab, December 2011, 69(12):3615-3627.
3. Giovannella L., M.Clark P., Chiovato L.,Duntas L., Elisei R., Feldt-Rasmussen U., Leenhardt L., Luster M., Schalin-Jantti C., Schott M., Seregin E., Rimmele H., Smith J., and A Verburg F. DIAGNOSIS OF ENDOCRINE DISEASE. Thyroglobulin measurement using highly sensitive assay in patients with differentiated thyroid cancer: a clinical position paper. European Journal of Endocrinology (2014) 171, R33-R46.
4. Spencer C., LoPresti J., and Fatemi S. How sensitive (second-generation) thyroglobulin measurement is changing paradigms for monitoring patients with differentiated thyroid cancer, in the absence or presence of thyroglobulin autoantibodies. Curr Opin Endocrinol Diabetes Obes 2014, 21:394-404.
5. Gunter E. Analyte stability and freeze-thaw information Specimen solution LCC.

DIESSE Diagnostica Senese
S.p.A.
 Via delle Rose 10
 53035 Monteriggioni (Siena)
 Italy

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F.DIAGNOSIS OF ENDOCRINE DISEASE.
Thyroglobulin measurement using highly sensitive assay in patients with differentiated thyroid cancer: a clinical position paper. European Journal of Endocrinology (2014) 171, R33-R46.

4. Spencer C., LoPresti J., and Fatemi S. How sensitive (second-generation) thyroglobulin measurement is changing paradigms for monitoring patients with differentiated thyroid cancer, in the absence or presence of thyroglobulin autoantibodies. Curr Opin Endocrinol Diabetes Obes 2014, 21:394-404.
5. Gunter E. Analyte stability and freeze-thaw information Specimen solution LCC.

DIESSE Diagnostica Senese
S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

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17. BIBLIOGRAFIE

1. Indrasena B.S.H. Use of thyroglobulin as a tumor marker. World J Biol Chem 2017 February 26; 8(1) ; 8(1): 81-85.
2. Spencer C. A. Clinical Utility of Thyroglobulin Antibody (TgAb) Measurement for Patients with Differentiated Thyroid Cancers (DTC). J Clin Endocrinol Metab, December 2011, 69(12):3615-3627.
3. Giovannella L., M.Clark P., Chiovato L., Duntas L., Elisei R., Feldt-Rasmussen U., Leenhardt L., Luster M., Schalin-Jantti C., Schott M., Seregini E., Rimmele H., Smith J., and A Verburg F. DIAGNOSIS OF ENDOCRINE DISEASE. Thyroglobulin measurement using highly sensitive assay in patients with differentiated thyroid cancer: a clinical position paper. European Journal of Endocrinology (2014) 171, R33-R46.
4. Spencer C., LoPresti J., and Fatemi S. How sensitive (second-generation) thyroglobulin measurement is changing paradigms for monitoring patients with differentiated thyroid cancer, in the absence or presence of thyroglobulin autoantibodies. Curr Opin Endocrinol Diabetes Obes 2014, 21:394-404.

5. Gunter E. Analyte stability and freeze-thaw information
Specimen solution LCC.

DIESSE Diagnostica Senese
S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

	EN Date of manufacture ES Fecha de fabricación IT Data di fabbricazione	FR Date de fabrication GR Ημερομηνία Παραγωγής PT Data de fabrico
	EN Use By ES Fecha de caducidad IT Utilizzare entro	FR Utiliser jusqu'à GR Ημερομηνία λήξης PT Prazo de validade
	EN Do not reuse ES No reutilizar IT Non riutilizzare	FR Ne pas réutiliser GR Μην κάνετε επαναληπτική χρήση PT Não reutilizar
	EN Caution, consult accompanying documents ES Atención, ver instrucciones de uso IT Attenzione, vedere le istruzioni per l'uso	FR Attention voir notice d'instructions GR Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα PT Atenção, consulte a documentação incluída
	EN Manufacturer ES Fabricante IT Fabbicante	FR Fabricant GR Κατασκευαστής PT Fabricante
	EN Contains sufficient for <n> tests ES Contenido suficiente para <n> ensayos IT Contenuto sufficiente per "n" saggi	FR Contenu suffisant pour "n" tests GR Περιεχόμενο επαρκές για «n» εξετάσεις PT Conteúdo suficiente para "n" ensaios
	EN Temperature limitation ES Límite de temperatura IT Limiti di temperatura	FR Limites de température GR Περιορισμοί θερμοκρασίας PT Limites de temperatura
	EN Consult Instructions for Use ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso	FR Consulter les instructions d'utilisation GR Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização
	EN Biological risks ES Riesgo biológico IT Rischio biologico	FR Risques biologiques GR Βιολογικοί κίνδυνοι PT Risco biológico
	EN Catalogue number ES Número de catálogo IT Numero di catalogo	FR Référence du catalogue GR Αριθμός καταλόγου PT Referência de catálogo
	EN In Vitro Diagnostic Medical Device ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro	FR Dispositif médical de diagnostic in vitro GR In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro
	EN Batch code ES Código de lote IT Codice del lotto	FR Code du lot GR Αριθμός Παρτίδας PT Código do lote