CHORUS

DENGUE IgG

DIESSE
REF 81154

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16. BIBLIOGRAFIA

1. World Health Organization. (1997). Dengue Haemorrhagic Fever: diagnosis, treatment, prevention and control, p34-36. 2nd edition. World Health Organisation Office of Publications Geneva, Switzerland.
2. Guzman MG, Halstead SB, Artsob H, et al., Dengue: a continuing global threat, in Nat. Rev. Microbiol., vol. 8, 12 Suppl, dicembre 2010
3. Chatchen S, Sabchareon A, Sirivichayakul C. Serodiagnosis of asymptomatic dengue infection. Asian Pac J Trop Med. 2017 Jan;10(1):11-14.

INSTRUCTIONS FOR USE

CHORUS DENGUE IgG

**For the qualitative determination of anti-Dengue
IgG antibodies**

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the qualitative determination of IgG class antibodies against Dengue in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Dengue fever (DF) is caused by four antigenically different serotypes of the Dengue virus (Denv-1, Denv-2, Denv-3 e Denv-4). DF is a mosquito-borne viral disease which causes over 50-100 million infections every year in the tropical and subtropical regions of the world.

The infection is associated with high fever, intense headache, myalgia, arthralgia and rash. Severe complications of Dengue infection include haemorrhagic fever, especially during the second infection.

Primary infection (caused by any virus serotype) is characterized by a rise in virus specific IgM antibodies, 4-5 days after the appearance of symptoms. While these are detectable in serum only for 30-90 days, the IgG antibodies persist for a long time, often for all the life.

During a secondary infection, the level of Dengue virus specific IgM antibodies in serum is generally low but IgG levels rise rapidly to much higher levels than that observed during the primary dengue infection.

Detection of IgG can be used to confirm the diagnosis also in the late phases of the infection. Detection of IgM, in a symptomatic person, is considered diagnostic.

3. PRINCIPLE OF THE METHOD

The Chorus Dengue IgG device is ready to use for the detection of IgG antibodies against Dengue, in the Chorus/Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human serum.

After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgG monoclonal antibodies conjugated to horseradish peroxidase. The unbound conjugate is eliminated, and the peroxidase substrate added. The blue colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.

The results are expressed in Index (OD sample/OD cut-off)

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.

The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the

table published on Diesse website
(<http://www.diesse.it/en/Support/Download/strumento:39/>)


5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with Dengue (Type 2) antigen

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing phenol 0.05%, Bronidox 0.02% and indicator to reveal the presence of the serum.

Position 2: CONJUGATE

Contents: anti-human IgG monoclonal antibodies labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

In which the operator must place the undiluted serum

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.175 ml

Contents: Diluted human serum containing anti-Dengue IgG antibodies and preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL

1 x 0.425 ml

Contents: Diluted human serum containing anti-Dengue IgG antibodies and preservative. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration

must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in Index (OD sample/OD cut-off).

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is > 1.1
NEGATIVE: when the result is < 0.9
DOUBTFUL/EQUIVOCAL: for all values between 0.9 and 1.1

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. ANALYTICAL SPECIFICITY

4 samples (2 Negative and 2 Positive) were spiked with the following potentially interfering factors and then tested :

Rheumatoid Factor (44 – 220 IU/ml)
Bilirubin (4.5 mg/dl – 45 mg/dl)
Triglycerides (10 mg/dl – 250 mg/dl)
Hemoglobin (5 mg/ml – 30 mg/ml)

The presence in the serum sample of the interfering substances described above does not affect the test result.

13. CROSS-REACTIONS

25 samples, positive to Zika, Parvovirus IgG, Herpes Simplex virus, Influenza, Epstein-Barr virus and Varicella were tested. No significant cross-reactions were found except for Zika virus.

14. METHOD COMPARISON

In an experimentation 159 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	34	0	34
	-	1	124	125
	Total	35	124	159

Percent Positive Agreement (~Diagnostic Sensitivity):

97.1 % CI_{95%}: 85.5 - 99.4

Percent Negative Agreement: (~Diagnostic Specificity):

100.0 % CI_{95%}: 97.0 - 100.0

The agreement between the two methods is excellent with a Cohen's Kappa of 0.97.

15. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision	
	Mean (Index)	CV%	Mean (Index)	CV%
1	0.1	50.0*	0.2	30.0*
2	0.3	-	0.3	-
3	1.0	8.0	1.0	5.0
4	4.1	8.8	4.0	4.3
5	3.6	4.2	3.7	4.6

* Artifact caused by the known fault of Variation Coefficient which becomes extremely sensitive to even very small changes in the mean when the mean value is near zero.

Sample	Precision between batches	
	Mean (Index)	CV%
1	0.2	-
2	0.4	15.0
3	1.2	5.0
4	4.5	4.7
5	4.0	11.8

16. REFERENCES

1. World Health Organization. (1997). Dengue Haemorrhagic Fever: diagnosis, treatment, prevention and control, p34-36. 2nd edition. World Health Organisation Office of Publications Geneva, Switzerland.
2. Guzman MG, Halstead SB, Artsob H, et al., Dengue: a continuing global threat, in Nat. Rev. Microbiol., vol. 8, 12 Suppl, dicembre 2010
3. Chatchen S, Sabchareon A, Sirivichayakul C. Serodiagnosis of asymptomatic dengue infection. Asian Pac J Trop Med. 2017 Jan;10(1):11-14.

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