**ENZYWELL**

**HERPES SIMPLEX 1+2 IgM**

**REF 91021 (96 tests)**

Prodotto da/Manufactured by/ Hersteller:  
 Κατασκευάζεται από/ Fabricado por/ Fabriqué par:  
 DIESSE Diagnostica Senese  
 Strada dei Laghi 39  
 53035 Monteriggioni (Siena) - Italy

|                                                                                                                                                                                                                                                                                                                                                      | Capitolo<br>Section<br>Capítulo |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------|
| Modifiche introdotte nella revisione corrente<br>Changes introduced in the current revision<br>Eingeführte Änderungen in der vorliegenden Fassung<br>Τροποποιήσεις που εισήχθησαν στην τρέχουσα διόρθωση<br>Cambios introducidos en la revisión actual<br>Modifications apportées à la révision courante<br>Alterações introduzidas na revisão atual | 4                               |

**INDICE / INDEX / INHALT/ ΠΕΡΙΕΧΟΜΕΝΑ / INDICE / SOMMAIRE / CONTEÚDO**

1. UTILIZZAZIONE / INTENDED USE / VERWENDUNGSZWECK/ ΥΠΟΔΕΙΚΝΥΟΜΕΝΗ ΧΡΗΣΗ / INDICACIONES DE USO/ UTILISATION/UTILIZAÇÃO PRETENDIDA
2. INTRODUZIONE / SUMMARY AND EXPLANATION OF TEST / PHYSIOLOGIE/ ΕΙΣΑΓΩΓΗ / RESUMEN Y EXPLICACIÓN DEL TEST / INTERET CLINIQUE / SUMÁRIO E EXPLICAÇÃO DO TESTE
3. PRINCIPIO DEL METODO / PRINCIPLE OF THE TEST / TESTPRINZIP/ ΑΡΧΗ ΤΗΣ ΜΕΘΟΔΟΥ / PRINCIPIO DEL MÉTODO / PRINCIPE/ PRINCÍPIOS DO TESTE
4. COMPOSIZIONE DEL KIT E PREPARAZIONE DEI REAGENTI / KIT CONTENTS AND REAGENT PREPARATION / INHALT DES KITS UND VORBEREITUNG DER REAGENZIEN / ΣΥΝΘΕΣΗ ΤΟΥ ΚΙΤ ΚΑΙ ΠΡΟΕΤΟΙΜΑΣΙΑ ΤΩΝ ΑΝΤΙΑΠΑΣΤΗΡΙΩΝ / COMPONENTES DEL KIT Y PREPARACIÓN DEL REACTIVO / CONTENU DU KIT ET PREPARATION DES REACTIFS / CONTEÚDA DO KIT E PREPARAÇÃO DO REAGENTE /
5. MODALITA' DI CONSERVAZIONE E STABILITA' DEI REAGENTI / STORAGE AND STABILITY OF REAGENTS / LAGERUNG UND STABILITÄT DER REAGENZIEN / ΤΡΟΠΟΣ ΑΠΟΘΗΚΕΥΣΗΣ ΚΑΙ ΣΤΑΘΕΡΟΤΗΤΑ ΤΩΝ ΑΝΤΙΑΠΑΣΤΗΡΙΩΝ/ CONSERVACIÓN Y ESTABILIDAD DE LOS REACTIVOS / MODALITES DE CONSERVATION ET STABILITE DES REACTIFS / ARMAZENAMENTO E ESTABILIDADE
6. PRECAUZIONI / PRECAUTIONS / HINWEISE UND VORSICHTSMASSNAHMEN/ ΠΡΟΦΥΛΑΞΕΙΣ / PRECAUCIONES DE USO / PRECAUTIONS D'UTILISATION / PRECAUÇÕES
7. TIPO DI CAMPIONE E CONSERVAZIONE / TYPE AND STORAGE OF SAMPLE / ART UND LAGERUNG DER PROBEN / ΕΙΔΟΣ ΔΕΙΓΜΑΤΟΣ ΚΑΙ ΑΠΟΘΗΚΕΥΣΗ / TIPO DE MUESTRA Y CONSERVACION / TYPE D'ECHANTILLON ET CONSERVATION / TIPO E ARMAZENAMENTO DAS AMOSTRAS
8. PROCEDIMENTO / TEST PROCEDURE / TESTDURCHFÜHRUNG / ΔΙΑΔΙΚΑΣΙΑ / PROCEDIMIENTO DEL TEST / PROCEDURE / PROCEDIMENTO DE TESTE /
9. SCHEMA DEL SAGGIO / SCHEME OF TEST PROCEDURE / TESTSCHEMA / ΔΙΑΓΡΑΜΜΑ ΤΗΣ ΑΝΑΛΥΣΗΣ/ ESQUEMA DEL PROCEDIMIENTO DEL TEST / SCHEMA D'ANALYSE/ ESQUEMA DE PROCEDIMENTO DE TESTE /
10. VALIDAZIONE DEL TEST / TEST VALIDATION / TESTVALIDIERUNG / ΕΓΚΥΡΟΤΗΤΑ ΤΟΥ ΤΕΣΤ / VALIDACIÓN DEL TEST / VALIDATION DU TEST / VALIDAÇÃO DO TESTE /
11. INTERPRETAZIONE DEL TEST / INTERPRETATION OF RESULTS / INTERPRETATION DER ERGEBNISSE / ΕΡΜΗΝΕΙΑ ΤΩΝ ΑΠΟΤΕΛΕΣΜΑΤΩΝ / INTERPRETACIÓN DE RESULTADOS / INTERPRETATION DU TEST / INTERPRETAÇÃO DOS RESULTADOS
12. LIMITAZIONI DELLA PROCEDURA / LIMITATIONS OF THE PROCEDURE / GRENZEN DES VERFAHRENS / ΠΕΡΙΟΡΙΣΜΟΙ ΣΤΗΝ ΜΕΘΟΔΟ / LIMITACIONES DEL TEST / LIMITES DU TEST/ LIMITAÇÕES DO PROCEDIMENTO /
13. SPECIFICITA' ANALITICA / ANALYTICAL SPECIFICITY / ANALYTISCHE SPEZIFITÄT / ΑΝΑΛΥΤΙΚΗ ΕΙΔΙΚΟΤΗΤΑ/ ESPECIFICIDAD ANALÍTICA / SPECIFICITE ANALYTIQUE / ESPECIFICIDADE ANALÍTICA /
14. SENSIBILITA' E SPECIFICITA' DIAGNOSTICA / DIAGNOSTIC SENSITIVITY AND SPECIFICITY / DIAGNOSTISCHE SENSITIVITÄT UND SPEZIFITÄT / ΕΥΑΙΣΘΗΣΙΑ ΚΑΙ ΔΙΑΓΝΩΣΤΙΚΗ ΕΙΔΙΚΟΤΗΤΑ / SENSIBILIDAD Y ESPECIFICIDAD DIAGNOSTICAS / SENSIBLITE ET SPECIFICITE DIAGNOSTIQUES / ESPECIFICIDADE E SENSIBILIDADE DO DIAGNÓSTICO
15. PRECISIONE / PRECISION / PRÄZISION / ΑΚΡΙΒΕΙΑ / PRECISIÓN / PRECISION / PRECISÃO
16. GUIDA DEI PROBLEMI DI UTILIZZO / "TROUBLE SHOOTING" / "TROUBLE SHOOTING" (PROBLEMLÖSUNGEN) / ΟΔΗΓΟΣ ΕΠΙΛΥΣΗΣ ΠΡΟΒΛΗΜΑΤΩΝ / GUIA DE RESOLUCION DE PROBLEMAS / RESOLUTION DES PROBLEMES / RESOLUÇÃO DE PROBLEMAS
17. BIBLIOGRAFIA / REFERENCES / LITERATUR / ΒΙΒΛΙΟΓΡΑΦΙΑ / BIBLIOGRAFÍA / BIBLIOGRAPHIE / REFERÊNCIAS

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### **ENZYWELL HERPES SIMPLEX 1+2 IgM**

**REF** 91021

**(Italiano)**

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**17. BIBLIOGRAFIA**

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
3. B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
4. C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
5. M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
6. D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgM antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).
7. R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).
8. J.E. Kuhn et al.: Analysis of the IgM and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).

## INSTRUCTIONS FOR USE

### **ENZYWELL HERPES SIMPLEX 1+2 IgM**

**REF** 91021

(English)

#### **1. INTENDED USE**

**IMMUNOENZYMATIC CAPTURE METHOD FOR THE QUALITATIVE DETERMINATION OF IgM-CLASS ANTIBODIES TO HERPES SIMPLEX VIRUS (TYPES 1 AND 2) IN HUMAN SERUM, AS AN AID IN THE DIAGNOSIS OF HERPES SIMPLEX VIRUS INFECTION.**

#### **2. SUMMARY AND EXPLANATION OF THE TEST**

The Herpes simplex virus (HSV) is a member of the Herpesviridae family, of which two types are known: type 1 (HSV-1) and type 2 (HSV-2) which present slight antigenic differences. HSV-1 causes chiefly oral-facial lesions, while HSV-2 is mainly responsible for genital lesions, but this distinction is not binding, both types occasionally causing infection in either anatomical site. HSV may also cause a form of ocular cheratitis, and lesions of the central nervous system.

HSV can affect practically the whole population. The primary infection is often in a subclinical form and is rarely diagnosed. After a latency period of variable duration, reactivation may occur and viral replication may or may not give rise to clinical lesions. Infection contracted during birth is of particular interest, this being an important cause of morbidity and mortality. It is therefore important to determine the immunitary state of women during pregnancy in order to detect serum conversion. The assay of specific IgM is important for the diagnosis of neonatal infection and encephalitis caused by HSV. Moreover, the presence of specific IgM indicates viral activity in progress, although it is not possible to distinguish between primary infection and reactivation.

#### **3. PRINCIPLE OF THE TEST**

The test for the assay of HSV IgM is based on the principle of the capture of these immunoglobulins and the subsequent identification of those which are specific, making use of their ability to bind an antigen conjugated to peroxidase (1-8). The capture is performed using monoclonal antibodies bound to the solid phase (microtiter wells). The antigen is composed of purified, inactivated HSV type 1 and type 2. The conjugate is composed of specific anti-HSV monoclonal antibodies labelled with peroxidase.

#### **4. KIT CONTENTS AND REAGENT PREPARATION**

- Reagents are sufficient for 96 determinations.

**Bring to room temperature before use.**

**MT PLATE** MICROPLATE 12x8 wells coated with anti-human IgM monoclonal antibodies.

Use: open the package at the opposite end from the code (M followed by the lot number) which is useful for identification purposes, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

**CONTROL +** POSITIVE CONTROL 1 x 1.6 mL

Contents: Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour is proportional to the relative antibody titer.

**CONTROL CUT-OFF** CUT OFF CONTROL 1 x 2.5 mL

Contents: Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour is proportional to the relative antibody titer.

**Ag** ANTIGEN Freeze-dried powder x 6 vials

Contents: Purified Herpes Simplex Virus, in Phosphate buffer with lactose.

Preparation: reconstitute with the conjugate volume shown on the label, mixing by inversion.

**CONJ** CONJUGATE 1 x 18 mL

Contents: monoclonal antibodies labelled with peroxidase, in phosphate buffer with phenol 0.05% and Bronidox 0.02%.

Preparation: ready for use.

The immunocomplex should be prepared about 45 min. before use.

**CONTROL IgM -** IgM NEGATIVE CONTROL (PF93900) 1 x 1.6 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

**WASH BUF 10x** WASH BUFFER 10X (PF93603) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5%.

Preparation: dilute the required volume 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

**SAMP DIL 2** DILUENT 2 (PF93611) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

To be used to dilute samples.

Contents: Proteic solution in phosphate buffer with sodium azide 0,09% containing methyl orange as dye.

**SUBS TMB** SUBSTRATE (PF93619) 1x 12 mL Ready for use **INTERCHANGEABLE BETWEEN LOTS**

Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilised in citrate buffer 0.05 mol/L (pH 3.8).

**H<sub>2</sub>SO<sub>4</sub> 0.3 M** STOP SOLUTION (PF93602) 1x16 mL **INTERCHANGEABLE BETWEEN LOTS**

H<sub>2</sub>SO<sub>4</sub> 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

#### **MATERIALS REQUIRED BUT NOT PROVIDED**

- Incubator at 37°C
- Microplate reader (wavelength 450 or 450/620 nm, with linearity up to OD >= 2000)
- Microplate washer (preferable) able to dispense volumes in the range 225-375 µl
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

#### **5. STORAGE AND STABILITY OF REAGENTS**

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

#### **Reagents have a limited stability after opening and/or preparation**

| REAGENT               | CONDITIONS                                                           |
|-----------------------|----------------------------------------------------------------------|
| Microplate            | 5 weeks at 2/8°C, polythene bag                                      |
| Controls              | 5 weeks at 2/8°C                                                     |
| Conjugate             | 5 weeks at 2/8°C                                                     |
| Reconstituted antigen | must be used the same day; cannot be frozen                          |
| Substrate             | up to the expiry date at 2/8°C, 1 week at 15-30°C; store in the dark |
| Sample Diluent        | up to the expiry date at 2/8°C                                       |
| Wash Buffer           | 2 weeks at 2/8°C, 5 days at 15/30°C.                                 |
| Stop Solution         | up to the expiry date at 2/8°C                                       |

## 6. PRECAUTIONS

### FOR IN VITRO DIAGNOSTIC USE ONLY.

#### Caution:

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

**Waste disposal:** serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

#### Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
  - a) The Wash Buffer contains detergents
  - b) The conjugate contains phenol
  - c) The substrate is acid
  - d) The controls contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

#### Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
6. Do not allow wells to become dry during the assay procedure.
7. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
8. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
9. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO<sub>2</sub> incubators must not be used.
10. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
11. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.

12. Use of the kit with automatic instruments must be validated by the user.
13. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:
  - installation and particular requisites
  - operating principles, instructions, precautions and risks
  - manufacturer's specifications and instrument performance
  - servicing and maintenance.

## **7. TYPE AND STORAGE OF SAMPLE**

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be carefully shaken before the test. Heat inactivation can lead to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric or contaminated samples should be avoided. **The test cannot be performed on human plasma.**

## **8. TEST PROCEDURE**

### **Manual Technique**

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Wash Buffer 10x (100 mL + 900 mL H<sub>2</sub>O).
- Prepare the immunocomplex by reconstituting the freeze-dried powder with the conjugate (volume shown on the label).

Dilute samples 1:101 distributing 10 µL of serum into 1 mL of diluent. Leave one well for the blank, performed using 100 µL of the substrate mixture. Dispense 100 µL of each diluted sample per well (duplicate testing is recommended). Place UNDILUTED controls in a strip (100 µL in each well). The minimum requisite is 1 negative control, 2 cut-off and 1 positive control.

Wells are covered with protective film and incubated for 45 minutes at 37°C. After washing four times for 30 seconds (300 µl), add 100 µL of the immunocomplex (Antigen/monoclonal antibodies anti-HSV labelled with HRP) to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film. The plate is washed again 4 times, as described above. Finally, the substrate is distributed, 100 µL/well.

After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µL of Stop Solution.

The adsorbance (O.D.) is read at 450 nm or at 450/620 nm within 30 min.

|                                    |
|------------------------------------|
| <b>9. SCHEME OF TEST PROCEDURE</b> |
|------------------------------------|

### **Manual Technique**

STEP 1 Place 100 µL of diluted sample//controls to the wells of the strips.

-  
Incubate for 45 min. at 37°C

-  
Wash 4 times (300 µL)

STEP 2 Add 100 µL of immunocomplex to each well

-  
Incubate for 45 min. at 37°C

-  
Wash 4 times (300 µL)

STEP 3 Add 100 µL of Substrate to each well

-  
Incubate for 15 min. at R.T.

STEP 4 Add 100 µL of Stop Solution

-  
Read absorbance at 450 nm within 30 min

**10. TEST VALIDATION**

Subtract the value of the blank ( $\leq 0.150$ ) from all the other readings. The O.D. values of the control Cut-off serum must be within 25% of the mean value if tested in triplicate. Disregard any abnormal value and recalculate the mean. The Positive control must have an O.D. at least 1.5 times that of the Cut-Off serum. The ratio between Negative Control and Cut-off must be  $\leq 0.6$ . The OD cut-off must be  $\geq 0.2$  at 450 nm and  $\geq 0.16$  at 450/620 nm.

**11. INTERPRETATION OF THE RESULTS****Qualitative results**

If the adsorbance of the sample is higher than that of the Cut-Off, the sample is positive for the presence of specific IgM.

Calculate the ratio between the O.D. value of the sample and that of the Cut-off (INDEX). The sample is considered:

Positive: if the ratio is  $> 1.2$ .

Doubtful:  $\pm 20\%$  of the Cut-Off.

Negative: if the ratio is  $< 0.8$ .

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

**12. LIMITATIONS OF THE PROCEDURE**

As with many other serological tests, the results obtained serve only as an aid to diagnosis and should always be interpreted together with other findings.

Paired specimens are required for determining seroconversion. Collect the first specimen as soon as possible after onset of symptoms, and the second specimen 2-3 weeks later. If the first sample is taken too late during the course of infection, seroconversion may not be detected. If the test results obtained for paired samples are inconclusive, reassay the specimens. Results must always be interpreted together with an evaluation of the clinical picture and results of other diagnostic investigations.

**13. ANALYTICAL SPECIFICITY**

67 serum samples containing potential interfering substances were tested:

- Antinuclear antibodies (ANA) (n=12)
- Rheumatoid factor (n=8)
- Heterophile antibodies (n=5)
- Bilirubin (n=9)
- Triglycerides (n=10).
- Cytomegalovirus IgM (n=23)

No interference was noted in any case.

**14. DIAGNOSTIC SENSITIVITY AND SPECIFICITY**

In a clinical trial performed in a hospital laboratory, 313 samples were analysed with the Diesse kit in comparison with another commercial method. All discordant samples were analysed with a third commercial method. The results are reported in the following table.

|        |   | REFERENCE METHOD |     |
|--------|---|------------------|-----|
|        |   | +                | -   |
| DIESSE | + | 43               | 6   |
|        | - | 2                | 262 |

The Enzywell kit has a sensitivity of 95.6% and a specificity of 97.8%.

**15. PRECISION****"In run" Precision**

| Sample         | HSM 1<br>(Negative < Cut Off) | HSM 2<br>(Positive > Cut Off) | HSM 3<br>(Positive) | Cut Off | Positive Control |
|----------------|-------------------------------|-------------------------------|---------------------|---------|------------------|
| n (replicates) | 23                            | 24                            | 24                  | 12      | 12               |
| O.D.           | 0.202                         | 0.643                         | 0.816               | 0.638   | 1.404            |
| CV%            | 23%                           | 14%                           | 8%                  | 15%     | 6%               |

**“Between run” Precision**

|                 | INDEX |     |
|-----------------|-------|-----|
|                 | Media | CV% |
| <b>Pos.Con.</b> | 5.8   | 6   |
| <b>HSM1</b>     | 0.4   | 14  |
| <b>HSM2</b>     | 1.2   | 13  |
| <b>HSM3</b>     | 2.0   | 11  |

**Precision between lots**

| Sample      | INDEX      |            |            | Average | CV% |
|-------------|------------|------------|------------|---------|-----|
|             | Lot n. 128 | Lot n. 129 | Lot n. 130 |         |     |
| <b>HSM1</b> | 0.4        | 0.4        | 0.3        | 0.4     | 16  |
| <b>HSM2</b> | 1.3        | 1.3        | 1.1        | 1.2     | 9   |
| <b>HSM3</b> | 2.2        | 1.9        | 1.7        | 1.9     | 13  |

**16. TROUBLE SHOOTING GUIDE**

| PROBLEM                      | POSSIBLE SOURCE                                           | TEST OR ACTION                                                                                        |
|------------------------------|-----------------------------------------------------------|-------------------------------------------------------------------------------------------------------|
| Invalid run (all negative)   | One or more reagents not added or added in wrong sequence | Recheck procedure<br>Check for unused solutions. Repeat test.                                         |
|                              | Unreactive plate                                          | Check the code on the package containing the plate (see package insert paragraph 4 for correct code). |
|                              |                                                           | Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow). Repeat test          |
| Invalid run (all positive)   | Contamination of substrate                                | Take new aliquot of substrate.                                                                        |
|                              | Inadequate washing                                        | Ensure that wash apparatus works well                                                                 |
| Poor precision               | Incomplete washing of wells                               | Ensure that wash apparatus works well                                                                 |
|                              | Inadequate aspiration of wells                            | Ensure that wash apparatus works well                                                                 |
|                              | Pipetting error                                           | Check pipette function                                                                                |
|                              | Reagent addition too slow                                 | Avoid drying of the plate after washing step.<br>Add reagents immediately                             |
|                              | Presence of bubbles                                       | Avoid air bubbles during pipetting.                                                                   |
|                              | Optical pathway not clean                                 | Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.           |
| Inadequate Color development | Incorrect incubation times or temperature                 | Check for temperature control and time monitoring                                                     |
|                              |                                                           | Adhere to recommended instruction for use.                                                            |
|                              | Inadequate volume of substrate added to the plate         | Check pipette function.                                                                               |

**17. REFERENCES**

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
3. B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
4. C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
5. M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
6. D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgM antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).

7. R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. *J. Med. Virology* 16: 1247 (1985).
8. J.E. Kuhn et al.: Analysis of the IgM and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. *J. Medical Virology* 23: 135 (1987).

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**17. LITERATUR**

9. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
10. S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
11. B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
12. C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
13. M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
14. D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgM antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).
15. R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).
16. J.E. Kuhn et al.: Analysis of the IgM and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).

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- 1.G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
- 2.S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
- 3.B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
- 4.C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
- M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
- 5.D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgM antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).
- 6.R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).
- 7.J.E. Kuhn et al.: Analysis of the IgM and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).

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