ENZY-WELL

DIESSE

SARS-CoV-2 IgM

DIESSE Diagnostica Senese
S.p.A.
Via delle Rose, 10
53035 Monteriggioni (Siena)
Italy

REF 91401 (96 tests)

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ENZY-WELL SARS-CoV-2 IgM

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18. BIBLIOGRAFIA

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
2. Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
5. Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

DIESSE Diagnostica Senese S.p.A.
 Via delle Rose 10
 53035 Monteriggioni (Siena)
 Italy

INSTRUCTIONS FOR USE

ENZY-WELL SARS-CoV-2 IgM

For the qualitative determination of IgM class antibodies to SARS-CoV-2 virus

For In Vitro Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the qualitative determination of IgM class antibodies to SARS-CoV-2 virus in human serum.

2. INTRODUCTION

Coronaviruses are single positive-stranded RNA viruses belonging to *Coronaviridae* family, which are divided into four genera: *Alpha-*, *Beta-*, *Delta-* and *Gammacoronavirus*.

CoVs are commonly found in many species of animals. Occasionally, the animal CoVs can acquire genetic mutations by errors during genome replication or recombination mechanism, which can further expand their tropism to humans. A total of six human CoV types were identified to be responsible for causing human respiratory ailments, which include two alpha CoVs and four beta CoVs (the most recent two are SARS-CoV and MERS-CoV). Typically, these CoVs cause asymptomatic infection or severe acute respiratory illness, including fever, cough, and shortness of breath. However, other symptoms such as gastroenteritis and neurological diseases of varying severity have also been reported. Since December 2019, the seventh novel human betacoronavirus (2019-nCoV) causing pneumonia outbreak was reported by the World Health Organization. This novel virus is linked with an outbreak of febrile respiratory illness in the city of Wuhan, Hubei Province, China. Some patients with 2019-nCoV have developed severe pneumonia, pulmonary oedema, ARDS, or multiple organ failure and have died. At present, information regarding the epidemiology and clinical features of pneumonia caused by 2019-nCoV is scarce. The monitoring of antibodies' titer can serve as an indication of the status of infection: during the first days, IgA and IgM appear; as the disease progresses, IgG becomes predominant.

3. PRINCIPLE OF THE TEST

The test is based on the ELISA technique (Enzyme Linked ImmunoSorbent Assay). The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with dilute human serum. After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgM monoclonal antibodies conjugated to horseradish peroxidase. The unbound conjugate is eliminated, and the peroxidase substrate added. The blue color which develops is proportional to the concentration of specific antibodies present in the serum sample. When the enzymatic reaction is interrupted by the addition of a sulphuric acid solution, a yellow coloring develops. The color, proportional to the quantity of specific antibodies present in the serum, can be easily read using a microplate reader.

4. WARNINGS AND PRECAUTIONS

This kit contains materials of human origin which have been tested and gave a negative response by methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens and performing the assay.
3. Wash hands thoroughly when finished.
4. Non-disposable apparatus should be sterilized after use, the preferred method is to autoclave for 1 h at 121°C. All the disposable apparatus should be eliminated according to law.
5. 2M hydrochloric acid, used for washing glassware, is corrosive and should be handled with appropriate care. In case of contact with skin or eyes, wash thoroughly with water.
6. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minutes exposure to 1% sodium hypochlorite may be sufficient to ensure effective decontamination.
7. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

1. The OD of Cut-Off, controls and samples can be slightly different among different plates. For such reason if during the same run strips from different plates are used, even if the lot is the same, it is necessary to repeat the determination of the Cut-Off.
2. Allow all reagents to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature.
3. Open the package containing the strips after 30 minutes at room temperature at least.
4. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
5. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.

6. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or deionized water.
7. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
8. Do not allow wells to become dry during the test.
9. Do not keep the samples in auto-defrosting freezers.
10. Care should be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the substrate and conjugate.
11. Care should be taken to avoid touching the rim of the well with conjugate.
12. Do not "blow-out" from microplates.
13. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at the temperature indicated in the procedure either in a water bath with a rack to support the plates, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate Operating Manual for further details. CO₂ incubators must not be used.
14. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
15. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
16. Care should be taken to avoid contaminating the microplate wells with disposable gloves powder.
17. The use of the kit with automated equipment has to be validated by the user.
18. For each instrument used, read the Operating Manual, in particular to obtain additional information on the following points:
 - installation and particular requisites
 - operating principles, instructions, precautions and risks
 - manufacturer's specifications and instrument performances
 - technical servicing and maintenance

5. KIT COMPOSITION AND REAGENT PREPARATION

Reagents are sufficient for 96 determinations.

- Bring to room temperature before use.

MT PLATE MICROPLATE 12x8 strip

Content: 1 microplate (96 wells) coated with inactivated native antigen.

CONJ CONJUGATE 1x16 mL

Content: anti-human IgM monoclonal antibodies, labeled with peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%. Liquid, ready for use.

CONTROL - NEGATIVE CONTROL 1x0.8ml

Content: Proteic solution not containing specific antibodies capable of binding the antigen present on the microplate, and preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1x0.8 mL

Content: Proteic solution containing specific antibodies capable of binding the antigen present on the microplate and preservative. Liquid, ready for use.

CONTROL CUT-OFF CUT-OFF CONTROL 1x0.8 mL

Content: Proteic solution containing specific antibodies capable of binding the antigen present on the microplate, and preservative. Liquid, ready for use.

SAMP DIL 2 DILUENT 2 1x100 mL

(PF93611)

INTERCHANGEABLE BETWEEN LOTS

Content: Proteic solution in phosphate buffer with Sodium Azide 0.09% and methyl orange as dye.

Use: To be used to dilute samples. Liquid, ready for use.

SAMPLE SORBENT 4 SAMPLE SORBENT 4 1x7 mL

(PF30401)

INTERCHANGEABLE BETWEEN LOTS

Content: Proteic solution in phosphate buffer with Sodium Azide 0.09%, methyl orange as dye and anti-human IgG antibodies.

Use: To be used to dilute samples (in the microwell). Liquid, ready for use.

WASH BUF 10x WASH BUFFER 10X 1x100 mL

(PF93603)

INTERCHANGEABLE BETWEEN LOTS

Content: Phosphate buffered saline (PBS), concentrated 10 times; contains Brij 0.5%.

Preparation: Dilute the required volume 1:10 with distilled or deionized water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SUBS TMB SUBSTRATE 1x12 mL

(PF93619)

INTERCHANGEABLE BETWEEN LOTS

Content: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in citrate buffer 0.05 mol/L (pH 3.8). Liquid, ready for use.

H₂SO₄ 0.3 M STOP SOLUTION 1x16 mL

(PF93602)

INTERCHANGEABLE BETWEEN LOTS

Content: H₂SO₄ 0.3 mol/L. Liquid, ready for use.

ADHESIVE FILMS (2).

MATERIALS REQUIRED BUT NOT PROVIDED

- Microplate reader (wavelength 450 or 450/620 nm, with linearity up to OD ≥ 2.000)
- Microplate washer (optional) able to dispense volumes in the range 225-375 µl
- Distilled or deionized water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl of solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

- Absorbent tissue

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2-8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation:

MICROPLATE	4 weeks at 2-8°C
CONJUGATE	8 weeks at 2-8°C
POSITIVE CONTROL	8 weeks at 2-8°C
NEGATIVE CONTROL	8 weeks at 2-8°C
CUT-OFF CONTROL	8 weeks at 2-8°C
DILUENT 2	up to the expiry date at 2-8°C
SAMPLE SORBENT 4	up to the expiry date at 2-8°C
SUBSTRATE	up to the expiry date at 2-8°C; 1 week at 15-30°C; store in the dark
WASHING BUFFER	2 weeks at 2-8°C; 5 days at 15-30°C
STOP SOLUTION	up to the expiry date at 2-8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C.

Avoid several freeze-thawing cycles.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

PREPARATION

Bring all the reagents to room temperature (18-30°C) before use.

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Wash Buffer 10x (100 mL + 900 mL H₂O).
- Dilute samples 1:51 distributing 10 µL of serum into 0.500 mL of diluent.

PROCEDURE

1. Distribution of the samples:
Dispense 100 µL of control (positive, negative or cut-off).
Dispense 50 µL of sample sorbent 4 + 50 µL diluted sample per well (duplicate testing is recommended).
The minimum requisite is 1 negative control, 1 positive control and 2 cut-off controls. Leave one well for the blank
2. Incubation:
Incubate the plate covered with an adhesive protective film for 60 ± 5 minutes at room temperature.

3. Washing:
Remove the adhesive film, aspirate the contents of all the wells and rinse 4 times, filling each well with 300 µL of washing solution. Wait 30 seconds between each wash.
4. Distribution of the conjugate:
Dispense 100 µL of conjugate in each well.
5. Conjugate Incubation:
Incubate the plate covered with an adhesive protective film for 60 ± 5 minutes at room temperature.
6. Washing:
Remove the adhesive film, aspirate the contents of all the wells and rinse 4 times, filling each well with 300 µL of washing solution. Wait 30 seconds between each wash.
7. Distribution of the substrate:
Dispense 100 µL of substrate in each well.
8. Substrate incubation:
Incubate the plate for 15 minutes at room temperature.
9. Interruption of the reaction:
Dispense 100 µL of stop solution, in the same order as that followed for point 4.
10. Reading:
Read the O.D. at 450 nm or 450/620 nm within 30 minutes. Repeat the reading at 405 nm in case of O.D. > 2.000.

9. SCHEME OF TEST PROCEDURE FOR SARS-CoV-2 IgM

- | | |
|--------|--|
| STEP 1 | Place 100 µL of controls (positive, negative and cut-off) and 50 µL of sample sorbent 4 + 50 µL diluted sample in the wells of the strips. Mix well. |
| | - |
| | Incubate for 60 ± 5 min. at R.T. |
| | - |
| | Wash 4 times (300 µL) |
| | - |
| STEP 2 | Add 100 µL of conjugate to each well |
| | - |
| | Incubate for 60 ± 5 min. at R.T. |
| | - |
| | Wash 4 times (300 µL) |
| | - |
| STEP 3 | Add 100 µL of Substrate to each well |
| | - |
| | Incubate for 15 min. at R.T. |
| | - |
| STEP 4 | Add 100 µL of Stop Solution |
| | - |
| | Read the O.D. at 450 nm or 450/620 nm within 30 min |

10. TEST VALIDATION

Blank: < 0.150

Negative control: the ratio between the negative control OD and cut-off OD must be < 0.6.

Positive control: the positive control must have an OD at least 1.5 times higher than the cut-off.

Cut-off control: the OD of the cut-off must be > 0.2.

Repeat the test, if one of the results of the control sera is not within the acceptability range.

If the problem persists contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

11. INTERPRETATION OF THE RESULTS

Calculate the ratio between the OD value of the sample and that of the Cut-off (Index).

The sample is considered:

POSITIVE: when the Index is > 1.1

DOUBTFUL: for all the values between 0.9 and 1.1

NEGATIVE: when the Index is < 0.9

12. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

13. ANALITICAL SPECIFICITY

4 samples (2 Negative and 2 Positive) were spiked with the following potentially interfering factors and then tested:

Rheumatoid factor (44-220 IU/ml)
 Bilirubin (4.5-45 mg/dl)
 Triglycerides (10-250 mg/dl)
 Hemoglobin (5-30 mg/ml)

Minimal interference could be observed with samples containing high concentrations of Rheumatoid Factor, however in no case the misinterpretation of the results (false positive or false negative) has been observed.

No interference was found with EDTA and heparin in 70 plasma samples.

14. CROSS-REACTIONS

128 samples positive to Mycoplasma p. (5), Cytomegalovirus (4), RSV (2), Influenza B (1), Influenza A (5), Adenovirus (3), Rheumatoid Factor (4), ANA (22), SARS-CoV (18), Human coronaviruses (21), Epstein-Barr Virus (13), HIV (10), Rubella (3), Chlamydomphila pneumoniae (5), Parainfluenza (5), Legionella pneumoniae (2) e Bordetella pertussis (5) were tested.

No significant cross-reactions were found except for a minimal cross-reactivity against Mycoplasma p. (1), Reumatoid Factor (1), ANA (3), Human coronaviruses (4), SARS-CoV (5) e Epstein-Barr Virus (3).

15. METHOD COMPARISON

In an experimentation 597 samples have been tested with Diesse kit and with an immunofluorescent test.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	57	10	67
	-	8	522	530
	Total	65	532	597

Percent Positive Agreement (~Diagnostic Sensitivity):

87.7% CI95%: 77.5-93.6.

Percent Negative Agreement (~Diagnostic Specificity):

98.1% CI95%: 96.6-99.0.

The analysis of the results of 29 samples with known date of symptoms onset showed that the sensitivity of the ELISA tests increases 10 days after the symptoms onset. 87.5% sensitivity is reached after 12 days.

Days for symptoms onset	Sensitivity
<10	58.3%
>10	70.6%
>12	87.5%

The analysis of the results of 282 samples with known pre-pandemic sampling date (before June 2019), showed that the specificity of the ELISA tests is 99%.

16. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision		Precision between batches	
	Mean Index	CV%	Mean Index	CV%	Mean Index	CV%
1	1.6	8.8	1.5	13.3	0.2	25.0*
2	0.5	10.0	2.1	11.9	2.0	7.5
3	1.9	2.6	0.8	11.3	0.5	8.0
4	0.6	10.0	0.7	10.0	1.8	10.6

* Artifact caused by the known fault of Variation Coefficient which becomes extremely sensitive to even very small changes in the mean when the mean value is near zero.

17. TROUBLESHOOTING

PROBLEM	POSSIBLE ERROR	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure. Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert).
		Check for moisture in unused plate. (Silica gel desiccant must be pale yellow). Repeat test.
Invalid run (all positive)	Contamination of substrate	Collect new aliquot of substrate.
	Inadequate washing	Ensure that plate washer works well
Poor precision	Inadequate aspiration of wells	Ensure that plate washer works well
	Pipetting error	Check pipette function
	Addition of reagents too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of air-bubbles	Avoid air-bubbles formation during pipetting.
	Optical pathway not clean	Check instrument light source and detector

		for dirt. Wipe bottom of plate with soft tissue.
Inadequate color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
	Inadequate volume of substrate added	Check pipette function

18. REFERENCES

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
2. Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
5. Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

DIESSE Diagnostica Senese S.p.A.
 Via delle Rose 10
 53035 Monteriggioni (Siena)
 Italy

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- Xintian et al. «Evolution of the novel coronavirus from the ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission». *Science China Life Sciences*; January 2020
- Chen et al. «Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study». *The Lancet*, Vol. 395 (507-513) February 15, 2020
- Huang et al. «Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China». *The Lancet*, Vol. 395 (497-506) February 15, 2020

- Yuxian et al. «Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus». *Journal of Clinical Microbiology* - Vol. 43 n8 (p.3718-37296); August 2005
- Woo et al. «Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia». *Journal of Clinical Microbiology* - Vol. 42 n5 (2306-2309); May 2004

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
18. BIBLIOGRAFIA

- Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; janeiro 2020
- Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) 15 de fevereiro 2020
- Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) 15 de fevereiro 2020
- Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); agosto 2005
- Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); maio 2004

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Itália

	EN ES IT DE	Use By Fecha de caducidad Utilizzare entro Verwendbar bis	FR GR PT CZ	Utiliser jusque Ημερομηνία λήξης Prazo de validade Použijte do
	EN ES IT DE	Caution, consult accompanying documents Atención, ver instrucciones de uso Attenzione, vedere le istruzioni per l'uso Achtung, die Gebrauchsanleitung lesen	FR GR PT CZ	Attention voir notice d'instructions Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα Atenção, consulte a documentação incluída Pozor, řiďte se návodem k použití
	EN ES IT DE	Manufacturer Fabricante Fabbicante Hersteller	FR GR PT CZ	Fabricant Κατασκευαστής Fabricante Výrobce
	EN ES IT DE	Contains sufficient for <n> tests Contenido suficiente para <n> ensayos Contenuto sufficiente per "n" saggi Inhalt reicht für „n“ Tests	FR GR PT CZ	Contenu suffisant pour "n" tests Περιεχόμενο επαρκές για «n» εξετάσεις Conteúdo suficiente para "n" ensaios Obsah stačí na < n > testů
	EN ES IT DE	Temperature limitation Límite de temperatura Limiti di temperatura Temperaturgrenzwerte	FR GR PT CZ	Limites de température Περιορισμοί θερμοκρασίας Limites de temperatura Teplotní omezení
	EN ES IT DE	Consult Instructions for Use Consulte las instrucciones de uso Consultare le istruzioni per l'uso Die Gebrauchsanleitung lesen	FR GR PT CZ	Consulter les instructions d'utilisation Συμβουλευτείτε τις οδηγίες χρήσης Consulte as instruções de utilização Čtěte návod k použití
	EN ES IT DE	Biological risks Riesgo biológico Rischio biológico Biologisches Risiko	FR GR PT CZ	Risques biologiques Βιολογικοί κίνδυνοι Risco biológico Biologická rizika
	EN ES IT DE	Batch code Código de lote Codice del lotto Chargennummer	FR GR PT CZ	Code du lot Αριθμός Παρτίδας Código do lote Kód šarže