DIESS
DIESS

CHORUS

SARS-CoV-2 IgG

REF 81400

DIESS Diagnostica Senese
S.p.A.
Via delle Rose, 10
53035 Monteriggioni (Siena)
Italy

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CHORUS SARS-CoV-2 IgG

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16. BIBLIOGRAFIE

- Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
- Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
- Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
- Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory

Syndrome-Associated Coronavirus” Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005

5. Woo at al. “Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia” Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

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Percent Positive Agreement (~Sensibilitatea diagnosticului):
 100.0% CI95%: 98.2 – 100.0

Percent Negative Agreement: (~Specificitatea diagnosticului):
 99.6% CI95%: 97.9 – 99.9

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16. BIBLIOGRAFIE

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein

- for risk of human transmission” Science China Life sciences; January 2020
2. Chen et al “Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study”. The Lancet, Vol 395 (507-513) February 15, 2020
 3. Huang et al “Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China”. The Lancet, Vol 395 (497-506) February 15, 2020
 4. Yuxian et al. “Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus” Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
 5. Woo at al. “Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia” Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

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
16. LITERATUR

1. Xintian et al. „Evolution of the novel coronavirus from the ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission“, Science China Life Sciences, Januar 2020
2. Chen et al „Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study“. The Lancet, Vol 395 (507-513), 15. Februar 2020
3. Huang et al „Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China“. The Lancet, Vol 395 (497–506), 15. Februar 2020
4. Yuxian et al. „Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus“ Journal of clinical microbiology - Vol 43 Nr. 8 (S. 3718–37296); August 2005
5. Woo at al. „Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia“ Journal of clinical microbiology - Vol 42 Nr. 5 (2306–2309); Mai 2004

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italien

| | | |
|---|--|--|
|  | EN Date of manufacture ES Fecha de fabricación IT Data di fabbricazione CS Datum výroby RO Data fabricatie | FR Date de fabrication GR Ημερομηνία Παραγωγής PT Data de fabrico DE Herstellungsdatum |
|  | EN Use By ES Fecha de caducidad IT Utilizzare entro CS Použitelné do RO A se folosi pana la | FR Utiliser jusque GR Ημερομηνία λήξης PT Prazo de validade DE Verwendbar bis |
|  | EN Do not reuse ES No reutilizar IT Non riutilizzare CS Nepoužívejte opakovaně RO A nu se refolosi | FR Ne pas réutiliser GR Μην κάνετε επαναληπτική χρήση PT Não reutilizar DE Nicht wieder verwenden |
|  | EN Caution, consult accompanying documents ES Atención, ver instrucciones de uso IT Attenzione, vedere le istruzioni per l'uso CS Pozor, čtěte přiložené dokumenty RO Atentie, consultati documentele insoitoare | FR Attention voir notice d'instructions GR Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα PT Atenção, consulte a documentação incluída DE Achtung, die Gebrauchsanleitung lesen |
|  | EN Manufacturer ES Fabricante IT Fabbicante CS Výrobce RO Productator | FR Fabricant GR Κατασκευαστής PT Fabricante DE Hersteller |
|  | EN Contains sufficient for <n> tests ES Contenido suficiente para <n> ensayos IT Contenuto sufficiente per "n" saggi CS Obsah stačí na <n> testů RO Continut sufficient pt <n> teste | FR Contenu suffisant pour "n" tests GR Περιεχόμενο επαρκές για «n» εξετάσεις PT Conteúdo suficiente para "n" ensaios DE Inhalt reicht für „n“ Tests |
|  | EN Temperature limitation ES Límite de temperatura IT Limiti di temperatura CS Teplotní omezení RO Limita da temperatura | FR Limites de température GR Περιορισμοί θερμοκρασίας PT Limites de temperatura DE Temperaturgrenzwerte |
|  | EN Consult Instructions for Use ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso CS Čtěte návod k použití RO Pentru utilizare consultati instructiunile | FR Consulter les instructions d'utilisation GR Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização DE Die Gebrauchsanleitung lesen |
|  | EN Biological risks ES Riesgo biológico IT Rischio biologico CS Biologická rizika RO Risk biologic | FR Risques biologiques GR Βιολογικοί κίνδυνοι PT Risco biológico DE Biologisches Risiko |
|  | EN Catalogue number ES Número de catálogo IT Numero di catalogo CS Katalogové číslo RO Numar de catalog | FR Référence du catalogue GR Αριθμός καταλόγου PT Referência de catálogo DE Katalognummer |
|  | EN In Vitro Diagnostic Medical Device ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro CS Lékařské vybavení pro diagnostiku <i>in vitro</i> RO Dizpositiv medical pentru diagnosticare <i>in vitro</i> | FR Dispositif médical de diagnostic in vitro GR In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro DE Medizinisches In-vitro-Diagnostikum |
|  | EN Batch code ES Código de lote IT Codice del lotto CS Kód šarže RO Lot | FR Code du lot GR Αριθμός Παρτίδας PT Código do lote DE Chargennummer |