CHORUS
Brucella IgG

DIESSE

REF 81164

DIESSE Diagnostica Senese
S.p.A.
Via delle Rose, 10
53035 Monteriggioni (Siena)
Italy

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CHORUS Brucella IgG

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15. BIBLIOGRAFIA

- Brucella Subgroup of the Northern Ireland Regional Zoonoses Group (2004). *Diagnosis and Management of Acute Brucellosis in Primary Care*.
- Franco, María Pía; Mulder, Maximilian; Gilman, Robert H; Smits, Henk L (2007). *Human brucellosis*. The Lancet Infectious Diseases. 7 (12): 775–786.
- Al Dahouk, Sascha; Nöckler, Karsten (July 2011). *Implications of laboratory diagnosis on brucellosis therapy*. Expert Review of Anti-infective Therapy. 9 (7): 833–845.
- Mantur, B.; Parande, A.; Amarnath, S.; Patil, G.; Walvekar, R.; Desai, A.; Parande, M.; Shinde, R.; Chandrashekar, M.; Patil, S. (2010). *ELISA versus Conventional Methods of Diagnosing Endemic Brucellosis*. American Journal of Tropical Medicine and Hygiene. 83 (2): 314–318.
- Annual Epidemiological Report for 2016 Brucellosis – Surveillance Report*. European Centre for Disease Prevention and Control.
- Paoli C.; Campisi E.; Gargani G. (1988). *Determinazione semiquantitativa delle IgM nella Brucellosi umana*. Atti XVII° Congresso Nazionale A.M.C.L.I. (Rimini 1988).
- M.J. Corbel (principal author) (2006). *Brucellosis in humans and animals*. WHO Library Cataloguing-in-Publication Data.
- D.T. Berman; B.L. Wilson; E. Moreno; R.D. Angus; L.M. Jones (Apr. 1980). *Characterization of Brucella abortus soluble antigen employed in immunoassay*. Journal of Clinical Microbiology; 355-362.
- P.M. Muñoz; C.M. Marín; D. Monreal; D. Gonzalez; B. Garin-Bastuji; R. Díaz; R.C. Mainar-Jaime; I. Moriyón; J.M. Blasco (Jan. 2005). *Efficacy of several serological tests and antigens for diagnosis of bovine Brucellosis in the presence of false-positive serological results due to Yersinia enterocolitica O:9*. Clinical and Diagnostic Laboratory Immunology; 141-151.
- B. Bonfini; G. Chiarenza; V. Paci; F. Sacchini; R. Salini; G. Vesco; S. Villari; K. Zilli; M. Tottarelli (2018). *Cross-reactivity in serological tests for brucellosis: a comparison for immune response of Escherichia coli O157:H7 and Yersinia Enterocolitica O:9 vs Brucella spp.* Veterinaria italiana 2018, 54 (2), 107-114.

11. J.S. García del Pozo; S.L. Ortuño; E. Navarro; J. Solera (2014). *Detection of IgM Anti-brucella Antibody in the Absence of IgGs: a challenge for the clinical interpretation of Brucella serology*. PLOS neglected Tropical Diseases. Volume 8, Issue 12, e3390
12. P.R. Murray; K.S. Rosenthal; M.A. Pfaller (2009). *Microbiologia Medica Settima Edizione, Brucella*, 318-321

INSTRUCTIONS FOR USE

CHORUS Brucella IgG

For the qualitative determination of anti-Brucella IgG antibodies

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the qualitative determination of IgG class antibodies against Brucella in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Brucellosis is an infectious disease caused by small, gram-negative, rod-shaped bacteria: there are 7 distinct species but only four species infect humans: *B. abortus*, *B. melitensis*, *B. suis* and *B. canis*.

People are infected by contact with infected animals or by eating meat or drinking unpasteurized milk; infected humans are not contagious. Brucellosis, for which endangered people are butchers, farmers, owners of pets, veterinaries and laboratory staff, appears in prevalence during winter and spring and the incubation period is between one and three weeks, but may be as long as two months. It is a systemic disease, which can damage different joints and organs and is characterized by not specific symptoms such as fever, sweating; sometimes it can be lymph-adenomegaly, hepatomegaly and splenomegaly.

The monitoring of antibodies' titer can serve as an indication of the status of infection: during the first days, IgM is the only immunoglobulin that appears; as the disease progresses, IgM recedes quantitatively and IgG becomes predominant. In chronic brucellosis, IgG may be produced for an extended period.

3. PRINCIPLE OF THE METHOD

The Chorus Brucella IgG device is ready to use for the detection of IgG antibodies against Brucella, in the Chorus/Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human serum.

After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgG monoclonal antibodies conjugated to horseradish peroxidase. The unbound conjugate is eliminated, and the peroxidase substrate added. The blue colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments. The results are expressed in Index (OD sample/OD cut-off).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minutes exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.

The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website

(<http://www.diesse.it/en/Support/Download/strumento:39/>)

5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.


7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use strongly hemolyzed, lipemic, icteric samples, serum not completely coagulated or samples presenting microbial contamination and heat-inactivated.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with B.abortus antigen

Position 5: UNCOATED MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing phenol 0.05%, Bronidox 0.02% and indicator to reveal the presence of the serum.

Position 2: CONJUGATE

Contents: anti-human IgG monoclonal antibodies labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

in which the sample is transferred.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.175 ml

Contents: Diluted serum containing antibodies anti-Brucella and preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.425 ml

Contents: Diluted serum containing antibodies anti-Brucella and preservative. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Chorus/Chorus TRIO Instrument
- Distilled or deionized water

- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in

Index (OD sample/OD cut-off).

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is > 1.2

NEGATIVE: when the result is < 0.8

DOUBTFUL/EQUIVOCAL: for all values between 0.8 and 1.2

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. CROSS-REACTIONS

It is not possible to exclude cross-reactions of antibodies directed against *Escherichia coli* O157, *Francisella tularensis*, *Yersinia enterocolitica*, *Vibrio cholerae*, *Salmonella* and among different biovars of *Brucella* spp. except for *Brucella canis*.

13. METHOD COMPARISON

In an experimentation 101 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	20	1	21
	-	0	80	80
	Total	20	81	101

Percent Positive Agreement (~Diagnostic Sensitivity):

100% CI_{95%}: 83.9-99.8

Percent Negative Agreement: (~Diagnostic Specificity):

98.8% CI_{95%}: 93.3-99.8

The agreement between the two methods is excellent with a Cohen's Kappa of 0.97.

14. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision	
	Mean Index	CV%	Mean Index	CV%
1	2.6	3.5	2.5	8.8
2	2.7	3.0	2.5	4.8
3	0.6	6.7	0.6	6.7

Sample	Precision between batches		Precision between instruments	
	Mean Index	CV%	Mean Index	CV%
1	2.6	5.8	2.6	9.3
2	2.7	1.9	2.7	8.9
3	0.6	6.7	0.6	11.9

15. REFERENCES

1. Brucella Subgroup of the Northern Ireland Regional Zoonoses Group (2004). *Diagnosis and Management of Acute Brucellosis in Primary Care*.
2. Franco, María Pía; Mulder, Maximilian; Gilman, Robert H; Smits, Henk L (2007). *Human brucellosis*. The Lancet Infectious Diseases. 7 (12): 775–786.
3. Al Dahouk, Sascha; Nöckler, Karsten (July 2011). *Implications of laboratory diagnosis on brucellosis therapy*. Expert Review of Anti-infective Therapy. 9 (7): 833–845.
4. Mantur, B.; Parande, A.; Amarnath, S.; Patil, G.; Walvekar, R.; Desai, A.; Parande, M.; Shinde, R.; Chandrashekar, M.; Patil, S. (2010). *ELISA versus Conventional Methods of Diagnosing Endemic Brucellosis*. American Journal of Tropical Medicine and Hygiene. 83 (2): 314–318.
5. *Annual Epidemiological Report for 2016 Brucellosis – Surveillance Report*. European Centre for Disease Prevention and Control.
6. Paoli C.; Campisi E.; Gargani G. (1988). *Determinazione semiquantitativa delle IgM nella Brucellosi umana*. Atti XVII° Congresso Nazionale A.M.C.L.I. (Rimini 1988).
7. M.J. Corbel (principal author) (2006). *Brucellosis in humans and animals*. WHO Library Cataloguing-in-Publication Data.
8. D.T. Berman; B.L. Wilson; E. Moreno; R.D. Angus; L.M. Jones (Apr. 1980). *Characterization of Brucella abortus soluble antigen employed in immunoassay*. Journal of Clinical Microbiology; 355-362.
9. P.M. Muñoz; C.M. Marín; D. Monreal; D. Gonzalez; B. Garin-Bastuji; R. Díaz; R.C. Mainar-Jaime; I. Moriyón; J.M. Blasco (Jan. 2005). *Efficacy of several serological tests and antigens for diagnosis of bovine Brucellosis in the presence of false-positive serological results due to Yersinia enterocolitica O:9*. Clinical and Diagnostic Laboratory Immunology; 141-151.
10. B. Bonfini; G. Chiarenza; V. Paci; F. Sacchini; R. Salini; G. Vesco; S. Villari; K. Zilli; M. Tottarelli (2018). *Cross-reactivity in serological tests for brucellosis: a comparison for immune response of Escherichia coli O157:H7 and Yersinia Enterocolitica O:9 vs Brucella spp.* Veterinaria italiana 2018, 54 (2), 107-114.
11. J.S. García del Pozo; S.L. Ortuño; E. Navarro; J. Solera (2014). *Detection of IgM Anti-brucella Antibody in the Absence of IgGs: a challenge for the clinical interpretation of Brucella serology*. PLOS neglected Tropical Diseases. Volume 8, Issue 12, e3390.
12. P.R. Murray; K.S. Rosenthal; M.A. Pfaller (2009). *Microbiologia Medica Settima Edizione, Brucella*, 318-321

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15. BIBLIOGRAFÍA

1. Brucella Subgroup of the Northern Ireland Regional Zoonoses Group (2004). *Diagnosis and Management of Acute Brucellosis in Primary Care*.
2. Franco, María Pia; Mulder, Maximilian; Gilman, Robert H; Smits, Henk L (2007). *Human brucellosis*. The Lancet Infectious Diseases. 7 (12): 775–786.
3. Al Dahouk, Sascha; Nöckler, Karsten (July 2011). *Implications of laboratory diagnosis on brucellosis therapy*. Expert Review of Anti-infective Therapy. 9 (7): 833–845.
4. Mantur, B.; Parande, A.; Amarnath, S.; Patil, G.; Walvekar, R.; Desai, A.; Parande, M.; Shinde, R.; Chandrashekar, M.; Patil, S. (2010). *ELISA versus Conventional Methods of Diagnosing Endemic Brucellosis*. American Journal of Tropical Medicine and Hygiene. 83 (2): 314–318.
5. *Annual Epidemiological Report for 2016 Brucellosis – Surveillance Report*. European Centre for Disease Prevention and Control.
6. Paoli C.; Campisi E.; Gargani G. (1988). *Determinazione semiquantitativa delle IgM nella Brucellosi umana*. Atti XVII° Congresso Nazionale A.M.C.L.I. (Rimini 1988).
7. M.J. Corbel (principal author) (2006). *Brucellosis in humans and animals*. WHO Library Cataloguing-in-Publication Data.
8. D.T. Berman; B.L. Wilson; E. Moreno; R.D. Angus; L.M. Jones (Apr. 1980). *Characterization of Brucella abortus soluble antigen employed in immunoassay*. Journal of Clinical Microbiology; 355-362.
9. P.M. Muñoz; C.M. Marín; D. Monreal; D. Gonzalez; B. Garin-Bastuji; R. Díaz; R.C. Mainar-Jaime; I. Moriyón; J.M. Blasco (Jan. 2005). *Efficacy of several serological tests and antigens for diagnosis of bovine Brucellosis in the presence of false-positive serological results due*

Referencia		
+	-	Total

- to *Yersinia enterocolitica* O:9. Clinical and Diagnostic Laboratory Immunology; 141-151.
10. B. Bonfini; G. Chiarenza; V. Paci; F. Sacchini; R. Salini; G. Vesco; S. Villari; K. Zilli; M. Tottarelli (2018). *Cross-reactivity in serological tests for brucellosis: a comparison for immune response of Escherichia coli O157:H7 and Yersinia Enterocolitica O:9 vs Brucella spp.* Veterinaria italiana 2018, 54 (2), 107-114.
 11. J.S. García del Pozo; S.L. Ortuño; E. Navarro; J. Solera (2014). *Detection of IgM Anti-brucella Antibody in the Absence of IgGs: a challenge for the clinical interpretation of Brucella serology.* PLOS neglected Tropical Diseases. Volume 8, Issue 12, e3390
 12. P.R. Murray; K.S. Rosenthal; M.A. Pfaller (2009). *Microbiologia Medica Settima Edizione, Brucella*, 318-321

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15. BIBLIOGRAPHIE

- Brucella Subgroup of the Northern Ireland Regional Zoonoses Group (2004). *Diagnosis and Management of Acute Brucellosis in Primary Care*.
- Franco, María Pía; Mulder, Maximilian; Gilman, Robert H; Smits, Henk L (2007). *Human brucellosis*. The Lancet Infectious Diseases. 7 (12): 775–786.
- Al Dahouk, Sascha; Nöckler, Karsten (July 2011). *Implications of laboratory diagnosis on brucellosis therapy*. Expert Review of Anti-infective Therapy. 9 (7): 833–845.
- Mantur, B.; Parande, A.; Amarnath, S.; Patil, G.; Walvekar, R.; Desai, A.; Parande, M.; Shinde, R.; Chandrashekar, M.; Patil, S. (2010). *ELISA versus Conventional Methods of Diagnosing Endemic Brucellosis*. American Journal of Tropical Medicine and Hygiene. 83 (2): 314–318.
- Annual Epidemiological Report for 2016 Brucellosis – Surveillance Report*. European Centre for Disease Prevention and Control.
- Paoli C.; Campisi E.; Gargani G. (1988). *Determinazione semiquantitativa delle IgM nella Brucellosi umana*. Atti XVII° Congresso Nazionale A.M.C.L.I. (Rimini 1988).
- M.J. Corbel (principal author) (2006). *Brucellosis in humans and animals*. WHO Library Cataloguing-in-Publication Data.
- D.T. Berman; B.L. Wilson; E. Moreno; R.D. Angus; L.M. Jones (Apr. 1980). *Characterization of Brucella abortus*

- soluble antigen employed in immunoassay. Journal of Clinical Microbiology*; 355-362.
9. P.M. Muñoz; C.M. Marín; D. Monreal; D. Gonzalez; B. Garin-Bastuji; R. Díaz; R.C. Mainar-Jaime; I. Moriyón; J.M. Blasco (Jan. 2005). *Efficacy of several serological tests and antigens for diagnosis of bovine Brucellosis in the presence of false-positive serological results due to Yersinia enterocolitica O:9*. *Clinical and Diagnostic Laboratory Immunology*; 141-151.
 10. B. Bonfini; G. Chiarenza; V. Paci; F. Sacchini; R. Salini; G. Vesco; S. Villari; K. Zilli; M. Tottarelli (2018). *Cross-reactivity in serological tests for brucellosis: a comparison for immune response of Escherichia coli O157:H7 and Yersinia Enterocolitica O:9 vs Brucella spp.* *Veterinaria italiana* 2018, 54 (2), 107-114.
 11. J.S. García del Pozo; S.L. Ortuño; E. Navarro; J. Solera (2014). *Detection of IgM Anti-brucella Antibody in the Absence of IgGs: a challenge for the clinical interpretation of Brucella serology*. *PLOS neglected Tropical Diseases*. Volume 8, Issue 12, e3390
 12. P.R. Murray; K.S. Rosenthal; M.A. Pfaller (2009). *Microbiologia Medica Settima Edizione, Brucella*, 318-321

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15. BIBLIOGRAFIE

1. Brucella Subgroup of the Northern Ireland Regional Zoonoses Group (2004). *Diagnosis and Management of Acute Brucellosis in Primary Care*.
2. Franco, María Pía; Mulder, Maximilian; Gilman, Robert H; Smits, Henk L (2007). *Human brucellosis*. The Lancet Infectious Diseases. 7 (12): 775–786.
3. Al Dahouk, Sascha; Nöckler, Karsten (July 2011). *Implications of laboratory diagnosis on brucellosis therapy*. Expert Review of Anti-infective Therapy. 9 (7): 833–845.
4. Mantur, B.; Parande, A.; Amarnath, S.; Patil, G.; Walvekar, R.; Desai, A.; Parande, M.; Shinde, R.; Chandrashekar, M.; Patil, S. (2010). *ELISA versus Conventional Methods of Diagnosing Endemic Brucellosis*. American Journal of Tropical Medicine and Hygiene. 83 (2): 314–318.
5. *Annual Epidemiological Report for 2016 Brucellosis – Surveillance Report*. European Centre for Disease Prevention and Control.
6. Paoli C.; Campisi E.; Gargani G. (1988). *Determinazione semiquantitativa delle IgM nella Brucellosi umana*. Atti XVII° Congresso Nazionale A.M.C.L.I. (Rimini 1988).
7. M.J. Corbel (principal author) (2006). *Brucellosis in humans and animals*. WHO Library Cataloguing-in-Publication Data.
8. D.T. Berman; B.L. Wilson; E. Moreno; R.D. Angus; L.M. Jones (Apr. 1980). *Characterization of Brucella abortus soluble antigen employed in immunoassay*. Journal of Clinical Microbiology; 355-362.
9. P.M. Muñoz; C.M. Marín; D. Monreal; D. Gonzalez; B. Garin-Bastuji; R. Díaz; R.C. Mainar-Jaime; I. Moriyón; J.M. Blasco (Jan. 2005). *Efficacy of several serological tests and antigens for diagnosis of bovine Brucellosis in the presence of false-positive serological results due to Yersinia enterocolitica O:9*. Clinical and Diagnostic Laboratory Immunology; 141-151.
10. B. Bonfini; G. Chiarenza; V. Paci; F. Sacchini; R. Salini; G. Vesco; S. Villari; K. Zilli; M. Tottarelli (2018). *Cross-reactivity in serological tests for brucellosis: a comparison for immune response of Escherichia coli O157:H7 and Yersinia Enterocolitica O:9 vs Brucella spp.* Veterinaria italiana 2018, 54 (2), 107-114.
11. J.S. García del Pozo; S.L. Ortuño; E. Navarro; J. Solera (2014). *Detection of IgM Anti-brucella Antibody in the Absence of IgGs: a challenge for the clinical interpretation of Brucella serology*. PLOS neglected Tropical Diseases. Volume 8, Issue 12, e3390
12. P.R. Murray; K.S. Rosenthal; M.A. Pfaller (2009). *Microbiologia Medica Settima Edizione, Brucella*, 318-321

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15. BIBLIOGRAFIJA

1. Brucella Subgroup of the Northern Ireland Regional Zoonoses Group (2004). *Diagnosis and Management of Acute Brucellosis in Primary Care*.
2. Franco, María Pía; Mulder, Maximilian; Gilman, Robert H; Smits, Henk L (2007). *Human brucellosis*. The Lancet Infectious Diseases. 7 (12): 775–786.
3. Al Dahouk, Sascha; Nöckler, Karsten (July 2011). *Implications of laboratory diagnosis on brucellosis therapy*. Expert Review of Anti-infective Therapy. 9 (7): 833–845.
4. Mantur, B.; Parande, A.; Amarnath, S.; Patil, G.; Walvekar, R.; Desai, A.; Parande, M.; Shinde, R.; Chandrashekar, M.; Patil, S. (2010). *ELISA versus Conventional Methods of Diagnosing Endemic Brucellosis*. American Journal of Tropical Medicine and Hygiene. 83 (2): 314–318.
5. *Annual Epidemiological Report for 2016 Brucellosis – Surveillance Report*. European Centre for Disease Prevention and Control.
6. Paoli C.; Campisi E.; Gargani G. (1988). *Determinazione semiquantitativa delle IgM nella Brucellosi umana*. Atti XVII° Congresso Nazionale A.M.C.L.I. (Rimini 1988).
7. M.J. Corbel (principal author) (2006). *Brucellosis in humans and animals*. WHO Library Cataloguing-in-Publication Data.
8. D.T. Berman; B.L. Wilson; E. Moreno; R.D. Angus; L.M. Jones (Apr. 1980). *Characterization of Brucella abortus soluble antigen employed in immunoassay*. Journal of Clinical Microbiology; 355-362.
9. P.M. Muñoz; C.M. Marín; D. Monreal; D. Gonzalez; B. Garin-Bastuji; R. Díaz; R.C. Mainar-Jaime; I. Moriyón; J.M. Blasco (Jan. 2005). *Efficacy of several serological tests and antigens for diagnosis of bovine Brucellosis in the presence of false-positive serological results due to Yersinia enterocolitica O:9*. Clinical and Diagnostic Laboratory Immunology; 141-151.
10. B. Bonfini; G. Chiarenza; V. Paci; F. Sacchini; R. Salini; G. Vesco; S. Villari; K. Zilli; M. Tottarelli (2018). *Cross-reactivity in serological tests for brucellosis: a comparison for immune response of Escherichia coli O157:H7 and Yersinia Enterocolitica O:9 vs Brucella spp*. Veterinaria italiana 2018, 54 (2), 107-114.
11. J.S. García del Pozo; S.L. Ortuño; E. Navarro; J. Solera (2014). *Detection of IgM Anti-brucella Antibody in the Absence of IgGs: a challenge for the clinical interpretation of Brucella serology*. PLOS neglected Tropical Diseases. Volume 8, Issue 12, e3390
12. P.R. Murray; K.S. Rosenthal; M.A. Pfaller (2009). *Microbiologia Medica Settima Edizione, Brucella*, 318-321

	EN Date of manufacture ES Fecha de fabricación IT Data di fabbricazione CS Datum výroby	FR Date de fabrication EL Ημερομηνία Παραγωγής PT Data de fabrico HR Datum proizvodnje
	EN Use By ES Fecha de caducidad IT Utilizzare entro CS Použití do	FR Utiliser jusque EL Ημερομηνία λήξης PT Prazo de validade HR Upotrijebiti do
	EN Do not reuse ES No reutilizar IT Non riutilizzare CS Nepoužívejte opakovaně	FR Ne pas réutiliser EL Μην κάνετε επαναληπτική χρήση PT Não reutilizar HR Ne koristiti ponovno
	EN Caution, consult accompanying documents ES Atención, ver instrucciones de uso IT Attenzione, vedere le istruzioni per l'uso CS Pozor, nahlédněte do průvodních dokumentů	FR Attention voir notice d'instructions EL Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα PT Atenção, consulte a documentação incluída HR Pozor, vidi upute za uporabu
	EN Manufacturer ES Fabricante IT Fabbicante CS Výrobce	FR Fabricant EL Κατασκευαστής PT Fabricante HR Proizvodač
	EN Contains sufficient for <n> tests ES Contenido suficiente para <n> ensayos IT Contenuto sufficiente per "n" saggi CS Obsahuje dostatečné množství pro <n> testů	FR Contenu suffisant pour "n" tests EL Περιεχόμενο επαρκές για «ν» εξετάσεις PT Conteúdo suficiente para "n" ensaios HR Sadržaj dovoljan za „tot.“ testova
	EN Temperature limitation ES Límite de temperatura IT Limiti di temperatura CS Teplotní omezení	FR Limites de température EL Περιορισμοί θερμοκρασίας PT Limites de temperatura HR Temperaturne granice
	EN Consult Instructions for Use ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso CS Přečtěte si návod k použití	FR Consulter les instructions d'utilisation EL Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização HR Pogledati upute za uporabu
	EN Biological risks ES Riesgo biológico IT Rischio biologico CS Biologická rizika	FR Risques biologiques EL Βιολογικοί κίνδυνοι PT Risco biológico HR Biološka opasnost
	EN Catalogue number ES Número de catálogo IT Numero di catalogo CS Katalogové číslo	FR Référence du catalogue EL Αριθμός καταλόγου PT Referência de catálogo HR Kataloški broj
	EN In Vitro Diagnostic Medical Device ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro CS Diagnostický zdravotnický prostředek in vitro	FR Dispositif médical de diagnostic in vitro EL In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro HR In Medicinski proizvod za dijagnostiku in vitro
	EN Batch code ES Código de lote IT Codice del lotto CS Kód šarže	FR Code du lot EL Αριθμός Παρτίδας PT Código do lote HR Šifra serije
	EN CE marking of conformity ES Marcado CE de conformidad IT Marcatura CE di conformità CS Označení shody CE	FR Marquage de conformité CE EL Σημάνση συμμορφωσης CE PT Marcação CE de conformidade HR CE oznaka sukladnosti