CHORUS

Herpes Simplex 1 IgG Recombinant

DIESSE
REF 81023

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16. BIBLIOGRAFIA

1. Arvaja M., Lehtinen M., Koskela P., Lappalainen M., Paavonen J., Vesikari T.; "Serological evaluation of Herpes Simplex Virus type 1 and type 2 infections in pregnancy"; Sexual Transmission Infections , 1999; 75: 168-171
2. Ashley R.L. et al.; "Genital herpes: review of the epidemic and potential use of type-specific serology"; Clinical of Microbiology Review, 1999; 12: 1-8
3. Ashley R.L., et al.; "Comparison of Western blot (immunoblot) and glycoprotein G-specific immunoblot assay for detecting antibodies to Herpes Simplex Virus type 1 and 2 in human sera"; Journal of Clinical Microbiology, 1988; 26: 662-667
4. Ashley R.L., Wu L., Pickering J.W., Tu M.C., Schonorenberg L; " Premarket evaluation of commercial glycoprotein G-based enzyme immunoassay for Herpes Simplex Virus type-specific antibodies"; Journal of Clinical Microbiology, 1998; 36: 294-295
5. Hashido M., Lee F. K., Inouye S., Kawana T.; "Detection of Herpes Simplex Virus type-specific antibodies by an Enzyme-Linked Immunosorbent Assay based on Glycoprotein G"; Journal of Medical Virology, 1997; 53: 319-323
6. Whittington W.L., Celum C.L., Cent A., Ashley R.L.; "Use of a glycoprotein G-based type-specific assay to detect antibodies to Herpes Simplex Virus type 2 among persons attending sexually transmitted disease clinics"; Sexual Transmission Diseases, 2001; 28: 99-104

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INSTRUCTIONS FOR USE

CHORUS Herpes Simplex 1 IgG Recombinant

**For the qualitative determination of IgG
antibodies anti-Herpes Simplex Virus (type 1)**

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the qualitative determination of IgG class antibodies to Herpes Simplex Virus (type 1) in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

The Herpes simplex virus (HSV) is a member of the Herpesviridae family, of which two types are known: type 1 (HSV-1) and type 2 (HSV-2) which present slight antigenic differences. HSV-1 causes chiefly oral-facial lesions, while HSV-2 is mainly responsible for genital lesions, but this distinction is not binding, both types occasionally causing infection in either anatomical site. HSV may also cause a form of ocular ch keratitis, and lesions of the central nervous system. HSV can affect practically the whole population. The primary infection is often in a subclinical form and is rarely diagnosed. After a latency period of variable duration, reactivation may occur and viral replication may or may not give rise to clinical lesions. Infection contracted during birth is of particular interest, this being an important cause of morbidity and mortality. It is therefore important to determine the immunitary state of women during pregnancy in order to detect serum conversion. The assay of specific IgG is important to establish the serological state of the patient.

3. PRINCIPLE OF THE METHOD

The Chorus Herpes Simplex 1 IgG Recombinant device is ready to use for the detection of IgG antibodies against Herpes Simplex Virus (type 1) in the Chorus/Chorus TRIO instruments. The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). The recombinant antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human serum. After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgG monoclonal antibodies conjugated to horseradish peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added. The blue colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments. The results are expressed in Index (OD sample/OD cut-off).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.


4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.
5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with recombinant HSV1 antigen, composed of recombinant G1 glycoprotein specific for Herpes Simplex Virus Type 1.

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution containing phenol 0.05%, Bronidox 0.02%, indicator to reveal the presence of the serum and cells extracts

Position 2: CONJUGATE

Contents: monoclonal antibodies anti-human IgG labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

In which the operator must place the undiluted serum

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.175 ml

Contents: Diluted human serum containing IgG antibodies anti-Herpes Simplex Virus (type 1) and preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.425 ml

Contents: Diluted human serum containing IgG antibodies anti-Herpes Simplex Virus (type 1) and preservative. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- SAMPLE DILUENT/NEGATIVE CONTROL **REF** 83607
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in Index (OD sample/OD cut-off).

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is > 1.1

NEGATIVE: when the result is < 0.9

DOUBTFUL/EQUIVOCAL: for all values between 0.9 and 1.1

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. ANALYTICAL SPECIFICITY

5 samples (2 Negative, 1 Cut-Off and 2 Positive) were spiked with the following potentially interfering factors and then tested:

Rheumatoid Factor (44 – 220 UI/ml)
Bilirubin (4.5 mg/dl – 45 mg/dl)
Triglycerides (10 mg/dl – 250 mg/dl)
Hemoglobin (5 mg/ml – 30 mg/ml)

The presence in the serum sample of the interfering substances described above except for Rheumatoid Factor does not affect the test result.

13. CROSS-REACTIONS

56 samples positive to Rubella, Epstein-Barr Virus, Varicella, Cytomegalovirus, HSV 2 and Toxoplasma were tested.

No significant cross-reactions were found.

14. METHOD COMPARISON

In an experimentation 153 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	80	0	80
	-	6	67	73
	Total	86	67	153

Percent Positive Agreement (~Diagnostic Sensitivity):

93.0% CI_{95%}: 85.4-96.7

Percent Negative Agreement: (~Diagnostic Specificity):

100.0% CI_{95%}: 94.6-99.9

15. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision	
	Mean Index	CV%	Mean Index	CV%
1	0.1	40.0*	0.1	-
2	0.3	-	0.3	-
3	0.5	-	0.5	12.0
4	0.8	3.8	0.8	10.0
5	1.3	5.4	1.4	7.9
6	1.8	4.4	1.8	8.3
7	3.2	4.7	3.0	7.3
8	4.5	5.1	4.2	7.4
9	6.4	3.6	6.6	10.6
10	6.9	4.5	7.2	12.5

* Artifact caused by the known fault of Variation Coefficient which becomes extremely sensitive to even very small changes in the mean when the mean value is near zero

Sample	Precision between batches		Precision between instruments	
	Mean Index	CV%	Mean Index	CV%
1	0.1	-	0.1	-
2	0.3	-	0.3	-
3	0.5	-	0.5	-
4	0.8	7.5	0.8	-
5	1.5	4.0	1.5	8.0
6	1.8	5.6	1.8	-
7	3.2	4.7	3.2	4.7
8	4.5	1.3	4.5	1.3
9	7.4	9.5	7.4	7.6
10	8.2	5.4	8.2	4.4

16. REFERENCES

1. Arvaja M., Lehtinen M., Koskela P., Lappalainen M., Paavonen J., Vesikari T.; "Serological evaluation of Herpes Simplex Virus type 1 and type 2 infections in pregnancy"; Sexual Transmission Infections, 1999; 75: 168–171
2. Ashley R.L. et al.; "Genital herpes: review of the epidemic and potential use of type-specific serology"; Clinical of Microbiology Review, 1999; 12: 1-8
3. Ashley R.L., et al.; "Comparison of Western blot (immunoblot) and glycoprotein G-specific immunoblot

- assay for detecting antibodies to Herpes Simplex Virus type 1 and 2 in human sera"; Journal of Clinical Microbiology, 1988; 26: 662-667
4. Ashley R.L., Wu L., Pickering J.W., Tu M.C., Schonorenberg L; " Premarket evaluation of commercial glycoprotein G-based enzyme immunoassay for Herpes Simplex Virus type-specific antibodies"; Journal of Clinical Microbiology, 1998; 36: 294-295
 5. Hashido M., Lee F. K., Inouye S., Kawana T.; "Detection of Herpes Simplex Virus type-specific antibodies by an Enzyme-Linked Immunosorbent Assay based on Glycoprotein G"; Journal of Medical Virology, 1997; 53: 319-323
 6. Whittington W.L., Celum C.L., Cent A., Ashley R.L.; "Use of a glycoprotein G-based type-specific assay to detect antibodies to Herpes Simplex Virus type 2 among persons attending sexually transmitted disease clinics"; Sexual Transmission Diseases, 2001; 28: 99-104

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3. Ashley R.L., et al.; "Comparison of Western blot (immunoblot) and glycoprotein G-specific immunoblot assay for detecting antibodies to Herpes Simplex Virus type 1 and 2 in human sera"; Journal of Clinical Microbiology, 1988; 26: 662-667
4. Ashley R.L., Wu L., Pickering J.W., Tu M.C., Schonorenberg L.; " Premarket evaluation of commercial glycoprotein G-based enzyme immunoassay for Herpes Simplex Virus type-specific antibodies"; Journal of Clinical Microbiology, 1998; 36: 294-295
5. Hashido M., Lee F. K., Inouye S., Kawana T.; "Detection of Herpes Simplex Virus type-specific antibodies by an Enzyme-Linked Immunosorbent Assay based on Glycoprotein G"; Journal of Medical Virology, 1997; 53: 319-323
6. Whittington W.L., Celum C.L., Cent A., Ashley R.L.; "Use of a glycoprotein G-based type-specific assay to detect antibodies to Herpes Simplex Virus type 2 among persons attending sexually transmitted disease clinics"; Sexual Transmission Diseases, 2001; 28: 99-104

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type 1 and 2 in human sera"; Journal of Clinical Microbiology, 1988; 26: 662-667

4. Ashley R.L., Wu L., Pickering J.W., Tu M.C., Schonorenberg L.; "Premarket evaluation of commercial glycoprotein G-based enzyme immunoassay for Herpes Simplex Virus type-specific antibodies"; Journal of Clinical Microbiology, 1998; 36: 294-295
5. Hashido M., Lee F. K., Inouye S., Kawana T.; "Detection of Herpes Simplex Virus type-specific antibodies by an Enzyme-Linked Immunosorbent Assay based on Glycoprotein G"; Journal of Medical Virology, 1997; 53: 319-323
6. Whittington W.L., Celum C.L., Cent A., Ashley R.L.; "Use of a glycoprotein G-based type-specific assay to detect antibodies to Herpes Simplex Virus type 2 among persons attending sexually transmitted disease clinics"; Sexual Transmission Diseases, 2001; 28: 99-104

DIESE Diagnostica Senese S.p.A.
Strada dei Laghi 39
53035 Monteriggioni (SIENA)
Italien

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pregnancy"; Sexual Transmission Infections , 1999; 75: 168-171

- Ashley R.L. et al.; "Genital herpes: review of the epidemic and potential use of type-specific serology"; Clinical of Microbiology Review, 1999; 12: 1-8
- Ashley R.L., et al.; "Comparison of Western blot (immunoblot) and glycoprotein G-specific immunoblot assay for detecting antibodies to Herpes Simplex Virus type 1 and 2 in human sera"; Journal of Clinical Microbiology, 1988; 26: 662-667
- Ashley R.L., Wu L., Pickering J.W., Tu M.C., Schonorenberg L.; " Premarket evaluation of commercial glycoprotein G-based enzyme immunoassay for Herpes Simplex Virus type-specific antibodies"; Journal of Clinical Microbiology, 1998; 36: 294-295
- Hashido M., Lee F. K., Inouye S., Kawana T.; "Detection of Herpes Simplex Virus type-specific antibodies by an Enzyme-Linked Immunosorbent Assay based on Glycoprotein G"; Journal of Medical Virology, 1997; 53: 319-323
- Whittington W.L., Celum C.L., Cent A., Ashley R.L.; "Use of a glycoprotein G-based type-specific assay to detect antibodies to Herpes Simplex Virus type 2 among persons attending sexually transmitted disease clinics"; Sexual Transmission Diseases, 2001; 28: 99-104

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16. BIBLIOGRAFIA

- Arvaja M., Lehtinen M., Koskela P., Lappalainen M., Paavonen J., Vesikari T.; "Serological evaluation of Herpes Simplex Virus type 1 and type 2 infections in pregnancy"; Sexual Transmission Infections , 1999; 75: 168-171
- Ashley R.L. et al.; "Genital herpes: review of the epidemic and potential use of type-specific serology"; Clinical of Microbiology Review, 1999; 12: 1-8
- Ashley R.L., et al.; "Comparison of Western blot (immunoblot) and glycoprotein G-specific immunoblot assay for detecting antibodies to Herpes Simplex Virus type 1 and 2 in human sera"; Journal of Clinical Microbiology, 1988; 26: 662-667
- Ashley R.L., Wu L., Pickering J.W., Tu M.C., Schonorenberg L.; "Premarket evaluation of commercial glycoprotein G-based enzyme immunoassay for Herpes Simplex Virus type-specific antibodies"; Journal of Clinical Microbiology, 1998; 36: 294-295
- Hashido M., Lee F. K., Inouye S., Kawana T.; "Detection of Herpes Simplex Virus type-specific antibodies by an Enzyme-Linked Immunosorbent Assay based on Glycoprotein G"; Journal of Medical Virology, 1997; 53: 319-323
- Whittington W.L., Celum C.L., Cent A., Ashley R.L.; "Use of a glycoprotein G-based type-specific assay to detect antibodies to Herpes Simplex Virus type 2 among persons attending sexually transmitted disease clinics"; Sexual Transmission Diseases, 2001; 28: 99-104

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