CHORUS

TRAb

DIESSE

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17. BIBLIOGRAFIA

1. Menconi F, Marcocci C, Marinò M. Diagnosis and classification of Graves' disease. Autoimmun Rev. 2014 Apr-May;13(4-5):398-402.
2. Tozzoli R, Bagnasco M, Giavarina D, Bizzaro N. TSH receptor autoantibody immunoassay in patients with Graves' disease: improvement of diagnostic accuracy over different generations of methods. Systematic review and meta-analysis. Autoimmun Rev. 2012 Dec;12(2):107-13.
3. Barbesino G, Tomer Y. Clinical review: Clinical utility of TSH receptor antibodies. J Clin Endocrinol Metab. 2013 Jun;98(6):2247-55.
4. Soh SB, Aw TC. Laboratory Testing in Thyroid Conditions - Pitfalls and Clinical Utility. Ann Lab Med. 2019 Jan;39(1):3-14.
5. Smith BR, Bolton J, Young S, Collyer A, Weeden A, Bradbury J, Weightman D, Perros P, Sanders J, Furmaniak J. A new assay for thyrotropin receptor autoantibodies. Thyroid. 2004 Oct;14(10):830-5.

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INSTRUCTIONS FOR USE

CHORUS TRAb

**For the quantitative determination of antibodies
anti-TSH receptor**

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the quantitative determination of antibodies anti-TSH receptor in human serum, using a disposable device applied on the Chorus TRIO instrument.

2. INTRODUCTION

Anti-TSH receptor antibodies are IgG immunoglobulins which, by binding to the thyrotropin receptor (Thyroid-stimulating hormone receptor TSHR), cause its stimulation or inhibition. The link created between TRAb and the TSH receptor prevents endogenous TSH from carrying out its physiological function. The dosage of anti-TSH receptor antibodies, with highly sensitive methods, assumes importance for the diagnosis of Graves' disease. The TRAb family can be divided into two categories: TRAb which by binding to TSHR mimic the biological action of TSH, or the production of thyroid hormones. These autoantibodies are present in most patients (TSAb). TRAb which by binding to TSHR block it, the resulting result is the absence of thyroid hormone production (TBAb). The reason why the TSH receptor behaves in two ways, stimulating or inhibiting the production of T3 and T4, is due to its particular structure.

The receptor is a 764 AA polypeptide, consisting of three portions such as extracellular, transmembrane and intracellular. In the extracellular portion there is the binding site of TSH and therefore of TRAb (TSAb and TBAb). The regions that are recognized as epitopes by TRAb appear to be very close, so much so that some amino acids appear to be shared. This discovery was possible thanks to the isolation of the human monoclonal autoantibody and its expansion on an animal model: mouse and rat. Several monoclonal antibodies (MAbs) have been isolated, including the most important: the stimulating Mabs M22 and the inhibiting Mabs 5C9.

3. PRINCIPLE OF THE METHOD

The Chorus TRAb device is ready to use for the detection of TSH receptor autoantibodies (TRAb) in Chorus TRIO instruments version 4.0.6 [65] and later. The test is based on an indirect competitive two-step ELISA method. The autoantibodies present in the serum and subsequently the biotinylated M22 monoclonal antibody compete with the binding of the TSH receptor (recombinant antigen absorbed in the solid phase).

The biotinylated M22 antibody binds streptavidin peroxidase in a second step: adding the specific substrate, TMB (tetramethylbenzidine), will cause the development of a blue

color. Subsequently the blue becomes yellow, adding sulfuric acid which blocks the enzymatic reaction; staining is inversely proportional to the concentration of TRAb autoantibodies present in the sample.

The disposable devices contain all the reagents to perform the test in the Chorus TRIO instrument.

The results are expressed in U/L.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by CE marked methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.

Waste disposal: serum samples, calibrators and strips, once used, must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30-minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged. Do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.

4. The devices are for use with the Chorus TRIO instrument; the Instructions for Use and the Instrument Operating Manual must be carefully followed.

The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website

(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)


5. Check that the Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 9: Space for application of bar code label

Position 8: CONJUGATE

Contents: M22 anti-TRAb monoclonal antibody labeled with Biotin, in protein solution containing preservative.

Position 7: MICROPLATE WELL

Coated with TSH receptor

Position 6: MICROPLATE WELL

Uncoated

Position 5: SAMPLE DILUENT

Contents: saline protein solution containing surfactant and preservative

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: STOP SOLUTION

Contents: 0.3 M sulfuric acid solution

Position 2: STREPTAVIDIN-POD

Contents: Streptavidin conjugated with horseradish peroxidase, in phosphate buffer solution containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

In which the sample is transferred.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.350 ml

Contents: Protein solution containing preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL

1 x 0.650 ml

Contents: Protein solution containing preservative. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	4 weeks at 2/8°C
CALIBRATOR	4 weeks at 2/8°C
POSITIVE CONTROL	4 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 7 days at 2/8°C, or frozen for longer periods at -20°C.

Avoid repeated freezing and thawing.

Do not keep the samples in auto-defrosting freezers.

Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.

2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 100 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the instrument Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus TRIO instrument expresses the result in U/L calculated on the basis of a lot-dependent graph stored in the instrument.

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is ≥ 0.4 U/L
NEGATIVE: when the result is < 0.4 U/L

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. CALIBRATION RANGE

Calibration range 0.40-30.0 U/L

For samples > 30 U/L repeat the test after diluting the sample in physiological solution. **DO NOT USE** Negative Control/Sample Diluent (REF. 83607).

13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, were between 0.0 and 0.24 U/L.

14. ANALYTICAL SPECIFICITY

3 samples (1 Negative and 2 Positive) were spiked with the following potentially interfering factors and then tested:

Rheumatoid factor (99 IU/ml)
Bilirubin (0.2 mg/ml)
Hemoglobin (2.4 mg/ml)

The presence in the serum sample of the interfering substances described above does not affect the test result.

15. METHOD COMPARISON

In an experimentation 176 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

Correlation	r	95%CI
Pearson	0.97	0.96-0.98
Spearman	0.91	0.88-0.93

The correlation between the two methods is very high.

16. PRECISION AND REPEATABILITY

Sample	Within run		Between run	
	Mean (U/L)	CV%	Mean (U/L)	CV%
1	0.5	14.0	1.1	13.6
2	2.0	11.5	3.1	6.8
3	10.8	10.2	8.1	8.1
4	26.9	1.4	15.0	11.0

Sample	Between lots		Between Instruments	
	Mean (U/L)	CV%	Mean (U/L)	CV%
1	2.3	13.0	1.1	13.6
2	7.6	12.9	3.1	14.2
3	14.4	10.8	7.9	6.6
4	-	-	14.2	14.1

17. REFERENCES

1. Menconi F, Marcocci C, Marinò M. Diagnosis and classification of Graves' disease. *Autoimmun Rev.* 2014 Apr-May;13(4-5):398-402.
2. Tozzoli R, Bagnasco M, Giavarina D, Bizzaro N. TSH receptor autoantibody immunoassay in patients with Graves' disease: improvement of diagnostic accuracy over different generations of methods. *Systematic review and meta-analysis.* *Autoimmun Rev.* 2012 Dec;12(2):107-13.
3. Barbesino G, Tomer Y. Clinical review: Clinical utility of TSH receptor antibodies. *J Clin Endocrinol Metab.* 2013 Jun;98(6):2247-55.
4. Soh SB, Aw TC. Laboratory Testing in Thyroid Conditions - Pitfalls and Clinical Utility. *Ann Lab Med.* 2019 Jan;39(1):3-14.
5. Smith BR, Bolton J, Young S, Collyer A, Weeden A, Bradbury J, Weightman D, Perros P, Sanders J, Furmaniak J. A new assay for thyrotropin receptor autoantibodies. *Thyroid.* 2004 Oct;14(10):830-5.

DIESE Diagnostica Senese S.p.A.
Strada dei Laghi 39
53035 Monteriggioni (SIENA)
Italy

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- disease: improvement of diagnostic accuracy over different generations of methods. Systematic review and meta-analysis. *Autoimmun Rev.* 2012 Dec; 12(2):107-13.
3. Barbesino G, Tomer Y. Clinical review: Clinical utility of TSH receptor antibodies. *J Clin Endocrinol Metab.* 2013 Jun; 98(6):2247-55.
 4. Soh SB, Aw TC. Laboratory Testing in Thyroid Conditions - Pitfalls and Clinical Utility. *Ann Lab Med.* 2019 Jan; 39(1):3-14.
 5. Smith BR, Bolton J, Young S, Collyer A, Weeden A, Bradbury J, Weightman D, Perros P, Sanders J, Furmaniak J. A new assay for thyrotropin receptor autoantibodies. *Thyroid.* 2004 Oct; 14(10):830-5.

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17. BIBLIOGRAFIA

1. Menconi F, Marcocci C, Marinò M. Diagnosis and classification of Graves' disease. *Autoimmun Rev.* 2014 Apr-May;13(4-5):398-402.
2. Tozzoli R, Bagnasco M, Giavarina D, Bizzaro N. TSH receptor autoantibody immunoassay in patients with Graves' disease: improvement of diagnostic accuracy over different generations of methods. *Systematic review and meta-analysis.* *Autoimmun Rev.* 2012 Dec;12(2):107-13.
3. Barbésino G, Tomer Y. Clinical review: Clinical utility of TSH receptor antibodies. *J Clin Endocrinol Metab.* 2013 Jun;98(6):2247-55.
4. Soh SB, Aw TC. Laboratory Testing in Thyroid Conditions - Pitfalls and Clinical Utility. *Ann Lab Med.* 2019 Jan;39(1):3-14.
5. Smith BR, Bolton J, Young S, Collyer A, Weeden A, Bradbury J, Weightman D, Perros P, Sanders J,

Furmaniak J. A new assay for thyrotropin receptor autoantibodies. Thyroid. 2004 Oct;14(10):830-5.

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17. BIBLIOGRAFIE

- Menconi F, Marcocci C, Marinò M. Diagnosis and classification of Graves' disease. Autoimmun Rev. 2014 Apr-May;13(4-5):398-402.
- Tozzoli R, Bagnasco M, Giavarina D, Bizzaro N. TSH receptor autoantibody immunoassay in patients with Graves' disease: improvement of diagnostic accuracy over different generations of methods. Systematic review and meta-analysis. Autoimmun Rev. 2012 Dec;12(2):107-13.
- Barbesino G, Tomer Y. Clinical review: Clinical utility of TSH receptor antibodies. J Clin Endocrinol Metab. 2013 Jun;98(6):2247-55.
- Soh SB, Aw TC. Laboratory Testing in Thyroid Conditions - Pitfalls and Clinical Utility. Ann Lab Med. 2019 Jan;39(1):3-14.
- Smith BR, Bolton J, Young S, Collyer A, Weeden A, Bradbury J, Weightman D, Perros P, Sanders J, Furmaniak J. A new assay for thyrotropin receptor autoantibodies. Thyroid. 2004 Oct;14(10):830-5.

DIESE Diagnostica Senese
S.p.A.
Strada dei Laghi 39
53035 Monteriggioni (SIENA)
Italy

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