# CHORUS

## SARS-CoV-2 Antigen

# DIESSE

**REF**      **81410**

DIESSE Diagnostica Senese  
S.p.A.  
Via delle Rose, 10  
53035 Monteriggioni (Siena)  
Italy

	Capitolo Section Κεφάλαιο Kapitel Capítulo Capítulo Chapitre Rozdział Capitol
Modifiche introdotte nella revisione corrente Changes introduced in the current revision Αλλαγές που εισήχθησαν στην τρέχουσα αναθεώρηση Änderungen in der aktuellen Revision eingeführt Cambios introducidos en la revisión actual Alterações introduzidas na revisão atual Modifications introduites dans la révision actuelle Zmiany wprowadzone w aktualnej rewizji Modificări introduse în actuala revizuire	<b>6,7,12,13</b>

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### **CHORUS SARS-CoV-2 Antigen**

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**17. BIBLIOGRAFIA**

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
2. Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
5. Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

**DIESSE Diagnostica Senese S.p.A.**  
Via delle Rose 10  
53035 Monteriggioni (SI)  
Italy

## INSTRUCTIONS FOR USE

### **CHORUS SARS-CoV-2 Antigen**

#### **For the semi-quantitative determination of SARS-CoV-2 virus in nasopharyngeal swab**

#### **For *In Vitro* Diagnostic Use Only**

##### **1. INTENDED USE**

Immunoenzymatic method for the semiquantitative determination of SARS-CoV-2 virus in nasopharyngeal swab, using a disposable device applied on the Chorus TRIO instruments.

##### **2. INTRODUCTION**

Coronaviruses are single positive-stranded RNA viruses belonging to *Coronaviridae* family, which are divided into four genera: *Alpha-*, *Beta-*, *Delta-* and *Gammacoronavirus*. CoVs are commonly found in many species of animals. Occasionally, the animal CoVs can acquire genetic mutations by errors during genome replication or recombination mechanism, which can further expand their tropism to humans. A total of six human CoV types were identified to be responsible for causing human respiratory ailments, which include two alpha CoVs and four beta CoVs (the most recent two are SARS-CoV and MERS-CoV). Typically, these CoVs cause asymptomatic infection or severe acute respiratory illness, including fever, cough, and shortness of breath. However, other symptoms such as gastroenteritis and neurological diseases of varying severity have also been reported. Since December 2019, the seventh novel human betacoronavirus (2019-nCoV) causing pneumonia outbreak was reported by the World Health Organization. This novel virus is linked with an outbreak of febrile respiratory illness in the city of Wuhan, Hubei Province, China. Some patients with 2019-nCoV have developed severe pneumonia, pulmonary oedema, ARDS, or multiple organ failure and have died. At present, information regarding the epidemiology and clinical features of pneumonia caused by 2019-nCoV is scarce. The monitoring of antibodies' titer can serve as an indication of the status of infection.

##### **3. PRINCIPLE OF THE METHOD**

The Chorus SARS-CoV-2 Antigen device is ready to use for the detection of SARS-CoV-2 virus antigen in nasopharyngeal swab, using the Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme Linked Immunosorbent Assay).

The anti-SARS-CoV-2 monoclonal capture antibodies are bound to the solid phase. The specific antigen binds to solid-phase antibodies following incubation with the samples. After washing to remove the proteins that have not reacted, the conjugate consisting of anti-SARS-CoV-2 monoclonal antibodies conjugated with horseradish peroxidase is added: the specific antigen binds to the conjugate following incubation.

After washings to eliminate the proteins which have not reacted, the peroxidase substrate is added.

The enzymatic reaction is stopped by the adding of the Stop Solution, which makes the solution turn to yellow. The colour which develops is proportional to the concentration of specific antigens present in the sample.

The disposable devices contain all the reagents to perform the test in the Chorus TRIO instruments.

The results are expressed in Arbitrary Units (AU/ml).

#### **4. WARNINGS AND PRECAUTIONS**

##### **FOR *IN VITRO* DIAGNOSTIC USE ONLY**

**This kit contains materials of human origin which have been tested and gave a negative response by CE marked methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.**

**Waste disposal: samples, calibrators and strips, once used, must be treated as infectious residuals and eliminated according to law.**

##### **Health and Safety Information**

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30-minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

##### **Analytical Precautions**

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged. Do not use devices which

are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.

- The devices are for use with the Chorus TRIO instrument; the Instructions for Use and the Instrument Operating Manual must be carefully followed.

**The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website**

(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)


- Check that the Chorus TRIO instrument is set up correctly (see Operating Manual).
- Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
- Avoid using self-defrosting freezers for the storage of the samples.
- Defective barcodes can be inserted manually in the instrument (see Operating Manual).
- Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
- Do not use the device after the expiry date.
- Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

## 5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

**DD** DEVICES 6 packages each containing 6 devices

Description:

**Position 8:** Space for application of bar code label

**Position 7:** Empty

**Position 6:** MICROPLATE WELL

Coated with anti-SARS-CoV-2 antibody

**Position 5:** Uncoated MICROPLATE WELL

**Position 4:** TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/mL and H<sub>2</sub>O<sub>2</sub> 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

**Position 3:** STOP SOLUTION

Contents: Sulfuric acid solution 0.3 M

**Position 2:** CONJUGATE

Contents: anti-SARS-CoV-2 monoclonal antibodies labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

**Position 1:** EMPTY WELL

in which the sample is transferred.

**Use:** equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

**CALIBRATOR** CALIBRATOR

1 x 0.475 ml

Contents Protein solution containing recombinant antigen (Nucleocapsid) and preservative. Liquid, ready for use.

**CONTROL +** POSITIVE CONTROL 1 x 0.925 ml

Contents: Protein solution containing recombinant antigen (Nucleocapsid) and preservative. Liquid, ready for use.

## MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-1000 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

## 6. STORAGE AND STABILITY OF REAGENTS

**Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the positive control (see section 9, Test validation).**

The expiry date is printed on each component and on the kit label.

**Reagents are stable after opening and/or preparation up to the expiration date, when stored correctly at 2-8°C**

## 7. SPECIMEN COLLECTION AND STORAGE

The sample is represented by viral transport medium of nasopharyngeal swab.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh sample should be tested within 24 hours of collection and should be kept at 2/8°C.

Heat-inactivation or freezing at -80°C can lead to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

## 8. ASSAY PROCEDURE

**The test Chorus SARS-CoV-2 Antigen (REF. 81410) cannot be performed simultaneously with the following tests:**

CHORUS BORDETELLA PERTUSSIS TOXIN IgA	REF 81274
CHORUS BORDETELLA PERTUSSIS TOXIN IgG	REF 81272
CHORUS BORRELLIA IgG CSF	REF 81048
CHORUS BORRELLIA IgG	REF 81046
CHORUS BORRELLIA IgM	REF 81047
CHORUS CALPROTECTIN	REF 86600
CHORUS CALPROTECTIN K	REF 86601
CHORUS CHIKUNGUNYA IgG	REF 81133
CHORUS CHIKUNGUNYA IgM	REF 81134
CHORUS DIPHtheria IgG	REF 81266
CHORUS LEGIONELLA URINARY ANTIGEN	REF 81301
CHORUS POLIOVIRUS IgG	REF 81264
CHORUS Q FEVER PHASE II IgG	REF 81171
CHORUS Q FEVER PHASE II IgM	REF 81172
CHORUS TETANUS IgG	REF 81270
CHORUS TICK-BORNE ENCEPHALITIS VIRUS IgG	REF 81276
CHORUS TICK-BORNE ENCEPHALITIS VIRUS IgM	REF 81277
CHORUS TOXOCARA IgG	REF 81280

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 150 µl of undiluted test sample in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

### 9. TEST VALIDATION

Use the Positive Control to check the validity of the obtained results. It should be used as reported in the instrument Operating Manual. If the instrument signals that the control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected. If the result of the control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554  
 Fax: 0039 0577 366605  
 email: scientificsupport@diesse.it

### 10. INTERPRETATION OF THE RESULTS

The Chorus TRIO instrument expresses the result in Arbitrary Units (AU/ml) calculated on the basis of a lot-dependent graph stored in the instrument.

The test on the examined sample can be interpreted as follows:

POSITIVE: when the result is > 16.5 AU/ml  
 NEGATIVE: when the result is < 13.5 AU/ml  
 EQUIVOCAL: for all values between 13.5 and 16.5 AU/ml

If the result is equivocal, repeat the test.

### 11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient. The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

### 12. CALIBRATION RANGE

Calibration range 3.0-100.0 AU/ml.  
 For samples > 100.0 AU/ml retest the sample diluted in the solution contained in a fresh nasopharyngeal swab.

### 13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 samples from healthy donors, were between 3.0 and 11.8 AU/ml.

### 14. CROSS-REACTIONS

24 samples, positive for *Acinetobacter baumannii*, Adenovirus (2-7), *Bordetella pertussis*, *Branhamella catarrhalis*, *Candida albicans*, *Candida glabrata*, Coxsackie virus (type B2-B3-B6), Echovirus (25 -7-31-30-3- 11-14), *Escherichia coli*, Influenza type A, Influenza type B, *Klebsiella pneumoniae*, *Legionella* (ser 1-2-3-4-5-6), *Mycoplasma pneumoniae*, *Neisseria gonorrhoeae*, Parainfluenza virus (type 1-3), Poliovirus (1-2-3), Rhino Virus,

*Staphylococcus epidermidis*, *Streptococcus* gr. A, *Streptococcus* gr. B, *Streptococcus* gr. C, *Streptococcus* gr. G, *Streptococcus pneumoniae*, were tested. No significant cross-reactions were found.

Furthermore, the kit does not cross-react with MERS Coronavirus and Human Coronavirus (NL63, 229E, OC43).

Cross-reaction with SARS-CoV virus is possible.

### 15. METHOD COMPARISON

In an experimentation the following samples were tested with the Diesse kit:

- 73 PCR positive samples, expressed in CT for various titers (high: CT <25, medium: CT between 25 and 30, low: CT between 30 and 35)
- 168 PCR negative samples

The diagnostic sensitivity is shown in the following table.

CT Titer	Sensitivity	CI <sub>95%</sub>
CT<25	100.0%	83.9 - 99.8
CT<30	66.7%	50.3 - 79.7
CT<35	41.1%	30.5 - 52.5

No false positives were found among the negative samples; therefore, the specificity is 100% (CI<sub>95%</sub> = 97.8 - 100.0).

### 16. PRECISION AND REPEATABILITY

Sample	Within run		Between run	
	Mean AU/ml	CV%	Mean AU/ml	CV%
1	50.0	7.5	35.2	10.3
2	22.5	14.9	12.7	5.2
3	7.5	10.1	5.8	10.3

Sample	Between lots		Between Instruments	
	Mean AU/ml	CV%	Mean AU/ml	CV%
1	82.1	14.9	32.3	15.0
2	21.3	12.0	13.3	12.4
3	6.1	12.3	6.3	14.9

### 17. REFERENCES

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" *Science China Life sciences*; January 2020
2. Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". *The Lancet*, Vol 395 (507-513) February 15, 2020
3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". *The Lancet*, Vol 395 (497-506) February 15, 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" *Journal of clinical microbiology-vol 43 n8 (p.3718-37296)*; august 2005
5. Woo et al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" *Journal of clinical microbiology -vol 42 n5 (2306-2309)*; May 2004

DIESE Diagnostica Senese S.p.A.  
Via delle Rose 10  
53035 Monteriggioni (Siena)  
Italy

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### 17. BIBLIOGRAFIE

- Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
- Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
- Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020

4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
5. Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

**DIESSE Diagnostica Senese S.p.A.**  
Via delle Rose 10  
53035 Monteriggioni (SI)  
Italy

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## 17. LITERATUR

1. Xintian et al. „Evolution of the novel coronavirus from the ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission“, Science China Life Sciences, Januar 2020
2. Chen et al „Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study“. The Lancet, Vol 395 (507–513), 15. Februar 2020
3. Huang et al „Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China“. The Lancet, Vol 395 (497–506), 15. Februar 2020
4. Yuxian et al. „Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus“ Journal of clinical microbiology - Vol 43 Nr. 8 (S. 3718–37296); August 2005
5. Woo at al. „Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia“ Journal of clinical microbiology - Vol 42 Nr. 5 (2306–2309); Mai 2004

**DIESE Diagnostica Senese S.p.A.**  
Via delle Rose 10  
53035 Monteriggioni (SI)  
Italy

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1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
2. Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); August 2005
5. Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

**DIESE Diagnostica Senese S.p.A.**  
Via delle Rose 10  
53035 Monteriggioni (SI)  
Italy

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CHORUS BORDETELLA PERTUSSIS TOXIN IgA	REF. 81274
CHORUS BORDETELLA PERTUSSIS TOXIN IgG	REF. 81272
CHORUS BORRELLIA IgG CSF	REF. 81048
CHORUS BORRELLIA IgG	REF. 81046
CHORUS BORRELLIA IgM	REF. 81047
CHORUS CALPROTECTIN	REF. 86600
CHORUS CALPROTECTIN K	REF. 86601
CHORUS CHIKUNGUNYA IgG	REF. 81133
CHORUS CHIKUNGUNYA IgM	REF. 81134
CHORUS DIPHTHERIA IgG	REF. 81266
CHORUS LEGIONELLA URINARY ANTIGEN	REF. 81301
CHORUS POLIOVIRUS IgG	REF. 81264
CHORUS Q FEVER PHASE II IgG	REF. 81171
CHORUS Q FEVER PHASE II IgM	REF. 81172
CHORUS TETANUS IgG	REF. 81270
CHORUS TICK-BORNE ENCEPHALITIS VIRUS IgG	REF. 81276
CHORUS TICK-BORNE ENCEPHALITIS VIRUS IgM	REF. 81277
CHORUS TOXOCARA IgG	REF. 81280

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1. Xintian et al. «Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission» Science China Life sciences; January 2020
2. Chen et al «Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study». The Lancet, Vol 395 (507-513) February 15, 2020
3. Huang et al «Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China». The Lancet, Vol 395 (497-506) February 15, 2020
4. Yuxian et al. «Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus» Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
5. Woo at al. «Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia» Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

**DIESSE Diagnostica Senese S.p.A.**  
 Via delle Rose 10  
 53035 Monteriggioni (SI)  
 Italy

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## 17. BIBLIOGRAFIA

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; janeiro 2020
2. Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) 15 de fevereiro 2020
3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) 15 de fevereiro 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); agosto 2005
5. Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); maio 2004

**DIESSE Diagnostica Senese S.p.A.**  
 Via delle Rose 10  
 53035 Monteriggioni (SI)  
 Italy

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**17. BIBLIOGRAFIE**

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
2. Chen et al „Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
3. Huang et al „Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
5. Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

**DIESSE Diagnostica Senese S.p.A.**  
Via delle Rose 10  
53035 Monteriggioni (SI)  
Italy

	IT	Data di fabbricazione	ES	Fecha de fabricación
	EN	Date of manufacture	FR	Date de fabrication
	CZ	Datum výroby	PT	Data de fabrico
	DE	Herstellungsdatum	RO	Data fabricatiei
	GR	Ημερομηνία Παραγωγής		
	IT	Utilizzare entro	ES	Fecha de caducidad
	EN	Use By	FR	Utiliser jusque
	CZ	Použitelné do	PT	Prazo de validade
	DE	Verwendbar bis	RO	A se folosi pana la
	GR	Ημερομηνία λήξης		
	IT	Non riutilizzare	ES	No reutilizar
	EN	Do not reuse	FR	Ne pas réutiliser
	CZ	Nepoužívejte opakovaně	PT	Não reutilizar
	DE	Nicht wieder verwenden	RO	A nu se refolosi
	GR	Μην κάνετε επαναληπτική χρήση		
	IT	Attenzione, vedere le istruzioni per l'uso	ES	Atención, ver instrucciones de uso
	EN	Caution, consult accompanying documents	FR	Attention voir notice d'instructions
	CZ	Pozor, čtěte příložené dokumenty	PT	Atenção, consulte a documentação incluída
	DE	Achtung, die Gebrauchsanleitung lesen	RO	Atentie, consultați documentele insoitoare
	GR	Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα		
	IT	Fabbricante	ES	Fabricante
	EN	Manufacturer	FR	Fabricant
	CZ	Výrobce	PT	Fabricante
	DE	Hersteller	RO	Productator
	GR	Κατασκευαστής		
	IT	Contenuto sufficiente per "n" saggi	ES	Contenido suficiente para <n> ensayos
	EN	Contains sufficient for <n> tests	FR	Contenu suffisant pour "n" tests
	CZ	Obsah stačí na <n> testů	PT	Conteúdo suficiente para "n" ensaios
	DE	Inhalt reicht für „n“ Tests	RO	Continut sufficient pt <n> teste
	GR	Περιεχόμενο επαρκές για «n» εξετάσεις		
	IT	Limiti di temperatura	ES	Límite de temperatura
	EN	Temperature limitation	FR	Limites de température
	CZ	Teplotní omezení	PT	Limites de temperatura
	DE	Temperaturgrenzwerte	RO	Limita da temperatura
	GR	Περιορισμοί θερμοκρασίας		
	IT	Consultare le istruzioni per l'uso	ES	Consulte las instrucciones de uso
	EN	Consult Instructions for Use	FR	Consulter les instructions d'utilisation
	CZ	Čtěte návod k použití	PT	Consulte as instruções de utilização
	DE	Die Gebrauchsanleitung lesen	RO	Pentru utilizarea consultați instrucțiunile
	GR	Συμβουλευτείτε τις οδηγίες χρήσης		
	IT	Rischio biologico	ES	Riesgo biológico
	EN	Biological risks	FR	Risques biologiques
	CZ	Biologická rizika	PT	Risco biológico
	DE	Biologisches Risiko	RO	Risk biologic
	GR	Βιολογικοί κίνδυνοι		
<b>REF</b>	IT	Numero di catalogo	ES	Número de catálogo
	EN	Catalogue number	FR	Référence du catalogue
	CZ	Katalogové číslo	PT	Referência de catálogo
	DE	Katalognummer	RO	Numar de catalog
	GR	Αριθμός καταλόγου		
<b>IVD</b>	IT	Dispositivo medico-diagnostico in vitro	ES	Producto sanitario para diagnóstico in vitro
	EN	In Vitro Diagnostic Medical Device	FR	Dispositif médical de diagnostic in vitro
	CZ	Lékařské vybavení pro diagnostiku in vitro	PT	Dispositivo médico para diagnóstico in vitro
	DE	Medizinisches In-vitro-Diagnostikum	RO	Dizpositiv medical pentru diagnosticare in vitro
	GR	In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν		
<b>LOT</b>	IT	Codice del lotto	ES	Código de lote
	EN	Batch code	FR	Code du lot
	CZ	Kód šarže	PT	Código do lote
	DE	Chargennummer	RO	Lot
	GR	Αριθμός Παρτίδας		