# CHORUS

## SARS-CoV-2 IgG

**REF 81400**

	Kapitola Capitol Kapitel
Změny zavedené v aktuální revizi Modificări introduse în actuala revizuire Änderungen in der aktuellen Revision eingeführt	12 – 13– 14

# DIESSE

DIESSE Diagnostica Senese  
S.p.A.  
Via delle Rose, 10  
53035 Monteriggioni (Siena)  
Italy

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### CHORUS SARS-CoV-2 IgG

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### 16. BIBLIOGRAFIE

- Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
- Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
- Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
- Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory

Syndrome-Associated Coronavirus” Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005

5. Woo at al. “Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia” Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

**DIESSE Diagnostica Senese S.p.A.**  
**Via delle Rose 10**  
**53035 Monteriggioni (Siena)**  
**Italy**

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Percent Positive Agreement (~Sensibilitatea diagnosticului):  
 98.2% CI95%: 95.8 – 99.2

Percent Negative Agreement (~Specificitatea diagnosticului):  
 99.7% CI95%: 98.2 – 99.9

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Zile de la debutul simptomelor	Sensibilitate
≤ 10	93.8%
> 10	98.4%
> 12	98.4%

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## 16. BIBLIOGRAFIE

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein

- for risk of human transmission” Science China Life sciences; January 2020
2. Chen et al “Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study”. The Lancet, Vol 395 (507-513) February 15, 2020
  3. Huang et al “Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China”. The Lancet, Vol 395 (497-506) February 15, 2020
  4. Yuxian et al. “Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus” Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
  5. Woo at al. “Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia” Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

**DIESSE Diagnostica Senese S.p.A.**  
Via delle Rose 10  
53035 Monteriggioni (Siena)  
Italy

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

### 16. LITERATUR

1. Xintian et al. „Evolution of the novel coronavirus from the ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission“, Science China Life Sciences, Januar 2020
2. Chen et al. „Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study“. The Lancet, Vol 395 (507-513), 15. Februar 2020
3. Huang et al. „Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China“. The Lancet, Vol 395 (497–506), 15. Februar 2020
4. Yuxian et al. „Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus“ Journal of clinical microbiology - Vol 43 Nr. 8 (S. 3718–37296); August 2005
5. Woo at al. „Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia“ Journal of clinical microbiology - Vol 42 Nr. 5 (2306–2309); Mai 2004

**DIESSE Diagnostica Senese S.p.A.**  
Via delle Rose 10  
53035 Monteriggioni (Siena)  
Italien

	EN Date of manufacture ES Fecha de fabricación IT Data di fabbricazione CS Datum výroby RO Data fabricatie	FR Date de fabrication GR Ημερομηνία Παραγωγής PT Data de fabrico DE Herstellungsdatum
	EN Use By ES Fecha de caducidad IT Utilizzare entro CS Použitelné do RO A se folosi pana la	FR Utiliser jusque GR Ημερομηνία λήξης PT Prazo de validade DE Verwendbar bis
	EN Do not reuse ES No reutilizar IT Non riutilizzare CS Nepoužívejte opakovaně RO A nu se refolosi	FR Ne pas réutiliser GR Μην κάνετε επαναληπτική χρήση PT Não reutilizar DE Nicht wieder verwenden
	EN Caution, consult accompanying documents ES Atención, ver instrucciones de uso IT Attenzione, vedere le istruzioni per l'uso CS Pozor, čtěte přiložené dokumenty RO Atentie, consultati documentele insoitoare	FR Attention voir notice d'instructions GR Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα PT Atenção, consulte a documentação incluída DE Achtung, die Gebrauchsanleitung lesen
	EN Manufacturer ES Fabricante IT Fabbicante CS Výrobce RO Productator	FR Fabricant GR Κατασκευαστής PT Fabricante DE Hersteller
	EN Contains sufficient for <n> tests ES Contenido suficiente para <n> ensayos IT Contenuto sufficiente per "n" saggi CS Obsah stačí na <n > testů RO Continut sufficient pt <n> teste	FR Contenu suffisant pour "n" tests GR Περιεχόμενο επαρκές για «n» εξετάσεις PT Conteúdo suficiente para "n" ensaios DE Inhalt reicht für „n“ Tests
	EN Temperature limitation ES Límite de temperatura IT Limiti di temperatura CS Teplotní omezení RO Limita da temperatura	FR Limites de température GR Περιορισμοί θερμοκρασίας PT Limites de temperatura DE Temperaturgrenzwerte
	EN Consult Instructions for Use ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso CS Čtěte návod k použití RO Pentru utilizare consultati instructiunile	FR Consulter les instructions d'utilisation GR Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização DE Die Gebrauchsanleitung lesen
	EN Biological risks ES Riesgo biológico IT Rischio biologico CS Biologická rizika RO Risk biologic	FR Risques biologiques GR Βιολογικοί κίνδυνοι PT Risco biológico DE Biologisches Risiko
	EN Catalogue number ES Número de catálogo IT Numero di catalogo CS Katalogové číslo RO Numar de catalog	FR Référence du catalogue GR Αριθμός καταλόγου PT Referência de catálogo DE Katalognummer
	EN In Vitro Diagnostic Medical Device ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro CS Lékařské vybavení pro diagnostiku <i>in vitro</i> RO Dizpositiv medical pentru diagnosticare <i>in vitro</i>	FR Dispositif médical de diagnostic in vitro GR In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro DE Medizinisches In-vitro-Diagnostikum
	EN Batch code ES Código de lote IT Codice del lotto CS Kód šarže RO Lot	FR Code du lot GR Αριθμός Παρτίδας PT Código do lote DE Chargennummer