CHORUS

SARS-CoV-2 IgA

DIESSE
DIESSE

REF 81402

DIESSE Diagnostica Senese
S.p.A.
Via delle Rose, 10
53035 Monteriggioni (Siena)
Italy

	Capitolo Section Κεφάλαιο Kapitel Capítulo Capítulo Chapitre Rozdział Capitol
Modifiche introdotte nella revisione corrente Changes introduced in the current revision Αλλαγές που εισήχθησαν στην τρέχουσα αναθεώρηση Änderungen in der aktuellen Revision eingeführt Cambios introducidos en la revisión actual Alterações introduzidas na revisão atual Modifications introduites dans la révision actuelle Zmiany wprowadzone w aktualnej rewizji Modificări introduse în actuala revizuire	12 – 13– 14

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CHORUS SARS-CoV-2 IgA

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Percent Positive Agreement (~Sensibilità Diagnostica):
 95.7% CI95%: 91.3-97.9.

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16. BIBLIOGRAFIA

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein

- for risk of human transmission” Science China Life sciences; January 2020
2. Chen et al “Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study”. The Lancet, Vol 395 (507-513) February 15, 2020
 3. Huang et al “Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China”. The Lancet, Vol 395 (497-506) February 15, 2020
 4. Yuxian et al. “Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus” Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
 5. Woo at al. “Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia” Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

INSTRUCTIONS FOR USE

CHORUS SARS-CoV-2 IgA

**For the qualitative determination of IgA
antibodies anti-SARS-CoV-2 virus**

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the qualitative determination of IgA class antibodies to SARS-CoV-2 virus in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Coronaviruses are single positive-stranded RNA viruses belonging to *Coronaviridae* family, which are divided into four genera: *Alpha-*, *Beta-*, *Delta-* and *Gammacoronavirus*. CoVs are commonly found in many species of animals. Occasionally, the animal CoVs can acquire genetic mutations by errors during genome replication or recombination mechanism, which can further expand their tropism to humans. A total of six human CoV types were identified to be responsible for causing human respiratory ailments, which include two alpha CoVs and four beta CoVs (the most recent two are SARS-CoV and MERS-CoV). Typically, these CoVs cause asymptomatic infection or severe acute respiratory illness, including fever, cough, and shortness of breath. However, other symptoms such as gastroenteritis and neurological diseases of varying severity have also been reported. Since December 2019, the seventh novel human betacoronavirus (2019-nCoV) causing pneumonia outbreak was reported by the World Health Organization. This novel virus is linked with an outbreak of febrile respiratory illness in the city of Wuhan, Hubei Province, China. Some patients with 2019-nCoV have developed severe pneumonia, pulmonary oedema, ARDS, or multiple organ failure and have died. At present, information regarding the epidemiology and clinical features of pneumonia caused by 2019-nCoV is scarce. The monitoring of antibodies' titer can serve as an indication of the status of infection: during the first days, IgA and IgM appear; as the disease progresses, IgG becomes predominant.

3. PRINCIPLE OF THE METHOD

The Chorus SARS-CoV-2 IgA device is ready to use for the detection of IgA antibodies against SARS-CoV-2 virus, using Chorus/Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human serum.

After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgA monoclonal antibodies conjugated to horseradish peroxidase. The unbound conjugate is eliminated, and the

peroxidase substrate added. The blue colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments. The results are expressed in Index (OD sample/OD cut-off).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully

followed and the Instrument Operating Manual must be consulted.

The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website

(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)


5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with inactivated native antigen

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic solution in phosphate buffer with Sodium Azide 0.09% and methyl orange as dye.

Position 2: CONJUGATE

Contents: anti-human IgA monoclonal antibodies, labeled with peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

in which the sample is transferred.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.175 ml

Contents: Protein solution containing specific antibodies capable of binding the antigen present on the microplate, and preservative Liquid, ready for use.

CONTROL + POSITIVE CONTROL

1 x 0.425 ml

Contents: Protein solution containing specific antibodies capable of binding the antigen present on the microplate, and preservative Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER [REF] 83606
- CLEANING SOLUTION 2000 [REF] 83609
- SANITIZING SOLUTION [REF] 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT [REF] 83607
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

The reagents have stability after opening and/or preparation until the expiration date, if stored correctly at 2-8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C.

Avoid repeated thawing-freezing cycles.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.

- Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in Index (OD sample/OD cut-off).

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is > 1.1
 NEGATIVE: when the result is < 0.9
 DOUBTFUL/EQUIVOCAL: for all values between 0.9 and 1.1

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. ANALYTICAL SPECIFICITY

4 samples (2 Negative, 1 Cut-Off and 1 Positive) were spiked with the following potentially interfering factors and then tested:

Rheumatoid Factor - 110 IU/mL
 Bilirubin - 18 mg/dL
 Triglycerides - 150 mg/dL
 Hemoglobin - 15 mg/mL

The presence in the serum sample of the interfering substances described above does not affect the test result.

13. CROSS-REACTIONS

95 samples, positive to Cytomegalovirus (4), Mycoplasma pneumoniae (5), Influenza B (1), Influenza A (5), RSV (1), Adenovirus (3), ANA (16), Rheumatoid Factor (2), Rubella Virus (3), Chlamydia pneumoniae (5), Parainfluenza Virus (5), Bordetella pertussis (5), Legionella pneumoniae (2), Human Coronavirus 229E Nucleoprotein IgG (18), Human Coronavirus HKU1 (5), Human Coronavirus OC43 (5), Human Coronavirus NL63 (5) and Human Coronavirus 229E (5) were tested.

No significant cross-reactions were found except for Influenza A (1 sample), ANA (4 samples), Rheumatoid Factor (1 sample),

Parainfluenza Virus (3 samples) and Human Coronavirus 229E Nucleoprotein IgG (3 samples).

14. METHOD COMPARISON

In an experimentation 400 samples have been tested with Diesse kit and with another commercial kit.

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	155	1	156
	-	7	237	244
	Total	162	238	400

Percent Positive Agreement (~Diagnostic Sensitivity):
 95.7% CI95%: 91.3-97.9.

Percent Negative Agreement (~Diagnostic Specificity):
 99.6% CI95%: 97.6-99.9.

The analysis of the results of 83 samples with known date pre pandemic (before June 2017) showed that the sensitivity of the ELISA tests is 99%.

The analysis of the results of 81 samples with known date of symptoms onset showed the following results:

Days from symptoms onset	Sensitivity
≤ 10	77.8%
> 10	96.2%
> 12	95.9%

15. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision	
	Mean (Index)	CV%	Mean (Index)	CV%
1	3.9	13.3	4.2	10.5
2	1.5	12.0	1.6	8.8
3	0.7	14.3	0.9	6.7

Sample	Precision between batches		Precision between instruments	
	Mean (Index)	CV%	Mean (Index)	CV%
1	3.3	13.0	4.2	14.3
2	1.6	10.6	1.5	13.3
3	0.8	12.5	0.8	13.8

16. REFERENCES

- Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
- Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
- Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
- Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
- Woo et al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus

Nucleocapsid protein for serodiagnosis of SARS
Coronavirus Pneumonia" Journal of clinical microbiology -
vol 42 n5 (2306-2309); May 2004

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

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16. BIBLIOGRAFIE

- Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
- Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020

3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
5. Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

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16. LITERATUR

1. Xintian et al. „Evolution of the novel coronavirus from the ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission“, Science China Life Sciences, Januar 2020
2. Chen et al. „Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study“. The Lancet, Vol 395 (507-513), 15. Februar 2020
3. Huang et al. „Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China“. The Lancet, Vol 395 (497–506), 15. Februar 2020
4. Yuxian et al. „Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus“ Journal of clinical microbiology - Vol 43 Nr. 8 (S. 3718–37296); August 2005
5. Woo et al. „Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia“ Journal of clinical microbiology - Vol 42 Nr. 5 (2306–2309); Mai 2004

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italien

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- Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
- Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
- Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
- Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
- Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

DIESE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

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16. BIBLIOGRAFÍA

1. Xintian et al. «Evolution of the novel coronavirus from the ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission». Science China Life Sciences; January 2020
2. Chen et al. «Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study». The Lancet, Vol. 395 (507-513) February 15, 2020
3. Huang et al. «Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China». The Lancet, Vol. 395 (497-506) February 15, 2020
4. Yuxian et al. «Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus». Journal of Clinical Microbiology - Vol. 43 n8 (p.3718-37296); August 2005
5. Woo et al. «Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia». Journal of Clinical Microbiology - Vol. 42 n5 (2306-2309); May 2004

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16. BIBLIOGRAPHIE

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
2. Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
5. Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italie

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16. BIBLIOGRAFIA

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
2. Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) 15 de fevereiro 2020
3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) 15 de fevereiro 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); agosto 2005
5. Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); maio 2004

DIESSE Diagnostica Senese S.p.A.
 Via delle Rose 10
 53035 Monteriggioni (Siena)
 Itália

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Percent Positive Agreement (~Sensibilitatea diagnosticului): 95.7% CI95%: 91.3 – 97.9.

Percent Negative Agreement (~Specificitatea diagnosticului): 99.6% CI95%: 97.6 – 99.9.

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Zile de debut al simptomelor
 Sensibilitate

Zile de debut al simptomelor	Sensibilitate
≤ 10	77.8%
> 10	96.2%
> 12	95.9%

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
16. BIBLIOGRAFIE

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" *Science China Life sciences*; January 2020
2. Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". *The Lancet*, Vol 395 (507-513) February 15, 2020
3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". *The Lancet*, Vol 395 (497-506) February 15, 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" *Journal of clinical microbiology-vol 43 n8 (p.3718-37296)*; august 2005
5. Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" *Journal of clinical microbiology - vol 42 n5 (2306-2309)*; May 2004

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

	IT Data di fabbricazione EN Date of manufacture CZ Datum výroby DE Herstellungsdatum GR Ημερομηνία Παραγωγής	ES Fecha de fabricación FR Date de fabrication PT Data de fabrico RO Data fabricatiei
	IT Utilizzare entro EN Use By CZ Použitelné do DE Verwendbar bis GR Ημερομηνία λήξης	ES Fecha de caducidad FR Utiliser jusque PT Prazo de validade RO A se folosi pana la
	IT Non riutilizzare EN Do not reuse CZ Nepoužívejte opakovaně DE Nicht wieder verwenden GR Μην κάνετε επαναληπτική χρήση	ES No reutilizar FR Ne pas réutiliser PT Não reutilizar RO A nu se refolosi
	IT Attenzione, vedere le istruzioni per l'uso EN Caution, consult accompanying documents CZ Pozor, čtěte příložené dokumenty DE Achtung, die Gebrauchsanleitung lesen GR Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα	ES Atención, ver instrucciones de uso FR Attention voir notice d'instructions PT Atenção, consulte a documentação incluída RO Atentie, consultați documentele însoțitoare
	IT Fabricante EN Manufacturer CZ Výrobce DE Hersteller GR Κατασκευαστής	ES Fabricante FR Fabricant PT Fabricante RO Productator
	IT Contenuto sufficiente per "n" saggi EN Contains sufficient for <n> tests CZ Obsah stačí na <n> testů DE Inhalt reicht für „n“ Tests GR Περιεχόμενο επαρκές για «n» εξετάσεις	ES Contenido suficiente para <n> ensayos FR Contenu suffisant pour "n" tests PT Conteúdo suficiente para "n" ensaios RO Continut sufficient pt <n> teste
	IT Limiti di temperatura EN Temperature limitation CZ Teplotní omezení DE Temperaturgrenzwerte GR Περιορισμοί θερμοκρασίας	ES Límite de temperatura FR Limites de température PT Limites de temperatura RO Limita da temperatura
	IT Consultare le istruzioni per l'uso EN Consult Instructions for Use CZ Čtěte návod k použití DE Die Gebrauchsanleitung lesen GR Συμβουλευτείτε τις οδηγίες χρήσης	ES Consulte las instrucciones de uso FR Consulter les instructions d'utilisation PT Consulte as instruções de utilização RO Pentru a utiliza consultați instrucțiunile
	IT Rischio biologico EN Biological risks CZ Biologická rizika DE Biologisches Risiko GR Βιολογικοί κίνδυνοι	ES Riesgo biológico FR Risques biologiques PT Risco biológico RO Risk biologic
	IT Numero di catalogo EN Catalogue number CZ Katalogové číslo DE Katalognummer GR Αριθμός καταλόγου	ES Número de catálogo FR Référence du catalogue PT Referência de catálogo RO Numar de catalog
	IT Dispositivo medico-diagnostico in vitro EN In Vitro Diagnostic Medical Device CZ Lékařské vybavení pro diagnostiku in vitro DE Medizinisches In-vitro-Diagnostikum GR In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν	ES Producto sanitario para diagnóstico in vitro FR Dispositif médical de diagnostic in vitro PT Dispositivo médico para diagnóstico in vitro RO Dizpozitiv medical pentru diagnosticare in vitro
	IT Codice del lotto EN Batch code CZ Kód šarže DE Chargennummer GR Αριθμός Παρτίδας	ES Código de lote FR Code du lot PT Código do lote RO Lot