CHORUS

DIESSE

SARS-CoV-2 “NEUTRALIZING” Ab

DIESSE Diagnostica Senese
S.p.A.

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18. BIBLIOGRAFIA

- Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
- Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020

- Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
- Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
- Woo et al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004
- Andreano et al. "SARS-CoV-2 escape in vitro from a highly neutralizing COVID-19 convalescent plasma" bioRxiv. Preprint. 2020 Dec 28
- <https://doi.org/10.1016/j.jinf.2020.12.024> (variante inglese)
- <https://doi.org/10.1016/j.jinf.2021.01.025> (variante brasiliana)
- <https://doi.org/10.1016/j.jinf.2021.01.007> (variante sudafricana)

INSTRUCTIONS FOR USE

CHORUS SARS-CoV-2 "NEUTRALIZING" Ab

For the quantitative determination of total antibodies anti-S1 SARS-CoV-2 virus as an aid in the determination of the immune reaction to SARS-CoV-2 or to the vaccine

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for quantitative determination of total antibodies anti-S1 SARS-CoV-2 virus in human serum, using a disposable device applied on the Chorus TRIO instrument. The test is intended as an aid in the determination of the immune reaction to SARS-CoV-2 or to the vaccine.

2. INTRODUCTION

Coronaviruses are positive single-stranded RNA viruses belonging to the Coronaviridae family, which are divided into four genera: Alpha, Beta, Delta and Gammacoronavirus. CoVs are commonly found in many animal species. Occasionally, animal CoVs may acquire genetic mutations by errors during genome replication, which can further expand their tropism to humans. A total of six human CoV types have been identified as responsible for human respiratory disorders, which include two CoV alpha and four CoV beta (the two most recent being SARS-CoV and MERS-CoV). Typically, these CoVs cause asymptomatic infection or severe acute respiratory illness, including fever, cough, and shortness of breath. However, other symptoms such as gastroenteritis and neurological diseases of varying severity have also been reported. Since December 2019, the World Health Organization has reported the seventh human betacoronavirus (2019-nCoV) that caused a pneumonia epidemic. This new virus is linked to an outbreak of febrile respiratory disease in Wuhan city, Hubei province, China. Some patients with 2019-nCoV developed severe pneumonia, pulmonary edema, ARDS, or multiple organ failure and died. At present, information regarding the epidemiology and clinical features of pneumonia caused by 2019-nCoV is scarce. Monitoring of the total anti-S1 antibody titer can serve as an indication of immune protection for post-COVID-19 patients or for vaccinated individuals.

This serological assay helps to identify individuals exposed to the SARS-CoV-2 virus and its British, Brazilian, South African variants; it will be useful to decide the extent of containment measures in COVID-19 disease.

3. PRINCIPLE OF THE METHOD

The Chorus SARS-CoV-2 "NEUTRALIZING" Ab device is ready to use for the detection of total antibodies anti-S1 SARS-CoV-2 virus, in the Chorus TRIO instrument.

The SARS-CoV-2 immunoglobulin (total Ig) anti-S1 SARS-CoV-2 virus assay is based on the principle of analysis by competition. The SARS-CoV-2 anti-S1 antibodies present in the test sample compete with the tracer (peroxidase-conjugated SARS-CoV-2 anti-S1 therapeutic monoclonal antibody) to occupy the binding sites, available in limited numbers, of the antigen fixed on the solid phase. The monoclonal tracer is a human therapeutic antibody directed against the RBD of the S1 subunit of the Spike protein of the SARS-CoV-2 virus. It shows a strong neutralizing activity against the virus and its variants (English, South African and Brazilian). The higher the concentration of antibodies in the test sample, the lower the amount of conjugated antibodies that bind. The unbound components are eliminated by washing and the bound enzymatic activity is evaluated colorimetrically by transformation of a chromogenic substrate. The intensity of the color developed is inversely proportional to the concentration of antibodies in the sample under examination.

The disposable devices contain all the reagents to perform the test when applied on the Chorus TRIO instrument.

The results are expressed in Binding Antibody Units (BAU/ml) calculated with reference to the first international standard WHO 20/136 for anti-SARS-CoV-2.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by CE-marked methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30-minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially

biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.


1. **Discard devices which show the substrate (well 5) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.
The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website
(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)
5. Check that the Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 9: Space for application of bar code label

Position 8: SAMPLE DILUENT SORBENT

Contents: Proteic colution containing phenol 0.05% and Bronidox 0.02%.

Position 7: MICROPLATE WELL

Coated with inactivated native SARS-CoV-2 antigen.

Position 6: MICROPLATE WELL

Coated with inactivated native SARS-CoV-2 antigen.

Position 5: TMB SUBSTRATE

Contents: Tetramethylbenzidine and H₂O₂ stabilized in buffer.

Position 4: STOP SOLUTION

Contents: 0.3 M sulfuric acid solution.

Position 3: SAMPLE DILUENT

Contents: Protein solution in phosphate buffer containing 0.09% sodium azide and dye.

Position 2: CONJUGATE

Contents: Human monoclonal antibodies (SARS-CoV-2 anti-S1 therapeutic monoclonal antibody) labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

in which the sample is transferred.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 1 ml

Contents: Protein solution containing specific antibodies capable of binding the antigen present on the microplate, and preservative Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 1 ml

Contents: Protein solution containing specific antibodies capable of binding the antigen present on the microplate, and preservative Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 – 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents are stable after opening and/or preparation up to the expiration date, when stored correctly at 2-8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.
Heat-inactivation can rise to erroneous results.
The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

The test Chorus SARS-CoV-2 “NEUTRALIZING” Ab (REF. 81408) cannot be performed simultaneously with the following tests:

CHORUS BORDETELLA PERTUSSIS TOXIN IgA	REF 81274
CHORUS BORDETELLA PERTUSSIS TOXIN IgG	REF 81272
CHORUS BORRELIA IgG	REF 81046
CHORUS BORRELIA IgG CSF	REF 81048
CHORUS BORRELIA IgM	REF 81047
CHORUS CALPROTECTIN	REF 86600
CHORUS CALPROTECTIN K	REF 86601
CHORUS CHIKUNGUNYA IgG	REF 81133
CHORUS CHIKUNGUNYA IgM	REF 81134
CHORUS DIPHTHERIA IgG	REF 81266
CHORUS HSV 1 SCREEN	REF 81018
CHORUS HSV 2 SCREEN	REF 81019
CHORUS LEGIONELLA URINARY ANTIGEN	REF 81301
CHORUS POLIOVIRUS IgG	REF 81264
CHORUS PSA	REF 86550
CHORUS Q FEVER PHASE II IgG	REF 81171
CHORUS Q FEVER PHASE II IgM	REF 81172
CHORUS SARS-CoV-2 SCREEN SERUM	REF 81406
CHORUS SYPHILIS SCREEN RECOMBINANT	REF 81100
CHORUS TETANUS IgG	REF 81270
CHORUS TICK-BORNE ENCEPHALITIS VIRUS IgG	REF 81276
CHORUS TICK-BORNE ENCEPHALITIS VIRUS IgM	REF 81277
CHORUS TOXOCARA IgG	REF 81280

Carry out the calibration and validation of the test at each analysis session, using one device for the Calibrator and one device for the Positive Control.

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 100 µl of undiluted test serum in well no. 1 of each device
4. Place the devices in the Chorus TRIO instrument. Perform the calibration and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the obtained results. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected. If the result of the positive control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus TRIO instrument expresses the result in Binding Antibody Units, BAU/ml, (first international standard WHO 20/136 for anti-SARS-CoV-2), calculated on the basis of a batch-dependent curve stored in the instrument.

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is > 50.0 BAU/ml

NEGATIVE: when the result is < 20.0 BAU/ml

EQUIVOCAL: for all values between 20.0-50.0 BAU/ml

If the result is equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. CALIBRATION RANGE

Calibration range 20-1500 BAU/ml

For samples > 1500 BAU/ml retest the diluted sample in the Negative Control/Sample Diluent (REF 83607 - not supplied with the kit).

13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, are negative, with values < 20.0 BAU/ml.

14. CROSS-REACTIONS

102 samples, positive to *Mycoplasma pneumoniae* (1), *Chlamydia pneumoniae* (3), Cytomegalovirus (4), Adenovirus (4), Influenza B (1), Influenza A (1), RSV (2), Rheumatoid Factor (4), Parainfluenza Virus (4), Bordetella pertussis (5), ANA (10), Epstein-Barr Virus (5), Influenza A H1N1 (5), Dengue Virus (5), Coronavirus HKU1 (5), Coronavirus OC43 (5), Coronavirus NL63 (5), Coronavirus 229E (5) e HIV were tested.

No significant cross-reactions were found.

15. ANALYTICAL SPECIFICITY

4 samples (2 Negative and 2 Positive) were spiked with the following potentially interfering factors and then tested:

Rheumatoid Factor – 110 IU/mL
Bilirubin – 18 mg/dL
Triglycerides – 150 mg/dL
Hemoglobin – 15 mg/mL

The presence in the serum sample of hemoglobin or triglycerides can affect the test result.

16. METHOD COMPARISON

In an experimentation 1106 samples have been tested with Diesse kit and with PRNT seroneutralization:

- 651 pre-pandemic samples
- 455 pandemic samples

Data are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	242	2	244
	-	1	861	862
	Total	243	863	1106

Percent Positive Agreement (~Diagnostic Sensitivity):

99.6% CI_{95%}: 97.7 - 99.9

Percent Negative Agreement: (~Diagnostic Specificity):

99.8% CI_{95%}: 99.2 - 99.9

The agreement between the two methods is excellent with a Cohen's Kappa of 1.0.

17. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run Precision	
	Mean (BAU/ml)	CV%	Mean (BAU/ml)	CV%
1	311.6	6.3	295.3	5.8
2	59.7	14.9	68.2	14.1
3	544.5	10.6	410.7	12.1
4	542.5	14.4	581.7	6.9

Sample	Precision between batches		Precision between instruments	
	Mean (BAU/ml)	CV%	Mean (BAU/ml)	CV%
1	260.9	6.3	296.5	5.3
2	59.7	12.7	577.4	14.7
3	141.3	3.9	426.8	13.4
4	107.3	13.6	577.4	14.7

18. BIBLIOGRAPHY

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
2. Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
5. Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004
6. Andreano et al. "SARS-CoV-2 escape in vitro from a highly neutralizing COVID-19 convalescent plasma" bioRxiv. Preprint. 2020 Dec 28
7. <https://doi.org/10.1016/j.jinf.2020.12.024> (British variant)
8. <https://doi.org/10.1016/j.jinf.2021.01.025> (Brazilian variant)
9. <https://doi.org/10.1016/j.jinf.2021.01.007> (South African variant)

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18. BIBLIOGRAPHY

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
2. Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
5. Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004
6. Andreano et al. "SARS-CoV-2 escape in vitro from a highly neutralizing COVID-19 convalescent plasma" bioRxiv. Preprint. 2020 Dec 28
7. <https://doi.org/10.1016/j.jinf.2020.12.024> (British variant)
8. <https://doi.org/10.1016/j.jinf.2021.01.025>(Brazilian variant)
9. <https://doi.org/10.1016/j.jinf.2021.01.007>(South African variant)

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18. BIBLIOGRAPHY

- Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
- Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
- Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
- Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
- Woo et al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004
- Andreano et al. "SARS-CoV-2 escape in vitro from a highly neutralizing COVID-19 convalescent plasma" bioRxiv. Preprint. 2020 Dec 28
- <https://doi.org/10.1016/j.jinf.2020.12.024> (British variant)
- <https://doi.org/10.1016/j.jinf.2021.01.025> (Brazilian variant)
- <https://doi.org/10.1016/j.jinf.2021.01.007> (South African variant)

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18. BIBLIOGRAPHY

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
2. Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
5. Woo et al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004
6. Andreano et al. "SARS-CoV-2 escape in vitro from a highly neutralizing COVID-19 convalescent plasma" bioRxiv. Preprint. 2020 Dec 28
7. <https://doi.org/10.1016/j.jinf.2020.12.024> (British variant)

8. <https://doi.org/10.1016/j.jinf.2021.01.025>(Brazilian variant)
9. <https://doi.org/10.1016/j.jinf.2021.01.007>(South African variant)

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CHORUS BORDETELLA PERTUSSIS TOXIN IgA	REF. 81274
CHORUS BORDETELLA PERTUSSIS TOXIN IgG	REF. 81272
CHORUS BORRELIA IgG	REF. 81046
CHORUS BORRELIA IgG CSF	REF. 81048
CHORUS BORRELIA IgM	REF. 81047
CHORUS CALPROTECTIN	REF. 86600
CHORUS CALPROTECTIN K	REF. 86601
CHORUS CHIKUNGUNYA IgG	REF. 81133
CHORUS CHIKUNGUNYA IgM	REF. 81134
CHORUS DIPHTHERIA IgG	REF. 81266
CHORUS HSV 1 SCREEN	REF. 81018
CHORUS HSV 2 SCREEN	REF. 81019
CHORUS LEGIONELLA URINARY ANTIGEN	REF. 81301
CHORUS POLIOVIRUS IgG	REF. 81264
CHORUS PSA	REF. 86550
CHORUS Q FEVER PHASE II IgG	REF. 81171
CHORUS Q FEVER PHASE II IgM	REF. 81172

CHORUS SARS-CoV-2 SCREEN SERUM	REF. 81406
CHORUS SYPHILIS SCREEN RECOMBINANT	REF. 81100
CHORUS TETANUS IgG	REF. 81270
CHORUS TICK-BORNE ENCEPHALITIS VIRUS IgG	REF. 81276
CHORUS TICK-BORNE ENCEPHALITIS VIRUS IgM	REF. 81277
CHORUS TOXOCARA IgG	REF. 81280

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18. BIBLIOGRAPHY

1. Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
2. Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020
3. Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
4. Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
5. Woo et al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004
6. Andreano et al. "SARS-CoV-2 escape in vitro from a highly neutralizing COVID-19 convalescent plasma" bioRxiv. Preprint. 2020 Dec 28
7. <https://doi.org/10.1016/j.jinf.2020.12.024> (British variant)
8. <https://doi.org/10.1016/j.jinf.2021.01.025> (Brazilian variant)
9. <https://doi.org/10.1016/j.jinf.2021.01.007> (South African variant)

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18. BIBLIOGRAPHY

- Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" *Science China Life sciences*; January 2020
- Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". *The Lancet*, Vol 395 (507-513) February 15, 2020
- Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". *The Lancet*, Vol 395 (497-506) February 15, 2020
- Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" *Journal of clinical microbiology-vol 43 n8 (p.3718-37296)*; august 2005
- Woo et al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" *Journal of clinical microbiology - vol 42 n5 (2306-2309)*; May 2004
- Andreano et al. "SARS-CoV-2 escape in vitro from a highly neutralizing COVID-19 convalescent plasma" *bioRxiv*. Preprint. 2020 Dec 28
- <https://doi.org/10.1016/j.jinf.2020.12.024> (British variant)
- <https://doi.org/10.1016/j.jinf.2021.01.025> (Brazilian variant)
- <https://doi.org/10.1016/j.jinf.2021.01.007> (South African variant)

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18. BIBLIOGRAPHY

- Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" *Science China Life sciences*; January 2020
- Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study". *The Lancet*, Vol 395 (507-513) February 15, 2020
- Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". *The Lancet*, Vol 395 (497-506) February 15, 2020
- Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" *Journal of clinical microbiology*-vol 43 n8 (p.3718-37296); august 2005
- Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" *Journal of clinical microbiology* - vol 42 n5 (2306-2309); May 2004
- Andreano et al. "SARS-CoV-2 escape in vitro from a highly neutralizing COVID-19 convalescent plasma" *bioRxiv*. Preprint. 2020 Dec 28
- <https://doi.org/10.1016/j.jinf.2020.12.024> (British variant)
- <https://doi.org/10.1016/j.jinf.2021.01.025> (Brazilian variant)
- <https://doi.org/10.1016/j.jinf.2021.01.007> (South African variant)

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Percent Positive Agreement (~Sensibilitatea diagnosticului):

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18. BIBLIOGRAPHY

- Xintian et al. "Evolution of the novel coronavirus from ongoing Wuhan outbreak and modeling of its spike protein for risk of human transmission" Science China Life sciences; January 2020
- Chen et al "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan,

China: a descriptive study". The Lancet, Vol 395 (507-513) February 15, 2020

- Huang et al "Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China". The Lancet, Vol 395 (497-506) February 15, 2020
- Yuxian et al. "Identification of immunodominant epitopes on the membrane protein of Severe Acute Respiratory Syndrome-Associated Coronavirus" Journal of clinical microbiology-vol 43 n8 (p.3718-37296); august 2005
- Woo at al. "Detection of specific antibodies to Severe Acute Respiratory Syndrome (SARS) Coronavirus Nucleocapsid protein for serodiagnosis of SARS Coronavirus Pneumonia" Journal of clinical microbiology - vol 42 n5 (2306-2309); May 2004
- Andreano et al. "SARS-CoV-2 escape in vitro from a highly neutralizing COVID-19 convalescent plasma" bioRxiv. Preprint. 2020 Dec 28
- <https://doi.org/10.1016/j.jinf.2020.12.024> (British variant)
- <https://doi.org/10.1016/j.jinf.2021.01.025> (Brazilian variant)
- <https://doi.org/10.1016/j.jinf.2021.01.007> (South African variant)

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