CHORUS

LKM-1

REF 86080

REF 86080/12

DIESSE Diagnostica Senese
S.p.A.

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INSTRUCTIONS FOR USE

CHORUS LKM-1

For the semi-quantitative determination of anti-liver-kidney microsomal IgG antibodies (LKM)

Only for *in vitro* diagnostic use

1. USE

Enzyme immunoassay for the semi-quantitative determination of anti-liver-kidney microsomal IgG class antibodies (LKM) in human serum with a single-use device applied to Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Autoimmune hepatitis (AIH) is a progressive liver disease of unknown origin, which responds well to immunosuppressive therapy, but which, if not treated, has a negative prognosis. An accurate and early diagnosis therefore plays a fundamental role. AIH is characterised by the histological picture of periportal hepatitis in the absence of viral markers, by hypergammaglobulinemia and in most patients, by the presence of autoantibodies in the serum. Anti-nuclear antibodies (ANA), smooth muscle antibodies (SMA), anti-liver-kidney microsomal (LKM, liver-kidney microsomes) antibodies, and antibodies against soluble liver antigen (SLA) are marker autoantibodies for AIH. 52% of AIH patients are ANA and/or SMA positive, 20% SLA positive and 3% have anti-LKM-1 antibodies. All these antibodies are very important for diagnosing an AIH, however, so far only Anti-Soluble-Liver-Antigens have been highly specific for AIH. ANAs/SMA are found in 10-15% of patients with viral hepatitis or other immune-mediated diseases. LKM-1s are also found in patients with hepatitis C.

Anti-LKM antibodies can be divided into three types in relation to their respective target antigens. Anti-LKM-1 antibodies are directed against cytochrome p450 IID6, a 50 kDa cytoplasmic protein found in hepatocytes and renal proximal tubular cells. LKM-2 antibodies are associated with ticrynafen-induced hepatitis. The target antigen of LKM-2 is the cytochrome p450 IIC9, a cytochrome p450 isoenzyme, which catalyses the metabolic oxidation of the drug. Anti-LKM-3 antibodies are associated with chronic hepatitis D and their target antigen is UDP-glucuronosyltransferase.

AIH associated with anti-LKM-1s occurs predominantly in childhood, especially in girls aged 2 to 14 years. The determination of anti-LKM-1 antibodies is therefore of paramount importance in paediatrics.

3. PRINCIPLE OF THE METHOD

The Chorus LKM-1 device is ready for use to determine anti-liver-kidney microsomal IgG class antibodies (LKM) in Chorus/Chorus TRIO instruments.

The test is based on the ELISA (enzyme-linked immunosorbent assay) principle. The antigen is bound to the solid phase. Specific immunoglobulins bind to the antigen following incubation with diluted human serum. After washing to remove unreacted proteins, incubation is performed with the conjugated compound consisting of human anti-immunoglobulin antibodies conjugated to horseradish peroxidase. The unbound conjugate is removed and the substrate for the peroxidase is added. The colour that develops is proportional to the concentration of specific antibodies in the serum under examination.

The single-use devices contain all the reagents to perform the test in the Chorus/Chorus TRIO instruments.

The results are expressed in Arbitrary Units (AU/ml).

4. PRECAUTIONS

ONLY FOR *IN VITRO* DIAGNOSTIC USE.

This kit contains materials of human origin which have been tested and found negative for both HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. Since no diagnostic test can offer a complete guarantee on the absence of infectious agents, any material of human origin must be considered potentially infected. All reagents and samples must be handled in accordance with the safety standards normally adopted in the laboratory.

Disposal of residues: the serum samples, calibrators and the strips used must be treated as infected residues, therefore disposed of in accordance with the provisions of current laws.

Personal safety warnings

1. Do not pipette with your mouth.
2. Use single-use gloves and eye protection when handling samples.
3. Wash your hands thoroughly after inserting the devices into the Chorus/Chorus TRIO instrument.
4. Refer to the Safety Data Sheet (available on request) for the safety characteristics of the reagents contained in the kit.
5. Neutralised acids and other liquid waste must be disinfected by adding sodium hypochlorite in a volume sufficient to obtain a final concentration of at least 1%. Exposure to 1% hypochlorite sodium for 30 minutes should be sufficient to ensure effective disinfection.
6. Any spillage of potentially infected materials must be removed immediately with absorbent paper and the contaminated area must be decontaminated, for example with 1% hypochlorite sodium, before continuing work. If an acid is present, sodium hypochlorite must not be used before the area has been dried. All materials used to decontaminate any accidental spillage, including gloves, must be discarded as potentially infected waste. Do not autoclave materials containing hypochlorite sodium.

Analytical warnings

Before use, bring the devices to be used to room temperature (18-30°C) and use within 60 minutes.

1. **Discard the devices with blue coloured substrate (well 4).**
2. When adding the sample to the well, check that it is perfectly distributed at the bottom.
3. Check the actual presence of reagents in the device and the integrity of the device itself. Do not use devices that lack any reagents and/or that have any foreign bodies in the reaction well when visually inspected.
4. The devices must be used in conjunction with the Chorus/Chorus TRIO instrument, strictly following the Instructions for Use and the User Manual of the instrument.
5. Check that the Chorus/Chorus TRIO instrument is set correctly (see User Manual).
6. Do not alter the barcode placed on the handle of the device in order to allow it to be read correctly by the instrument.
7. Avoid the use of self-defrosting freezers for sample storage.
8. Defective barcodes can be manually entered into the instrument (see User Manual).
9. Do not expose the devices to strong lighting or hypochlorite vapours during storage and use.
10. The use of heavily haemolysed, lipemic, icteric, non-fully coagulated serum samples or of samples with microbial contamination may be a source of errors.
11. Do not use the device after its expiry date
12. **Check that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004)**

5. COMPOSITION OF THE KIT AND PREPARATION OF REAGENTS

The kit is sufficient for 36 tests (REF 86080).


The kit is sufficient for 12 tests (REF 86080/12).

DD DEVICES

6 packs of 6 devices each (REF 86080).

2 packs of 6 devices each (REF 86080/12).

Description:

Position 8: Space available for barcode label

Position 7: Empty

Position 6: MICROPLATE WELL

Sensitised with recombinant cytochrome p450 IID6.

Position 5: MICROPLATE WELL

Not sensitised.

Position 4: TMB SUBSTRATE

Content: Tetramethylbenzidine stabilised in citrate buffer

Position 3: DILUENT FOR SAMPLES

Content: saline protein solution containing Proclin (0.1%).

Position 2: CONJUGATE

Content: Peroxidase-labelled monoclonal anti-human IgG antibodies.

Position 1: EMPTY WELL

Where the user must dispense the undiluted serum.

Usage: bring a bag to room temperature, open the bag, pick up the necessary devices; store the others in the bag containing the silica gel, let the air out and **seal** by pressing on the closure. Store at 2/8°C.

CALIBRATOR CALIBRATOR 1 x 0.175 mL

Content: Protein solution containing specific antibodies capable of binding the microplate and preservative. Liquid, ready to use.

CONTROL + POSITIVE CONTROL 1 x 0.425 mL

Content: Protein solution containing specific antibodies capable of binding the microplate and preservative. Liquid, ready to use.

OTHER MATERIAL REQUIRED, BUT NOT SUPPLIED

- WASHING BUFFER AUTOIMMUNITY (REF 86004)
- CLEANING SOLUTION 2000 (REF 83609)
- SANITISING SOLUTION (REF 83604 - 83608)
- CHORUS NEGATIVE CONTROL/DILUENT SAMPLE (REF 83607)
- Chorus/Chorus TRIO instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test tubes, etc.
- Micropipettes capable of accurately withdrawing volumes of 50-200 µl.
- Single-use gloves
- 5% sodium hypochlorite solution
- Containers for the collection of potentially infected materials

6. METHOD OF STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the event of an incorrect storage temperature, the calibration must be repeated and the correctness of the result checked using the control serum (see chapter 9: Test validation).

The expiry date is printed on each component and on the outer label of the package.

Reagents have limited stability after opening and/or preparation:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SAMPLE TYPE AND STORAGE

The sample type is serum obtained from blood collected by normal venous blood sampling and handled as required by standard laboratory procedures.

The consequences of using other body fluids are unknown.

Fresh serum can be kept for 4 days at 2/8°C; for longer storage periods, freeze at -20°C.

The sample can undergo up to a maximum of 3 thaws.
Avoid the use of self-defrosting freezers for sample storage. After thawing, carefully shake the sample before dosing.
Heat inactivation may give erroneous results. Sample quality can be seriously affected by microbial contamination which can lead to erroneous results.

8. PROCEDURE

1. Open the bag (side containing the snap closure), take the number of devices needed to perform the tests and store the others by closing the bag after letting the air out.
2. Visually check the status of the device according to the instructions in chapter 4 Analytical Warnings.
3. Add 50µl of undiluted serum to be analysed to well no. 1, at each batch change use a calibrator device.
4. Introduce the devices on the Chorus/Chorus TRIO instrument. Perform the calibration (if required) and the test as reported in the Instruction Manual of the instrument.

9. TEST VALIDATION

Use the positive control serum to verify the correctness of the result obtained, processing it as indicated in the User Manual of the instrument. If the instrument indicates that the control serum has a value falling outside the acceptability limit, the calibration must be performed again. Previous results will be corrected automatically.

If the control serum result continues to be outside the acceptability range, contact Scientific Support.

Phone: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE TEST

The Chorus/Chorus TRIO instrument provides the result in Arbitrary Units (AU/ml) calculated on the basis of a batch-dependent graph stored in the instrument.

The test on the serum under examination can be interpreted as follows:

POSITIVE: when the result is > 18 AU/ml

NEGATIVE: when the result is <12 AU/ml

DOUBTFUL/EQUIVOCAL: when the result is between 12 and 18 AU/ml

In case of doubtful/equivocal result, repeat the test. If the result remains doubtful/equivocal, repeat the sampling.

11. TEST LIMITATIONS

All values obtained require a careful interpretation, without ever disregarding other indicators relating to the same patient.

The test, in fact, cannot be used on its own for a clinical diagnosis and the result of the test must be evaluated together with data from the patient's medical history and/or from other diagnostic investigations.

12. CALIBRATION RANGE

Calibration range 3.0-300 AU/ml.

For samples > 300 AU/ml repeat the test, pre-diluting the sample with the Negative Control/Sample Diluent (PF83607- not supplied with the kit).

13. ANALYTICAL SPECIFICITY

2 samples were tested (1 Negative, 1 Positive) to which the following interfering agents were added:

Rheumatoid Factor (44 IU/ml - 220 IU/ml)

Bilirubin (4.5 mg/dl - 45 mg/dl)

Triglycerides (10 mg/dl - 250 mg/dl)

Haemoglobin (5 mg/ml - 30 mg/ml)

The presence of the above interfering substances in the test serum does not alter the test result.

14. CROSS-REACTIVITY

15 samples, positive for AMA-M2, SMA and HCV were tested.

No significant cross-reactions were detected.

15. COMPARISON STUDIES

In one trial, 263 samples were analysed with a Diesse kit and another commercially available kit.

The trial data are outlined below:

		Reference		
		+	-	Total
Diesse	+	4	0	4
	-	0	259	259
	Total	4	259	263

Positive Percent Agreement (~Diagnostic Sensitivity):

100.0% CI_{95%}: 51.0-99.4

Negative Percent Agreement: (~Diagnostic Specificity):

100.0% CI_{95%}: 98.5-100.0

The degree of agreement between the two methods is excellent with a value of K (Cohen's coefficient) of 1.

16. ACCURACY AND REPEATABILITY

Sample	Within the session		Between sessions	
	Mean (AU/ml)	CV%	Mean (AU/ml)	CV%
1	3.8	12.4	3.0	14.3
2	37.1	10.8	33.9	8.7
3	59.5	11.2	61.3	13.5

Sample	Between batches		Between instruments	
	Mean (AU/ml)	CV%	Mean (AU/ml)	CV%
1	5.0	3.4	3.4	14.1
2	39.7	8.0	30.2	13.9
3	72.0	11.1	59.3	14.2

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