**CHORUS
ASPERGILLUS
GALACTOMANNAN
Ag**

REF 81178

DIESSE

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INSTRUCTIONS FOR USE

CHORUS ASPERGILLUS GALACTOMANNAN Ag

**For the qualitative determination of the
Aspergillus Galactomannan antigen**

Only for in vitro diagnostic use

1. USE

Enzyme immunoassay method for the qualitative determination of the Galactomannan antigen of *Aspergillus* in human serum and bronchial alveolar lavage fluid (BAL) with single-use device applied to Chorus TRIO instruments.

2. INTRODUCTION

Aspergillosis is an opportunistic infection caused by the inhalation of *Aspergillus* fungal spores, commonly present in the environment.

Invasive forms, on the rise over the past 10 years, represent the most serious infections.

Aspergillus tends to infect open spaces, such as lung cavities caused by previous diseases (e.g. bronchiectasis, cancer, tuberculosis), sinuses, or external auditory canals. In particular, *A. fumigatus* is the most frequent cause of invasive lung disease.

Traditional methods of diagnosis often rely on non-specific diagnostic or radiological tests (clinical symptoms, CT scans, chest X-rays, etc.). Since *Aspergillus* is widespread in the environment, the positivity of sputum cultures may be due to environmental contamination.

The test for galactomannan antigen in serum and bronchial alveolar lavage fluid (BAL) appears to be a method capable of supporting the diagnosis of invasive Aspergillosis. In invasive pulmonary Aspergillosis, the Galactomannan test in BAL is often the only option for patients with thrombocytopenia, for whom biopsy is contraindicated.

3. PRINCIPLE OF THE METHOD

The test is based on the ELISA (enzyme-linked immunosorbent assay) principle.

The monoclonal antibody specific for the antigen is bound to the solid phase. During incubation, the antigen present in the sample binds to the monoclonal antibody present in the solid phase and to the one conjugated with horseradish peroxidase, forming a "sandwich".

After washing to remove unreacted proteins, the peroxidase substrate is added. The enzymatic reaction is subsequently stopped by adding the Stop Solution which causes the solution to turn yellow.

The single-use devices contain all the reagents to perform the test when applied to the Chorus TRIO instruments.

The results are expressed in Index (OD sample/OD cut-off)

4. PRECAUTIONS

ONLY FOR IN VITRO DIAGNOSTIC USE.

This kit contains materials of human origin which have been tested and found negative for both HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. Since no diagnostic test can offer a complete guarantee on the absence of infectious agents, any material of human origin must be considered potentially infected. All reagents and samples must be handled in accordance with the safety standards normally adopted in the laboratory.

Disposal of residues: the samples, calibrators and the strips used must be treated as infected residues, therefore disposed of in accordance with the provisions of current laws.

Personal safety warnings

1. Do not pipette with your mouth.
 2. Use single-use gloves and eye protection when handling samples.
 3. Wash your hands thoroughly after inserting the devices into the Chorus TRIO instrument.
 4. Refer to the Safety Data Sheet (available on request) for the safety characteristics of the reagents contained in the kit.
 5. Neutralised acids and other liquid waste must be disinfected by adding sodium hypochlorite in a volume sufficient to obtain a final concentration of at least 1%. Exposure to 1% hypochlorite sodium for 30 minutes should be sufficient to ensure effective disinfection.
 6. Any spillage of potentially infected materials must be removed immediately with absorbent paper and the contaminated area must be decontaminated, for example with 1% hypochlorite sodium, before continuing work. If an acid is present, sodium hypochlorite must not be used before the area has been dried.
- All materials used to decontaminate any accidental spillage, including gloves, must be discarded as potentially infected waste.
- Do not autoclave materials containing hypochlorite sodium.

Analytical warnings

Before use, bring the bag containing the devices to room temperature (18-30°C) for at least 30 minutes.

Use the devices within 60 minutes.

1. **Discard the devices with blue coloured substrate (well 4).**
2. When adding the sample to the well, check that it is perfectly distributed at the bottom.
3. Check the actual presence of reagents in the device and the integrity of the device itself. Do not use devices that lack any reagents and/or that have any foreign bodies in the reaction well when visually inspected.
4. The devices must be used in conjunction with the Chorus TRIO instrument, strictly following the Instructions for Use and the User Manual of the instrument.

Use of the kit is only possible with CHORUS TRIO with an updated version of software. Make sure that the

software installed in the instrument matches or has a Release (date) above the one shown in the table published on the Diesse website

(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)


5. Check that the Chorus TRIO instrument is set correctly (see User Manual).
6. Do not alter the barcode placed on the handle of the device in order to allow it to be read correctly by the instrument.
7. Avoid the use of self-defrosting freezers for sample storage.
8. Defective barcodes can be manually entered into the instrument (see User Manual).
9. Do not expose the devices to strong lighting or hypochlorite vapours during storage and use.
10. The use of heavily haemolysed, lipemic, icteric, non-fully coagulated serum samples or of samples with microbial contamination may be a source of errors.
11. Do not use the device after its expiry date
12. **Check that the instrument is connected to the Washing Buffer (Ref. 83606)**

5. COMPOSITION OF THE KIT AND PREPARATION OF REAGENTS

The kit is sufficient for 36 tests

DD DEVICES 6 packs of 6 devices each

Description:

Position 8: Space available for barcode label

Position 7: Empty

Position 6: MICROPLATE WELL

Sensitised with anti-galactomannan monoclonal antibody

Position 5: MICROPLATE WELL

Not sensitised.

Position 4: TMB SUBSTRATE

Content: Tetramethylbenzidine stabilised in citrate buffer

Position 3: STOP SOLUTION

Content: Sulphuric acid 0.3 M solution

Position 2: CONJUGATE

Content: Peroxidase-labelled anti-galactomannan monoclonal antibody, in phosphate buffer solution containing 0.05% phenol and 0.02% Bronidox

Position 1: EMPTY WELL

Where the sample is transferred.

Usage: bring a bag to room temperature, open the bag, pick up the necessary devices; store the others in the bag containing the silica gel, let the air out and **seal** by pressing on the closure. Store at 2/8°C.

CALIBRATOR CALIBRATOR 1 x 0.950 ml

Content: Human serum containing galactomannan and preservative. Liquid, to be treated as the test samples.

CONTROL + POSITIVE CONTROL 2 x 1.250 ml

Content: Human serum containing galactomannan and preservative. Liquid, to be treated as the test samples.

SAMPLE DILUENT SAMPLE DILUENT 1 x 7 ml

Content: EDTA solution. Liquid, ready to use

OTHER MATERIAL REQUIRED, BUT NOT SUPPLIED:

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITISING SOLUTION **REF** 83604
- Chorus TRIO instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test tubes, etc.
- Micropipettes capable of accurately withdrawing volumes of 50-200µl.
- Single-use gloves
- 5% sodium hypochlorite solution
- Containers for the collection of potentially infected materials
- Vortex
- 1.5 ml centrifuge tubes
- Block heater up to 100°C or water bath at 100°C.
- Bench-top centrifuge for 1.5 ml tubes able to obtain 10,000xg

6. METHOD OF STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the event of an incorrect storage temperature, the calibration must be repeated and the correctness of the result checked using the positive control (see chapter 9: Test validation).

The expiry date is printed on each component and on the outer label of the package.

Reagents have limited stability after opening and/or preparation:

	After opening	After treatment
DEVICES	8 weeks at 2/8°C	-
CALIBRATOR	8 weeks at 2/8°C	5 days
POSITIVE CONTROL	8 weeks at 2/8°C	5 days
SAMPLE DILUENT	8 weeks at 2/8°C	-

7. SAMPLE TYPE AND STORAGE

The sample type is serum obtained from blood collected by normal venous blood sampling and handled as required by standard laboratory procedures, or bronchoalveolar lavage fluid (BAL) which must be collected in a sterile saline solution.

Fresh serum can be kept for 4 days at 2/8°C; for longer storage periods (at least up to 3 months), freeze at -20°C. The sample can undergo up to a maximum of 3 thaws.

Avoid the use of self-defrosting freezers for sample storage. After thawing, carefully shake the sample before dosing.

Sample quality can be seriously affected by microbial contamination or fungal spores, which can lead to erroneous results.

8. PROCEDURE

The CHORUS Aspergillus Galactomannan Ag test (REF 81178) can be performed simultaneously ONLY with the following tests:

CHORUS CYTOMEGALOVIRUS IgG AVIDITY	REF 81092
CHORUS RUBELLA IgG AVIDITY	REF 81095
CHORUS TOXOPLASMA IgG AVIDITY	REF 81098

SAMPLE PREPARATION

The Positive Control and the Calibrator must be processed as the serum/fluid BAL samples:

- Pipette 300 µl of Calibrator, Positive Control and each sample into 1.5 ml centrifuge tubes (not supplied)
- Add 100 µl of Sample Diluent to each tube.
- Mix well
- Heat the tubes for 3 minutes in a block heater at 100°C or 3 minutes in a water bath at 100°C.
- Centrifuge the tubes at 10,000xg for 10 minutes
- Take the supernatant and use it as a sample.

It is possible to store the treated sample for 48 hours at 2-8°C.

PROCEDURE

1. Open the bag (side containing the snap closure), take the number of devices needed to perform the tests and store the others by closing the bag after letting the air out.
2. Visually check the status of the device according to the instructions in chapter 4 Analytical Warnings.
3. Add 75 µl of sample to be analysed to well no. 1. At each batch change use a calibrator device.
4. Introduce the devices on the Chorus TRIO instrument. Perform the calibration (if required) and the test as reported in the User Manual of the instrument.

9. TEST VALIDATION

Use the positive control to verify the correctness of the result obtained, processing it as indicated in the User Manual of the instrument. If the instrument indicates that the positive control has a value falling outside the acceptability limit, the calibration must be performed again. Previous results will be corrected automatically.

If the positive control result continues to be outside the acceptability range, contact Scientific Support.

Phone: 0039 0577 319554

Fax: 0039 0577 366605

email: scientificsupport@diesse.it

10. INTERPRETATION OF THE TEST

The Chorus TRIO instrument provides the result in Index (sample OD/cut-off OD).

The test on the sample under examination can be interpreted as follows:

POSITIVE: when the result is > 1.1

NEGATIVE: when the result is < 0.9

DOUBTFUL/EQUIVOCAL: when the result is between 0.9 and 1.1

In case of doubtful/equivocal result, repeat the test. If the result remains doubtful/equivocal, repeat the sampling.

11. TEST LIMITATIONS

All values obtained require a careful interpretation, without ever disregarding other indicators relating to the same patient.

The test, in fact, cannot be used on its own for a clinical diagnosis and the result obtained must always be evaluated together with data from the patient's medical history and/or from other diagnostic investigations.

12. ANALYTICAL SPECIFICITY

4 samples were tested (1 Negative, 1 Cut-Off and 2 Positive ones) to which the following interfering agents were added:

Rheumatoid factor (44 - 110 IU/ml)
 Bilirubin (4.5 mg/dl - 18 mg/dl)
 Triglycerides (10 mg/dl - 150 mg/dl)
 Haemoglobin (5 mg/ml - 15 mg/ml)
 Ibuprofen (2.19 mg/dl)
 Acetaminophen (1.56 mg/dl)
 Oxymetazoline (1.25 mg/dl)

The presence in the serum under examination of the interfering substances listed above (with the exception of haemoglobin at the level of 15 mg/ml) does not alter the test result.

13. CROSS-REACTIVITY

81 samples, positive for Cytomegalovirus (11), Toxoplasma (10), Treponema (10), Rubella (11), Mycoplasma (4), SARS-CoV-2 (10), ANA (25) were tested.

In addition, bacterial suspensions in saline solution of *Candida krusei*, *Candida albicans*, *Candida tropicalis*, *Candida glabrata*, *Saccharomyces cerevisiae*, *Cryptococcus neoformans* and SARS-CoV-2 viral antigen were also tested.

No cross-reactions were detected.

14. COMPARISON STUDIES

In one trial, 305 samples were analysed with Diesse kit and with another commercially available kit.

The trial data are outlined below:

		Reference		
		+	-	Total
Diesse	+	85	11	96
	-	6	203	209
	Total	91	214	305

Positive Percent Agreement (~Diagnostic Sensitivity):

93.4% CI_{95%}: 86.3 - 96.9

Negative Percent Agreement (~Diagnostic Specificity):

94.4% CI_{95%}: 91.0 - 97.1

Positive Predictive Value: 88.5 CI_{95%}: 84.9 - 92.1

Negative Predictive Value: 97.1 CI_{95%}: 95.2 - 99.0

The degree of agreement between the two methods is excellent with a value of K (Cohen's coefficient) of 0.86

15. PRECISION AND REPEATABILITY

Sample	Within the session		Between sessions	
	Media Index	CV%	Media Index	CV%
1	0.2	5.0	0.3	-
2	0.3	3.3	0.4	12.5
3	2.0	2.5	1.8	2.8
4	2.8	3.6	2.9	4.1

5	0.8	3.8	0.8	5.0
---	-----	-----	-----	-----

Sample	Between batches		Between instruments	
	Media Index	CV%	Media Index	CV%
1	0.3	13.3	0.3	13.3
2	0.3	13.3	0.4	12.5
3	1.8	6.7	1.7	7.1
4	2.8	11.1	2.8	3.9
5	0.8	-	0.8	6.3

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