CHORUS Promonitor Adalimumab

REF 86702

DIESSE

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18. BIBLIOGRAFIA

1. Chen, D. Y. et al "Significant associations of antidrug antibody levels with serum drug trough levels and therapeutic response of adalimumab and etanercept treatment in rheumatoid arthritis" Ann Rheum Dis 2014
2. Bartelds, G. M., et al "High levels of human anti-human antibodies to adalimumab in a patient not responding to adalimumab treatment" Ann.Rheum.Dis. 2006; 65:1249-1250.
3. Bartelds, G. M., et al "Clinical response to adalimumab: relationship to anti-adalimumab antibodies and serum adalimumab concentrations in rheumatoid arthritis" Ann.Rheum.Dis. 2007; 66:921-926.
4. Bartelds, G. M., et al "Development of antidrug antibodies against adalimumab and association with disease activity and treatment failure during long-term follow-up". JAMA 2011; 305:1460-1468.
5. Lecluse, L. L., et al "Extent and clinical consequences of antibody formation against adalimumab in patients with plaque psoriasis" Arch.Dermatol. 2010;146:127-132.

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INSTRUCTIONS FOR USE

CHORUS Promonitor Adalimumab

For the quantitative detection of Adalimumab (ADL) in human serum

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Chorus Promonitor Adalimumab is an immunoassay kit for automated quantitative detection of Adalimumab (ADL) in human serum using a disposable device applied on the Chorus TRIO instrument.

The product, used in conjunction with other clinical and laboratory findings, is useful as an aid in the management of patients treated with Adalimumab (ADL), like patients with rheumatic and gastrointestinal diseases.

It must be used by professional laboratory users only.

2. INTRODUCTION

Adalimumab (ADL) is a recombinant human monoclonal antibody expressed in Chinese Hamster Ovary (CHO) cells used in the treatment of inflammatory diseases. This includes rheumatoid arthritis (RA), polyarticular juvenile idiopathic arthritis (JIA), psoriatic arthritis (PA), ankylosing spondylitis (AS), Crohn's disease (CD) and Psoriasis (PsO). ADL is an immunomodulator biological drug which binds TNF α , a cytokine that plays an important role in the pathogenesis of these diseases, inhibits binding with its receptor, and prevents triggering the inflammatory cascade. As for other biological drugs, ADL can lead to the development of anti-drug antibodies (ADAs). Studies confirm a significant inverse correlation between ADL bioavailability and ADAs presence. The clinical response to ADL closely follows the trough drug levels and the presence of ADAs. It has been demonstrated that immunogenicity of ADL is associated with loss of clinical response in RA, CD and PsO. Furthermore, patients who develop anti-ADL antibodies are more likely not to respond to the therapy and those with therapeutic drug levels are more likely to have an adequate clinical response to ADL therapy in RA and IBD (Inflammatory Bowel Disease).

3. PRINCIPLE OF THE METHOD

The Chorus Promonitor Adalimumab device is ready to use for the detection of Adalimumab (ADL) in human serum, in the Chorus TRIO instrument.

The assay is based on a "sandwich" ELISA method.

The microwell strips are pre-coated with an anti-ADL human monoclonal antibody. Calibrator, positive control and samples are added to separate device.

ADL is able to bind to preimmobilized anti-ADL following incubation with the sample. After washing to remove the unbound sample, a HRP-labelled anti-ADL antibody (conjugate) is added. A second incubation step allows the conjugate to bind to the ADL that has become attached to the microwells. Monoclonal antibodies to ADL do not compete with TNF α for the

binding to ADL. After washing to remove the excess of unbound conjugate, the peroxidase substrate is added. The enzymatic reaction is stopped by adding Stop Solution, which makes the solution turn to yellow. The color which develops is proportional to the amount of drug in the patient sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus TRIO instrument.

The results are expressed in $\mu\text{g/ml}$ calculated in reference to 1st International Standard for Adalimumab (NIBSC code: 16/236).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30-minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the bags with the devices to room temperature (18-30°C) for at least 30 minutes before use.

Use the devices within 60 min.

1. **Discard devices which show the substrate (well 2) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which

are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.

4. The devices are for use with the Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.

The use of the kit is only possible with CHORUS TRIO with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website

(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)


5. Check that the Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 9: Space for application of bar code label

Position 8: DILUTION BUFFER

Contents: Buffer solution containing preservative

Position 7: MICROPLATE WELL

Coated with an anti-ADL human monoclonal antibody

Position 6: MICROPLATE WELL

Coated with an anti-ADL human monoclonal antibody

Position 5: STOP SOLUTION

Contents: 0.3 M sulfuric acid solution

Position 4: CONJUGATE

HRP-labelled anti-ADL antibody and buffer containing preservative

Position 3: DILUTION BUFFER

Contents: Buffer solution containing preservative

Position 2: TMB SUBSTRATE

Contents: Tetramethylbenzidine stabilized in citrate buffer

Position 1: EMPTY WELL

in which the sample is transferred.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CALIBRATOR	CALIBRATOR	1 x 0.400 ml
Contents: Human serum containing ADL, and preservative. Liquid, ready for use.		
CONTROL +	POSITIVE CONTROL	1 x 0.800 ml
Contents: Human serum containing ADL, and preservative. Liquid, ready for use.		

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Chorus TRIO instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

It is strongly recommended to avoid any harsh treatment of the samples.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

The test CHORUS Promonitor Adalimumab (REF. 86702) can be performed simultaneously ONLY with the other CHORUS Promonitor tests.

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.

3. Dispense the following volumes in well no. 1 of each device:

SAMPLE	50 µl/device
CALIBRATOR	130 µl/device
POSITIVE CONTROL	130 µl/device

4. Place the devices in the Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the obtained results. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected. If the result of the positive control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus TRIO instrument expresses the result in µg/ml, calculated on the basis of a batch-dependent curve stored in the instrument.

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is ≥ 0.30 µg/ml

NEGATIVE: when the result is < 0.30 µg/ml

11. LIMITATIONS

1. Trough ADL levels should be measured. The optimum sample collection time is immediately before the ADL injection.
2. This assay measures free ADL. Although it could detect TNF α -bound ADL, this contribution is neglectable, since the concentration of circulating TNF α is in the range of picograms per milliliter.
3. This assay does not measure ADL bound to anti-ADL antibodies (immunocomplexes) which may have been developed by the patient due to drug immunogenicity.
4. This assay is not intended for use as an aid in diagnosis, neither as a companion diagnostic test. The clinician should interpret test results in the context of other variables like sampling window, methodology, anti-ADL antibodies and other clinical findings like disease activity.
5. Although the assay has been confirmed to work with plasma, the assay performance characteristics have not been established for matrices other than serum.
6. This product should only be used by appropriately trained laboratory professional personnel.

12. CALIBRATION RANGE

Calibration range: 0.30-55.0 µg/ml

For samples > 55 µg/ml retest the diluted sample in the Negative Control/Sample Diluent (REF 83607 - not supplied with the kit).

13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, are negative, with values < 0.30 µg/ml.

14. ANALYTICAL SPECIFICITY

3 samples were tested (1 Negative, 1 Low Positive and 1 High Positive) to which the following interfering agents were added:

Bilirubin (4.5 mg/dl - 45 mg/dl)
Triglycerides (250 mg/dl - 1500 mg/dl)
Hemoglobin (2.5 mg/ml - 10 mg/ml)

The presence in the serum sample of the interfering substances listed above (except for hemoglobin) does not affect the test result.

15. CROSS-REACTIONS

15 samples, not containing ADL, but containing various concentrations of IFX 5 samples (11.2 - 48.4 µg/ml), ETN 5 samples (10.0 µg/ml) and HAMA 5 samples (38 - 86 ng/ml) were tested.

No cross-reactions were found.

16. METHOD COMPARISON

In an experimentation 153 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Tot.
Diesse	+	47	0	47
	-	0	106	106
	Tot.	47	106	153

Percent Positive Agreement (~Diagnostic Sensitivity):

100.0% CI_{95%}: 92.4-99.9

Percent Negative Agreement (~Diagnostic Specificity):

100.0% CI_{95%}: 96.5-100.0

Positive Predictive Value (PPV): 100.0% CI_{95%}: 100.0-100.0

Negative Predictive Value (NPV): 100.0% CI_{95%}: 100.0-100.0

The agreement between the two methods is excellent with a Cohen's Kappa of 1.

Out of 153, for 47 positive samples the correlation between Diesse kit and a competitor kit was tested.

Data are summarized in the following table:

Correlation	r	95%CI
Pearson	0.97	0.94-0.98
Spearman	0.97	0.95-0.99

The correlation between the two methods is very high.

17. PRECISION AND REPEATABILITY

Sample	Within run		Between run	
	Mean (µg/ml)	CV%	Mean (µg/ml)	CV%
1	16.1	7.0	16.7	6.5
2	9.4	3.7	9.1	3.8
3	3.9	5.4	3.9	7.4
4	0.47	12.8	0.47	10.6
5	0.30	-	0.30	-
6	0.30	-	0.30	-

Sample	Between lots		Between instruments	
	Mean (µg/ml)	CV%	Mean (µg/ml)	CV%
1	19.7	3.8	17.5	13.9
2	9.3	6.8	8.8	12.8
3	5.1	6.7	4.0	12.8
4	0.50	6.0	0.46	13.0
5	0.30	-	0.30	-
6	0.30	-	0.30	-

18. BIBLIOGRAPHY

1. Chen, D. Y. et al "Significant associations of antidrug antibody levels with serum drug trough levels and

therapeutic response of adalimumab and etanercept treatment in rheumatoid arthritis" Ann Rheum Dis 2014

2. Bartelds, G. M., et al "High levels of human anti-human antibodies to adalimumab in a patient not responding to adalimumab treatment" Ann.Rheum.Dis. 2006; 65:1249-1250.
3. Bartelds, G. M., et al "Clinical response to adalimumab: relationship to anti-adalimumab antibodies and serum adalimumab concentrations in rheumatoid arthritis" Ann.Rheum.Dis. 2007; 66:921-926.
4. Bartelds, G. M., et al "Development of antidrug antibodies against adalimumab and association with disease activity and treatment failure during long-term follow-up". JAMA 2011; 305:1460-1468.
5. Lecluse, L. L., et al "Extent and clinical consequences of antibody formation against adalimumab in patients with plaque psoriasis" Arch.Dermatol. 2010,146:127-132.

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18. BIBLIOGRAFIE

1. Chen, D. Y. et al "Significant associations of antidrug antibody levels with serum drug trough levels and therapeutic response of adalimumab and etanercept treatment in rheumatoid arthritis" Ann Rheum Dis 2014
2. Bartelds, G. M., et al "High levels of human anti-human antibodies to adalimumab in a patient not responding to adalimumab treatment" Ann.Rheum.Dis. 2006; 65:1249-1250.
3. Bartelds, G. M., et al "Clinical response to adalimumab: relationship to anti-adalimumab antibodies and serum adalimumab concentrations in rheumatoid arthritis" Ann.Rheum.Dis. 2007; 66:921-926.
4. Bartelds, G. M., et al "Development of antidrug antibodies against adalimumab and association with disease activity and treatment failure during long-term follow-up". JAMA 2011; 305:1460-1468.
5. Lecluse, L. L., et al "Extent and clinical consequences of antibody formation against adalimumab in patients with plaque psoriasis" Arch.Dermatol. 2010,146:127-132.

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18. LITERATUR

- Chen, D. Y. et al "Significant associations of antidrug antibody levels with serum drug trough levels and therapeutic response of adalimumab and etanercept treatment in rheumatoid arthritis" Ann Rheum Dis 2014
- Bartelds, G. M., et al "High levels of human anti-human antibodies to adalimumab in a patient not responding to adalimumab treatment" Ann.Rheum.Dis. 2006; 65:1249-1250.
- Bartelds, G. M., et al "Clinical response to adalimumab: relationship to anti-adalimumab antibodies and serum adalimumab concentrations in rheumatoid arthritis" Ann.Rheum.Dis. 2007; 66:921-926.
- Bartelds, G. M., et al "Development of antidrug antibodies against adalimumab and association with disease activity and treatment failure during long-term follow-up". JAMA 2011; 305:1460-1468.
- Lecluse, L. L., et al "Extent and clinical consequences of antibody formation against adalimumab in patients with plaque psoriasis" Arch.Dermatol. 2010;146:127-132.

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- Chen, D. Y. et al "Significant associations of antidrug antibody levels with serum drug trough levels and therapeutic response of adalimumab and etanercept treatment in rheumatoid arthritis" Ann Rheum Dis 2014
- Bartelds, G. M., et al "High levels of human anti-human antibodies to adalimumab in a patient not responding to adalimumab treatment" Ann.Rheum.Dis. 2006; 65:1249-1250.
- Bartelds, G. M., et al "Clinical response to adalimumab: relationship to anti-adalimumab antibodies and serum adalimumab concentrations in rheumatoid arthritis" Ann.Rheum.Dis. 2007; 66:921-926.
- Bartelds, G. M., et al "Development of antidrug antibodies against adalimumab and association with disease activity and treatment failure during long-term follow-up". JAMA 2011; 305:1460-1468.
- Lecluse, L. L., et al "Extent and clinical consequences of antibody formation against adalimumab in patients with plaque psoriasis" Arch.Dermatol. 2010,146:127-132.

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18. BIBLIOGRAFIA

1. Chen, D. Y. et al "Significant associations of antidrug antibody levels with serum drug trough levels and therapeutic response of adalimumab and etanercept treatment in rheumatoid arthritis" Ann Rheum Dis 2014

2. Bartelds, G. M., et al "High levels of human anti-human antibodies to adalimumab in a patient not responding to adalimumab treatment" Ann.Rheum.Dis. 2006; 65:1249-1250.
3. Bartelds, G. M., et al "Clinical response to adalimumab: relationship to anti-adalimumab antibodies and serum adalimumab concentrations in rheumatoid arthritis" Ann.Rheum.Dis. 2007; 66:921-926.
4. Bartelds, G. M., et al "Development of antidrug antibodies against adalimumab and association with disease activity and treatment failure during long-term follow-up". JAMA 2011; 305:1460-1468.
5. Lecluse, L. L., et al "Extent and clinical consequences of antibody formation against adalimumab in patients with plaque psoriasis" Arch.Dermatol. 2010;146:127-132.

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DIESE Diagnostica Senese S.p.A.