W.W.-DAT**DIESSE****REF 26010**

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INSTRUCTIONS FOR USE

W.W.-DAT

For the determination of anti-Salmonella and anti-Brucella specific agglutinating antibodies

For In vitro Diagnostic Use Only

1. INTENDED USE

Method for the determination of specific agglutinating antibodies in human serum using colored Salmonella and Brucella suspensions.

2. INTRODUCTION

Serological studies of patients with persistent fever can be useful in detecting typhoid fever and brucellosis.

The test consists in the observation of a direct agglutination reaction between the bacterial suspension and the serum sample which is suspected of containing specific antibodies. The present method, known as the Widal-Wright reaction, allows a semiquantitative titration of the serum agglutinins.

3. PRINCIPLE OF THE METHOD

The serum sample is mixed with a bacterial suspension at the optimal concentration for the semi-automated test on the AUTO-DAT instrument and/or in slide.

The titer is attributed with a software processing of the image of the agglutinate, which is obtained in microplate wells, and stored in the instrument.

The bacterial suspensions can also be used to perform the test in slide and / or tube.

4. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

All material of human origin must be handled as potentially infectious. This kit contains materials of human origin which have been tested and gave a negative response for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens and performing the assay.
3. Wash hands thoroughly when finished.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use.

1. After use, return reagents to 2-8°C.
2. Do not use strongly hemolyzed, lipemic, icteric samples, serum not completely coagulated or samples presenting microbial contamination.
3. Do not use the reagents after expiry date.
4. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
5. Do not modify the Test Procedure.

5. KIT COMPOSITION AND REAGENT PREPARATION

Reagents are sufficient to perform 100 tests in automation with AUTO-DAT (700 tests).

Reagents are sufficient to analyze 60 samples in slide.

BACTERIAL SUSPENSION BACTERIAL SUSPENSION 1 x 2.5 mL

of the following suspensions:

- S. THYPHI O
- S. THYPHI H
- S. PARATHYPHI A O
- S. PARATHYPHI A H
- S. PARATHYPHI B O
- S. PARATHYPHI B H
- BRUCELLA

Contents: Suspension of stained, killed bacteria at an optimized concentration for slide test performance.

Contains sodium azide 0.09% as preservative.

CONTROL + POSITIVE CONTROL 1 x 1 mL

Contents: Polyvalent positive control serum of animal origin, containing sodium azide 0.09% as preservative. Contains antibodies which agglutinate all the bacterial suspensions of the kit.

CONTROL - NEGATIVE CONTROL

1 x 1 mL

Contents: Negative control serum containing Proclin and Gentamicin as preservative.

BRING REAGENTS TO ROOM TEMPERATURE BEFORE USE**OTHER MATERIALS:**

- Slides (3)
- Mixing sticks (100)

MATERIALS REQUIRED BUT NOT PROVIDED:

- Pipettes with disposable needles
- Physiological saline
- 96-well microplates with flat bottom, Auto-DAT Microplate [REF] 26005
- Single use tubes for tube test.
- AUTO-DAT instrument ([REF] 26000)

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

The stability of the reagents does not change after opening the bottle, provided that the user is careful to keep the product safe from possible microbial contamination.

7. TYPE AND STORAGE OF SAMPLE

Fresh serum may be stored at 2/8°C for 7 days. For longer storage, freeze at -20°C, and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Use immediately after thawing.

8. PROCEDURE

Bring the reagent to room temperature and mix the vial of bacterial suspension carefully but not vigorously.

INSTRUMENTAL APPLICATION WITH AUTO-DAT

Pipette 25 µL of sample and 25 µL of reagent into a microplate well ([REF] 26005). Insert the microplate into the instrument and press START. The instrument, at the end of the reaction (which lasts about 15 minutes) will print the result in titer.

SCREENING TEST

If the test on AUTO-DAT cannot be performed, reagents may be used for a screening test.

Distribute 20 µL of test serum into the slide circle. Add 40 µL of suspension and mix; rotate the slide and read after exactly one minute. Positive samples will give agglutination while negative samples will remain homogeneous.

TUBE TRITATION

Dilute the sample 1:20 with physiological saline (0.1 mL serum + 1.9 mL saline). Take 1 mL from the first tube and make serial dilutions with 1 mL of saline, as shown in the

scheme. Add 40 µL of suspension to each tube, mix and incubate overnight at 37°C (see the scheme shown below). The blank is performed by mixing the suspension with saline.

Tube	1	2	3	4	5	6	7
Saline Solution ml	1.9	1	1	1	1	1	1
Serum ml	0.1						
Mix and transfer ml	1	1	1	1	1	1	1
Dilution/ Titer	1/20	1/40	1/80	1/160	1/320	1/640	1/1280

9. TEST VALIDATION ON AUTO-DAT

Using 25 µL of **Positive Control** serum, perform the test on AUTO-DAT and compare the result with the titers reported on the certificate of analysis. +/- one titer dilution deviations are accepted.

If the result is different from the expected, contact the Scientific Support:

Tel: 0039 0577 319554

Fax: 0039 0577 366605

email: scientificsupport@diesse.it

10. INTERPRETATION OF THE TEST**Positive reaction**

Instrumental application with AUTO-DAT: the titer is printed.

Screening test: presence of agglutinate compared to negative control.

Tube titration: presence of granular (antigen O) or flocculant (antigen H) agglutinate with clear supernatant in comparison with the blank.

Negative reaction

Instrumental application with AUTO-DAT: the printing reports that the sample is negative.

Screening test: absence of agglutinate compared to negative control. Homogeneous suspension.

Tube titration: absence of agglutinate compared to negative control. Homogeneous suspension.

The titer is given by the highest dilution able to agglutinate the present bacteria.

Titers 1:80 are considered suspect; higher titers are diagnostically significant.

11. LIMITATIONS OF THE TEST

The prozone has been studied with the Brucella suspension which is more subject to this phenomenon. No problems have been found with high-titer samples (1/640).

Since the phenomenon cannot be absolutely excluded for suspected samples, repeat the test by diluting the sample at least 1:4. A negative reaction with the Brucella suspension can occur in presence of incomplete antibodies in chronic infections. The evaluation of the presence of incomplete antibodies can be performed adding 10 µl of positive control (REF 26007) to the reaction mix that resulted negative. After

new incubation, the absence of agglutination indicates the presence of incomplete antibodies in the tested serum.

In the slide and automated test, the agglutination of the H antigens suspension can also be caused by antibodies that react with the lipopolysaccharide. Differentiation is best obtained by examining the flaky agglutination in the tube.

12. METHOD COMPARISON – DIAGNOSTIC SENSITIVITY AND SPECIFICITY

An experimentation performed on human sera obtained from hospital laboratories, the titers of which were taken as reference, allowed the evaluation of the sensitivity of the suspensions *S. typhi* O, *S. typhi* H, *S. paratyphi* BO, *S. paratyphi* BH. When *S. typhi* antigen H was determined on 30 positive samples, a sensitivity of 96.8% was found, while the suspension *S. typhi* antigen O tested on 51 positive samples gave 2 false negative results and a sensitivity of 96.2%. The *S. paratyphi* B antigen H and antigen O suspensions were tested with 12 and 24 positive samples, respectively. A sensitivity of 100% and 91.7% were obtained respectively.

The suspensions *S. paratyphi* A antigen O and antigen H were tested in parallel with another commercial product. 12 positive and 3 doubtful results were found. The Diesse reagent gave 1 false negative and therefore the sensitivity was 95.0% considering the doubtful results as positive. With the same criteria, 45 positive samples for Brucellosis were tested. The Diesse reagent gave 2 false negative results and therefore the sensitivity was 95.7%.

Human sera found negative in another commercial test were also tested, with the following results:

Suspension	No. of samples	Specificity %
<i>S. typhi</i> antigene O	93	97.4
<i>S. typhi</i> antigene H	58	95.1
<i>S. paratyphi</i> B antigene O	96	97.0
<i>S. paratyphi</i> B antigene H	93	95.9
<i>S. paratyphi</i> A antigene O	82	98.8
<i>S. paratyphi</i> A antigene H	80	97.6
Brucella	83	97.6

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