RPR-DAT

REF 26030

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INSTRUCTIONS FOR USE

RPR-DAT

Non-treponemal test for the qualitative and semi-quantitative determination of reagins in serum to be applied on the AUTO-DAT instrument.

For In Vitro Diagnostic Use Only

1. INTENDED USE

Non-treponemal test for the qualitative and semi-quantitative determination of reagins in serum.

2. INTRODUCTION

Infection with *Treponema pallidum* gives rise to a complex antibody response. Non treponemal tests are based on the detection of reagins, an antibody class present in syphilis and occasionally in other acute and chronic conditions. Reagins can be detected in the serum about 4-6 weeks after infection, or 1-3 weeks after appearance of the chancre. Purified beef heart extract (cardiolipin) fortified with lecithin and cholesterol, is used as antigen for reagin tests. Non treponemal tests are routinely used in syphilis serology because of their practical advantages and reproducibility, even though they are not always 100% specific.

3. PRINCIPLE OF THE METHOD

The serum sample reacts with the colloidal suspension of cardiolipin, lecithin and cholesterol mixed with microparticles to perform the semi-automated test on the AUTO-DAT instrument.

The analysis is carried out through a software processing of the agglutinate image, which is obtained in microplate wells, and which is stored in the instrument.

In case of malfunction of the instrument the reagent can be used with manual procedure.

4. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY

All material of biological origin must be considered potentially infectious. All reagents and samples must be handled in accordance with the safety precautions adopted in the laboratory.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens and performing the assay.
3. Wash hands thoroughly when finished.

4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the device to room temperature (18-30°C) before use.

1. After use, and return reagents to 2-8°C.
2. Do not use the reagents after the expiry date.
3. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
4. Do not modify the Test Procedure. Do not substitute reagents using reagents from other suppliers or other lots, unless it is specifically indicated that the reagent is interchangeable between lots.
5. Do not use lipemic sera

5. KIT CONTENTS AND REAGENT PREPARATION

Reagents are sufficient for 500 determinations in automation with AUTO-DAT.

ANTIGEN

ANTIGEN

2 x 4 mL

Contents: the antigen is composed of a stabilised mixture of cardiolipin, cholesterol, lecithin and carbon microparticles.

Attention: Store in an upright position when not in use. Bring to room temperature before use. Do not shake too vigorously.

CONTROL +

POSITIVE CONTROL

1 x 0.5 mL

Contents: Human serum, reactive for syphilis, diluted in a protein solution with sodium azide 0.09%.

Liquid, ready for use.

CONTROL -

NEGATIVE CONTROL/SAMPLE DILUENT

1 x 4 mL

Contents: Human serum containing Proclin and Gentamicin as preservatives.

Liquid, ready for use.

BRING THE REAGENTS TO ROOM TEMPERATURE BEFORE USE.

MT PLATE

6 FLAT BOTTOM MICROPLATES (8x12)

MATERIALS REQUIRED BUT NOT PROVIDED

- AUTO-DAT Instrument (REF 26000)
- Rotating mixer
- Micropipette with variable volume
- Microplate mixer

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C

After the first use, the microplate must be stored at room temperature.

The expiry date is printed on each component and on the kit label.

The stability of the reagents does not change after opening the bottle, provided that the user is careful to keep the product safe from possible microbial contamination.

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of fresh or defrosted serum. Fresh serum may be stored at 2/8°C for 5 days.

For longer storage, freeze at -20°C; it can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Use immediately after thawing.

8. ASSAY PROCEDURE**INSTRUMENTAL APPLICATION WITH AUTO-DAT****SCREENING TEST**

1. Bring reagents and samples to room temperature. Resuspend the reagent.
2. Draw up 25 µl of serum and drop on the bottom of the microplate well. Take care not to transfer any cellular elements. **Avoid the formation of bubbles.**
3. Repeat these steps for each sample to be tested.
4. Distribute 15 µl of antigen in each well. **Avoid the formation of bubbles.**
5. Verify that the sample-antigen mix is **evenly distributed on the bottom of the well**. If not, tap the plate on the side and / or gently tap it on the work surface.
6. Place the plates in the AUTO-DAT instrument and follow the instructions on the tablet.
7. At the end of the analysis, extract the plates from the AUTO-DAT instrument and, using the adhesive film, cover the wells in which the reaction occurred. Wells not previously used and not covered can be used for further analysis.

NOTE: Leave one strip microwell for Negative Control and one for Positive Control, to be tested exactly as the analyzed samples.

SEMIQUANTITATIVE TEST

To be carried out on the samples reactive to the screening test.

1. Add 25 µl of Negative Control/Sample Diluent to each well of the plate.

2. Add to the first well 25 µl of the sample that gave positive result.
3. Make doubling dilutions, mixing thoroughly before transferring to each step, as shown in the following scheme. **Avoid the formation of bubbles.**

Strip 1x8	A	B	C	D	E	F	G
Sample Diluent (µl)	25	25	25	25	25	25	25
Serum (µl)	25						
Mix and transfer (µl)							
Dilution/Titer	1/2	1/4	1/8	1/16	1/32	1/64	1/128

4. Follow the procedure described for the screening test from step 4 to step 7.

The titer will be given by the last dilution found reactive with the instrument.

MANUAL PROCEDURE

1. Prepare samples as described for the screening test (points 1-5).
2. Mix reagent and sample for 1 minute with a microplate mixer at maximum speed.
3. Place the plates on a rotating mixer for 10 minutes at 100 ± 4 rpm.
4. Analyze the obtained result.

9. TEST VALIDATION ON AUTO-DAT

Use the controls provided at each session.

Proceed as described in the paragraph "INSTRUMENTAL APPLICATION WITH AUTO-DAT".

If the result is different from the expected, contact the Scientific Support:

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE TEST**INSTRUMENTAL APPLICATION WITH AUTO-DAT**

At the end of the test the instrument will print the results as non-reactive (N), doubtful (D), reactive (P).

It is advisable to visually verify the obtained results for confirmation.

MANUAL PROCEDURE

Once the test is finished, the sample will be:

REACTIVE: the presence of clearly visible black aggregates is typical of a reactive sample.

The Positive Control must give evident agglutination.

NON-REACTIVE: No aggregate formation.

The Negative Control must not give aggregate formation.

NOTE: Any sample giving a reactive result should be tested in the semi-quantitative method (as described above).

ATTENTION: failure to add the antigen to the sample causes false positive results.

11. LIMITATIONS

Biological false-positive reactions occur occasionally with cardiolipin antigens, e.g. in cases of drug abuse, or in certain diseases such as lupus erythematosus, mononucleosis, malaria, leprosy.

A pro-zone reaction is sometimes found: complete or partial inhibition of reactivity occurs with undiluted sample and maximum reactivity is found with diluted sample. This can lead to a weakly positive result in the qualitative test, whereas a strong reaction will be found in the quantitative test after dilution. For this reason, all specimens that give doubtful or reactive results should be also tested in the semi-quantitative procedure.

On human sera the phenomenon has not been experienced up to titers of 1/64.

12. METHOD COMPARISON

In an experimentation, 159 samples were analyzed with Diesse kit and with another commercial kit. The experimental data are summarized below:

		Reference		
		+	-	Total
		48	1	49
Diesse	+	48	1	49
	-	1	109	110
	Totale	49	110	159

Percent Positive Agreement (Diagnostic Sensitivity):

98.0% CI95%: 89.2-99.6

Percent Negative Agreement: (Diagnostic Specificity): 99.1%

CI95%: 95.0-99.8

The agreement between the two methods is excellent with a with a Cohen's Kappa of 0.98.

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