CHORUS

BORDETELLA PERTUSSIS TOXIN IgG

REF 81272

DIESSE

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Correlazione	r	95%CI
Pearson	0.97	0.95-0.98
Spearman	0.97	0.95-0.98

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17. BIBLIOGRAFIA

- Riffelmann M., Krefeld H.K., Wirsing von König C.H., Xing D., He Q. (2012). Guidance and protocol for the serological diagnosis of human infection with Bordetella pertussis. European Centre for Disease Prevention and Control.
- Guiso N., Berbers G., Fry N.K., He Q., Riffelmann M., Wirsing von König C.H. (2011). What to do and what not to do in serological diagnosis of pertussis: recommendations from EU reference laboratories. Eur J Clin Microbiol Infect Dis (2011) 30:307-312.
- Riffelmann M., Thiel K., Schmetz J., Wirsing von Koenig C.H. (2010). Performance of Commercial Enzyme-Linked Immunosorbent Assays for detection of Antibodies to Bordetella Pertussis. Journal of Clinical Microbiology p. 4459-4463

INSTRUCTIONS FOR USE

CHORUS BORDETELLA PERTUSSIS TOXIN IgG

For the quantitative determination of anti-Bordetella pertussis toxin IgG antibodies

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the quantitative determination of IgG class antibodies against Bordetella pertussis toxin in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Whooping cough is an acute and very contagious infection of the respiratory tract. It causes inflammation and necrosis of the respiratory tract ciliated epithelium and it produces immunological and pathophysiological reactions.

Pertussis toxin (PT) is a protein produced by the Bordetella pertussis Gram-negative bacterium and it is one of the main virulent factors of the disease.

It is typically characterized by an incubation phase of 15-20 days, followed by a catarrhal stage of about 2 weeks, and in the end by a paroxysmal stage lasting 2 or 3 weeks. In young children more serious complications are constituted by further infections that can lead to otitis, pneumonia, bronchitis, or even neurological diseases.

The serological test is based on the determination of specific antibodies through the ELISA principle. IgA antibodies, specific for PT, can be detected after 1-2 weeks and persist for 6-24 months depending on the age. The IgA production is slower in children. IgA antibodies, in suckling infants, can be detected after 6-8 weeks.

IgG antibodies are not detectable before 2-3 weeks from the beginning of the disease and reach their maximum after 6-8 weeks.

3. PRINCIPLE OF THE METHOD

The Chorus Bordetella pertussis toxin IgG device is ready to use for the detection of IgG antibodies against Bordetella pertussis toxin, in the Chorus/Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted sample.

After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human IgG monoclonal antibodies conjugated to horseradish peroxidase. The unbound conjugate is eliminated, and the peroxidase substrate added. The blue colour which develops is proportional to the concentration of specific antibodies present in the examined sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments. The results are expressed in International Units (IU/ml) calculated in reference to WHO 1st IS NIBSC 06/140.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.

The use of the kit is only possible with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website

(<https://www.diesse.it/en/downloads/downloads/strum ento:39/>)


5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
11. Do not use the device after the expiry date.
12. Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with highly purified Bordetella pertussis toxin antigen

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine 0.26 mg/ml and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT

Contents: Proteic buffer solution.

Position 2: CONJUGATE

Contents: anti-human IgG antibodies labelled with horseradish peroxidase.

Position 1: EMPTY WELL

in which the sample is transferred.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.450 mL

Contents: Diluted human serum containing anti- Bordetella pertussis toxin IgG antibodies and preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL

1 x 0.900 mL

Contents: Diluted human serum containing anti- Bordetella pertussis toxin IgG antibodies and preservative. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh sample may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense in well no. 1 of each device the following:

SAMPLE	50 µl/device
CALIBRATOR	130 µl/device
POSITIVE CONTROL	130 µl/device

4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in International Units (IU/ml) (WHO 1st IS NIBSC 06/140), calculated on the basis of a batch-dependent curve stored in the instrument.

The test on the examined sample can be interpreted as follows:

POSITIVE: when the result is > 100 IU/ml

NEGATIVE: when the result is < 40 IU/ml

DOUBTFUL/EQUIVOCAL: for all values between 40 and 100 IU/ml

If the result is doubtful/equivocal, it is recommended to perform the test of IgA antibodies too.

Eventually, collect a new sample after 3-4 weeks to observe the increase or decrease of IgG antibodies concentration.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. CALIBRATION RANGE

Calibration range 1-350 IU/ml.

For samples > 350 IU/ml repeat the test pre-diluting the sample with Negative Control/Sample Diluent (PF83607 – not supplied with the kit).

13. ANALYTICAL SPECIFICITY

2 samples (1 Negative and 1 Positive) were spiked with the following potentially interfering factors and then tested:

Bilirubin (0.4 mg/ml)
Triglycerides (20 mg/ml)
Hemoglobin (5 mg/ml)

The presence in the serum sample of the interfering substances described above does not affect the test result.

14. CROSS-REACTIONS

127 samples, positive to *Chlamydomydia pneumoniae*, *Chlamydia trachomatis*, Influenza A, Influenza B, Parainfluenza, *Mycoplasma pneumoniae*, RSV, Adenovirus, SARS-CoV-2, RF, ANA and ASCA were tested.

No significant cross-reactions were found.

15. METHOD COMPARISON

In an experimentation 69 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

Correlation	r	95%CI
Pearson	0.97	0.95-0.98
Spearman	0.97	0.95-0.98

The correlation between the two methods is very high.

16. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision	
	Mean (IU/ml)	CV%	Mean (IU/ml)	CV%
1	109.0	6.6	114.7	5.8
2	69.7	6.0	72.9	6.9
3	34.4	6.0	34.9	5.9

Sample	Precision between batches		Precision between instruments	
	Mean (IU/ml)	CV%	Mean (IU/ml)	CV%
1	104.0	10.9	107.1	9.4
2	71.0	9.2	68.3	14.7
3	32.3	11.6	32.1	14.7

17. REFERENCES

1. Riffelmann M., Krefeld H.K., Wirsing von König C.H., Xing D., He Q. (2012). Guidance and protocol for the serological diagnosis of human infection with *Bordetella pertussis*. European Centre for Disease Prevention and Control.
2. Guiso N., Berbers G., Fry N.K., He Q., Riffelmann M., Wirsing von König C.H. (2011). What to do and what not to do in serological diagnosis of pertussis: recommendations from EU reference laboratories. *Eur J Clin Microbiol Infect Dis* (2011) 30:307-312.
3. Riffelmann M., Thiel K., Schmetz J., Wirsing von König C.H. (2010). Performance of Commercial Enzyme-Linked Immunosorbent Assays for detection of Antibodies to *Bordetella Pertussis*. *Journal of Clinical Microbiology* p. 4459-4463

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1. Riffelmann M., Krefeld H.K., Wirsing von König C.H., Xing D., He Q. (2012). Guidance and protocol for the serological diagnosis of human infection with *Bordetella pertussis*. European Centre for Disease Prevention and Control.
2. Guiso N., Berbers G., Fry N.K., He Q., Riffelmann M., Wirsing von König C.H. (2011). What to do and what not to do in serological diagnosis of pertussis: recommendations from EU reference laboratories. *Eur J Clin Microbiol Infect Dis* (2011) 30:307-312.
3. Riffelmann M., Thiel K., Schmetz J., Wirsing von Koenig C.H. (2010). Performance of Commercial Enzyme-Linked Immunosorbent Assays for detection of Antibodies to *Bordetella Pertussis*. *Journal of Clinical Microbiology* p. 4459-4463

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17. BIBLIOGRAPHIE

- Riffelmann M., Krefeld H.K., Wirsing von König C.H., Xing D., He Q. (2012). Guidance and protocol for the serological diagnosis of human infection with *Bordetella pertussis*. European Centre for Disease Prevention and Control.
- Guiso N., Berbers G., Fry N.K., He Q., Riffelmann M., Wirsing von König C.H. (2011). What to do and what not to do in serological diagnosis of pertussis: recommendations from EU reference laboratories. *Eur J Clin Microbiol Infect Dis* (2011) 30:307-312.
- Riffelmann M., Thiel K., Schmetz J., Wirsing von König C.H. (2010). Performance of Commercial Enzyme-Linked Immunosorbent Assays for detection of Antibodies to *Bordetella Pertussis*. *Journal of Clinical Microbiology* p. 4459-4463

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17. BIBLIOGRAFIA

1. Riffelmann M., Krefeld H.K., Wirsing von König C.H., Xing D., He Q. (2012). Guidance and protocol for the serological diagnosis of human infection with Bordetella pertussis. European Centre for Disease Prevention and Control.
2. Guiso N., Berbers G., Fry N.K., He Q., Riffelmann M., Wirsing von König C.H. (2011). What to do and what not to do in serological diagnosis of pertussis: recommendations from EU reference laboratories. Eur J Clin Microbiol Infect Dis (2011) 30:307-312.
3. Riffelmann M., Thiel K., Schmetz J., Wirsing von Koenig C.H. (2010). Performance of Commercial Enzyme-Linked Immunosorbent Assays for detection of Antibodies to Bordetella Pertussis. Journal of Clinical Microbiology p. 4459-4463

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17. BIBLIOGRAFIE

- Riffelmann M., Krefeld H.K., Wirsing von König C.H., Xing D., He Q. (2012). Guidance and protocol for the serological diagnosis of human infection with Bordetella pertussis. European Centre for Disease Prevention and Control.
- Guiso N., Berbers G., Fry N.K., He Q., Riffelmann M., Wirsing von König C.H. (2011). What to do and what not to do in serological diagnosis of pertussis: recommendations from EU reference laboratories. Eur J Clin Microbiol Infect Dis (2011) 30:307-312.
- Riffelmann M., Thiel K., Schmetz J., Wirsing von König C.H. (2010). Performance of Commercial Enzyme-Linked Immunosorbent Assays for detection of Antibodies to Bordetella Pertussis. Journal of Clinical Microbiology p. 4459-4463

	IT	Data di fabbricazione	ES	Fecha de fabricación
	EN	Date of manufacture	FR	Date de fabrication
	CZ	Datum výroby	PT	Data de fabrico
	DE	Herstellungsdatum	RO	Data fabricatiei
	GR	Ημερομηνία Παραγωγής		
	IT	Utilizzare entro	ES	Fecha de caducidad
	EN	Use By	FR	Utiliser jusqu'à
	CZ	Použitelné do	PT	Prazo de validade
	DE	Verwendbar bis	RO	A se folosi pana la
	GR	Ημερομηνία λήξης		
	IT	Non riutilizzare	ES	No reutilizar
	EN	Do not reuse	FR	Ne pas réutiliser
	CZ	Nepoužívejte opakovaně	PT	Não reutilizar
	DE	Nicht wieder verwenden	RO	A nu se refolosi
	GR	Μην κάνετε επαναληπτική χρήση		
	IT	Attenzione, vedere le istruzioni per l'uso	ES	Atención, ver instrucciones de uso
	EN	Caution, consult accompanying documents	FR	Attention voir notice d'instructions
	CZ	Pozor, čtěte příložené dokumenty	PT	Atenção, consulte a documentação incluída
	DE	Achtung, die Gebrauchsanleitung lesen	RO	Atentie, consultați documentele însoțitoare
	GR	Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα		
	IT	Fabbricante	ES	Fabricante
	EN	Manufacturer	FR	Fabricant
	CZ	Výrobce	PT	Fabricante
	DE	Hersteller	RO	Productator
	GR	Κατασκευαστής		
	IT	Contenuto sufficiente per "n" saggi	ES	Contenido suficiente para <n> ensayos
	EN	Contains sufficient for <n> tests	FR	Contenu suffisant pour "n" tests
	CZ	Obsah stačí na <n> testů	PT	Conteúdo suficiente para "n" ensaios
	DE	Inhalt reicht für „n“ Tests	RO	Continut suficient pt <n> teste
	GR	Περιεχόμενο επαρκές για «n» εξετάσεις		
	IT	Limiti di temperatura	ES	Límite de temperatura
	EN	Temperature limitation	FR	Limites de température
	CZ	Teplotní omezení	PT	Limites de temperatura
	DE	Temperaturgrenzwerte	RO	Limita da temperatura
	GR	Περιορισμοί θερμοκρασίας		
	IT	Consultare le istruzioni per l'uso	ES	Consulte las instrucciones de uso
	EN	Consult Instructions for Use	FR	Consulter les instructions d'utilisation
	CZ	Čtěte návod k použití	PT	Consulte as instruções de utilização
	DE	Die Gebrauchsanleitung lesen	RO	Pentru a utiliza consultați instrucțiunile
	GR	Συμβουλευτείτε τις οδηγίες χρήσης		
	IT	Rischio biologico	ES	Riesgo biológico
	EN	Biological risks	FR	Risques biologiques
	CZ	Biologická rizika	PT	Risco biológico
	DE	Biologisches Risiko	RO	Risk biologic
	GR	Βιολογικοί κίνδυνοι		
REF	IT	Numero di catalogo	ES	Número de catálogo
	EN	Catalogue number	FR	Référence du catalogue
	CZ	Katalogové číslo	PT	Referência de catálogo
	DE	Katalognummer	RO	Numar de catalog
	GR	Αριθμός καταλόγου		
IVD	IT	Dispositivo medico-diagnostico in vitro	ES	Producto sanitario para diagnóstico in vitro
	EN	In Vitro Diagnostic Medical Device	FR	Dispositif médical de diagnostic in vitro
	CZ	Lékařské vybavení pro diagnostiku in vitro	PT	Dispositivo médico para diagnóstico in vitro
	DE	Medizinisches In-vitro-Diagnostikum	RO	Dizpositiv medical pentru diagnosticare in vitro
	GR	Ιn Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν		
LOT	IT	Codice del lotto	ES	Código de lote
	EN	Batch code	FR	Code du lot
	CS	Kód šarže	PT	Código do lote
	DE	Chargennummer	RO	Lot
	GR	Αριθμός Παρτίδας		
CE	IT	Marcatura CE di conformità	ES	Marcado CE de conformidad
	EN	CE marking of conformity	FR	Marquage de conformité CE
	CZ	Označení shody CE	PT	Marcação CE de conformidade
	DE	CE-Konformität Skennzeichnung	RO	Marcajul de conformitate CE
	GR	Σημάνση συμμορφωσ CE		