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19. BIBLIOGRAFIA

1. Vaughn B.P. "A Practical Guide to Therapeutic Drug Monitoring of Biologic Medications for Inflammatory Bowel Disease" J Clin Med (2021) Oct 27;10(21):4990.
2. Gómez-Arango C. et All. "Cost-Effectiveness of Therapeutic Drug Monitoring-Guided Adalimumab Therapy in Rheumatic Diseases: A Prospective, Pragmatic Trial" Rheumatol Ther (2021) Sep;8(3):1323-1339
3. Sam M.J. et All. "Comparative Evaluation of 4 Commercially Available ELISA Kits for Measuring Adalimumab and Anti-adalimumab Antibodies" Ther Drug Monit (2020) Dec;42(6):821-828.
4. Mytych D. et All. "Evaluation of marketed kits for measurement of ABP 501, the first approved Adalimumab biosimilar, drug concentration and anti-drug antibody levels in patient serum" Ann Rheum Dis (2019) Abstract FRI0645.

INSTRUCTIONS FOR USE

CHORUS Promonitor ANTI-ADALIMUMAB

For the semiquantitative detection of anti-Adalimumab antibodies in human serum

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

CHORUS Promonitor ANTI-ADALIMUMAB is an immunoassay kit for automated semiquantitative detection of anti-Adalimumab antibodies in human serum using a disposable device applied on the CHORUS TRIO instrument. The product, used in conjunction with other clinical and laboratory findings, is useful as an aid in the management of patients treated with Adalimumab (ADL), like patients with rheumatic and gastrointestinal diseases. It must be used by professional laboratory users only.

2. INTRODUCTION

Adalimumab (ADL) is a recombinant human monoclonal antibody expressed in Chinese Hamster Ovary (CHO) cells used in the treatment of inflammatory diseases. This includes rheumatoid arthritis (RA), polyarticular juvenile idiopathic arthritis (JIA), psoriatic arthritis (PA), ankylosing spondylitis (AS), Crohn's disease (CD) and Psoriasis (PsO).

ADL is an immunomodulator biological drug which binds TNF α , a cytokine that plays an important role in the pathogenesis of these diseases, inhibits binding with its receptor, and prevents triggering the inflammatory cascade. As for other biological drugs, ADL can lead to the development of anti-drug antibodies (ADAs). Studies confirm a significant inverse correlation between ADL bioavailability and ADAs presence. The clinical response to ADL closely follows the trough drug levels and the presence of ADAs. It has been demonstrated that immunogenicity of ADL is associated with loss of clinical response in RA, CD and PsO. Furthermore, patients who develop anti-ADL antibodies are more likely not to respond to the therapy and those with therapeutic drug levels are more likely to have an adequate clinical response to ADL therapy in RA and IBD (Inflammatory Bowel Disease).

3. PRINCIPLE OF THE METHOD

The CHORUS Promonitor ANTI-ADALIMUMAB device is ready to use for the detection of anti-Adalimumab antibodies in human serum, in the CHORUS TRIO instrument.

The assay is based on a bridging ELISA method.

The bridging antibody ELISA method is widely used for anti-drug antibody screening and is well accepted as a standard for immunogenicity screening.

The microwell strips are pre-coated with ADL. Calibrator, positive control and samples are added to separate devices.

Anti-ADL antibodies bind to pre-immobilized ADL following incubation with the sample. The following addition and incubation of HRP-labelled ADL (conjugate) allows the conjugate to bind to the anti-ADL antibodies that have become attached to the microwells. After washing to remove the unbound conjugate, the

peroxidase substrate is added. The enzymatic reaction is stopped by adding Stop Solution, which makes the solution turn to yellow. The color which develops is proportional to the amount of anti-ADL binding antibodies in the patient sample.

The disposable devices contain all the reagents to perform the test when applied on the CHORUS TRIO instrument.

The results are expressed in Arbitrary Units (AU/ml).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the CHORUS TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30-minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the bags with the devices to room temperature (18-30°C) for at least 30 minutes before use.

Use the devices within 60 min.

1. **Discard devices which show the substrate (well 2) blue colored.**
2. If the sample is added into the well manually, verify that it is perfectly distributed.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.

4. The devices are for use with the CHORUS TRIO instrument; the Instructions for Use must be carefully followed, and the Instrument Operating Manual must be consulted.

The use of the kit is only possible with CHORUS TRIO with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website

(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)

5. Check that the CHORUS TRIO instrument is set up correctly (see CHORUS TRIO Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see CHORUS TRIO Operating Manual).
9. When CHORUS Autosampler (REF 81210) is used, refer to the related Operating Manual
10. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
11. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
12. Do not use the device after the expiry date.
13. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 9: Space for application of bar code label

Position 8: EMPTY

Position 7: MICROPLATE WELL

Coated with ADL

Position 6: MICROPLATE WELL

Coated with ADL

Position 5: STOP SOLUTION

Contents: 0.3 M sulfuric acid solution

Position 4: CONJUGATE

Contents: HRP-labelled ADL antibody and buffer containing preservative

Position 3: DILUTION BUFFER

Contents: Buffer solution containing preservative

Position 2: TMB SUBSTRATE

Contents: Tetramethylbenzidine stabilized in citrate buffer

Position 1: EMPTY WELL

in which the sample is transferred.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.650 ml

Contents: Human serum containing anti-ADL antibodies, and preservative.

Liquid, ready for use.

CONTROL + POSITIVE CONTROL

1 x 1.100 ml

Contents: Human serum containing anti-ADL antibodies, and preservative.

Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- CHORUS TRIO instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated, and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 5 days at 2/8°C, or frozen for longer periods at -20°C and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be mixed carefully before use.

It is strongly recommended to avoid any harsh treatment of the samples.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

The test CHORUS Promonitor ANTI-ADALIMUMAB (REF. 86703) can be performed simultaneously **ONLY** with the other CHORUS Promonitor tests.

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense the following volumes in well no. 1 of each device:

SAMPLE	130 µl/device
CALIBRATOR	130 µl/device
POSITIVE CONTROL	130 µl/device

At each change of batch, use a device for the calibrator.

4. Place the devices in the CHORUS TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the obtained results. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected. If the result of the positive control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The CHORUS TRIO instrument expresses the result in AU/ml, calculated on the basis of a batch-dependent curve stored in the instrument.

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is >10.0 AU/ml
NEGATIVE: when the result is ≤ 10.0 AU/ml

11. LIMITATIONS

1. Trough antibody levels should be measured. The optimum sample collection time is immediately before the ADL infusion.
2. This assay measures binding antibodies to ADL. It does not measure immunocomplexes formed between ADL and anti-ADL antibodies.
3. Bridging ELISA does not detect IgG4, since this isotype is bispecific due to Fab arm dynamic exchange and therefore cannot be crosslinked in the assay. However, this should not be clinically relevant, since IgG4 will not be present without other immunoglobulins G, even in patients with a predominant IgG4 response.
4. This assay is not intended for use as an aid in diagnosis, neither as a companion diagnostic test. The clinician should interpret test results in the context of other variables like sampling window, methodology, ADL levels and other clinical findings like disease activity.
5. The assay performance characteristics have not been established for matrices other than serum.
6. This product should only be used by appropriately trained laboratory professional personnel.

12. CALIBRATION RANGE

Calibration range: 6.0-1650.0 AU/ml

For samples > 1650.0 AU/ml retest the diluted sample in the Negative Control/Sample Diluent (REF 83607 - not supplied with the kit).

13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, are negative, with values < 6.0 AU/ml.

14. ANALYTICAL SPECIFICITY

3 samples were tested (1 Negative, 1 Low Positive and 1 High Positive) to which the following interfering agents were added:

Bilirubin (0.12 mg/ml – 0.45 mg/ml)
Triglycerides (3.75 mg/ml – 15 mg/ml)
Hemoglobin (2.5 mg/ml - 10 mg/ml)

The presence in the serum sample of the interfering substances listed above (except for bilirubin at a level of 0.45 mg/ml) does not affect the test result.

15. CROSS-REACTIONS

4 samples, not containing anti-ADL antibodies, but containing various concentrations of anti-IFX antibodies (33.0–1100.0 AU/ml), were tested.

No cross-reactions were found.

16. METHOD COMPARISON

In an experimentation 168 samples have been tested with Diesse kit and with Promonitor Anti-ADL ELISA kit. Data are summarized in the following table:

		Reference		
		+	-	Tot.
Diesse	+	44	1	45
	-	0	123	123
	Tot.	44	124	168

Percent Positive Agreement (~Diagnostic Sensitivity):

100.0% IC_{95%}: 92.0-99.9

Percent Negative Agreement: (~Diagnostic Specificity):

99.2% IC_{95%}: 95.5-99.8

Positive Predictive Value (PPV): 97.8% IC_{95%}: 95.6-100.0

Negative Predictive Value (NPV): 100.0% IC_{95%}: 100.0-100.0

The agreement between the two methods is excellent with a Cohen's Kappa of 0.97.

17. PRECISION AND REPEATABILITY

Sample	Within run		Between run	
	Mean (AU/ml)	CV%	Mean (AU/ml)	CV%
1	124.7	6.5	131.3	9.0
2	27.2	5.3	33.7	12.8
3	1187.2	5.5	1194.1	11.0
4	66.9	8.1	78.2	15.0
5	6.0	-	6.0	-
6	6.0	-	6.0	-

Sample	Between lots		Between instruments	
	Mean (AU/ml)	CV%	Mean (AU/ml)	CV%
1	121.7	13.2	123.7	13.8
2	33.0	14.3	32.2	13.6
3	1203.6	8.0	1123.0	10.5
4	72.8	4.9	73.4	13.8
5	6.0	-	6.0	-
6	6.0	-	6.0	-

18. LINEARITY

The linearity of the test was evaluated with 3 samples with known concentration of anti-ADL antibodies. For each sample serial dilutions were done; the linearity was demonstrated in the range 6.8-1650.0 AU/ml.

19. BIBLIOGRAPHY

1. Vaughn B.P. "A Practical Guide to Therapeutic Drug Monitoring of Biologic Medications for Inflammatory Bowel Disease" J Clin Med (2021) Oct 27;10(21):4990.
2. Gómez-Arango C. et All. "Cost-Effectiveness of Therapeutic Drug Monitoring-Guided Adalimumab Therapy in Rheumatic Diseases: A Prospective, Pragmatic Trial" Rheumatol Ther (2021) Sep;8(3):1323-1339
3. Sam M.J. et All. "Comparative Evaluation of 4 Commercially Available ELISA Kits for Measuring Adalimumab and Anti-adalimumab Antibodies" Ther Drug Monit (2020) Dec;42(6):821-828.
4. Mytych D. et All. "Evaluation of marketed kits for measurement of ABP 501, the first approved Adalimumab biosimilar, drug concentration and anti-drug antibody levels in patient serum" Ann Rheum Dis (2019) Abstract FRI0645.

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CHORUS Promonitor ANTI-ADALIMUMAB

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