CHORUS

INTRINSIC FACTOR

REF 86090

REF 86090/12

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17. BIBLIOGRAFIA

1. Carmel, R. (1992). *Reassessment of the relative prevalences of antibodies to gastric parietal cell and to intrinsic factor in patients with pernicious anaemia: influence of patient age and race*. Clin Exp Immunol, 89: 74-77 (11): 2617-20.
2. Oh R, Brown DL (2003). *Vitamin B12 deficiency*. Am Fam Physician, 67: 979-986.
3. Toh Ban-Hock, Alderuccio, F. (2004). *Pernicious Anaemia*. Autoimmunity 37 : 357-361.
4. Toh Ban-Hock, van Driel Ian R. (1997). *Pernicious Anaemia*. NEJM vol. 337 : 1441-1448.
5. NCBI Bookshelf. A service of the National Library of Medicine, National Institutes of Health. StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing; 2022 Jan-
Pernicious Anemia. Authors: Nicole M. Rodriguez1; Karen Shackelford.
6. Optimal management of pernicious anemia; Emmanuel Andres1,2Khalid Serraj; Department of Internal Medicine, Diabetes, and Metabolic Disorders, Competence Center of Autoimmune, Cytopenias in Adults, University, Hospital of Strasbourg, France; 3 Department of Internal Medicine, University Hospital of Oujda, Oujda, Morocco; Journal of Blood Medicine, 7 september 2012.
7. Diagnosis and classification of pernicious anemia Nicola Bizzaro Antonio Antico Autoimmunity Reviews Volume 13, Issues 4–5, April–May 2014, Pages 565-568.

INSTRUCTIONS FOR USE

CHORUS INTRINSIC FACTOR

For the semi-quantitative determination of anti-Intrinsic Factor IgG antibodies

Only for *in vitro* diagnostic use

1. USE

Enzyme immunoassay for the semi-quantitative determination of anti-Intrinsic Factor IgG class antibodies in human serum with a single-use device applied to Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Pernicious Anaemia represents the final stage of an autoimmune gastritis (type A gastritis) which leads to the destruction of the gastric mucosa. A small percentage of patients (10-15%) subject to autoimmune gastritis develops Pernicious Anaemia which in Western countries is the most widespread cause of vitamin B12 deficiency, due to the lack of the Intrinsic Factor, a glycoprotein necessary for the absorption of vitamin B12 in the gastrointestinal tract.

The Intrinsic Factor is produced by the cells of the gastric wall; Vitamin B12, in turn, is necessary for the formation of red blood cells.

Pernicious Anaemia does not generally appear before the age of 30 and is diagnosed at an average age of 60. In fact, a recent study revealed that almost 2% of individuals over 60 suffers from Pernicious Anaemia and that women are slightly more prone than men to suffer from it. The disease affects any ethnic group but is more common in populations with Scandinavian or Northern European origins.

Patients with Pernicious Anaemia develop autoantibodies to the parietal cells of the gastric mucosa and to the Intrinsic Factor. Antibodies against parietal cells have a sensitivity of 80-90% but can be detected up to 5% in healthy subjects; antibodies against Intrinsic Factor show a sensitivity of 50-70% with a 100% specificity in healthy subjects.

3. PRINCIPLE OF THE METHOD

The Chorus Intrinsic Factor device is ready to use for the determination of anti-Intrinsic Factor IgG antibodies, in Chorus/Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme Linked ImmunoSorbent Assay). The antigen is bound to the solid phase. Specific immunoglobulins bind to the antigen following incubation with diluted human serum. After washing to remove unreacted proteins, incubation is performed with the conjugated

compound consisting of human anti-immunoglobulin antibodies conjugated to horseradish peroxidase. The unbound conjugate is removed and the substrate for the peroxidase is added. The colour that develops is proportional to the concentration of specific antibodies in the serum under examination.

The single-use devices contain all the reagents to perform the test in the Chorus/Chorus TRIO instruments.

The results are expressed in Arbitrary Units (AU/ml).

4. WARNINGS AND PRECAUTIONS

ONLY FOR *IN VITRO* DIAGNOSTIC USE.

This kit contains materials of human origin which have been tested and found negative for both HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. Since no diagnostic test can offer a complete guarantee on the absence of infectious agents, any material of human origin must be considered potentially infected. All reagents and samples must be handled in accordance with the safety standards normally adopted in the laboratory.

Disposal of residues: the serum samples, calibrators and the strips used must be treated as infected residues, therefore disposed of in accordance with the provisions of current laws.

Personal safety warnings

1. Do not pipette with your mouth.
2. Use single-use gloves and eye protection when handling samples.
3. Wash your hands thoroughly after inserting the devices into the Chorus/Chorus TRIO instrument.
4. Refer to the Safety Data Sheet (available on request) for the safety characteristics of the reagents contained in the kit.
5. Neutralised acids and other liquid waste must be disinfected by adding sodium hypochlorite in a volume sufficient to obtain a final concentration of at least 1%. Exposure to 1% hypochlorite sodium for 30 minutes should be sufficient to ensure effective disinfection.
6. Any spillage of potentially infected materials must be removed immediately with absorbent paper and the contaminated area must be decontaminated, for example with 1% hypochlorite sodium, before continuing work. If an acid is present, sodium hypochlorite must not be used before the area has been dried. All materials used to decontaminate any accidental spillage, including gloves, must be discarded as potentially infected waste. Do not autoclave materials containing hypochlorite sodium.

Analytical warnings

Before use, bring the devices to be used to room temperature (18-30°C) and use within 60 minutes.

1. **Discard the devices with blue coloured substrate (well 4).**
2. When adding the sample to the well, check that it is perfectly distributed at the bottom.

3. Check the actual presence of reagents in the device and the integrity of the device itself. Do not use devices that lack any reagents and/or that have any foreign bodies in the reaction well when visually inspected.
4. The devices must be used in conjunction with the Chorus/Chorus TRIO instrument, strictly following the Instructions for Use and the User Manual of the instrument.
Use of the kit is only possible with an updated version of the software. Make sure that the software installed in the instrument matches or has a Release (date) above the one shown in the table published on the Diesse website
(<https://www.diesse.it/it/downloads/downloads/strumento:39/>)
5. Check that the Chorus/Chorus TRIO instrument is set correctly (see User Manual).
6. Do not alter the barcode placed on the handle of the device in order to allow it to be read correctly by the instrument.
7. Avoid the use of self-defrosting freezers for sample storage.
8. Defective barcodes can be manually entered into the instrument (see User Manual).
9. Do not expose the devices to strong lighting or hypochlorite vapours during storage and use.
10. The use of heavily haemolysed, lipemic, icteric, non-fully coagulated serum samples or of samples with microbial contamination may be a source of errors.
11. Do not use the device after its expiry date
12. **Check that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004)**

5. COMPOSITION OF THE KIT AND PREPARATION OF REAGENTS

The kit is sufficient for 36 tests (REF 86090).


The kit is sufficient for 12 tests (REF 86090/12).

DD DEVICES

6 packs of 6 devices each (REF 86090).

2 packs of 6 devices each (REF 86090/12).

Description:

Position 8: Space available for barcode label

Position 7: Empty

Position 6: MICROPLATE WELL

Sensitised with recombinant Intrinsic Factor.

Position 5: MICROPLATE WELL

Not sensitised.

Position 4: TMB SUBSTRATE

Content: Tetramethylbenzidine stabilised in citrate buffer

Position 3: DILUENT FOR SAMPLES

Content: Saline protein solution containing Proclin (0.1%)

Position 2: CONJUGATE

Content: Peroxidase-labelled anti-galactomannan monoclonal antibodies in phosphate buffer solution containing 0.05% phenol and 0.02% Bronidox.

Position 1: EMPTY WELL

Where the sample is transferred.

Usage: bring a bag to room temperature, open the bag, pick up the necessary devices; store the others in the bag containing the silica gel, let the air out and **seal** by pressing on the closure. Store at 2/8°C.

CALIBRATOR CALIBRATOR **1 x 0.175 ml**

Content: Diluted serum containing antibodies and preservative.

Liquid, ready to use

CONTROL + POSITIVE CONTROL **1 x 0.425 ml**

Content: Diluted serum containing antibodies and preservative.

Liquid, ready to use

OTHER MATERIAL REQUIRED, BUT NOT SUPPLIED:

- WASHING BUFFER AUTOIMMUNITY (REF 86004)
- CLEANING SOLUTION 2000 (REF 83609)
- SANITISING SOLUTION (REF 83604 - 83608)
- CHORUS NEGATIVE CONTROL/DILUENT SAMPLE (REF 83607)
- Chorus/Chorus TRIO instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test tubes, etc.
- Micropipettes capable of accurately withdrawing volumes of 50-200µl.
- Single-use gloves
- 5% sodium hypochlorite solution
- Containers for the collection of potentially infected materials

6. METHOD OF STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the event of an incorrect storage temperature, the calibration must be repeated and the correctness of the result checked using the control serum (see chapter 9: Test validation).

The expiry date is printed on each component and on the outer label of the package.

Reagents have limited stability after opening and/or preparation:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SAMPLE TYPE AND STORAGE

The sample type is serum obtained from blood collected by normal venous blood sampling and handled as required by standard laboratory procedures.

The consequences of using other body fluids are unknown.

Fresh serum can be kept for 4 days at 2/8°C; for longer storage periods, freeze at -20°C.

The sample can undergo up to a maximum of 3 thaws.
Avoid the use of self-defrosting freezers for sample storage. After thawing, carefully shake the sample before dosing.
Heat inactivation may give erroneous results. Sample quality can be seriously affected by microbial contamination which can lead to erroneous results.

8. PROCEDURE

1. Open the bag (side containing the snap closure), take the number of devices needed to perform the tests and store the others by closing the bag after letting the air out.
2. Visually check the status of the device according to the instructions in chapter 4 Analytical Warnings.
3. Add 50µl of undiluted serum to be analysed to well no. 1 of each device, at each batch change use a calibrator device.
4. Introduce the devices on the Chorus/Chorus TRIO instrument. Perform the calibration (if required) and the test as reported in the Instruction Manual of the instrument.

9. TEST VALIDATION

Use the positive control serum to verify the correctness of the result obtained, processing it as indicated in the User Manual of the instrument. If the instrument indicates that the control serum has a value falling outside the acceptability limit, the calibration must be performed again. Previous results will be corrected automatically.

If the control serum result continues to be outside the acceptability range, contact Scientific Support.

Phone: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE TEST

The Chorus/Chorus TRIO instrument provides the result in Arbitrary Units (AU/ml) calculated on the basis of a batch-dependent graph stored in the instrument.

The test on the serum under examination can be interpreted as follows:

POSITIVE: when the result is >18 AU/ml

NEGATIVE: when the result is <12 AU/ml

DOUBTFUL/EQUIVOCAL: when the result is between 12 and 18 AU/ml

In case of equivocal result, repeat the test. If the result remains equivocal, repeat the sampling.

11. TEST LIMITATIONS

All values obtained require a careful interpretation, without ever disregarding other indicators relating to the same patient.

The test, in fact, cannot be used on its own for a clinical diagnosis and the result of the test must be evaluated together with data from the patient's medical history and/or from other diagnostic investigations.

12. CALIBRATION RANGE

Calibration range 3.0-300.0 AU/ml.

For samples >300.0 AU/ml repeat the test, pre-diluting the sample with the Negative Control/Sample Diluent (REF 83607- not supplied with the kit).

13. ANALYTICAL SPECIFICITY

4 samples were tested (2 Negative, 2 Positive) to which the following interfering agents were added:

Rheumatoid factor (44 - 220 IU/ml)

Bilirubin (4.5 mg/dl - 45 mg/dl)

Triglycerides (10 mg/dl - 250 mg/dl)

Haemoglobin (5 mg/ml - 30 mg/ml)

The presence of the above interfering substances in the test serum does not alter the test result (with the exception of haemoglobin).

14. CROSS-REACTIVITY

39 samples, positive for Herpes virus (3), Rubella virus (4), Helicobacter pylori (4), Varicella zoster (3), Cytomegalovirus (4), Epstein-Barr virus (4), Toxoplasma gondii (6) and ANA (11) were tested.

No significant cross-reactions were detected.

15. COMPARISON STUDIES

In one trial, 309 samples were analysed with a Diesse kit and another commercially available kit

The trial data are outlined below:

		Reference		
		+	-	Total
Diesse	+	40	0	40
	-	0	269	269
	Total	40	269	309

Positive Agreement Rate: (~ Diagnostic Sensitivity):

100.0% CI_{95%}: 91.2-99.9

Negative Agreement Rate: (~ Diagnostic Specificity):

100.0% CI_{95%}: 98.6-100.0

Positive Predictive Value (PPV): 100.0% CI_{95%}: 100.0-100.0

Negative Predictive Value (NPV): 100.0% CI_{95%}: 100.0-100.0

The degree of agreement between the two methods is excellent with a value of K (Cohen's coefficient) of 1.

16. PRECISION

Sample	Within the session		Between sessions	
	Mean (AU/ml)	CV%	Mean (AU/ml)	CV%
1	66.6	4.5	55.2	13.0
2	14.2	3.4	12.1	11.2
3	3.4	11.5	3.5	12.0

Sample	Between batches		Between instruments	
	Mean (AU/ml)	CV%	Mean (AU/ml)	CV%
1	3.4	9.1	57.0	9.1
2	57.5	11.5	12.3	14.3
3	20.8	8.8	3.7	13.5

17. BIBLIOGRAPHY

1. Carmel, R. (1992). *Reassessment of the relative prevalences of antibodies to gastric parietal cell and to intrinsic factor in patients with pernicious anaemia: influence of patient age and race*. Clin Exp Immunol, 89: 74-77 (11): 2617-20.
2. Oh R, Brown DL (2003). *Vitamin B12 deficiency*. Am Fam Physician, 67: 979-986.
3. Toh Ban-Hock, Alderuccio, F. (2004). *Pernicious Anaemia*. Autoimmunity 37 : 357-361.
4. Toh Ban-Hock, van Driel Ian R. (1997). *Pernicious Anaemia*. NEJM vol. 337 : 1441-1448.
5. NCBI Bookshelf. A service of the National Library of Medicine, National Institutes of Health. StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing; 2022 Jan-
Pernicious Anemia. Authors: Nicole M. Rodriguez¹; Karen Shackelford.
6. Optimal management of pernicious anemia; Emmanuel Andres^{1,2}Khalid Serraj; Department of Internal Medicine, Diabetes, and Metabolic Disorders, Competence Center of Autoimmune, Cytopenias in Adults, University, Hospital of Strasbourg, France; 3 Department of Internal Medicine, University Hospital of Oujda, Oujda, Morocco; Journal of Blood Medicine, 7 september 2012.
7. Diagnosis and classification of pernicious anemia Nicola Bizzaro Antonio Antico Autoimmunity Reviews Volume 13, Issues 4–5, April–May 2014, Pages 565–568.

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17. BIBLIOGRAFIE

- Carmel, R. (1992). Reassessment of the relative prevalences of antibodies to gastric parietal cell and to intrinsic factor in patients with pernicious anaemia: influence of patient age and race. Clin Exp Immunol, 89: 74-77 (11): 2617-20.
- Oh R, Brown DL (2003). Vitamin B12 deficiency. Am Fam Physician, 67: 979-986.
- Toh Ban-Hock, Alderuccio, F. (2004). Pernicious Anaemia. Autoimmunity 37 : 357-361.
- Toh Ban-Hock, van Driel Ian R. (1997). Pernicious Anaemia. NEJM vol. 337 : 1441-1448.
- NCBI Bookshelf. A service of the National Library of Medicine, National Institutes of Health. StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing; 2022 Jan-
Pernicious Anemia. Authors: Nicole M. Rodriguez1; Karen Shackelford.
- Optimal management of pernicious anemia; Emmanuel Andres1,2Khalid Serraj; Department of Internal Medicine, Diabetes, and Metabolic Disorders, Competence Center of Autoimmune, Cytopenias in Adults, University, Hospital of Strasbourg, France; 3 Department of Internal Medicine, University Hospital of Oujda, Oujda, Morocco; Journal of Blood Medicine, 7 september 2012.
- Diagnosis and classification of pernicious anemia Nicola Bizzaro Antonio Antico Autoimmunity Reviews Volume 13, Issues 4–5, April–May 2014, Pages 565-568.

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17. LITERATUR

1. Carmel, R. (1992). Reassessment of the relative prevalences of antibodies to gastric parietal cell and to intrinsic factor in patients with pernicious anaemia: influence of patient age and race. Clin Exp Immunol, 89: 74-77 (11): 2617-20.
2. Oh R, Brown DL (2003). Vitamin B12 deficiency. Am Fam Physician, 67: 979-986.
3. Toh Ban-Hock, Alderuccio, F. (2004). Pernicious Anaemia. Autoimmunity 37 : 357-361.
4. Toh Ban-Hock, van Driel Ian R. (1997). Pernicious Anaemia. NEJM vol. 337 : 1441-1448.
5. NCBI Bookshelf. A service of the National Library of Medicine, National Institutes of Health. StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing; 2022 Jan-
Pernicious Anemia. Authors: Nicole M. Rodriguez¹; Karen Shackelford.
6. Optimal management of pernicious anemia; Emmanuel Andres^{1,2}Khalid Serraj; Department of Internal Medicine, Diabetes, and Metabolic Disorders, Competence Center of Autoimmune, Cytopenias in Adults, University, Hospital of Strasbourg, France; 3 Department of Internal Medicine, University Hospital of Oujda, Oujda, Morocco; Journal of Blood Medicine, 7 september 2012.
7. Diagnosis and classification of pernicious anemia Nicola Bizzaro Antonio Antico Autoimmunity Reviews Volume 13, Issues 4–5, April–May 2014, Pages 565-568.

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- Carmel, R. (1992). Reassessment of the relative prevalences of antibodies to gastric parietal cell and to intrinsic factor in patients with pernicious anaemia: influence of patient age and race. Clin Exp Immunol, 89: 74-77 (11): 2617-20.
- Oh R, Brown DL (2003). Vitamin B12 deficiency. Am Fam Physician, 67: 979-986.
- Toh Ban-Hock, Alderuccio, F. (2004). Pernicious Anaemia. Autoimmunity 37 : 357-361.
- Toh Ban-Hock, van Driel Ian R. (1997). Pernicious Anaemia. NEJM vol. 337 : 1441-1448.
- NCBI Bookshelf. A service of the National Library of Medicine, National Institutes of Health. StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing; 2022 Jan-Pernicious Anemia. Authors: Nicole M. Rodriguez1; Karen Shackelford.
- Optimal management of pernicious anemia; Emmanuel Andres1,2Khalid Serraj; Department of Internal Medicine, Diabetes, and Metabolic Disorders, Competence Center of Autoimmune, Cytopenias in Adults, University, Hospital of Strasbourg, France; 3 Department of Internal Medicine, University Hospital of Oujda, Oujda, Morocco; Journal of Blood Medicine, 7 september 2012.
- Diagnosis and classification of pernicious anemia Nicola Bizzaro Antonio Antico Autoimmunity Reviews Volume 13, Issues 4–5, April–May 2014, Pages 565-568.

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1. Carmel, R. (1992). Reassessment of the relative prevalences of antibodies to gastric parietal cell and to intrinsic factor in patients with pernicious anaemia: influence of patient age and race. Clin Exp Immunol, 89: 74-77 (11): 2617-20.
2. Oh R, Brown DL (2003). Vitamin B12 deficiency. Am Fam Physician, 67: 979-986.
3. Toh Ban-Hock, Alderuccio, F. (2004). Pernicious Anaemia. Autoimmunity 37 : 357-361.
4. Toh Ban-Hock, van Driel Ian R. (1997). Pernicious Anaemia. NEJM vol. 337 : 1441-1448.
5. NCBI Bookshelf. A service of the National Library of Medicine, National Institutes of Health. StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing; 2022 Jan-
Pernicious Anemia. Authors: Nicole M. Rodriguez¹; Karen Shackelford.
6. Optimal management of pernicious anemia; Emmanuel Andres^{1,2}Khalid Serraj; Department of Internal Medicine, Diabetes, and Metabolic Disorders, Competence Center of Autoimmune, Cytopenias in Adults, University, Hospital of Strasbourg, France; 3 Department of Internal Medicine, University Hospital of Oujda, Oujda, Morocco; Journal of Blood Medicine, 7 september 2012.

7. Diagnosis and classification of pernicious anemia Nicola Bizzaro Antonio Antico Autoimmunity Reviews Volume 13, Issues 4–5, April–May 2014, Pages 565-568.

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