CHORUS

Zika IgM Capture

REF 81161

DIESSE

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16. BIBLIOGRAFIA

1. Song B.H. et al.; "Zika virus: History, epidemiology, transmission, and clinical presentation"; *Journal of Neuroimmunology*, **2017**; 15(308):50-64
2. Javed F., Manzoor K.N., et al.; "Zika virus: what we need to know?"; *Journal of Basic Microbiology*, **2018**; 58(1):3-16
3. Pierson T.C., Diamond M.S.; "The emergence of Zika virus and its new clinical syndromes"; *Nature*, **2018**; 560(7720):573-581
4. Plourde A.R., Bloch E.M.; "A Literature Review of Zika Virus"; *Emerging Infectious Diseases*, **2016**; 22(7):1185-1192
5. Musso D., Gubler D.J.; "Zika Virus"; *Clinical Microbiology Reviews*, **2016**; 29(3):487-524
6. Calvet G. et al.; "Detection and sequencing of Zika virus from amniotic fluid of fetuses with microcephaly in Brazil: a case study"; *The Lancet Infectious Diseases*, **2016**; 16(6):653-660
7. Singh T., Lopez et al.; "Efficient transplacental IgG transfer in women infected with Zika virus during pregnancy"; *PLoS Neglected Tropical Diseases*, **2019**; 13(8):e0007648

INSTRUCTIONS FOR USE

CHORUS Zika IgM Capture

For the qualitative determination of anti-Zika virus IgM class antibodies

For in vitro diagnostic use only

1. USE

Enzyme immunoassay method for the qualitative determination of anti-Zika virus IgM class antibodies in human serum with a single-use device applied to CHORUS TRIO instruments.

2. INTRODUCTION

Zika virus is an Arbovirus belonging to the genus Flavivirus, of the Flaviviridae family, transmitted by Aedes mosquitoes, but unlike other arboviruses, can be transmitted during sexual intercourse and vertically from the mother to the foetus through the placenta and amniotic fluid, resulting in teratogenic effects on the embryo: the lack of prevention mechanisms against ZIKV makes this pathogen a re-emergent cause of congenital infections. An important epidemic is the one that took place in 2015 in Brazil: in 2016 the WHO International Health Regulation Emergency Committee, declared Zika virus infection as an international public health emergency. Zika virus disease is an acute viral disease with clinical manifestations similar to those of Dengue, which tends to self-resolve within 4-6 days. However, what makes the Zika virus a threat to human health are neurological complications, which are the main cause of disability and death induced by this pathogen. Direct detection of the virus or components of the virus (NS1) is possible only during viraemia (maximum 1 week); specific antibodies are clearly detectable a few days after onset of the symptoms. Detection of specific IgM antibodies or a significant increase in specific IgG titre is evidence of acute ZIKV infection. Possible cross-reactions with other flaviviruses such as Dengue, West Nile, Yellow Fever or TBE virus should be emphasised.

3. PRINCIPLE OF THE METHOD

The CHORUS Zika IgM Capture device is ready to use for the determination of anti-Zika virus **IgM** antibodies in CHORUS TRIO instruments.

The test is based on the capture ELISA principle (Enzyme Linked ImmunoSorbent Assay). Human anti-IgM monoclonal antibodies are bound to the solid phase. IgM immunoglobulins bind to anti-IgM antibodies following incubation with a diluted sample.

After washing to remove proteins that have not reacted, incubation is carried out with the Zika NS1 antigen bound to specific anti-Zika virus monoclonal antibodies conjugated with horseradish peroxidase.

The unbound conjugate is removed and the substrate for the peroxidase is added.

The enzymatic reaction is subsequently stopped by adding the Stop Solution which causes the solution to turn yellow.

The colour that develops is proportional to the concentration of specific antibodies in the serum under examination.

The single-use devices contain all the reagents to perform the test when applied to the CHORUS TRIO instruments.

The results are expressed in Index (sample OD/cut-off OD).

4. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

This kit contains materials of human origin which have been tested and found negative for both HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. Since no diagnostic test can offer a complete guarantee on the absence of infectious agents, any material of human origin must be considered potentially infected. All reagents and samples must be handled in accordance with the safety standards normally adopted in the laboratory.

Disposal of residues: the serum samples, calibrators and strips used must be treated as infected residues, therefore disposed of in accordance with the provisions of current laws.

Personal safety warnings

1. Do not pipette with your mouth.
2. Use single-use gloves and eye protection when handling samples.
3. Wash your hands thoroughly after inserting the devices into the CHORUS TRIO instrument.
4. Refer to the Safety Data Sheet (available on request) for the safety characteristics of the reagents contained in the kit.
5. Neutralised acids and other liquid waste must be disinfected by adding sodium hypochlorite in a volume sufficient to obtain a final concentration of at least 1%. Exposure to 1% sodium hypochlorite for 30 minutes should be sufficient to ensure effective disinfection.
6. Any spillage of potentially infected materials must be removed immediately with absorbent paper and the contaminated area must be decontaminated, for example with 1% sodium hypochlorite, before continuing work. If an acid is present, sodium hypochlorite must not be used before the area has been dried.

All materials used to decontaminate any accidental spillage, including gloves, must be discarded as potentially infected waste.

Do not autoclave materials containing sodium hypochlorite.

Analytical warnings

Before use, bring the devices to be used to room temperature (18-30°C) and use within 60 minutes.

1. **Discard the devices with blue coloured substrate (well 4).**
2. When adding the sample to the well, check that it is perfectly distributed at the bottom.
3. Check the actual presence of reagents in the device and the integrity of the device itself. Do not use devices that lack any reagents and/or that have any foreign bodies in the reaction well when visually inspected.

4. The devices must be used in conjunction with the CHORUS TRIO instrument, strictly following the Instructions for Use and the User Manual of the instrument.

Use of the kit is only possible with an updated version of the software. Make sure that the software installed in the instrument matches or has a release (Rel.) date later than the one shown in the table published on the Diesse website

(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)


5. Check that the CHORUS TRIO instrument is set correctly (see User Manual).
6. Do not alter the barcode placed on the handle of the device in order to allow it to be read correctly by the instrument.
7. Avoid the use of self-defrosting freezers for sample storage.
8. Defective barcodes can be manually entered into the instrument (see User Manual).
9. Do not expose the devices to strong lighting or hypochlorite vapours during storage and use.
10. The use of heavily haemolysed, lipemic, icteric, non-fully coagulated serum samples or of samples with microbial contamination may be a source of errors.
11. Do not use the device after its expiry date.
12. **Check that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests

DD DEVICES 6 packs of 6 devices each

Description:

Position 9: Space available for barcode label

Position 8: LYOPHILIC ANTIGEN

Content: Zika Virus NS1 recombinant protein

Position 7: MICROPLATE WELL

Sensitised with human anti-IgM monoclonal antibodies

Position 6: MICROPLATE WELL

Not sensitised

Position 5: STOP SOLUTION

0.3 M sulphuric acid solution

Position 4: TMB SUBSTRATE

Content: Tetramethylbenzidine stabilised in citrate buffer

Position 3: DILUENT FOR SAMPLES

Content: Protein solution in phosphate buffer with 0.09% sodium azide and colouring

Position 2: CONJUGATE

Content: anti-Zika virus monoclonal antibodies marked with peroxidase in phosphate buffer solution containing proteins and 0.1% phenol

Position 1: EMPTY WELL

Where the sample is transferred.

Usage: bring a bag to room temperature, open the bag, take out the necessary devices; store the others in the bag containing the silica gel, let the air out and seal by pressing on the closure. Store at 2-8°C.

CALIBRATOR

1 x 0.450 ml

Content: protein solution containing proteins able to specifically bind the solid phase and immunocomplex, and preservative.

Liquid, ready to use.

CONTROL +

POSITIVE CONTROL

1 x 0.900 ml

Content: protein solution containing proteins able to specifically bind the solid phase and immunocomplex, and preservative.

Liquid, ready to use.

OTHER MATERIAL REQUIRED, BUT NOT SUPPLIED:

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- CHORUS TRIO instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test tubes, etc.
- Micropipettes capable of accurately withdrawing volumes of 50-200µl.
- Single-use gloves
- 5% sodium hypochlorite solution
- Containers for the collection of potentially infected materials

6. STORAGE METHOD AND REAGENT STABILITY

Reagents must be stored at 2-8°C. In the event of an incorrect storage temperature, the calibration must be repeated and the correctness of the result checked using the control serum (see section 9: Test validation).

The expiry date is printed on each component and on the outer label of the package.

Reagents have limited stability after opening and/or preparation:

DEVICES	8 weeks at 2-8°C
CALIBRATOR	8 weeks at 2-8°C
POSITIVE CONTROL	8 weeks at 2-8°C

7. SAMPLE TYPE AND STORAGE

The sample type is serum obtained from blood collected by normal venous blood sampling and handled as required by standard laboratory procedures.

The consequences of using other body fluids are unknown.

Fresh serum can be kept for 5 days at 2-8°C; for longer storage periods, freeze at -20°C.

Avoid repeated thawing.

Avoid the use of self-defrosting freezers for sample storage. After thawing, carefully shake the sample before dosing.

Inactivation to heat might give erroneous results.

Sample quality can be seriously affected by microbial contamination which can lead to erroneous results.

8. PROCEDURE

The CHORUS Zika IgM Capture test (REF. 81161) can be carried out simultaneously **ONLY** with the following kits:

CHORUS CYTOMEGALOVIRUS IgM CAPTURE	REF 81012
CHORUS TOXOPLASMA IgA CAPTURE	REF 81044
CHORUS TOXOPLASMA IgM CAPTURE	REF 81042

CHORUS ZIKA IgG
CHORUS ZIKA NS1

REF 81158
REF 81159

1. Open the bag (side containing the snap closure), take the number of devices needed to perform the tests and store the others by closing the bag after letting the air out.
2. Visually check the status of the device according to the instructions in section 4 Analytical Warnings.
3. Dispense into well no. 1 of each device:

SAMPLE 50 µl/device
CALIBRATOR 130 µl/device
POSITIVE CONTROL 130 µl/device

At each batch change use a device for the calibrator.

4. Introduce the devices on the CHORUS TRIO instrument. Perform the calibration (if required) and the test as indicated in the User Manual of the instrument.

9. TEST VALIDATION

Use the positive control serum to verify the correctness of the result obtained, processing it as indicated in the User Manual of the instrument. If the instrument indicates that the control serum has a value outside the acceptability limit, the calibration must be performed again. Previous results will be corrected automatically.

If the control serum result continues to be outside the acceptability range, contact Scientific Support.

Phone: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE TEST

The CHORUS TRIO instrument provides the result in Index (sample OD/cut-off OD).

The test on the serum under examination can be interpreted as follows:

POSITIVE: when the result is >1.1

NEGATIVE: when the result is <0.9

AMBIGUOUS: when the result is between 0.9 and 1.1

In case of an ambiguous result, repeat the test. If the result remains ambiguous, repeat the sampling.

11. TEST LIMITATIONS

All values obtained require a careful interpretation, without ever disregarding other indicators relating to the same patient. The test, in fact, cannot be used on its own for a clinical diagnosis and the result obtained must always be evaluated together with data from the patient's medical history and/or from other diagnostic investigations.

12. ANALYTICAL SPECIFICITY

5 samples were tested (2 Negative, 1 at Cut-Off and 2 Positive) to which the following interfering agents were added:

Rheumatoid Factor (44-220 IU/mL)
Bilirubin (4.5-45 mg/dL)
Triglycerides (10-250 mg/dL)
Hemoglobin (5-30 mg/mL)

The presence of the above interfering substances in the test serum does not alter the test result.

13. CROSS-REACTIVITY

22 samples positive for Dengue (12) and West Nile (10) were tested.

Cross-reactions were detected with only 2 samples positive to West Nile Virus.

14. COMPARISON STUDIES

In one trial, 159 samples were analysed with a Diesse kit and with another commercially available kit.

The trial data are outlined below:

		Predicate		
		+	-	Total
Diesse	+	14	2	16
	-	1	142	143
	Total	15	144	159

Positive percentage agreement (~ diagnostic sensitivity): 93.3
CI_{95%}: 70.0 – 98.7

Negative percentage agreement (~ diagnostic specificity):
98.6 CI_{95%}: 95.1 – 99.6

Positive Predictive Value (PPV): 87.5 CI_{95%}: 82.4 - 92.6

Negative Predictive Value (NPV): 99.3 CI_{95%}: 98.0 - 100.0

The degree of agreement between the two methods is excellent with a K value (Cohen's coefficient) of 0.88

15. ACCURACY AND REPEATABILITY

Sample	Within the session		Between sessions	
	Mean (Index)	CV%	Mean (Index)	CV%
1	4.2	11.9	4.0	14.5
2	0.4	15.0	0.3	13.3
3	0.3	13.3	0.3	3.3
4	1.1	10.9	1.1	10.9
5	1.7	14.7	1.8	12.8

Sample	Between batches		Between instruments	
	Mean (Index)	CV%	Mean (Index)	CV%
1	4.2	13.9	4.1	14.6
2	0.4	13.4	0.4	12.5
3	0.4	12.4	0.4	15.0
4	1.0	7.4	1.1	11.8
5	1.8	13.2	1.9	12.6

16. BIBLIOGRAPHY

1. Song B.H. et al.; "Zika virus: History, epidemiology, transmission, and clinical presentation"; *Journal of Neuroimmunology*, **2017**; 15(308):50-64
2. Javed F., Manzoor K.N., et al.; "Zika virus: what we need to know?"; *Journal of Basic Microbiology*, **2018**; 58(1):3-16
3. Pierson T.C., Diamond M.S.; "The emergence of Zika virus and its new clinical syndromes"; *Nature*, **2018**; 560(7720):573-581

4. Plourde A.R., Bloch E.M.; "A Literature Review of Zika Virus"; *Emerging Infectious Diseases*, **2016**; 22(7):1185-1192
5. Musso D., Gubler D.J.; "Zika Virus"; *Clinical Microbiology Reviews*, **2016**; 29(3):487-524
6. Calvet G. et al.; "Detection and sequencing of Zika virus from amniotic fluid of fetuses with microcephaly in Brazil: a case study"; *The Lancet Infectious Diseases*, **2016**; 16(6):653-660
7. Singh T., Lopez et al.; "Efficient transplacental IgG transfer in women infected with Zika virus during pregnancy"; *PLoS Neglected Tropical Diseases*, **2019**; 13(8):e0007648

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3. Pierson T.C., Diamond M.S.; "The emergence of Zika virus and its new clinical syndromes"; *Nature*, **2018**; 560(7720):573-581
4. Plourde A.R., Bloch E.M.; "A Literature Review of Zika Virus"; *Emerging Infectious Diseases*, **2016**; 22(7):1185-1192
5. Musso D., Gubler D.J.; "Zika Virus"; *Clinical Microbiology Reviews*, **2016**; 29(3):487-524
6. Calvet G. et al.; "Detection and sequencing of Zika virus from amniotic fluid of fetuses with microcephaly in Brazil: a case study"; *The Lancet Infectious Diseases*, **2016**; 16(6):653-660
7. Singh T., Lopez et al.; "Efficient transplacental IgG transfer in women infected with Zika virus during pregnancy"; *PLoS Neglected Tropical Diseases*, **2019**; 13(8):e0007648

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