# VES-TEC VACU-TEC VACU-CODE

**DIESSE**

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#### *For In vitro Diagnostic Use Only*

#### 1. INTENDED USE

**VES-TEC:** Single-use test tube containing an anticoagulant solution (sodium citrate) to be used for whole blood collection for the determination of Erythrocyte Sedimentation Rate (ESR) in DIESSE VES instruments.

**VACU-TEC, VACU CODE:** Single-use vacuum tube containing anticoagulant (sodium citrate) to be used for manual sampling and whole blood collection for the determination of Erythrocyte Sedimentation Rate (ESR) in DIESSE VES instruments.

All the tubes must be used by professional laboratory users only.

#### 2. INTRODUCTION AND PRINCIPLE OF THE METHOD

VES-TEC, VACU-TEC e VACU-CODE tubes allow to collect the volume of blood ( $1.0 \pm 0.2$  ml) necessary for the following determination of Erythrocytes Sedimentation Rate (ESR), performed using the manual VES-RACK system or the instruments of the DIESSE VES-Matic Original line (reading system)

#### 3. WARNINGS AND PRECAUTIONS

##### **FOR IN VITRO DIAGNOSTIC USE ONLY**

**Waste disposal:** the tubes, once used, must be treated as infectious residuals and eliminated according to law.

#### Health and Safety Information

1. All material of human origin must be considered potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling samples.
2. The presence of the outer tube is intended as a protection for the operator from accidental breakages of the tube during the sample collection.
3. Wear disposable gloves, proper labcoat and eye protection while handling specimens.
4. Wash hands thoroughly after finishing the blood collection.

5. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

#### Analytical Precautions

1. Store the tubes in their bag.
2. Handle the bag with care: it contains FRAGILE material.
3. **Before use, carefully check the tube for its integrity.**
4. **For gamma irradiated vacuum tubes, proper filling ensures sterility is maintained until the expiration date**
5. **Discard tubes showing cracks.**
6. **It is advisable not to use tubes if the bag is damaged.**
7. **Avoid the tubes from coming into contact with solvents such as ethyl alcohol.**
8. Check that the shape of the blood-collection holder is compatible with the DIESSE evacuated tubes.
9. To avoid hemolysis and opening of the cap during blood sampling, use 20G needles.
10. Avoid hook-shaped needles.
11. To be certain of perforating the rubber cap correctly, the holder needle must not be less than 15 mm in length.
12. Do not use the device after the expiry date.

#### 4. COMPOSITION

**VES-TEC:** Polystyrene tubes with polyethylene caps inserted in an outer polyethylene tube onto which the label is applied.

**VACU-TEC:** Evacuated metacrylate tubes with a silicone rubber cap, inserted in an outer polyethylene tube onto which the label is applied.

**VACU-CODE:** Evacuated metacrylate tubes with a T-cross-section, with silicone rubber caps and bar-code label.

All the tubes contain  $0.37 \pm 0.03$  ml of 0.105 M sodium citrate and are packaged in triplex-layer bags each containing 100 tubes or multiples of 100.

The sterile tubes are irradiated with gamma rays.

#### **MATERIALS REQUIRED BUT NOT PROVIDED**

- NEEDLES FOR BLOOD COLLECTION
- HOLDER FOR VACUUM-COLLECTION.

#### 5. STORAGE AND STABILITY OF REAGENTS

**The tubes should be stored between 0-28°C, and are stable up to the expiry date.**

**After opening, the tubes are stable for 2 months.**

## 6. SPECIMEN COLLECTION AND STORAGE

The sample is composed of whole blood, collected in accordance with the International Committee for Standardization in Hematology (ICSH) recommendations:

- the blood must be drawn in 30 seconds to avoid venostasis;
- it is recommended to use a suitable size needle for blood collection (Ø 0.7 mm).

Use the sample within 2 hours of collection, or 12 hours if stored at 2/8°C. Bring to room temperature before performing the test.

## 7. ASSAY PROCEDURE

### Method to be used with VES-TEC tube

The sample is drawn using needle and syringe.

- Fill the tube with blood up to the mark
- Mix carefully with the anticoagulant by repeatedly turning upside down
- Transfer the tube to the reading system

### Method to be used with VACU-TEC and VACU-CODE tubes

The sample is drawn using a holder.

- Insert the needle in the patient's arm
- Place the tube in the holder and perforate the center of the cap by pressing the bottom of the tube. It will fill automatically.
- Mix the blood with the anticoagulant by repeatedly turning upside down.  
Safety collar tube: after mixing the sample, take the tube out of the outer tube holding on the top part and untwisting the bottom part.
- Transfer to the reading system.

During the whole process, all the conditions of maximum protection of the operator must be maintained.

Contact the Scientific Support, in case of problems:

Tel: 0039 0577 319554  
Fax: 0039 0577 366605  
email: [scientificsupport@diesse.it](mailto:scientificsupport@diesse.it)

## 8. INCIDENT REPORTING

If a serious incident has occurred in relation to this device in the European Union market territory, please report it to the manufacturer and the competent authority of its Member State without delay.

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