CHORUS suPAR Extended

REF 81414

DIESSE

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17 BIBLIOGRAFIA

1. Hartmann Rasmussen L. et All. Soluble urokinase plasminogen activator receptor (suPAR) as a biomarker of systemic chronic inflammation. Front Immunol. (2021)
2. Lafon T. et All. Prognostic performance of endothelial biomarkers to early predict clinical deterioration of patients with suspected bacterial infection and sepsis admitted to the emergency department. Ann Intensive Care. (2020)
3. Desmedt S et al. The intriguing role of soluble urokinase receptor in inflammatory diseases. Crit Rev Clin Lab Sci. 2017 Mar;54(2):117-133
4. Rasmussen LJH et al. Combining National Early Warning Score with Soluble Urokinase Plasminogen Activator Receptor (suPAR) Improves Risk Prediction in Acute Medical Patients: A Registry-Based Cohort Study. Crit Care Med. 2018 Dec;46:1961-8.
5. Rasmussen LJH, et al. Soluble urokinase plasminogen activator receptor (suPAR) in acute care: a strong marker of disease presence and severity, readmission and mortality. A retrospective cohort study. Emerg Med J. 2016 Nov;33:769-75.
6. Smith HW, Marshall CJ. Regulation of cell signalling by uPAR. Nat Rev Mol Cell Biol. 2010;11:23-36.
7. Enocsson H, Lukic T, Ziegelsch M, Kastbom A. Serum levels of the soluble urokinase plasminogen activator receptor (suPAR) correlates with disease activity in early rheumatoid arthritis and reflects joint damage over time. Transl Res. 2021 Jun;232:142-149.

INSTRUCTIONS FOR USE

CHORUS Extended suPAR

For the quantitative determination of suPAR

Only for *in vitro* diagnostic use

1. USE

CHORUS Extended suPAR is an enzyme immunoassay kit for the quantitative determination of suPAR. The test is conducted on human serum and plasma using a single-use device applied to the CHORUS TRIO instruments.

suPAR is a non-specific parameter of an inflammatory state, used as an aid in monitoring the patient's physiological or pathological state.

It must only be used by professional laboratory personnel.

2. INTRODUCTION

suPAR is the soluble form of the urokinase plasminogen activator receptor (uPAR) expressed on the cell membrane of immune cells, endothelial cells and smooth muscle cells and is upregulated in inflammation sites. uPAR interacts and cooperates with many ligands and receptors, mainly integrin, to facilitate intracellular signalling, cellular migration, cellular adhesion and tissue remodelling.

suPAR is released during immune inflammation or activation and although it is not specific for a disease, the circulating levels of suPAR reflect the severity and prognostic outcome of many infectious, inflammatory and autoimmune diseases.

One of the most interesting applications of the marker is that it can predict early clinical deterioration due to suspected infections and evaluate the severity of chronic inflammation.

3. METHOD PRINCIPLE

The CHORUS Extended suPAR test is based on the principle of the ELISA sandwich method.

The suPAR in the test sample binds to the specific antibodies present in the solid phase and in the conjugate. The unbound components are eliminated by washing and the bound enzyme activity is evaluated colorimetrically through transformation of a chromogenic substrate. The intensity of the colour that develops is directly proportional to the concentration of suPAR in the sample under test.

The single-use devices contain all the reagents to perform the test when applied to the CHORUS TRIO instruments.

The results are expressed in ng/ml.

4. PRECAUTIONS

ONLY FOR *IN VITRO* DIAGNOSTIC USE

Since no diagnostic test can offer a complete guarantee on the absence of infectious agents, any material of human origin must be considered potentially infected. All reagents and samples must be handled in accordance with the safety standards normally adopted in the laboratory.

Disposal of residues: the samples, calibrators and strips used must be treated as infected residues, therefore disposed of in accordance with the provisions of current laws.

Personal safety warnings

1. Do not pipette with your mouth.
2. Use single-use gloves and eye protection when handling samples.
3. Wash your hands thoroughly after inserting the devices into the CHORUS TRIO instrument.
4. Refer to the Safety Data Sheet (available on request) for the safety characteristics of the reagents contained in the kit.
5. Neutralised acids and other liquid waste must be disinfected by adding sodium hypochlorite in a volume sufficient to obtain a final concentration of at least 1%. Exposure to 1% sodium hypochlorite for 30 minutes should be sufficient to ensure effective disinfection.
6. Any spillage of potentially infected materials must be removed immediately with absorbent paper and the contaminated area must be decontaminated, for example with 1% sodium hypochlorite, before continuing work. If an acid is present, sodium hypochlorite must not be used before the area has been dried.

All materials used to decontaminate any accidental spillage, including gloves, must be discarded as potentially infected waste.

Do not autoclave materials containing sodium hypochlorite.

Analytical warnings

Before use, bring the devices to be used to room temperature (18-30°C) and use within 60 minutes.

1. **Discard the devices with blue coloured substrate (well 4).**
2. When adding the sample to the well, check that it is perfectly distributed at the bottom.
3. Check the actual presence of the reagents in the device and the integrity of the device itself. Do not use devices that lack any reagents and/or that have any foreign bodies in the reaction well when visually inspected.
4. The devices must be used in conjunction with the CHORUS TRIO instrument, strictly following the Instructions for Use and the User Manual of the instrument.

Use of the kit is only possible with an updated version of the software. Make sure that the software installed in the instrument matches or has a release (Rel.) date later than the one shown in the table published on the Diesse website.

(<https://www.diesse.it/en/downloads/downloads/strumento:39/>)


5. Check that the CHORUS TRIO instrument is set correctly (see User Manual).
6. Do not alter the barcode placed on the handle of the device in order to allow it to be read correctly by the instrument.
7. Avoid the use of self-defrosting freezers for sample storage.
8. Defective barcodes can be manually entered into the instrument (see User Manual).
9. Do not expose the devices to strong lighting or hypochlorite vapours during storage and use.
10. Do not use haemolysed, lipemic or icteric samples with a higher interfering concentration than that tested (according to the indications given in the "Analytical Specificity" chapter).
11. Do not use the device after its expiry date.
12. Check that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packs of 6 devices each

Description:

Position 9: Space available for barcode label.

Position 8: EMPTY

Position 7: MICROPLATE WELL

Sensitised with anti-suPAR monoclonal antibody.

Position 6: NON-SENSITISED WELL

Position 5: STOP SOLUTION

Content: 0.3 M sulphuric acid solution.

Position 4: TMB SUBSTRATE

Content: Buffer stabilized tetramethylbenzidine and H₂O₂.

Position 3: DILUENT FOR SAMPLES

Content: Protein solution in phosphate buffer containing 0.05% phenol, 0.02% Bronidox and ascitic fluid (0.35 mg/ml).

Position 2: CONJUGATE

Content: Peroxidase-labelled anti-suPAR monoclonal antibody, in phosphate buffer solution containing 0.05% phenol and 0.02% Bronidox.

Position 1: EMPTY WELL

Where the sample is transferred.

Usage: bring a bag to room temperature, open the bag, take out the necessary devices; store the others in the bag containing the silica gel, let the air out and **seal** by pressing on the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.65 ml

Content: Buffer containing suPAR and 0.05% Bronidox. Liquid, ready to use.

CONTROL + POSITIVE CONTROL

1 x 1.1 ml

Content: Buffer containing suPAR and 0.05% Bronidox. Liquid, ready to use.

OTHER MATERIAL REQUIRED, BUT NOT SUPPLIED:

- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- CHORUS TRIO instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test tubes, etc.
- Micropipettes capable of accurately withdrawing volumes of 50-200 µl.
- Single-use gloves
- 5% sodium hypochlorite solution
- Containers for the collection of potentially infected materials

6. STORAGE METHOD AND REAGENT STABILITY

The reagents must be stored at 2-8°C. In the event of an incorrect storage temperature, the calibration must be repeated and the correctness of the result checked using the positive control (see chapter 9: Test validation).

The expiry date is printed on each component and on the outer label of the package.

The reagents have limited stability after opening and/or preparation:

DEVICES	8 weeks at 2-8°C
CALIBRATOR	8 weeks at 2-8°C
POSITIVE CONTROL	8 weeks at 2-8°C

7. SAMPLE TYPE AND STORAGE

The sample type is plasma and serum (anticoagulants: heparin and EDTA) obtained from blood collected through venous blood sampling and handled as required by standard laboratory procedures.

The consequences of using other body fluids are unknown.

The fresh sample can be kept for 5 days at 2-8°C; for longer storage periods, freeze at ≤ -20°C.

The sample can undergo up to a maximum of 3 thaws. Avoid the use of self-defrosting freezers for sample storage. After thawing, carefully shake the sample before dosing.

Heat inactivation of the sample could give erroneous results.

Sample quality can be seriously affected by microbial contamination which can lead to erroneous results.

8. PROCEDURE

The CHORUS test on the EXTENDED suPAR (REF 81414) cannot be simultaneously performed with the following tests:

CHORUS TG	REF 86076
CHORUS TG	REF 86076\12
CHORUS 25OH VITAMIN D TOTAL	REF 86900
CHORUS PHOSPHOLIPID SCREEN-G	REF 86106
CHORUS PHOSPHOLIPID SCREEN-M	REF 86108
CHORUS PROMONITOR INFLIXIMAB	REF 86700
CHORUS PROMONITOR ANTI-INFLIXIMAB	REF 86701
CHORUS PROMONITOR ADALIMUMAB	REF 86702
CHORUS PROMONITOR ANTI-ADALIMUMAB	REF 86703
CHORUS ADENOVIRUS	REF 85210
CHORUS BORDETELLA PERTUSSIS	REF 85355
CHORUS BORRELIA	REF 85358
CHORUS BRUCELLA	REF 85373
CHORUS CAMPYLOBACTER JEJUNI	REF 85361
CHORUS CHLAMYDIA	REF 85342
CHORUS COXSACKIE VIRUS AT MIX	REF 85313
CHORUS COXSACKIE VIRUS B MIX	REF 85251
CHORUS ECHINOCOCCUS	REF 85328
CHORUS ECHO VIRUS N MIX	REF 85275
CHORUS ECHO VIRUS P MIX	REF 85221
CHORUS INFLUENZA A VIRUS	REF 85222
CHORUS INFLUENZA B VIRUS	REF 85224
CHORUS LEGIONELLA PNEUMOPHILA	REF 85352
CHORUS LEPTOSPIRA MIX	REF 85367
CHORUS LISTERIA MONOCYTOGENES	REF 85344
CHORUS MYCOPLASMA PNEUMONIAE	REF 85340
CHORUS PARAINFLUENZA MIX	REF 85225
CHORUS POLIOVIRUS MIX	REF 85283
CHORUS Q-FEVER	REF 85330
CHORUS REOVIRUS	REF 85234
CHORUS RESPIRATORY SYNCYTIAL VIRUS	REF 85236
CHORUS SHIGELLA FLEXNERI	REF 85370
CHORUS YERSINIA ENTEROCOLITICA	REF 85375

1. Open the bag (side containing the snap closure), take the number of devices needed to perform the tests and store the others by closing the bag after letting the air out.
2. Visually check the status of the device according to the instructions in section 4 Analytical Warnings.
3. Dispense 130 µl of undiluted sample to be analysed into well no. 1 of each device. At each batch change use a device for the calibrator.
4. Place the devices on the CHORUS TRIO instrument. Perform the calibration (if required) and the test as indicated in the User Manual of the instrument.

9. TEST VALIDATION

Use the positive control to verify the correctness of the result obtained, processing it as indicated in the User Manual of the instrument. If the instrument indicates that the value of the control is outside the acceptance limit, the calibration must be performed again. Previous results will be corrected automatically.

If the result of the control continues to be outside the acceptance range, contact Scientific Support.

Phone: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The CHORUS TRIO instrument provides the result in ng/ml calculated on the basis of a batch-dependent graph stored in the instrument.

The test on the sample being examined can be interpreted as follows:

- < 4 ng/ml - Physiological value
- 4-6 ng/ml - Deviation from the physiological value, indication of a likely ongoing inflammatory process
- > 6 ng/ml - Pathological value, strong indication of an ongoing inflammatory process (acute or chronic) and likely deterioration of the clinical condition

11. TEST LIMITATIONS

All values obtained require a careful interpretation, without ever disregarding other indicators relating to the same patient.

The result of the test must be evaluated together with data from the patient's medical history and/or from other diagnostic investigations.

12. CALIBRATION RANGE

Calibration range 0.5-25.0 ng/ml.

For samples > 25.0 ng/ml repeat the test, pre-diluting the sample with the Negative Control/Sample Diluent (PF83607-not supplied with the kit).

13. REFERENCE RANGES

The expected values in the normal population, determined by examining 120 samples from healthy donors, were < 4 ng/ml.

14. ANALYTICAL SPECIFICITY

5 samples were tested to which the following interfering agents were added:

Bilirubin (4.5 - 45 mg/dl)
Triglycerides (250 - 1500 mg/dl)
Haemoglobin (5 - 30 mg/ml)

The presence of high concentrations of haemoglobins (30 mg/ml) may alter the test result.

15. COMPARISON STUDIES

In one trial, 206 samples were analysed with a Diesse kit and with another reference kit on the market.

The trial data are outlined below:

		Reference		
		+	-	Total
Diesse	+	31	3	34
	-	0	172	172
	Total	31	175	206

Percent Positive Agreement (~Diagnostic Sensitivity):

100.0% CI_{95%}: 89.0 - 99.9

Percent Negative Agreement (~Diagnostic Specificity):

98.3% CI_{95%}: 95.1 - 99.4

The degree of agreement between the two methods is excellent with a K value (Cohen's coefficient) of 0.96.

The results obtained are reported below:

Correlation	r	95%CI
Pearson	0.96	0.95-0.97
Spearman	0.90	0.87-0.92

16. ACCURACY AND REPEATABILITY

Sample	Within the session		Between sessions	
	Mean (ng/ml)	CV%	Mean (ng/ml)	CV%
1	1.7	7.1	1.6	5.0
2	5.7	2.6	5.8	4.1
3	7.8	2.2	7.5	6.4
4	20.7	4.1	20.1	3.2

Sample	Between batches		Between instruments	
	Mean (ng/ml)	CV%	Mean (ng/ml)	CV%
1	1.8	7.7	1.8	12.2
2	5.6	5.3	5.7	7.4
3	7.7	6.8	7.8	9.1
4	18.7	4.5	19.9	5.3

17 BIBLIOGRAPHY

1. Hartmann Rasmussen L. et All. Soluble urokinase plasminogen activator receptor (suPAR) as a biomarker of systemic chronic inflammation. *Front Immunol.* (2021)
2. Lafon T. et All. Prognostic performance of endothelial biomarkers to early predict clinical deterioration of patients with suspected bacterial infection and sepsis admitted to the emergency department. *Ann Intensive Care.* (2020)
3. Desmedt S et al. The intriguing role of soluble urokinase receptor in inflammatory diseases. *Crit Rev Clin Lab Sci.* 2017 Mar;54(2):117-133
4. Rasmussen LJH et al. Combining National Early Warning Score with Soluble Urokinase Plasminogen Activator Receptor (suPAR) Improves Risk Prediction in Acute Medical Patients: A Registry-Based Cohort Study. *Crit Care Med.* 2018 Dec;46:1961-8.
5. Rasmussen LJH, et al. Soluble urokinase plasminogen activator receptor (suPAR) in acute care: a strong marker of disease presence and severity, readmission and mortality. A retrospective cohort study. *Emerg Med J.* 2016 Nov;33:769-75.
6. Smith HW, Marshall CJ. Regulation of cell signalling by uPAR. *Nat Rev Mol Cell Biol.* 2010;11:23-36.
7. Enocsson H, Lukic T, Ziegelsch M, Kastbom A. Serum levels of the soluble urokinase plasminogen activator receptor (suPAR) correlates with disease activity in early rheumatoid arthritis and reflects joint damage over time. *Transl Res.* 2021 Jun;232:142-149.

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17 BIBLIOGRAFIE

1. Hartmann Rasmussen L. et All. Soluble urokinase plasminogen activator receptor (suPAR) as a biomarker of systemic chronic inflammation. Front Immunol. (2021)
2. Lafon T. et All. Prognostic performance of endothelial biomarkers to early predict clinical deterioration of patients with suspected bacterial infection and sepsis admitted to the emergency department. Ann Intensive Care. (2020)
3. Desmedt S et al. The intriguing role of soluble urokinase receptor in inflammatory diseases. Crit Rev Clin Lab Sci. 2017 Mar;54(2):117-133
4. Rasmussen LJH et al. Combining National Early Warning Score with Soluble Urokinase Plasminogen Activator Receptor (suPAR) Improves Risk Prediction in Acute Medical Patients: A Registry-Based Cohort Study. Crit Care Med. 2018 Dec;46:1961-8.
5. Rasmussen LJH, et al. Soluble urokinase plasminogen activator receptor (suPAR) in acute care: a strong marker of disease presence and severity, readmission and mortality. A retrospective cohort study. Emerg Med J. 2016 Nov;33:769-75.
6. Smith HW, Marshall CJ. Regulation of cell signalling by uPAR. Nat Rev Mol Cell Biol. 2010;11:23-36.
7. Enocsson H, Lukic T, Ziegler M, Kastbom A. Serum levels of the soluble urokinase plasminogen activator receptor (suPAR) correlates with disease activity in early rheumatoid arthritis and reflects joint damage over time. Transl Res. 2021 Jun;232:142-149.

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17 LITERATUR

1. Hartmann Rasmussen L. et All. Soluble urokinase plasminogen activator receptor (suPAR) as a biomarker of systemic chronic inflammation. Front Immunol. (2021)
2. Lafon T. et All. Prognostic performance of endothelial biomarkers to early predict clinical deterioration of patients with suspected bacterial infection and sepsis admitted to the emergency department. Ann Intensive Care. (2020)
3. Desmedt S et al. The intriguing role of soluble urokinase receptor in inflammatory diseases. Crit Rev Clin Lab Sci. 2017 Mar;54(2):117-133
4. Rasmussen LJH et al. Combining National Early Warning Score with Soluble Urokinase Plasminogen Activator Receptor (suPAR) Improves Risk Prediction in Acute Medical Patients: A Registry-Based Cohort Study. Crit Care Med. 2018 Dec;46:1961-8.
5. Rasmussen LJH, et al. Soluble urokinase plasminogen activator receptor (suPAR) in acute care: a strong marker of disease presence and severity, readmission and mortality. A retrospective cohort study. Emerg Med J. 2016 Nov;33:769-75.
6. Smith HW, Marshall CJ. Regulation of cell signalling by uPAR. Nat Rev Mol Cell Biol. 2010;11:23-36.

7. Enocsson H, Lukic T, Ziegelasch M, Kastbom A. Serum levels of the soluble urokinase plasminogen activator receptor (suPAR) correlates with disease activity in early rheumatoid arthritis and reflects joint damage over time. Transl Res. 2021 Jun;232:142-149.

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- Hartmann Rasmussen L. et All. Soluble urokinase plasminogen activator receptor (suPAR) as a biomarker of systemic chronic inflammation. Front Immunol. (2021)
- Lafon T. et All. Prognostic performance of endothelial biomarkers to early predict clinical deterioration of patients with suspected bacterial infection and sepsis admitted to the emergency department. Ann Intensive Care. (2020)
- Desmedt S et al. The intriguing role of soluble urokinase receptor in inflammatory diseases. Crit Rev Clin Lab Sci. 2017 Mar;54(2):117-133
- Rasmussen LJH et al. Combining National Early Warning Score with Soluble Urokinase Plasminogen Activator Receptor (suPAR) Improves Risk Prediction in Acute Medical Patients: A Registry-Based Cohort Study. Crit Care Med. 2018 Dec;46:1961-8.
- Rasmussen LJH, et al. Soluble urokinase plasminogen activator receptor (suPAR) in acute care: a strong marker of disease presence and severity, readmission and mortality. A retrospective cohort study. Emerg Med J. 2016 Nov;33:769-75.
- Smith HW, Marshall CJ. Regulation of cell signalling by uPAR. Nat Rev Mol Cell Biol. 2010;11:23-36.

- Enocsson H, Lukic T, Ziegelsch M, Kastbom A. Serum levels of the soluble urokinase plasminogen activator receptor (suPAR) correlates with disease activity in early rheumatoid arthritis and reflects joint damage over time. Transl Res. 2021 Jun;232:142-149.

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CHORUS TG	REF. 86076
CHORUS TG	REF. 86076\12
CHORUS 25OH VITAMIN D TOTAL	REF. 86900
CHORUS PHOSPHOLIPID SCREEN-G	REF. 86106
CHORUS PHOSPHOLIPID SCREEN-M	REF. 86108
CHORUS PROMONITOR INFILIXIMAB	REF. 86700
CHORUS PROMONITOR ANTI-INFILIXIMAB	REF. 86701
CHORUS PROMONITOR ADALIMUMAB	REF. 86702
CHORUS PROMONITOR ANTI-ADALIMUMAB	REF. 86703
CHORUS ADENOVIRUS	REF. 85210
CHORUS BORDETELLA PERTUSSIS	REF. 85355
CHORUS BORRELIA	REF. 85358
CHORUS BRUCELLA	REF. 85373
CHORUS CAMPYLOBACTER JEJUNI	REF. 85361
CHORUS CHLAMYDIA	REF. 85342
CHORUS COXSACKIE VIRUS A MIX	REF. 85313
CHORUS COXSACKIE VIRUS B MIX	REF. 85251
CHORUS ECHINOCOCCUS	REF. 85328
CHORUS ECHO VIRUS N MIX	REF. 85275
CHORUS ECHO VIRUS P MIX	REF. 85221
CHORUS INFLUENZA A VIRUS	REF. 85222
CHORUS INFLUENZA B VIRUS	REF. 85224
CHORUS LEGIONELLA PNEUMOPHILA	REF. 85352
CHORUS LEPTOSPIRA MIX	REF. 85367
CHORUS LISTERIA MONOCYTOGENES	REF. 85344
CHORUS MYCOPLASMA PNEUMONIAE	REF. 85340
CHORUS PARAINFLUENZA MIX	REF. 85225
CHORUS POLIOVIRUS MIX	REF. 85283
CHORUS Q-FEVER	REF. 85330
CHORUS REOVIRUS	REF. 85234
CHORUS RESPIRATORY SYNCYTIAL VIRUS	REF. 85236
CHORUS SHIGELLA FLEXNERI	REF. 85370
CHORUS YERSINIA ENTEROCOLITICA	REF. 85375

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- Hartmann Rasmussen L. et All. Soluble urokinase plasminogen activator receptor (suPAR) as a biomarker of systemic chronic inflammation. Front Immunol. (2021)
- Lafon T. et All. Prognostic performance of endothelial biomarkers to early predict clinical deterioration of patients with suspected bacterial infection and sepsis admitted to the emergency department. Ann Intensive Care. (2020)
- Desmedt S et al. The intriguing role of soluble urokinase receptor in inflammatory diseases. Crit Rev Clin Lab Sci. 2017 Mar;54(2):117-133
- Rasmussen LJH et al. Combining National Early Warning Score with Soluble Urokinase Plasminogen Activator Receptor (suPAR) Improves Risk Prediction in Acute Medical Patients: A Registry-Based Cohort Study. Crit Care Med. 2018 Dec;46:1961-8.
- Rasmussen LJH, et al. Soluble urokinase plasminogen activator receptor (suPAR) in acute care: a strong marker of disease presence and severity, readmission and mortality. A retrospective cohort study. Emerg Med J. 2016 Nov;33:769-75.
- Smith HW, Marshall CJ. Regulation of cell signalling by uPAR. Nat Rev Mol Cell Biol. 2010;11:23-36.
- Enocsson H, Lukic T, Ziegelsch M, Kastbom A. Serum levels of the soluble urokinase plasminogen activator receptor (suPAR) correlates with disease activity in early rheumatoid arthritis and reflects joint damage over time. Transl Res. 2021 Jun;232:142-149.

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CHORUS TG	APT. 86076
CHORUS TG	APT. 86076/12
CHORUS 25OH VITAMIN D TOTAL	APT. 86900
CHORUS PHOSPHOLIPID SCREEN-G	APT. 86106
CHORUS PHOSPHOLIPID SCREEN-M	APT. 86108
CHORUS PROMONITOR INFILXIMAB	APT. 86700
CHORUS PROMONITOR ANTI-INFILXIMAB	APT. 86701
CHORUS PROMONITOR ADALIMUMAB	APT. 86702
CHORUS PROMONITOR ANTI-ADALIMUMAB	APT. 86703
CHORUS ADENOVIRUS	APT. 85210
CHORUS BORDETELLA PERTUSSIS	APT. 85355
CHORUS BORRELLIA	APT. 85358
CHORUS BRUCELLA	APT. 85373
CHORUS CAMPYLOBACTER JEJUNI	APT. 85361
CHORUS CHLAMYDIA	APT. 85342
CHORUS COXSACKIE VIRUS A MIX	APT. 85313
CHORUS COXSACKIE VIRUS B MIX	APT. 85251
CHORUS ECHINOCOCCUS	APT. 85328
CHORUS ECHO VIRUS N MIX	APT. 85275
CHORUS ECHO VIRUS P MIX	APT. 85221
CHORUS INFLUENZA A VIRUS	APT. 85222
CHORUS INFLUENZA B VIRUS	APT. 85224
CHORUS LEGIONELLA PNEUMOPHILA	APT. 85352
CHORUS LEPTOSPIRA MIX	APT. 85367
CHORUS LISTERIA MONOCYTOGENES	APT. 85344
CHORUS MYCOPLASMA PNEUMONIAE	APT. 85340
CHORUS PARAINFLUENZA MIX	APT. 85225
CHORUS POLIOVIRUS MIX	APT. 85283
CHORUS Q-FEVER	APT. 85330
CHORUS REOVIRUS	APT. 85234
CHORUS RESPIRATORY SYNCYTIAL VIRUS	APT. 85236
CHORUS SHIGELLA FLEXNERI	APT. 85370
CHORUS YERSINIA ENTEROCOLITICA	APT. 85375

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