CHORUS

**WEST NILE VIRUS
IgM**

REF 81176

CE

DIESSE

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INSTRUCTIONS FOR USE

CHORUS WEST NILE VIRUS IgM

For the qualitative determination of anti-West Nile Virus IgM antibodies

For in vitro diagnostic use only

1. USE

Enzyme immunoassay method for the qualitative determination of anti-West Nile Virus (WNV) IgM class antibodies in human serum with a single-use device applied to Chorus and Chorus TRIO instruments.

2. INTRODUCTION

The West Nile virus belongs to the Flaviviridae family. The virus reservoirs are wild birds and mosquitoes (most frequently type Culex), whose bites are the main means of transmission to humans.

The incubation period from the time of the bite of the infected mosquito varies between 2 and 14 days, but may also be 21 days in those with immune deficiency.

West Nile infections are asymptomatic in 80% of people but may cause fever, headache and nausea; in some cases meningitis and/or encephalitis may occur, which can cause serious brain damage or death. There is no specific therapy for West Nile fever. In most cases, symptoms go away on their own after a few days or may last for a few weeks. In the most severe cases, hospitalization is required, where treatments given include intravenous fluids and assisted breathing.

Diagnosis is mainly performed by serological testing for the determination of WNV-specific IgG and IgM antibodies. Alternatively, it may also be performed by PCR or viral culture on serum and cerebrospinal fluid specimens.

3. METHOD PRINCIPLE

The Chorus West Nile Virus IgM device is ready to use for the determination of anti-WNV IgG antibodies in Chorus/Chorus TRIO instruments.

The test is based on the ELISA (Enzyme-Linked Immunosorbent Assay) principle. The antigen is bound to the solid phase. Specific immunoglobulins bind to the antigen following incubation with diluted human serum. After washing to remove proteins that have not reacted, incubation is carried out with the conjugate consisting of human anti-IgM monoclonal antibodies conjugated with horseradish peroxidase. The unbound conjugate is removed and the substrate for the peroxidase is added. The blue colour that develops is proportional to the concentration of specific antibodies in the serum under examination.

The single-use devices contain all the reagents to perform the test when applied to the Chorus/Chorus TRIO instruments. The results are expressed in Index (sample OD/cut-off OD).

4. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

This kit contains materials of human origin which have been tested with FDA-approved tests and found negative for both HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. Since no diagnostic test can offer a complete guarantee on the absence of infectious agents, any material of human origin must be considered potentially infected. All reagents and samples must be handled in accordance with the safety standards normally adopted in the laboratory.

Disposal of residues: the serum samples, calibrators and strips used must be treated as infected residues, therefore disposed of in accordance with the provisions of current laws.

Personal safety warnings

1. Do not pipette with your mouth.
 2. Use single-use gloves and eye protection when handling samples.
 3. Wash your hands thoroughly after inserting the devices into the Chorus/Chorus TRIO instrument.
 4. Refer to the Safety Data Sheet (available on request) for the safety characteristics of the reagents contained in the kit.
 5. Neutralised acids and other liquid waste must be disinfected by adding sodium hypochlorite in a volume sufficient to obtain a final concentration of at least 1%. Exposure to 1% sodium hypochlorite for 30 minutes should be sufficient to ensure effective disinfection.
 6. Any spillage of potentially infected materials must be removed immediately with absorbent paper and the contaminated area must be decontaminated, for example with 1% sodium hypochlorite, before continuing work. If an acid is present, sodium hypochlorite must not be used before the area has been dried.
- All materials used to decontaminate any accidental spillage, including gloves, must be discarded as potentially infected waste.
- Do not autoclave materials containing sodium hypochlorite.

Analytical warnings

Before use, bring the devices to be used to room temperature (18-30°C) and use within 60 minutes.

1. **Discard the devices with blue coloured substrate (well 4).**
 2. When adding the sample to the well, check that it is perfectly distributed at the bottom.
 3. Check the actual presence of the reagents in the device and the integrity of the device itself. Do not use devices that lack any reagents and/or that have any foreign bodies in the reaction well when visually inspected.
 4. The devices must be used in conjunction with the Chorus/Chorus TRIO instrument, strictly following the Instructions for Use and the User Manual of the instrument.
- Use of the kit is only possible with an updated version of the software. Make sure that the software installed in the instrument matches or has a release (Rel.) date later than the one shown in the table published on the Diesse website**

(<http://www.diesse.it/it/Support/Download/strumento:39/>)


5. Check that the Chorus/Chorus TRIO instrument is set correctly (see User Manual).
6. Do not alter the barcode placed on the handle of the device in order to allow it to be read correctly by the instrument.
7. Avoid the use of self-defrosting freezers for sample storage.
8. Defective barcodes can be manually entered into the instrument (see User Manual).
9. Do not expose the devices to strong lighting or hypochlorite vapours during storage and use.
10. The use of heavily haemolysed, lipemic, icteric, non-fully coagulated serum samples or of samples with microbial contamination may be a source of errors.
11. Do not use the device after its expiry date
12. **Check that the instrument is connected to the Washing Buffer (Ref. 83606)**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests

DD DEVICES packs of 6 devices each

Description:

Position 8: Space available for barcode label

Position 7: Empty

Position 6: MICROPLATE WELL

Sensitized with WNV NS1 recombinant protein antigen **Position**

5: MICROPLATE WELL

Not sensitized.

Position 4: TMB SUBSTRATE

Content: 0.26 mg/mL tetramethylbenzidine and 0.01% H₂O₂ stabilized in 0.05 mol/L citrate buffer (pH 3.8).

Position 3: DILUENT FOR SAMPLES

Content: Protein solution containing 0.05% phenol, 0.02% Bronidox and an indicator to detect the presence of serum.

Position 2: CONJUGATE

Content: Peroxidase-labelled human anti-IgM monoclonal antibodies in phosphate buffer solution containing 0.05% phenol and 0.02% Bronidox.

Position 1: EMPTY WELL

Where the sample is transferred.

Usage: bring a bag to room temperature, open the bag, take out the necessary devices; store the others in the bag containing the silica gel, let the air out and **seal** by pressing on the closure. Store at 2-8°C.

CALIBRATOR 1 x 0.175 ml

Content: Buffer containing a preservative and specific antibodies capable of binding the antigen present on the micro-plate. Liquid, ready to use.

CONTROL + POSITIVE CONTROL 1 x 0.425 ml

Content: Buffer containing a preservative and specific antibodies capable of binding the antigen present on the micro-plate. Liquid, ready to use.

OTHER MATERIAL REQUIRED, BUT NOT SUPPLIED:

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITISING SOLUTION **REF** 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Chorus/Chorus TRIO instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test tubes, etc.
- Micropipettes capable of accurately withdrawing volumes of 50-200µl.
- Single-use gloves
- 5% sodium hypochlorite solution
- Containers for the collection of potentially infected materials

6. STORAGE METHOD AND REAGENT STABILITY

The reagents must be stored at 2-8°C. In the event of an incorrect storage temperature, the calibration must be repeated and the correctness of the result checked using the control serum (see section 9: Test validation).

The expiry date is printed on each component and on the outer label of the package.

The reagents have limited stability after opening and/or preparation:

DEVICES	8 weeks at 2-8°C
CALIBRATOR	8 weeks at 2-8°C
POSITIVE CONTROL	8 weeks at 2-8°C

7. SAMPLE TYPE AND STORAGE

The sample type is serum obtained from blood collected by normal venous blood sampling and handled as required by standard laboratory procedures.

The consequences of using other body fluids are unknown.

Fresh serum can be kept for 5 days at 2-8°C; for longer storage periods, freeze at -20°C.

The sample can undergo up to a maximum of 3 thaws.

Avoid the use of self-defrosting freezers for sample storage. After thawing, carefully shake the sample before dosing.

Heat inactivation may give erroneous results. Sample quality can be seriously affected by microbial contamination which can lead to erroneous results.

8. PROCEDURE

1. Open the bag (side containing the snap closure), take the number of devices needed to perform the tests and store the others by closing the bag after letting the air out.
2. Visually check the status of the device according to the instructions in section 4 Analytical Warnings.
3. Dispense into well no. 1 of each device:
 - SAMPLE: 50 µl/device
 - CALIBRATOR: 50 µl/device
 - POSITIVE CONTROL: 50 µl/device

At each batch change use a device for the calibrator.

4. Introduce the devices on the Chorus/Chorus TRIO instrument. Perform the calibration (if required) and the test as indicated in the User Manual of the instrument.

9. TEST VALIDATION

Use the positive control serum to verify the correctness of the result obtained, processing it as indicated in the User Manual of

the instrument. If the instrument indicates that the control serum has a value outside the acceptability limit, the calibration must be performed again. Previous results will be corrected automatically.

If the control serum result continues to be outside the acceptability range, contact Scientific Support.

Phone: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE TEST

The Chorus/Chorus TRIO instrument provides the result in Index (sample OD/cut-off OD).

The test on the serum under examination can be interpreted as follows:

POSITIVE: when the result is >1.1

NEGATIVE: when the result is <0.9

DOUBTFUL/UNCERTAIN: when the result is between 0.9 and 1.1

In case of a doubtful/uncertain result, repeat the test. If the result remains doubtful/uncertain, repeat the sampling.

11. TEST LIMITATIONS

All values obtained require a careful interpretation, without ever disregarding other indicators relating to the same patient.

The test, in fact, cannot be used on its own for a clinical diagnosis and the result obtained must always be evaluated together with data from the patient's medical history and/or from other diagnostic investigations.

12. ANALYTICAL SPECIFICITY

3 samples were tested (1 Negative, 1 at Cut-Off and 1 Positive) to which the following interfering agents were added:

Rheumatoid factor (110 IU/ml - 220 IU/ml)
Bilirubin (18 mg/dl - 45 mg/dl)
Lipoprotein (750 mg/dl - 1500 mg/dl)
Haemoglobin (15 mg/ml - 30 mg/ml)

The presence of the above interfering substances in the test serum does not alter the test result.

13. CROSS-REACTIVITY

22 samples, positive for Dengue (6), Zika (8), Rheumatoid Factor (6), Rheumatoid Factor/Cyclic Citrullinated Peptide (2) were tested.

Cross-reactions with Dengue were detected. Cross-reactions with Zika cannot be excluded.

14. COMPARISON STUDIES

In one trial, 133 samples were analysed with a Diesse kit and with another commercially available kit.

The trial data are outlined below:

		Reference		Total
		+	-	
Diesse	+	21	0	21
	-	1	111	112
	Total	22	111	133

Percent Positive Agreement (~Diagnostic Sensitivity):

95.5% CI_{95%}: 78.1-99.0

Percent Negative Agreement (~Diagnostic Specificity):

100.0% CI_{95%}: 96.6-100.0

The degree of agreement between the two methods is excellent with a K value (Cohen's coefficient) of 0.96.

15. ACCURACY AND REPEATABILITY

Sample	Within the session		Between sessions	
	Mean (Index)	CV%	Mean (Index)	CV%
1	2.2	6.4	2.0	9.5
2	1.8	5.6	1.5	8.0
3	0.9	10.0	0.8	12.5
4	0.1	11.4	0.1	10.0
5	0.3	13.3	0.2	10.0

Sample	Between batches		Between instruments	
	Mean (Index)	CV%	Mean (Index)	CV%
1	2.5	12.4	2.1	10.5
2	1.5	6.0	1.5	10.0
3	0.9	6.7	0.8	13.8
4	0.1	-	0.1	-
5	0.2	5.0	0.2	15.0

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