CHORUS

ANTI-CCP

REF 86094

REF 86094/12

CE
0123

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For In Vitro Diagnostic Use Only

1. INTENDED PURPOSE

The product Chorus ANTI-CCP (REF 86094-86094/12) is an immunoassay kit for automated semi-quantitative detection of IgG class antibodies against cyclic citrullinated peptide (anti-CCP).

Since the anti-CCP antibodies are highly specific for Rheumatoid Arthritis (RA) and their detection is one of the ACR (American College of Rheumatology) criteria for classification of RA, the kit is used as an aid to the relative diagnosis.

The test, performed in human serum using a disposable device applied on the Chorus and Chorus TRIO instruments, must be used by professional laboratory users only.

2. INTRODUCTION

Rheumatoid arthritis (RA) is a chronic inflammatory autoimmune disease characterized by joint swelling, joint tenderness and destruction of synovial joints, leading to functional limitations, severe disability with poor quality of life and premature mortality. High Rheumatoid Factor (RF) titers correlate with the risk of developing RA. However, RF positivity alone is not sufficient for diagnosis. Anti-citrullinated protein antibodies (ACPA, Anti-CCP), found in 63-76% of patients with confirmed RA, have a superb diagnostic specificity of up to 90%.

The standard and accepted means of defining RA is by use of classification criteria. The international classification criteria, set by American College of Rheumatology (ACR) and European League Against Rheumatism (EULAR) foresee 4 categories, which include joint involvement, serology, acute-phase reactants and duration of symptoms. Among the serology, one of the criteria to be determined is the level of antibodies against cyclic citrullinated peptide (Anti-CCP).

Since the product Chorus Anti-CCP is intended to determine the presence in patient serum of specific anti-CCP IgG antibodies, the device is used as an aid in the diagnosis of RA.

3. PRINCIPLE OF THE METHOD

The CHORUS anti-CCP device is ready to use for the detection of antibodies against CCP in the CHORUS instruments.

The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay).

The antigen is bound to the solid phase.

The specific immunoglobulins are bound to the antigen through incubation with diluted human serum.

After washing to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of anti-human immunoglobulins conjugated to horse radish peroxidase. The unbound conjugate is eliminated and the TMB substrate is added.

The colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test in the Chorus instruments.

The results are expressed in Arbitrary Unit (AU/ml).

4. WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips, once used, must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on DIESSE website: www.diesse.it) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the bags with the devices to room temperature (18-30°C) for at least 30 minutes before use.

Use the device within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.
5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Chorus Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use hemolyzed, lipemic, icteric samples having a interfering concentration higher than that tested (according to the indications reported in chapter "Analytical specificity")
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity REF 86004.**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 86094).


The kit is sufficient for 12 tests. (REF 86094/12).

DD DEVICES

6 packages each containing 6 devices (REF 86094).

2 packages each containing 6 devices (REF 86094/12).

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL coated with highly purified synthetic peptide containing citrulline residues (peptide concentration of 0.002 mg/ml)

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine in citrate buffer

Position 3: SAMPLE DILUENT

Contents: saline proteic solution with Proclin (0.1%)

Position 2: CONJUGATE

Contents: anti-human IgG monoclonal antibodies labeled with horse radish peroxidase (maximum concentration of 0.25 µg/ml), in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

in which undiluted serum must be added

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CALIBRATOR

1 x 0.175 ml

Contents: Diluted human serum (maximum concentration of 30%), containing anti-CCP IgG antibodies, Tween-20 0.2%, Proclin 300 0.1% and Methylorange: 0.0008%.

Liquid, ready for use.

CONTROL + POSITIVE CONTROL

1 x 0.425 ml

Contents: Diluted human serum (maximum concentration of 30%), containing anti-CCP IgG antibodies, Tween-20 0.2%, Proclin 300 0.1% and Methylorange 0.0008%.

Liquid, ready for use.

The device is a semi-quantitative test: Calibrator and Positive control cannot provide traceability to International Standard Units. Confidence in measurements is established with traceability to measurement standards as follows.

The relative exact range concentration is lot-dependent and is assigned during the releasing Quality control phase using a series of Working Calibrators.

The Working Calibrators are prepared and characterized, checking the consensus with a reference product.

MATERIALS REQUIRED BUT NOT PROVIDED

- Chorus Instruments
- WASHING BUFFER AUTOIMMUNITY (REF 86004)
- CLEANING SOLUTION 2000 (REF 83609)
- SANITIZING SOLUTION (REF 83604 - 83608)
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT (REF 83607)
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum from whole blood collected by venipuncture and handled following standard laboratory techniques.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored frozen for long periods at -20°C (for at least 61 months) and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers.

Heat-inactivation can rise to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. At each change of batch, perform the calibration.
4. Dispense the following volumes in well no. 1 of each device:

SAMPLE	50 µl/device
CALIBRATOR	50 µl/device
POSITIVE CONTROL	50 µl/device

5. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the control serum as reported in the Instrument Operating Manual to check the validity of the obtained results.

The acceptability range is reported on the enclosed Quality control certificate.

If the result of the positive control is out of the acceptability range, the calibration must be repeated.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in AU/ml calculated on the basis of a lot-dependent graph stored in the instrument.

The results can be interpreted as follows:

POSITIVE: when the result is > 18 AU/ml

NEGATIVE: when the result is < 12 AU/ml.

DOUBTFUL/EQUIVOCAL: for all values between 12 and 18 AU/ml

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new serum sample.

11. LIMITATIONS

The product must be used by professional laboratory users only. The test is not applicable to sample types different from serum.

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. CALIBRATION RANGE

Calibration range of 3.0-300 AU/ml.

For samples > 300 AU/ml retest the diluted sample in the Negative Control/Sample Diluent (PF83607-not supplied with the kit).

13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, were between 3.0 and 7.5 AU/ml.

Among the affected population the expected values, which have been determined by examining 72 sera from patients with a previous diagnosis of RA or previously characterized as positive, were between 26.4 and 300.0 AU/ml.

14. ANALYTICAL SPECIFICITY

5 samples were spiked with the following potentially interfering factors and then tested:

Rheumatoid factor (44-220 IU/ml)

Bilirubin (4.5-45 mg/dl)

Triglycerides (10-250 mg/dl)

Hemoglobin (5-30 mg/ml)

The presence in the serum sample of interfering substances described above do not affect the test result.

15. CROSS REACTION

232 samples from the sera collection were tested, their clinical characteristics are listed in the following table:

Samples Group	n° TOT	Neg < 12AU/ml	D/E 12-18 AU/ml	Pos > 18AU/ml
Healthy Donors	90	88	0	2
Routine samples	58	58	0	0
ANA	14	13	0	1
Thyroiditis	20	20	0	0
Rubella IgM Positive	5	5	0	0
Legionella IgM Positive	4	4	0	0

Chlamydia Positive	5	5	0	0
Salmonella Positive	3	3	0	0
Toxoplasma g. IgM Positive	6	6	0	0
Mycoplasma p. IgM Positive	9	9	0	0
Anti-streptolysin Positive	3	3	0	0
EBV IgM Positive	4	4	0	0
Parvovirus IgM Positive	5	5	0	0
LUE IgM Reactive	6	5	1	0

The clinical specificity is:

Healthy Donors 88/90 = 97.8% (CI_{95%}: 92.3-99.4%)

Routine Samples 58/58 = 100% (CI_{95%}: 93.8-100%)

Autoimmune Disease 33/34 = 97.1% (CI_{95%}: 85.1-99.5%)

Infectious Disease 50/50 = 100% (CI_{95%}: 92.7-100%)

16. METHOD COMPARISON

174 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Tot.
Diesse	+	72	2	74
	-	9	91	100
	Tot.	81	93	174

Percent Positive Agreement (~Diagnostic Sensitivity): 88.9%
CI_{95%}: 80.2 – 94.0

Percent Negative Agreement: (~Diagnostic Specificity): 97.8%
CI_{95%}: 92.5 – 99.4

Positive Predicted Value: 97.3% CI_{95%}: 94.8-99.7

Negative Predicted Value: 91.0% CI_{95%}: 86.7-96.6

The agreement between the two methods is high with a Cohen's Kappa of 0.88.

17. PRECISION AND REPEATABILITY

Sample	Within run		Between run	
	Mean (AU/ml)	CV%	Mean (AU/ml)	CV%
1	64.9	2.8	65.1	4.1
2	26.3	2.5	26.5	13.4
3	27.0	3.6	26.7	12.4
4	15.9	7.5	15.6	13.3
5	16.8	5.8	15.8	11.5
6	8.1	5.6	7.5	12.8

Sample	Between lots		Between Instruments	
	Mean (AU/ml)	CV%	Mean (AU/ml)	CV%
1	63.3	3.3	63.3	3.5
2	25.2	7.6	25.2	7.1
3	25.5	7.7	25.5	7.4
4	15.7	7.1	15.7	8.2
5	15.8	6.3	15.9	6.4
6	7.2	6.3	7.2	7.1

18. REFERENCES

1. Aletaha D. et All., 2010 Rheumatoid Arthritis Classification Criteria. An American College of Rheumatology/European League Against Rheumatism Collaborative Initiative. ARTHRITIS & RHEUMATISM Vol. 62, No. 9 (2010): 2569–2581
2. Atzeni F. et All., Biomarkers in Rheumatoid Arthritis. Isr Med Assoc J. (2017); 19(8): 512-516
3. Derksen V.F.A.M. et All., The role of autoantibodies in the pathophysiology of rheumatoid arthritis. Semin Immunopathol (2017) 39:437–446
4. Wu C.Y. et All., Anti-Citrullinated Protein Antibodies in Patients with Rheumatoid Arthritis: Biological Effects and Mechanisms of Immunopathogenesis. Int. J. Mol. Sci. (2020) 21, 4015
5. Tenstad H.B. et All., Use and utility of serologic tests for rheumatoid arthritis in primary care. Dan Med J 67/2 (2020).

19. INCIDENT REPORTING

If any serious incident in relation to this device has occurred in the European Union market territory, please report without delay to the manufacturer and competent authority of your Member State.

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