CHORUS

TG

REF 86076

REF 86076/12

CE
0123

DIESSE

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17. BIBLIOGRAFIA

1. Esfandiari N.H., et All. Biochemical Testing in Thyroid Disorders. Endocrinol Metab Clin North Am. (2017) September; 46(3): 631-648
2. Buddhike Sri Harsha Indrasena. Use of thyroglobulin as a tumour marker. World J Biol Chem (2017) February 26; 8(1): 81-85
3. Prpic M., et All. Thyroglobulin as a tumor marker in differentiated thyroid cancer – clinical considerations. Acta Clin Croat (2018); 57:518-527
4. Giovannella L., et All. Thyroglobulin measurement using highly sensitive assays in patients with differentiated thyroid cancer: a clinical position paper. European Journal of Endocrinology (2014) 171, R33-R46
5. <https://emedicine.medscape.com/article/2089532-overview#a1>
6. NCCLS CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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For *In Vitro* Diagnostic Use Only

1. INTENDED PURPOSE

The product Chorus TG is an immunoassay kit for automated quantitative determination of Thyroglobulin (TG), performed in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

High TG level is reported in various Thyroid diseases like hyperthyroidism, non-toxic goitre, thyroiditis and differentiated Thyroid carcinoma. The main clinical application for the TG determination, however, is represented by the post-operative monitoring of differentiated thyroid carcinoma. For this reason, the kit is intended to be used as an aid in the diagnosis and monitoring.

It must be used by professional laboratory users only.

2. INTRODUCTION

Thyroglobulin (TG) is a 660 kDa dimeric glycoprotein, composed of two identical 330 kDa monomers. It is synthesized, with the stimulation of Thyroid-Stimulating Hormone (TSH), by thyrocytes and released into the lumen of thyroid follicles. It plays a crucial role in the synthesis of the peripheral thyroid hormones Triiodothyronine (T3) and Thyroxine (T4).

Since the TG serum levels are proportional to the volume of thyroid tissue, it can be elevated in any thyroid pathology. High level of TG are found in patients with goiter, hyperthyroidism, Graves' disease and inflammatory or physical injury to the thyroid gland (thyroiditis).

TG is, however, primarily used as a tumor marker to evaluate the effectiveness of treatment for Differentiated Thyroid Cancer (DTC) and to monitor for residual or recurrent disease. DTC is the most common endocrine malignancy and is characterized by low morbidity and mortality. Initial treatment consists of total thyroidectomy.

Since TG is only produced by thyroid follicular cells, it is expected to be undetectable in patients who underwent total or near-total thyroidectomy. For this reason, TG measurement is the reference method for clinical management of patients previously treated for DTC.

Since, the product Chorus TG is intended to determinate the concentration of Thyroglobulin (TG) in patient serum, it is intended to be used as an aid in the diagnosis and monitoring of DTC.

3. PRINCIPLE OF THE METHOD

The Chorus TG devices are ready to use for the assay of Thyroglobulin (TG), in automation in Chorus and Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay).

Monoclonal antibodies highly specific for TG are bound to the solid phase. The TG is bound to the well through incubation with the human serum examined.

After washings to eliminate the proteins which have not reacted, incubation is performed with the enzyme tracer conjugate, composed of anti-TG monoclonal antibodies conjugated to horse radish peroxidase.

The unbound tracer conjugate is eliminated, and the peroxidase substrate added.

The colour which develops is proportional to the concentration of TG present in the serum sample.

The disposable devices contain all the reagents to perform the test in the Chorus/Chorus TRIO instruments.

The results are expressed in ng/ml, calculated in reference to BCR® certified reference material, BCR-457.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.

Waste disposal: serum samples, calibrators and strips, once used, must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on DIESSE website: www.diesse.it) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the bags with the devices to room temperature (18-30°C) for at least 30 minutes before use.

Use the device within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged. Do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use and the Instrument Operating Manual must be carefully followed and the Instrument Operating Manual must be consulted.
5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use hemolyzed, lipemic, icteric samples having a interfering concentration higher than that tested (according to the indications reported in chapter "Analytical specificity")
11. Do not use the device after the expiry date.
12. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 86076).


The kit is sufficient for 12 tests (REF 86076/12).

DD DEVICES

6 packages each containing 6 devices (REF 86076).

2 packages each containing 6 devices (REF 86076/12).

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with anti-Tg monoclonal antibodies (monoclonal antibodies concentration: 20 µg/ml)

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine in citrate buffer

Position 3: SAMPLE DILUENT

Contents: Proteic solution in phosphate buffer, with sodium azide.

Position 2: CONJUGATE

Contents: anti-TG monoclonal antibodies labeled with peroxidase (maximum concentration: 1 µg/ml).

Position 1: EMPTY WELL in which the patient sample has to be added.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in

the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.500 ml

Contents: Purified human thyroglobulin (maximum concentration: 0.5 µg/ml) in proteinic stabilizing buffer containing Proclin 300 0.2% and Ethylene Glycol 25%. Liquid, ready for use.

CONTROL + POSITIVE CONTROL

1 x 0.500 ml

Contents: Purified human thyroglobulin (maximum concentration: 0.5 µg/ml) in proteinic stabilizing buffer containing Proclin 300 0.2% and Ethylene Glycol 25%. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY (REF 86004)
- CLEANING SOLUTION 2000 (REF 83609)
- SANITIZING SOLUTION (REF 83604 - 83608)
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT (REF 83607)
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum from whole blood collected by venipuncture and handled following standard laboratory techniques.

According to the guideline NCCLS CLSI H18-A3 Serum specimens to be tested should be clotted before centrifugation; spontaneous and complete clotting normally occurs within 30 to 60 minutes at 22°C-25°C.

It is recommended that serum be physically separated, by centrifugation, from contact with cells as soon as possible with a maximum time limit of 2 hours from the time of collection.

The fresh serum may be frozen for long periods at -20°C±3°C (for at least 37 months), and can be thawed a maximum of 3 times.

Possible consequences, in case of use of other biological liquids, are not known.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

Heat-inactivation can rise to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. At each change of batch, perform the calibration.
4. Dispense the following volumes in well no. 1 of each device:

SAMPLE	140 µl/device
CALIBRATOR	140 µl/device
POSITIVE CONTROL	140 µl/device

5. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the instrument Operating Manual.

9. TEST VALIDATION

Use the control serum as reported in the Instrument Operating Manual to check the validity of the obtained results. The acceptability range is reported on the enclosed Quality control certificate.

If the result of the positive control is out of the acceptability range, the calibration must be repeated.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in ng/ml (BCR® certified reference material, BCR-457), calculated on the basis of a lot-dependent graph stored in the instrument.

The test on the examined serum can be interpreted as follows:

Patient	ng/ml
Euthyroidic patients	50 (expected average superior limit of normality)
Patients after total thyroidectomy	<1

Attention: do not consider the qualitative report N; D/E; P on the instrument.

11. LIMITATIONS

The product must be used by professional laboratory users only. The test is not applicable to sample types different from serum. All the values obtained require a careful interpretation that must consider other indicators relative to the patient. The test, indeed, cannot be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

The presence of anti-TG antibodies can influence the test results; for this reason the parallel assay of such auto-antibodies is recommended.

12. CALIBRATION RANGE

Calibration range 0.4-100.0 ng/ml.

For samples > 100.0 ng/ml retest the diluted sample in the Negative Control/Sample Diluent (PF83607-not supplied with the kit).

Functional sensitivity: 0.4 ng/ml.

13. REFERENCE RANGE

The reference range in normal population is given according to the literature as 0.5-53.0 ng/ml.

14. ANALITICAL SPECIFICITY

5 samples (3 Negative and 2 Positive) were spiked with the following potentially interfering factors and then tested:

Rheumatoid Factor (1.25 – 20 IU/ml)
Bilirubin (5 mg/dl – 100 mg/dl)
Triglycerides (50 mg/dl – 3000 mg/dl)
Hemoglobin (0.6 mg/dl – 10 mg/dl)

The presence in the serum sample of the interfering substances described above does not affect the test result.

15. METHOD COMPARISON

Cut-off: 1 ng/ml

In an experimentation 67 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Tot.
Diesse	+	28	2	30
	-	1	36	37
	Tot.	29	38	67

Percent Positive Agreement (~Diagnostic Sensitivity):

96.6% CI_{95%}: 82.8-99.2

Percent Negative Agreement (~Diagnostic Specificity):

94.7% CI_{95%}: 82.7-98.4

Cut-off: 50 ng/ml

In an experimentation 151 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table:

		Reference		
		+	-	Tot.
Diesse	+	63	5	68
	-	2	81	83
	Tot.	65	86	151

Percent Positive Agreement (~Diagnostic Sensitivity):

96.9% CI_{95%}: 89.4-99.1

Percent Negative Agreement (~Diagnostic Specificity):

94.2% CI_{95%}: 87.1-97.4

16. PRECISION AND REPEATABILITY**Intra assay precision**

	ng/ml			Mean	SD	CV%
Tg1	533.3	500.0	516.7	516.7	16.65	3.2
Tg2	60.0	61.7	58.3	60.0	1.70	2.8
Tg3	2.3	2.3	2.3	2.3	0.00	-

Inter assay precision

	ng/ml			Mean	SD	CV%
	Day1	Day 2	Day 3			
Tg1	528.1	448.9	500.8	482.7	31.53	6.5
	503.0	449.8	463.9			
	517.0	445.6	487.2			
Tg2	60.1	57.6	57.5	58.5	1.62	2.8
	61.9	59.0	57.3			
	58.2	56.9	57.7			
Tg3	2.0	1.7	1.8	1.8	0.14	7.8
	2.0	1.7	1.6			
	1.9	1.9	1.8			

Precision between instruments

ng/ml							Mean	SD	CV %
1471			1717						
Day	Day	Day	Day	Day	Day				
Tg 1	528.1	448.9	500.8	454.8	464.8	515.2	484.3	27.09	5.6
	503.0	449.8	463.9	474.4	484.6	517.2			
	517.0	445.6	487.2	469.4	479.5	513.8			
Tg 2	60.1	57.6	57.5	55.7	56.8	68.6	59.6	4.19	7.0
	61.9	59.0	57.3	56.5	57.7	69.1			
	58.2	56.9	57.7	57.3	58.4	66.8			
	2.0	1.7	1.8	1.5	1.5	2.1	1.8	0.20	11.1

Tg 3	2.0	1.7	1.6	1.5	1.5	2.1			
	1.9	1.9	1.8	1.7	1.7	1.7			

17. REFERENCES

1. Esfandiari N.H., et All. Biochemical Testing in Thyroid Disorders. Endocrinol Metab Clin North Am. (2017) September; 46(3): 631–648
2. Buddhike Sri Harsha Indrasena. Use of thyroglobulin as a tumour marker. World J Biol Chem (2017) February 26; 8(1): 81-85
3. Prpic M., et All. Thyroglobulin as a tumor marker in differentiated thyroid cancer – clinical considerations. Acta Clin Croat (2018); 57:518-527
4. Giovannella L., et All. Thyroglobulin measurement using highly sensitive assays in patients with differentiated thyroid cancer: a clinical position paper. European Journal of Endocrinology (2014) 171, R33–R46
5. <https://emedicine.medscape.com/article/2089532-overview#a1>
6. NCCLS CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

18. INCIDENT REPORTING

If any serious incident in relation to this device has occurred in the European Union market territory, please report without delay to the manufacturer and competent authority of your Member State.

19. SUMMARY OF SAFETY AND PERFORMANCE

This document, available on Eudamed database, is part of Technical Documentation.

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- Esfandiari N.H., et All. Biochemical Testing in Thyroid Disorders. Endocrinol Metab Clin North Am. (2017) September; 46(3): 631–648
- Buddhike Sri Harsha Indrasena. Use of thyroglobulin as a tumour marker. World J Biol Chem (2017) February 26; 8(1): 81-85
- Prpic M., et All. Thyroglobulin as a tumor marker in differentiated thyroid cancer – clinical considerations. Acta Clin Croat (2018); 57:518-527
- Giovannella L., et All. Thyroglobulin measurement using highly sensitive assays in patients with differentiated thyroid cancer: a clinical position paper. European Journal of Endocrinology (2014) 171, R33–R46
- <https://emedicine.medscape.com/article/2089532-overview#a1>
- NCCLS CLSI H18-A3 Procedures for the Handling and Processing of Blood Specimens; Approved Guideline — Third Edition. Vol. 24 No. 38

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