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INSTRUCTIONS FOR USE

CHORUS Promonitor ANTI-USTEKINUMAB

For the semiquantitative detection of anti-Ustekinumab antibodies in human serum

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

CHORUS Promonitor ANTI-USTEKINUMAB is an immunoassay kit for automated semiquantitative detection of anti-Ustekinumab antibodies in human serum using a disposable device applied on the CHORUS TRIO instrument.

The product, used in conjunction with other clinical and laboratory findings, is useful as an aid in the management of patients treated with Ustekinumab (UTK), like patients with inflammatory bowel disease (IBD), psoriasis and rheumatic diseases.

It must be used by professional laboratory users only.

2. INTRODUCTION

Ustekinumab (UTK) is a fully human IgG1 kappa monoclonal antibody that binds with high affinity to the p40 subunit of human interleukin (IL) 12 and IL23 and has been approved for the treatment of patients with moderate to severe Crohn's disease (CD), moderate to severe ulcerative colitis (UC), psoriatic arthritis and plaque psoriasis. The drug prevents the bioactivity of IL12 and IL23 by binding and neutralizing the shared p40 subunit, preventing interaction with the cell surface receptor protein IL12Rbeta1. Clinical studies show that 2.9% to 12.4% of patients develop anti-ustekinumab antibodies (ATU) over time. A gradual decrease in efficacy over time following an initial response to biologics is common. ATU may increase drug clearance in treated patients or neutralize drug effect, thus potentially contributing to loss of response. Therefore, it is important to monitor trough UTK concentrations in serum to correlate drug levels with loss of response to therapy and also the occurrence of adverse events such as serum sickness and hypersensitivity reactions. Currently, Ustekinumab quantification is performed in conjunction with the immunogenicity assessment for ATU.

3. PRINCIPLE OF THE METHOD

The CHORUS Promonitor ANTI-USTEKINUMAB device is ready to use for the detection of anti-Ustekinumab antibodies in human serum, in the CHORUS TRIO instrument.

The assay is based on a bridging ELISA method.

The bridging antibody ELISA method is widely used for anti-drug antibody screening and is well accepted as a standard for immunogenicity screening.

The microwell strips are pre-coated with Ustekinumab (UTK). Calibrator, positive control and samples are added to separate device.

Anti-UTK antibodies are able to bind to preimmobilized UTK following incubation with the sample. After washing to remove the unbound sample, an HRP-labelled UTK (conjugate) is added. A second incubation step allows the conjugate to bind to the anti-UTK antibodies, previously attached to the microwells. After washing to remove the unbounded conjugate, the peroxidase substrate is added. The enzymatic reaction is stopped by adding Stop Solution, which makes the solution turns to yellow. The color which develops is proportional to the amount of anti-UTK antibodies in the patient sample. The disposable devices contain all the reagents to perform the test when applied on the CHORUS TRIO instrument. The results are expressed in Arbitrary Units (AU/ml).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the CHORUS TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30-minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the bags with the devices to room temperature (18-30°C) for at least 30 minutes before use.

Use the devices within 60 min.


1. **Discard devices which show the substrate (well 3) blue colored.**
2. If the sample is added into the well manually, verify that it is perfectly distributed.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the CHORUS TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.
The use of the kit is only possible with CHORUS TRIO with an updated version of software. Make sure that the software installed in the instrument corresponds or has a Release (Rel.) subsequent to the one reported in the table published on Diesse website (<https://www.diesse.it/en/downloads/downloads/strumento:39/>)
5. Check that the CHORUS TRIO instrument is set up correctly (see CHORUS TRIO Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see CHORUS TRIO Operating Manual).
9. When CHORUS Autosampler (REF 81210) is used, refer to the related Operating Manual
10. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
11. The use of strongly hemolyzed, lipemic, icteric samples, of serum not completely coagulated or of samples presenting microbial contamination may all constitute sources of error.
12. Do not use the device after the expiry date.
13. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity (Ref. 86004).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests.

DD DEVICES 6 packages each containing 6 devices

Description:

Position 9: Space for application of bar code label

Position 8: EMPTY

Position 7: MICROPLATE WELL

Coated with UTK

Position 6: MICROPLATE WELL

Coated with UTK

Position 5: STOP SOLUTION

Contents: 0.3 M sulfuric acid solution

Position 4: CONJUGATE

Contents: HRP-labelled UTK in buffer containing preservative

Position 3: TMB SUBSTRATE

Contents: Tetramethylbenzidine stabilized in citrate buffer

Position 2: DILUTION BUFFER

Contents: Buffer solution containing preservative

Position 1: EMPTY WELL

in which the sample is transferred.

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR

1 x 0.650 ml

Contents: Human serum containing anti-UTK antibodies and preservative.

Liquid, ready for use.

CONTROL + POSITIVE CONTROL

1 x 1.100 ml

Contents: Human serum containing anti-UTK antibodies and preservative.

Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- Washing Buffer AUTOIMMUNITY **REF** 86004
- Cleaning Solution 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- CHORUS TRIO instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice.

Possible consequences, in case of use of other biological liquids, are not known.

The fresh serum may be stored for 5 days at 2/8°C, or frozen for longer periods at -20°C and can be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be mixed carefully before use.

It is strongly recommended to avoid any harsh treatment of the samples.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

The test CHORUS Promonitor ANTI-USTEKINUMAB (REF. 86713) can be performed simultaneously ONLY with the other CHORUS Promonitor tests.

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense the following volumes in well no. 1 of each device:

SAMPLE	130 µl/device
CALIBRATOR	130 µl/device
POSITIVE CONTROL	130 µl/device

At each change of batch, use a device for the calibrator.

4. Place the devices in the CHORUS TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the obtained results. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the positive control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The CHORUS TRIO instrument expresses the result in AU/ml, calculated on the basis of a batch-dependent curve stored in the instrument.

The test on the examined serum can be interpreted as follows:

POSITIVE: when the result is ≥ 3.5 AU/ml

NEGATIVE: when the result is ≤ 1.5 AU/ml

DOUBTFUL/EQUIVOCAL: when the result is between 1.6 and 3.4 AU/ml.

Doubtful/equivocal results are obtained in case of values which ranges between Limit of detection (< 1.5 AU/ml) and Low Limit of Quantitation (3.5 AU/ml)

11. LIMITATIONS

1. Trough antibody levels should be measured. The optimum sample collection time is immediately before the UTK injection.

2. This assay measures binding antibodies to UTK. It does not measure immunocomplexes formed between UTK and anti-UTK antibodies.
3. This assay is not intended for use as an aid in diagnosis, neither as a companion diagnostic test. The clinician should interpret test results in the context of other variables like sampling window, methodology, UTK levels and other clinical findings like disease activity.
4. The assay performance characteristics have not been established for matrices other than serum.
5. This product should only be used by appropriately trained laboratory professional personnel.

12. CALIBRATION RANGE

Calibration range : 1.5-296.0 AU/ml

For samples > 296.0 AU/ml retest the diluted sample in the Negative Control/Sample Diluent (REF 83607 - not supplied with the kit).

13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, are negative, with values < 1.5 AU/ml.

14. ANALYTICAL SPECIFICITY

3 samples were tested (1 Negative, 1 Low Positive and 1 High Positive) to which the following interfering agents were added:

Bilirubin (4.5 mg/dl – 45 mg/dl)

Triglycerides (250 mg/dl – 1500 mg/dl)

Hemoglobin (2.5 mg/ml – 10 mg/ml)

The presence of the above interfering substances in the test serum does not alter the test result.

3 Rheumatoid Factor (RF) positive samples 100 IU/ml (1) and 200 IU/mL (2) were spiked with Anti-Ustekinumab:

High level – 150.0 AU/ml

Low level – 20.0 AU/ml

Negative – 1.5 AU/ml

The presence of RF in the test sera does not alter the test result.

15. CROSS-REACTIONS

15 samples, not containing anti-Ustekinumab antibodies, but containing anti-IFX antibodies (10.1 – 474.6 AU/ml) (5 samples), anti-ADL antibodies (21.1 – 1100.4 AU/ml) (5 samples) and anti-VDZ antibodies (3.7 – 337.0 AU/ml) (5 samples), were tested.

No cross-reactions were found.

16. METHOD COMPARISON

In an experimentation 123 samples have been tested with Diesse kit and with Promonitor-ANTI-UTK ELISA kit.

Data are summarized in the following table:

		Reference		
		+	-	Tot.
Diesse	+	39	0	39
	-	0	84	84

	Tot.	39	84	123
--	-------------	----	----	-----

Percent Positive Agreement (~Diagnostic Sensitivity):

100.0% IC_{95%}: 91.0-99.9

Percent Negative Agreement: (~Diagnostic Specificity):

100.0% IC_{95%}: 95.6-99.0

Positive Predictive Value (PPV): 100.0% IC_{95%}: 100.0-100.0

Negative Predictive Value (NPV): 100.0% IC_{95%}: 100.0-100.0

The agreement between the two methods is excellent with a Cohen's Kappa of 1.0.

17. PRECISION AND REPEATABILITY

Sample	Within run		Between run	
	Mean (AU/ml)	CV%	Mean (AU/ml)	CV%
1	253.8	9.4	247.2	7.9
2	64.7	12.8	60.6	10.5
3	27.1	6.8	21.5	13.4
4	3.7	8.3	2.9	14.3
5	<1.5	-	<1.5	-
6	<1.5	-	<1.5	-

Sample	Between lots		Between instruments	
	Mean (AU/ml)	CV%	Mean (AU/ml)	CV%
1	258.4	14.0	257.6	12.6
2	42.6	14.0	58.2	13.1
3	18.5	14.3	24.2	14.5
4	6.4	2.7	3.0	15.0
5	<1.5	-	<1.5	-
6	<1.5	-	<1.5	-

18. LINEARITY

The linearity of the test was evaluated with 3 samples with known concentration of anti-UTK antibodies. For each sample serial dilutions were done; the linearity was demonstrated in the range 1.8 - 279.7 AU/ml.

19. BIBLIOGRAPHY

1. Lucas Guillo, Ferdinando D'Amico, Silvio Danese, Laurent Peyrin-Biroulet "Ustekinumab for Extra-intestinal Manifestations of Inflammatory Bowel Disease: A Systematic Literature Review" J Crohns Colitis 2021 Jul 5;15(7):1236-1243.
2. Sands BE, Sandborn WJ, Panaccione R, O'Brien CD, Zhang H, Johanns J, Adedokun OJ, Li K, Peyrin-Biroulet L, Van Assche G, Danese S, Targan S, Abreu MT, Hisamatsu T, Szapary P, Marano C; "Ustekinumab as Induction and Maintenance Therapy for Ulcerative Colitis". N Engl J Med. 2019 Sep 26;381(13):1201-1214.
3. Restellini S, Khanna R, Afif W. "Therapeutic Drug Monitoring With Ustekinumab and Vedolizumab in Inflammatory Bowel Disease." Inflamm Bowel Dis. 2018 Sep 15;24(10):2165-2172.
4. Menting SP, van den Reek JM, Baerveldt EM, de Jong EM, Prens EP, Lecluse LL, Wolbink GJ, Van der Kleij D, Spuls PI, Rispens T. The correlation of clinical efficacy, serum trough levels and antidrug antibodies in ustekinumab-treated patients with psoriasis in a clinical-practice setting. Br J Dermatol. 2015 Sep;173(3):855-7.
5. Klenske E, Osaba L, Nagore D, Rath T, Neurath MF, Atreya R. Drug Levels in the Maternal Serum, Cord Blood and Breast Milk of a Ustekinumab-Treated Patient with Crohn's Disease. J Crohns Colitis. 2019 Feb 1;13(2):267-269. Martín-González S, Urigoitia-Ugalde P, Careaga J, Nagore D, Navarro R, Izu-Belloso R, Blanch L, González Del Tánago J, Martínez A, Ruiz Del Agua A. Optimal concentration range of ustekinumab in patients with plaque-type psoriasis. J Am Acad Dermatol. 2019 Jun;80(6):1782-1784.

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CHORUS Promonitor ANTI-USTEKINUMAB

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Inflammatory Bowel Disease." Inflamm Bowel Dis. 2018 Sep 15;24(10):2165-2172.

4. Menting SP, van den Reek JM, Baerveldt EM, de Jong EM, Prens EP, Lecluse LL, Wolbink GJ, Van der Kleij D, Spuls PI, Rispens T. The correlation of clinical efficacy, serum trough levels and antidrug antibodies in ustekinumab-treated patients with psoriasis in a clinical-practice setting. Br J Dermatol. 2015 Sep;173(3):855-7.
5. Klenske E, Osaba L, Nagore D, Rath T, Neurath MF, Atreya R. Drug Levels in the Maternal Serum, Cord Blood and Breast Milk of a Ustekinumab-Treated Patient with Crohn's Disease. J Crohns Colitis. 2019 Feb 1;13(2):267-269. Martín-González S, Urigoitia-Ugalde P, Careaga J, Nagore D, Navarro R, Izu-Belloso R, Blanch L, González Del Tánago J, Martínez A, Ruiz Del Agua A. Optimal concentration range of ustekinumab in patients with plaque-type psoriasis. J Am Acad Dermatol. 2019 Jun;80(6):1782-1784.

	EN IT DE ES CS RO	Date of manufacture Data di fabbricazione Herstellungsdatum Fecha de fabricación Datum výroby Data di fabbricazione	FR GR PT PL RO	Date de fabrication Ημερομηνία Παραγωγής Data de fabrico Data produkcji Data fabricației
	EN IT DE ES CS RO	Use By Utilizzare entro Verwendbar bis Fecha de caducidad Použijte do Utilizzare entro	FR GR PT PL RO	Utiliser jusqu'e Ημερομηνία λήξης Prazo de validade Zużyć do dnia A se folosi până la
	EN IT DE ES CS RO	Do not reuse Non riutilizzare Nicht wieder verwenden No reutilizar Nepoužívejte znovu Nu reutilizați	FR GR PT PL RO	Ne pas réutiliser Μην κάνετε επαναληπτική χρήση Não reutilizar Nie używać ponownie A nu se refolosi
	EN IT DE ES CS RO	Caution, consult accompanying documents Attenzione, vedere le istruzioni per l'uso Achtung, die Gebrauchsanleitung lesen Atención, ver instrucciones de uso Pozor, řídte se návodem k použití Attenzione, vedere le istruzioni per l'uso	FR GR PT PL RO	Attention voir notice d'instructions Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα Atenção, consulte a documentação incluída Uwaga, patrz instrukcja obsługi Atenție, consultați documentele însoțitoare
	EN IT DE ES CS RO	Manufacturer Fabbricante Hersteller Fabricante Výrobce Fabbricante	FR GR PT PL RO	Fabricant Κατασκευαστής Fabricante Producuent Producător
	EN IT DE ES CS RO	Contains sufficient for <n> tests Contenuto sufficiente per "n" saggi Inhalt reicht für „n“ Tests Contenido suficiente para <n> ensayos Dostatečný obsah pro „n“ testů Contenuto sufficiente per "n" saggi	FR GR PT PL RO	Contenu suffisant pour "n" tests Περιεχόμενο επαρκές για «n» εξετάσεις Conteúdo suficiente para “n” ensaios Wystarczająca ilość materiału do „(liczby)” badań Conținut suficient pt <n> teste
	EN IT DE ES CS RO	Temperature limitation Limiti di temperatura Temperaturgrenzwerte Limite de temperatur Limity teploty Limiti di temperatura	FR GR PT PL RO	Limites de température Περιορισμοί θερμοκρασίας Limites de temperatura Granice temperatury Limită de temperatură
	EN IT DE ES CS RO	Consult Instructions for Use Consultare le istruzioni per l'uso Die Gebrauchsanleitung lesen Consulte las instrucciones de uso Řídte se návodem k použití Consultare le istruzioni per l'uso	FR GR PT PL RO	Consulter les instructions d'utilisation Συμβουλευτείτε τις οδηγίες χρήσης Consulte as instruções de utilização Należy zapoznać się z instrukcją obsługi Pentru utilizare consultați instrucțiunile
	EN IT DE ES CS RO	Catalogue number Numero di catalogo Katalognummer Número de catálogo Číslo katalogu Numero di catalogo	FR GR PT PL RO	Référence du catalogue Αριθμός καταλόγου Referência de catálogo Numer katalogowy Număr de catalog
	EN IT DE ES CS RO	In Vitro Diagnostic Medical Device Dispositivo medico-diagnostico in vitro Medizinisches In-vitro-Diagnostikum Producto sanitario para diagnóstico in vitro Zdravotnický nástroj k diagnostice in vitro Dispozitiv medical pentru diagnosticare in vitro	FR GR PT PL RO	Dispositif médical de diagnostic in vitro In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν Dispositivo médico para diagnóstico in vitro Wyrób medyczny do diagnostyki in vitro Dispozitiv medical pentru diagnosticare in vitro
	EN IT DE ES CS RO	Batch code Codice del lotto Chargennummer Código de lote Kód šarže Codice del lotto	FR GR PT PL RO	Code du lot Αριθμός Παρτίδας Código do lote Kod partii Lot
	EN IT DE ES CS RO	CE marking of conformity Marcatura CE di conformità CE-Konformität Skennzeichnung Marcado CE de conformidad Označení shody CE Marcajul de conformitate CE	FR GR PT PL RO	Marquage de conformité CE Σημόνση συμμορφώση CE Marcação CE de conformidade Oznakowanie zgodności CE Marcajul de conformitate CE