ENZYWELL

TOSCANA VIRUS IgG/IgM

REF 91067 (96 tests)

Prodotto da/Manufactured by/Fabricado por/
Fabriqué par: DIESSE Diagnostica Senese
Strada dei Laghi, 39
53035 Monteriggioni (Siena) – Italy

INDICE / INDEX / INDICE / SOMMAIRE

1. UTILIZZAZIONE / INTENDED USE / INDICACIONES DE USO/ UTILISATION
2. INTRODUZIONE / SUMMARY AND EXPLANATION OF TEST / RESUMEN Y EXPLICACIÓN DEL TEST / INTERET CLINIQUE
3. PRINCIPIO DEL METODO / PRINCIPLE OF THE TEST / PRINCIPIO DEL MÉTODO / PRINCIPE
4. COMPOSIZIONE DEL KIT E PREPARAZIONE DEI REAGENTI / KIT CONTENTS AND REAGENT PREPARATION / COMPONENTES DEL KIT Y PREPARACIÓN DEL REACTIVO / CONTENU DU KIT ET PREPARATION DES REACTIFS
5. MODALITA' DI CONSERVAZIONE E STABILITA' DEI REAGENTI / STORAGE AND STABILITY OF REAGENTS / CONSERVACIÓN Y ESTABILIDAD DE LOS REACTIVOS / MODALITES DE CONSERVATION ET STABILITE DES REACTIFS
6. PRECAUZIONI / PRECAUTIONS / PRECAUCIONES DE USO / PRECAUTIONS D'UTILISATION
7. TIPO DI CAMPIONE E CONSERVAZIONE / TYPE AND STORAGE OF SAMPLE / TIPO DE MUESTRA Y CONSERVACION / TYPE D'ECHANTILLON ET CONSERVATION
8. PROCEDIMENTO / TEST PROCEDURE / PROCEDIMIENTO DEL TEST / PROCEDURE /
9. SCHEMA DEL SAGGIO / SCHEME OF TEST PROCEDURE / ESQUEMA DEL PROCEDIMIENTO DEL TEST / SCHEMA D'ANALYSE
10. VALIDAZIONE DEL TEST / TEST VALIDATION / VALIDACIÓN DEL TEST / VALIDATION DU TEST
11. INTERPRETAZIONE DEL TEST / INTERPRETATION OF RESULTS / INTERPRETACIÓN DE RESULTADOS / INTERPRETATION DU TEST / INTERPRETAÇÃO DOS RESULTADOS
12. LIMITAZIONI DELLA PROCEDURA / LIMITATIONS OF THE PROCEDURE / LIMITACIONES DEL TEST / LIMITES DU TEST
13. SENSIBILITA' E SPECIFICITA' DIAGNOSTICA / DIAGNOSTIC SENSITIVITY AND SPECIFICITY / SENSIBILIDAD Y ESPECIFICIDAD DIAGNOSTICAS / SENSIBILITE ET SPECIFICITE DIAGNOSTIQUES
14. PRECISIONE / PRECISION / PRECISIÓN / PRECISION
15. GUIDA DEI PROBLEMI DI UTILIZZO / "TROUBLE SHOOTING" / GUIA DE RESOLUCION DE PROBLEMAS / RESOLUTION DES PROBLEMES
16. BIBLIOGRAFIA / REFERENCES / BIBLIOGRAFÍA / BIBLIOGRAPHIE

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ENZYWELL TOSCANA VIRUS IgG/IgM

REF 91067

(Italiano)

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17. BIBLIOGRAFIA

1. Nicoletti L. et al. 1991. Central Nervous System involvement during infection by Phlebovirus Toscana of residents in natural foci in Central Italy (1977-1988). *Am. J. Trop. Med. Hyg.* 45, 429-434.
2. Schwarz et al. 1996. Immunoblot detection of antibodies to Toscana Virus. *J. Medical Virology* 49, 83-86.
3. Valassina M. et al. 1998. Diagnostic potential of Toscana virus in protein expressed in *Escherichia coli*. *J. Clin. Microbiol.* 36, 3170-3172.
4. Soldateschi D. et al. 1999. Laboratory diagnosis of Toscana virus infection by enzyme immunoassay with recombinant viral nucleoprotein. *J. Clin. Microbiol.* 37, 649-652.

INSTRUCTIONS FOR USE

ENZYWELL TOSCANA VIRUS IgG/IgM

REF 91067

(English)

1. INTENDED USE

IMMUNOENZYMATIC METHOD FOR THE QUALITATIVE DETERMINATION OF IgG AND IgM ANTIBODIES TO TOSCANA VIRUS IN HUMAN SERUM, AS AN AID IN THE DIAGNOSIS OF TOSCANA VIRUS INFECTION.

2. SUMMARY AND EXPLANATION OF THE TEST

The Toscana Virus, belonging to the Phlebovirus genus of the Bunyaviridae family, was isolated for the first time in Italy in 1971. Since then, its presence in the mediterranean area has been established and its role as an etiological agent in human disease has been demonstrated. The disease caused by Toscana Virus is characterized by fever, strong headache and meningitis or meningoencephalitis. The infection is transmitted by hematophagous insects of the Phlebotomus genus and is seasonal with a peak in the summer months, due to the life cycle of the vector. Epidemiological data indicate that in areas where the Toscana Virus is endemic, many, perhaps the majority of cases of aseptic meningitis during the summer period are caused by this virus. Until now, the diagnosis of Toscana virus infection has been mainly performed serologically. Specific antibodies of the IgM class are present during the acute phase, while the presence of specific IgG antibodies in the absence of IgM indicates a previous infection. All the serological methods adopted up to the present time require propagation of the virus in cell cultures; these methods can only be performed in specialized laboratories, while the infection cannot be diagnosed in the majority of serological laboratories. As it has been shown that the viral nucleoprotein is the chief target of the anti-Toscana virus antibodies, this has been expressed by recombinant DNA technology in E. coli and used to set up an immunoenzymatic kit for the detection of specific IgG and IgM, as an aid in the serological study of this virus.

3. PRINCIPLE OF THE TEST

The test is based on the ELISA technique (Enzyme linked Immunosorbent Assay).

Determination for the IgG is performed after incubation of the serum sample with a solid phase precoated with the highly purified recombinant nucleoprotein of the Toscana virus. Specific IgG bind to the solid phase and are detected by further incubation with a mouse anti-human IgG monoclonal antibody conjugated with peroxidase.

The IgM are detected by the capture method. The total IgM present in the serum bind to the solid phase precoated with a mouse anti-human IgM monoclonal antibody. After washing, a complex composed of the biotinylated recombinant nucleoprotein of Toscana virus and the Streptavidin-peroxidase conjugate.

In both cases, the whole complex (for the IgG: specific IgG anti-Toscana virus with anti-IgG labelled with POD; for the IgM: IgM-nucleoprotein-biotin/Streptavidin-peroxidase) is determined by the addition of the substrate which becomes blue coloured in the presence of peroxidase. A stable yellow colour is obtained by addition of a blocking solution.

4. KIT CONTENTS AND REAGENT PREPARATION

N.B. The present kit contains all the reagents necessary for the assay of both IgG and IgM. The two types of assay must be performed separately.

- Reagents are sufficient for 2 x 48 determinations.
- **Bring to room temperature before use.**

For the IgG

MT PLATE MICROPLATE 3 x 2 strips coated with the recombinant nucleoprotein of Toscana virus for the determination of IgG.

Use: open the package on the opposite side from the code (M followed by the lot number) which is useful for identification purposes, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

CONTROL + IgG POSITIVE CONTROL 1 x 1.6 mL

Contents: Diluted human serum containing a known amount of anti-Toscana virus IgG, in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour of the controls is proportional to the relative antibody titer.

CONTROL CUT-OFF IgG CUT-OFF CONTROL 1 x 2 mL

Contents: Diluted human serum containing a known amount of anti-Toscana virus IgG in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour of the controls is proportional to the relative antibody titer.

CONJ IgG CONJUGATE 2 x 5 mL

Contents: anti-human IgG monoclonal antibodies labelled with Peroxidase, in phosphate buffer solution containing phenol 0.05% and Bronidox 0.02%. Ready for use without further dilution.

For the IgM

MT PLATE MICROPLATE 3 x 2 strips coated with anti-human IgM monoclonal antibodies for the determination of IgM.

Use: open the package on the opposite side from the code (M followed by the lot number) which is useful for identification purposes, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

Ag ANTIGEN Freeze-dried x 6 vials

Contents: Recombinant antigen of the Toscana virus linked to biotin in phosphate buffer containing lactose.

Preparation: reconstitute with the volume of conjugate reported on the label.

CONJ IgM CONJUGATE 2 x 5 mL

Contents: streptavidin conjugated to peroxidase in buffer solution. Ready for use without further dilution.

Stability: The unopened product is stable up to the expiry date. The immunocomplex must be prepared 45 minutes before use.

Reagents used in both assays

CONTROL - IgG/ IgM NEGATIVE CONTROL (PF93900) 2 x 1.6 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

WASH BUF 10x WASH BUFFER 10X (PF93603) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5% .

Preparation: dilute the required volume 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SAMP DIL 10 DILUENT 10 (PF93621) 1 x 100 mL Ready for use **INTERCHANGEABLE BETWEEN LOTS**

To be used to dilute samples.

Contents: PBS containing proteins 10% w/v and sodium azide 0,09% containing methyl orange as dye.

SUBS TMB SUBSTRATE (PF93619) 1 x 12 mL Ready for use **INTERCHANGEABLE BETWEEN LOTS**

Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilised in citrate buffer 0.05 mol/L (pH 3.8).

H₂SO₄ 0.3 M STOP SOLUTION (PF93602) 1x16 mL **INTERCHANGEABLE BETWEEN LOTS**
H₂SO₄ 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

MATERIALS REQUIRED BUT NOT PROVIDED.

- Incubator at 37°C
- Microplate reader (wave length 450 or 450/620 nm, with linearity up to OD >= 2000)

- Microplate washer (preferable) able to dispense volumes in the range 225-375 µL
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µL solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
Microplate IgG and IgM	5 weeks at 2/8°C, polythene bag
Controls	5 weeks at 2/8°C
Conjugate IgG	5 weeks at 2/8°C
Conjugate IgM	5 weeks at 2/8°C
Reconstituted antigen	the solution cannot be stored and must be used immediately after reconstitution
Substrate	until the expiry date at 2/8°C, 1 week at 15-30°C in the dark
Sample Diluent	until the expiry date at 2/8°C
Wash Buffer	2 weeks at 2/8°C, 5 days at 15/30°C.
Stop Solution	until the expiry date at 2/8°C

6. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

Caution: This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The Wash Buffer contains detergents
 - b) The conjugate contains phenol
 - c) The substrate is acid
 - d) The controls contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
6. Do not allow wells to become dry during the assay procedure.
7. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
8. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
9. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO₂ incubators must not be used.
10. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
11. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.
12. Use of the kit with automatic instruments must be validated by the user.
13. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:
 - installation and particular requisites
 - operating principles, instructions, precautions and risks
 - manufacturer's specifications and instrument performance
 - servicing and maintenance

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be carefully mixed before performing the test. Heat inactivation can lead to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric or contaminated samples should be avoided.

The test is not applicable to human plasma.

8. TEST PROCEDURE

- Prepare the required number of wells.
- Prepare the washing buffer by diluting the Rinsing Buffer 10x (100 mL + 900 mL H₂O).
- Prepare the antigen by reconstituting the freeze-dried product with the conjugate volume shown on the label.

Dilute samples 1:101 distributing 10 µL of serum into 1 mL of diluent. Dispense 100 µL of each diluted sample per well (duplicate testing is recommended). Place UNDILUTED controls in a strip (100 µL in each well). The minimum requisite is 1 negative control, 2 cut-off, 1 positive control for IgG and 2 negative controls for IgM. Leave one well for the blank, performed using 100 µL of the substrate mixture.

Wells are covered with protective film and incubated for 45 minutes at 37°C. After washing four times for 30 seconds (300 µl), add 100 µL of conjugate (anti-IgG for the determination of IgG, or the antigen/conjugate complex for determination of the IgM) to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film. The plate is washed again 4 times, as described above. Finally, the substrate is distributed, 100 µL/well. After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µL of Stop Solution.

The adsorbance (O.D.) is read at 450 nm or 450/620 nm within 30 min.

9. SCHEME OF TEST PROCEDURE

STEP 1	Place 100 µL of diluted sample/controls to the wells Incubate for 45 min. at 37°C Wash 4 times (300 µl)
STEP 2	Add 100 µL of conjugate to each well (anti-IgG for determination of the IgG; the antigen/conjugate complex for determination of the IgM) Incubate for 45 min. at 37°C Wash 4 times (300 µl)
STEP 3	Add 100 µL of Substrate to each well Incubate for 15 min. at R.T.
STEP 4	Add 100 µL of Stop Solution Read absorbance at 450 nm within 30 min

10. VALIDATION

Subtract the value of the blank (≤ 0.150) from all the other readings. For the IgG assay the OD of the IgG Cut-Off must be ≥ 0.200 . The Positive control must have an O.D. at least 1.5 times that of the Cut-Off serum. The ratio between Negative Control and Cut-off must be ≤ 0.6 . For the assay of the IgM the OD of the Negative control must be < 0.150 .

11. INTERPRETATION OF THE RESULTS

Calculation of the IgM Cut-Off: add 0.250 to the OD of the Negative Control.

Qualitative results

If the adsorbance of the sample is higher than that of the Cut-Off, the sample is positive for the presence of specific IgG.

Calculate the ratio (Index) between the average O.D. value of the sample and that of the Cut-off. The sample is considered:

Positive: if the ratio is > 1.2 .

Doubtful: $\pm 20\%$ of the Cut-Off.

Negative: if the ratio is < 0.8 .

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

12. LIMITATIONS OF THE PROCEDURE

The test result should be used in conjunction with information available from the evaluation of other diagnostic procedures.

13. ANALYTICAL SPECIFICITY

Cross-reactoins cannot be excluded with other correlated serotypes of the Phlebovirus genus such as Sicily and Naples.

14. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

97 samples were tested from patients belonging to the following 3 groups:

1. 22 sera from patients with neurological symptoms of meningitis caused by Toscana Virus.
2. 22 sera from pediatric patients
3. 53 samples of patients from the routine practice-

The samples were tested with the kit under study and the results were compared with an immunoenzymatic system which utilises the native antigen. Discordant results (0 in the IgM test, 10 in the IgG test) were further studied in Western Blotting. The results obtained are reported below:

Toscana virus IgG :

	Other method	
	P	N
Diesse P	38	2
N	2	55

Diagnostic sensitivity 95.0% (IC at 95%: 83.5 - 98.6%)

Diagnostic specificity 96.5 % (IC at 95%: 88.1 - 99.0%)

Toscana virus IgM

	Other method	
	P	N
Diesse P	19	0
N	0	78

Diagnostic sensitivity 100.% (IC al 95%: 83.2-100%)

Diagnosotic specificity 100 % (IC al 95%: 95.3-100%)

15. PRECISIONToscana virus IgG

Table 1 Precision “In run” (12 replicates of the Cut Off in thre different lots)

O.D.	Lot. N. 36	37	38
Cut Off n=12	0.384	0.234	0.386
CV%	8	11	6

Table 2 Precision between runs (average of the Index in three different lots)

Sample	Lots			average	cv%
	036	037	038		
1	0.79	0.65	0.5	0.65	22
2	1.74	1.36	1	1.55	17
3	2.4	3.2	2.5	2.7	16

Toscana Virus IgM

Table 3 Precision “in run” (12 replicates of a borderline sample)

O.D.	Lot. N. 36	37	38
Cut Off n=12	0.235	0.245	0.387
CV%	5	5	9

Table 4 Precision between runs (average of the Index in three different lots)

Sample	Lots			average	cv%
	036	037	038		
1	0.68	0.59	0.35	0.54	32
2	1.46	2.08	1.32	1.62	25
3	2.95	5.1	3.23	3.76	32

16. TROUBLE SHOOTING GUIDE

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert point 4 for correct code).
		Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow).Repeat test
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.
	Inadequate washing	Ensure that wash apparatus works well
Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate Color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

17. REFERENCES

1. Nicoletti L. et al. 1991. Central Nervous System involvement during infection by Phlebovirus Toscana of residents in natural foci in Central Italy (1977-1988). Am. J. Trop. Med. Hyg. 45, 429-434.
2. Schwarz et al. 1996. Immunoblot detection of antibodies to Toscana Virus. J. Medical Virology 49, 83-86.
3. Valassina M. et al. 1998. Diagnostic potential of Toscana virus in protein expressed in Escherichia coli. J. Clin. Microbiol. 36, 3170-3172.
4. Soldateschi D. et al. 1999. Laboratory diagnosis of Toscana virus infection by enzyme immunoassay with recombinant viral nucleoprotein. J. Clin. Microbiol. 37, 649-652.

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17. BIBLIOGRAPHIE

1. Nicoletti L. et al. 1991. Central Nervous System involvement during infection by Phlebovirus Toscana of residents in natural foci in Central Italy (1977-1988). Am. J. Trop. Med. Hyg. 45, 429-434.
5. Schwarz et al. 1996. Immunoblot detection of antibodies to Toscana Virus. J. Medical Virology 49, 83-86.
6. Valassina M. et al. 1998. Diagnostic potential of Toscana virus in protein expressed in Escherichia coli. J. Clin. Microbiol. 36, 3170-3172.
7. Soldateschi D. et al. 1999. Laboratory diagnosis of Toscana virus infection by enzyme immunoassay with recombinant viral nucleoprotein. J. Clin. Microbiol. 37, 649-652.

	EN ES IT	Date of manufacture Fecha de fabricación Data di fabbricazione	FR GR PT	Date de fabrication Ημερομηνία Παραγωγής Data de fabrico
	EN ES IT	Use By Fecha de caducidad Utilizzare entro	FR GR PT	Utiliser jusque Ημερομηνία λήξης Prazo de validade
	EN ES IT	Caution, consult accompanying documents Atención, ver instrucciones de uso Attenzione, vedere le istruzioni per l'uso	FR GR PT	Attention voir notice d'instructions Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα Atenção, consulte a documentação incluída
	EN ES IT	Manufacturer Fabricante Fabbricante	FR GR PT	Fabricant Κατασκευαστής Fabricante
	EN ES IT	Contains sufficient for <n> tests Contenido suficiente para <n> ensayos Contenuto sufficiente per "n" saggi	FR GR PT	Contenu suffisant pour "n" tests Περιεχόμενο επαρκές για «n» εξετάσεις Conteúdo suficiente para "n" ensaios
	EN ES IT	Temperature limitation Límite de temperatura Limiti di temperatura	FR GR PT	Limites de température Περιορισμοί θερμοκρασίας Limites de temperatura
	EN ES IT	Consult Instructions for Use Consulte las instrucciones de uso Consultare le istruzioni per l'uso	FR GR PT	Consulter les instructions d'utilisation Συμβουλευτείτε τις οδηγίες χρήσης Consulte as instruções de utilização
	EN ES IT	Biological risks Riesgo biológico Rischio biologico	FR GR PT	Risques biologiques Βιολογικοί κίνδυνοι Risco biológico
	EN ES IT	Catalogue number Número de catálogo Numero di catalogo	FR GR PT	Référence du catalogue Αριθμός καταλόγου Referência de catálogo
	EN ES IT	In Vitro Diagnostic Medical Device Producto sanitario para diagnóstico in vitro Dispositivo medico-diagnostico in vitro	FR GR PT	Dispositif médical de diagnostic in vitro In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν Dispositivo médico para diagnóstico in vitro
	EN ES IT	Batch code Código de lote Codice del lotto	FR GR PT	Code du lot Αριθμός Παρτίδας Código do lote

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