ENZYWELL

CMV SCREEN

REF 91013 (96 tests)

REF 91173 (6 x 96 tests)

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(English)

1. INTENDED USE

IMMUNOENZYMATIC KIT FOR THE QUALITATIVE DETERMINATION OF ANTI-CYTOMEGALOVIRUS IMMUNOGLOBULINS IN HUMAN SERUM OR PLASMA.


2. SUMMARY AND EXPLANATION OF THE TEST

Cytomegalovirus is a herpes virus transmitted by close human contact. No symptoms of infection are apparent in the majority of cases. However, the virus is very dangerous and may be fatal in immunodepressed patients. Serum-negative female patients who become infected during pregnancy may transmit the disease to the fetus. In 95% of cases this occurs without symptoms, but some neonates may present jaundice, hepato-splenomegaly and retarded psycho-motorial development. For this reason it is of great importance to determine the immunitary state of the patient and check for serum conversion, particularly for transfusion purposes.

3. PRINCIPLE OF THE PROCEDURE

The test for the assay of anti-Cytomegalovirus (CMV) immunoglobulins (total Ig) is based on the principle of analysis by competition. The anti-CMV antibodies present in the sample compete with the marker (peroxydase/anti-CMV monoclonal antibody conjugate) to occupy the few binding sites available on the antigen coated on the solid phase. The higher the antibody concentration in the sample, the lower the number of conjugated antibodies which are bound. The unbound components are eliminated by washing, and the bound enzymatic activity is determined colorimetrically by transformation of a chromogen substrate (the blue colouring of the reaction turns yellow after the addition of the Stop reagent). The intensity of the colour which develops is inversely proportional to the antibody concentration in the sample.

This kit also offers the opportunity to monitor sample addition in order to reduce the risk of incorrect results due to omission of the sample. The conjugate which has a yellow-orange colour in the absence of the sample, turns dark red in the presence of serum, plasma or negative or positive controls: The colour change is easily recorded if the plate is read before incubation at 540 nm.

4. REAGENTS AND REAGENT PREPARATION

- Reagents are sufficient for 96 tests (**REF** 91013) or 6x96 tests (**REF** 91173).

- **Bring to room temperature before use.**

Code 91013 contains:

MT PLATE MICROPLATE 1 plate with 12x8 wells coated with Cytomegalovirus antigens.

Use: open the package at the opposite end from the code (G, followed by the lot number) which is useful for its identification, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expell the air and seal by pressing the closure.

CAL CALIBRATOR 1 x 1 mL

Contents: Newborn Calf Serum in phosphate buffer 0.01 mol/L with phenol 0.05% and Bronidox 0.02%, liquid, ready for use.

CONTROL + POSITIVE CONTROL SERUM 1 x 1 mL

Contents: Dilute human serum in phosphate buffer 0.01 mol/L with BSA 1%, Bronidox 0.02% and phenol 0.05%, liquid, ready for use without further dilution.

CONJ CONJUGATE 1 x 16 mL Ready for use

Contents: anti-CMV monoclonal antibodies labelled with Horse Radish Peroxidase, in phosphate buffer with phenol 0.05% and Bronidox 0.02%.

WASH BUF 10x WASH BUFFER 10X (PF93603) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5%.

Preparation: dilute the required volume 1:10 with distilled or deionized water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SUBS TMB SUBSTRATE (PF93619) 1 x 12 mL Ready for use **INTERCHANGEABLE BETWEEN LOTS**

Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilized in citrate buffer 0.05 mol/L (pH 3.8).

H₂SO₄ 0.3 M STOP SOLUTION (PF93602) 1x16 mL **INTERCHANGEABLE BETWEEN LOTS**

H₂SO₄ 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

Code 91173 contains:

MT PLATE MICROPLATE 6 plates with 12x8 wells coated with Cytomegalovirus antigens.

Use: open the package at the opposite end from the code (G, followed by the lot number), remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expell the air and seal by pressing the closure.

CAL CALIBRATOR 1 x 1 mL

Contents: Newborn Calf Serum in phosphate buffer 0.01 mol/L with phenol 0.05% and Bronidox 0.02%, liquid, ready for use.

CONTROL + POSITIVE CONTROL SERUM 1 x 1 mL

Contents: Dilute human serum in phosphate buffer 0.01 mol/L with BSA 1%, Bronidox 0.02% and phenol 0.05%, liquid, ready for use without further dilution.

CONJ CONJUGATE 6 x 16 mL Ready for use

Contents: anti-CMV monoclonal antibodies labelled with Horse Radish Peroxidase, in phosphate buffer with phenol 0.05% and Bronidox 0.02%.

WASH BUF 10x WASH BUFFER 10X (PF93603) 3 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5%.

Preparation: dilute the required volume 1:10 with distilled or deionized water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SUBS TMB SUBSTRATE (PF93619) 6 x 12 mL Ready for use **INTERCHANGEABLE BETWEEN LOTS**

Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilized in citrate buffer 0.05 mol/L (pH 3.8).

H₂SO₄ 0.3 M STOP SOLUTION (PF93602) 1x130 mL **INTERCHANGEABLE BETWEEN LOTS**

H₂SO₄ 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (6).

POLYTHENE BAG (1).

MATERIALS REQUIRED BUT NOT PROVIDED

- Incubator at 37°C
- Microplate reader, wavelength 450 or 450/620 nm, with OD linearity up to 2,000 (at least).
- Microplate washer (preferable) able to dispense volumes in the range 225-375 µl
- Distilled or deionized water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
MICROPLATE	4 weeks at 2/8°C in the polythene bag
CALIBRATOR	1 week at room temperature
CONTROL SERUM	1 week at room temperature
CONJUGATE	1 week at room temperature.
SUBSTRATE	until the expiry date at 2/8°C, 1 week at 15-30°C; store in the dark
WASH BUFFER	2 weeks at 2/8°C, 5 days at 15/30°C.
STOP SOLUTION	until the expiry date at 2/8°C

6. PRECAUTIONS AND WARNINGS***FOR IN VITRO DIAGNOSTIC USE. STORE AT 2-8°C.*****Caution:**

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin. Waste disposal: the used devices containing materials of human origin must be handled as infectious material and disposed of according to the full requirements of this issue.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The Rinsing Buffer contains detergents
 - b) The conjugate and controls contain phenol
 - c) The substrate is acid.
 If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid 2 M used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. The extinction value of the calibrator, and also of the samples, may differ slightly from plate to plate. Therefore, if strips taken from different plates, even if from the same batch, are to be used in the same run, the calibrator must be repeated.
2. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.**
3. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
4. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
5. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
6. Avoid to use the self-defrosting freezers for the storage of samples.
7. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
8. Do not allow wells to become dry during the assay procedure.
9. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for use with the substrate and for use with the conjugate.
10. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
11. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to

support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO₂ incubators must not be used.

12. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
13. Use of highly hemolyzed samples, incompletely clotted sera or samples with microbial contamination may give rise to erroneous results.
14. Avoid contaminating microwells with the dust from disposable gloves.
15. The use of the kit with automated instruments must be validated by the user.
15. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:
 - installation and particular requisites
 - operating principles, instructions, precautions and risks
 - manufacturer's specifications and instrument performance
 - servicing and maintenance.

7. TYPE OF SPECIMENS AND STORAGE:

No special preparation of the patient is necessary. Serum or plasma samples can be used. The use of anticoagulants such as citrate, heparin or EDTA does not interfere in the test. Samples can be stored for 7 days at 2-8°C and can be thawed a maximum of 3 times. For longer storage, freeze at -20°C. Do not keep the samples in auto-defrosting freezers. Freezing/thawing of the samples or heat inactivation for 30 min at 56°C does not influence the results.

Strongly lipemic, icteric or contaminated samples should be avoided.

8. TEST PROCEDURE

1. Distribution of the samples:
Dispense 30 µL of calibrator and of the positive control in wells on the plate (preferably in duplicate), and 30 µL of the samples being tested in the remaining wells. Add 100 µL of conjugate to each well. Mix carefully.
2. Incubation:
Incubate the plate covered with the adhesive foil at 36-39°C for 60 ± 5 minutes.
3. Rinsing:
Remove the adhesive film, aspirate the contents of all the wells and rinse 4 times, filling each well with about 0.3 mL of washing solution. Wait 30" between each wash.
4. Distribution of the substrate:
Dispense 100 µL of the substrate in each well.
5. Substrate incubation:
Incubate the plate for 15 ± 2 minutes at room temperature (18-30°C).
6. Interruption of the reaction:
Dispense 100 µL of blocking solution, in the same order as that followed for point 4.
7. Reading:
Take the photometric readings within 30 minutes at 450 nm or 450/620 nm.

9. SCHEME OF TEST PROCEDURE

- | | |
|--------|--|
| STEP 1 | Place 30 µL of calibrator, positive control and serum sample in the wells. Add 100 µL of conjugate to all wells. Mix well. |
| STEP 2 | Incubate for 60 ± 5 min. at 37-40°C

Wash 4 times (300 µL) |
| STEP 3 | Place 100 µL of Substrate in each well

Incubate for 15 ± 2 min. at room temperature |
| STEP 4 | Add 100 µL of Stop Solution

Read absorbance at 450 ± 5 nm within 30 min |

10. TEST VALIDATION

The test is valid if:

O.D. Calibrator at 450 nm ≥ 0.600 , or ≥ 0.560 at 450/620 nm.

O.D. Positive Control 450 nm/O.D. Calibrator 450 nm ≤ 0.50 .

If the results obtained do not meet the test validation specifications, contact your local technical assistant.

11. INTERPRETATION OF THE TEST

The Cut-Off is calculated by multiplying the O.D. of the Calibrator by the factor

$$\mathbf{0.6}$$

Example: Calibrator: O.D. = 1.100 Factor = 0.6

$$\text{Cut-Off} = 1.100 \times 0.6 = \mathbf{0.660}.$$

If the absorbance value of the sample is lower than that of the Cut-Off, the sample is positive for the presence of anti-Cytomegalovirus immunoglobulins.

The reproductibility data suggest the use of a “grey zone” of $\pm 10\%$ around the Cut-off O.D. In the example, samples with OD in the range 0.6 to 0.72 would therefore be doubtful.

12. LIMITATIONS

The test is not able to discriminate between the presence of IgG and that of IgM.

The test cannot be used on its own for clinical diagnostic purposes. A negative result does not exclude the possibility of the disease.

The Test Procedure and Interpretation of the Test must be followed. Non-repeatable reactive results may be obtained with any EIA procedure.

13. ANALYTICAL SPECIFICITY

Samples containing antibodies against Epstein Barr Virus (n=40), Herpes Simplex Virus (n=33) and anti-Nuclear (n=44), and samples containing high levels of bilirubin (n=40), proteins (n=20), lipids (n=24), hemoglobin (n=20) and rheumatoid factor (n=40) were tested; in no case was any disagreement found when the results were compared with those found with another commercial method.

14. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

704 human sera from healthy blood donors were analyzed. In parallel, these samples were tested for the presence of specific anti-CMV IgG using a commercial ELISA kit. 418 negative samples and 265 positive samples were found, in agreement between the two methods.

20 samples were not in agreement, and were subjected to further study using two other commercial methods based on different analytical principles from the immunoenzymatic test. The use of these alternative methods led to the conclusion that there were 2 false negative and 1 false positive result. Based on the same criteria, the comparative method also showed a lack of sensitivity with 2 samples. In conclusion, the data found in the experimentation showed that the sensitivity of the test is 99.3% and the specificity 99.8%.

15. REPRODUCIBILITY

Reproducibility was studied by testing negative and positive samples: All the runs (n=27) were valid. The CV% of 6.4 and 5.8% were found respectively with the negative and positive controls. Variations between runs were studied using a panel of serum samples which were positive, negative and borderline for the presence of anti-CMV antibodies. The samples were tested in 10 different runs. The CV% values ranged from 1.4 to 4.4%.

16. TROUBLE SHOOTING GUIDE

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert paragraph 4 for correct code).
		Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow).Repeat test
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.
	Inadequate washing	Ensure that wash apparatus works well
Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well

	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate Color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

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