ENZYWELL

HERPES SIMPLEX 1+2 IgG

REF 91020 (96 tests)

Prodotto da/Manufactured by/ Hersteller:
Κατασκευάζεται από/ Fabricado por/ Fabriqué par:
DIESSE Diagnostica Senese
Strada dei Laghi 39
53035 Monteriggioni (Siena) - Italy

INDICE / INDEX / INDICE / CONTEÚDO

INDICE / INDEX / INHALT/ ΠΕΡΙΕΧΟΜΕΝΑ / INDICE / SOMMAIRE /CONTEÚDO

1. UTILIZZAZIONE / INTENDED USE / VERWENDUNGSZWECK/ ΥΠΟΔΕΙΚΝΥΟΜΕΝΗ ΧΡΗΣΗ / INDICACIONES DE USO/ UTILISATION/UTILIZAÇÃO PRETENDIDA
2. INTRODUZIONE / SUMMARY AND EXPLANATION OF TEST / PHYSIOLOGIE/ ΕΙΣΑΓΩΓΗ / RESUMEN Y EXPLICACIÓN DEL TEST / INTERET CLINIQUE / SUMÁRIO E EXPLICAÇÃO DO TESTE
3. PRINCIPIO DEL METODO / PRINCIPLE OF THE TEST / TESTPRINZIP/ ΑΡΧΗ ΤΗΣ ΜΕΘΟΔΟΥ / PRINCIPIO DEL MÉTODO / PRINCIPE/ PRINCÍPIOS DO TESTE
4. COMPOSIZIONE DEL KIT E PREPARAZIONE DEI REAGENTI / KIT CONTENTS AND REAGENT PREPARATION / INHALT DES KITS UND VORBEREITUNG DER REAGENZIEN / ΣΥΝΘΕΣΗ ΤΟΥ ΚΙΤ ΚΑΙ ΠΡΟΕΤΟΙΜΑΣΙΑ ΤΩΝ ΑΝΤΙΑΠΑΣΤΗΡΙΩΝ / COMPONENTES DEL KIT Y PREPARACIÓN DEL REACTIVO / CONTENU DU KIT ET PREPARATION DES REACTIFS / CONTEÚDA DO KIT E PREPARAÇÃO DO REAGENTE /
5. MODALITA' DI CONSERVAZIONE E STABILITA' DEI REAGENTI / STORAGE AND STABILITY OF REAGENTS / LAGERUNG UND STABILITÄT DER REAGENZIEN / ΤΡΟΠΟΣ ΑΠΟΘΗΚΕΥΣΗΣ ΚΑΙ ΣΤΑΘΕΡΟΤΗΤΑ ΤΩΝ ΑΝΤΙΑΠΑΣΤΗΡΙΩΝ/ CONSERVACIÓN Y ESTABILIDAD DE LOS REACTIVOS / MODALITES DE CONSERVATION ET STABILITE DES REACTIFS / ARMAZENAMENTO E ESTABILIDADE
6. PRECAUZIONI / PRECAUTIONS / HINWEISE UND VORSICHTSMASSNAHMEN/ ΠΡΟΦΥΛΑΞΕΙΣ / PRECAUCIONES DE USO / PRECAUTIONS D'UTILISATION / PRECAUÇÕES
7. TIPO DI CAMPIONE E CONSERVAZIONE / TYPE AND STORAGE OF SAMPLE / ART UND LAGERUNG DER PROBEN / ΕΙΔΟΣ ΔΕΙΓΜΑΤΟΣ ΚΑΙ ΑΠΟΘΗΚΕΥΣΗ / TIPO DE MUESTRA Y CONSERVACION / TYPE D'ECHANTILLON ET CONSERVATION / TIPO E ARMAZENAMENTO DAS AMOSTRAS
8. PROCEDIMENTO / TEST PROCEDURE / TESTDURCHFÜHRUNG / ΔΙΑΔΙΚΑΣΙΑ / PROCEDIMIENTO DEL TEST / PROCEDURE / PROCEDIMENTO DE TESTE /
9. SCHEMA DEL SAGGIO / SCHEME OF TEST PROCEDURE / TESTSCHEMA /ΔΙΑΓΡΑΜΜΑ ΤΗΣ ΑΝΑΛΥΣΗΣ/ ESQUEMA DEL PROCEDIMIENTO DEL TEST / SCHEMA D'ANALYSE/ ESQUEMA DE PROCEDIMENTO DE TESTE /
10. VALIDAZIONE DEL TEST / TEST VALIDATION / TESTVALIDIERUNG / ΕΓΚΥΠΟΤΗΤΑ ΤΟΥ ΤΕΣΤ / VALIDACIÓN DEL TEST / VALIDATION DU TEST / VALIDAÇÃO DO TESTE /
11. INTERPRETAZIONE DEL TEST / INTERPRETATION OF RESULTS / INTERPRETATION DER ERGEBNISSE / ΕΡΜΗΝΕΙΑ ΤΩΝ ΑΠΟΤΕΛΕΣΜΑΤΩΝ / INTERPRETACIÓN DE RESULTADOS / INTERPRETATION DU TEST / INTERPRETAÇÃO DOS RESULTADOS
12. LIMITAZIONI DELLA PROCEDURA / LIMITATIONS OF THE PROCEDURE / GRENZEN DES VERFAHRENS / ΠΕΡΙΟΡΙΣΜΟΙ ΣΤΗΝ ΜΕΘΟΔΟ / LIMITACIONES DEL TEST / LIMITES DU TEST/ LIMITAÇÕES DO PROCEDIMENTO /
13. SPECIFICITA' ANALITICA / ANALYTICAL SPECIFICITY / ANALYTISCHE SPEZIFITÄT /ΑΝΑΛΥΤΙΚΗ ΕΙΔΙΚΟΤΗΤΑ/ ESPECIFICIDAD ANALÍTICA / SPECIFICITE ANALYTIQUE / ESPECIFICIDADE ANALÍTICA /
14. SENSIBILITA' E SPECIFICITA' DIAGNOSTICA / DIAGNOSTIC SENSITIVITY AND SPECIFICITY / DIAGNOSTISCHE SENSITIVITÄT UND SPEZIFITÄT /ΕΥΑΙΣΘΗΣΙΑ ΚΑΙ ΔΙΑΓΝΩΣΤΙΚΗ ΕΙΔΙΚΟΤΗΤΑ / SENSIBILIDAD Y ESPECIFICIDAD DIAGNOSTICAS / SENSIBLITE ET SPECIFICITE DIAGNOSTIQUES / ESPECIFICIDADE E SENSIBILIDADE DO DIAGNÓSTICO
15. PRECISIONE / PRECISION / PRÄZISION / ΑΚΡΙΒΕΙΑ / PRECISIÓN / PRECISION / PRECISÃO
16. GUIDA DEI PROBLEMI DI UTILIZZO / “TROUBLE SHOOTING” / “TROUBLE SHOOTING” (PROBLEMLÖSUNGEN) /ΟΔΗΓΟΣ ΕΠΙΛΥΣΗΣ ΠΡΟΒΛΗΜΑΤΩΝ /GUIA DE RESOLUCION DE PROBLEMAS / RESOLUTION DES PROBLEMES / RESOLUÇÃO DE PROBLEMAS
17. BIBLIOGRAFIA / REFERENCES / LITERATUR / ΒΙΒΛΙΟΓΡΑΦΙΑ / BIBLIOGRAFÍA / BIBLIOGRAPHIE / REFERÊNCIAS

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ENZYWELL HERPES SIMPLEX 1+2 IgG

REF 91020

(Italiano)

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17. BIBLIOGRAFIA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
3. B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
4. C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
5. M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
6. D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgM antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).
7. R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).
8. J.E. Kuhn et al.: Analysis of the IgM and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).

INSTRUCTIONS FOR USE

ENZYWELL HERPES SIMPLEX 1+2 IgG

REF 91020

(English)

1. INTENDED USE

IMMUNOENZYMATIC KIT FOR THE QUALITATIVE DETERMINATION OF IgG-CLASS ANTIBODIES TO HERPES SIMPLEX VIRUS (TYPES 1 AND 2) IN HUMAN SERUM, AS AN AID IN THE DIAGNOSIS OF HERPES SIMPLEX VIRUS INFECTION.

2. SUMMARY AND EXPLANATION OF THE TEST

The Herpes simplex virus (HSV) is a member of the Herpesviridae family, of which two types are known: type 1 (HSV-1) and type 2 (HSV-2) which present slight antigenic differences. HSV-1 causes chiefly oral-facial lesions, while HSV-2 is mainly responsible for genital lesions, but this distinction is not binding, both types occasionally causing infection in either anatomical site. HSV may also cause a form of ocular cheratitis, and lesions of the central nervous system.

HSV can affect practically the whole population. The primary infection is often in a subclinical form and is rarely diagnosed. After a latency period of variable duration, reactivation may occur and viral replication may or may not give rise to clinical lesions. Infection contracted during birth is of particular interest, this being an important cause of morbidity and mortality. It is therefore important to determine the immunitary state of women during pregnancy in order to detect serum conversion. The assay of specific IgM is important for the diagnosis of neonatal infection and encephalitis caused by HSV. Moreover, the presence of specific IgM indicates viral activity in progress, although it is not possible to distinguish between primary infection and reactivation.

3. PRINCIPLE OF THE TEST

The test is based on the ELISA technique (Enzyme linked Immunosorbent Assay) (1-8).

The antigen, composed of purified and inactivated Herpes Virus types 1 and 2, is bound to the solid phase (8-well strips). The specific immunoglobulins are bound to the antigen through incubation with dilute human serum.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of human IgG monoclonal antibodies labelled with peroxidase. The unbound conjugate is eliminated and the peroxidase substrate is added.

The colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

4. KIT CONTENTS AND REAGENT PREPARATION

- Reagents are sufficient for 96 determinations.

Bring to room temperature before use.

MT PLATE MICROPLATE 12x8 wells coated with herpes simplex virus (coloured in violet).

Use: open the package at the opposite end from the code (H followed by the lot number) which is useful for identification purposes, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

CONTROL + POSITIVE CONTROL 1 x 1.6 mL

Contents: Diluted human serum, at known concentration of anti-HSV IgG, in Phosphate buffer 0.01 mol/L containing BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour is proportional to the relative antibody titer.

CONTROL CUT-OFF CUT OFF CONTROL 1 x 2 mL

Contents: Diluted human serum, at known concentration of anti-HSV IgG, in Phosphate buffer 0.01 mol/L containing BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour is proportional to the relative antibody titer.

CONJ CONJUGATE 1x16 mL

Contents: monoclonal antibodies labelled with Peroxidase, in phosphate buffer solution containing phenol 0.05% and Bronidox 0.02%. Ready for use without further dilution.

CONTROL IgG - IgG NEGATIVE CONTROL (PF93910) 1 x 1.6 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Human serum in Phosphate buffer 0.01 mol/L, with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Stability: The product is stable up to the expiry date when stored unopened at 2/8°C.

WASH BUF 10x WASH BUFFER 10X (PF93603) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5% .

Preparation: dilute the required volume 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SAMP DIL 2 DILUENT 2 (PF93611) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

To be used to dilute samples.

Contents: Proteic solution in phosphate buffer with sodium azide 0,09% containing methyl orange as dye.

SUBS TMB SUBSTRATE (PF93619) 1x12 mL Ready for use. **INTERCHANGEABLE BETWEEN LOTS**

Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilised in citrate buffer 0.05 mol/L (pH 3.8).

H₂SO₄ 0.3 M STOP SOLUTION (PF93602) 1x16 mL **INTERCHANGEABLE BETWEEN LOTS**

H₂SO₄ 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

MATERIALS REQUIRED BUT NOT PROVIDED

- Incubator at 37-40°C
- Microplate reader, wavelength 450 or 450/620 nm, with OD linearity up to 2,000 (at least).
- Microplate washer (preferable) able to dispense volumes in the range 225-375 µl
- Distilled or deionized water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
Microplate	6 weeks at 2/8°C, polythene bag
Controls	6 weeks at 2/8°C
Conjugate	6 weeks at 2/8°C
Substrate	up to the expiry date at 2/8°C, 1 week at 15-30°C; store in the dark
Sample Diluent	up to the expiry date at 2/8°C
Wash Buffer	2 weeks at 2/8°C, 5 days at 15/30°C.
Stop Solution	up to the expiry date at 2/8°C

6. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

Attention:

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be

handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The Wash Buffer contains detergents
 - b) The conjugate contains phenol
 - c) The substrate is acid
 - d) Calibrators contain Sodium Azide 0.09% as preservative which may react with copper and lead in plumbing to form potentially explosive metal azides.
 If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
6. Do not allow wells to become dry during the assay procedure.
7. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
8. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
9. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO₂ incubators must not be used.
10. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
11. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.
12. Use of the kit with automatic instruments must be validated by the user.
13. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:
 - installation and particular requisites
 - operating principles, instructions, precautions and risks
 - manufacturer's specifications and instrument performance

-servicing and maintenance.

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be carefully mixed before the test. Heat inactivation can lead to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric or contaminated samples should be avoided.

The test cannot be performed on human plasma.

8. TEST PROCEDURE

Manual Technique

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Wash Buffer 10x (100 mL + 900 mL H₂O).

Dilute samples 1:101 distributing 10 µL of serum into 1 mL of diluent. Leave one well for the blank, performed using 100 µL of the substrate mixture. Dispense 100 µL of each diluted sample per well (duplicate testing is recommended). Place UNDILUTED controls in a strip (100 µL in each well). The minimum requisite is 1 negative control, 2 cut-off and 1 positive control.

Wells are covered with protective film and incubated for 45 minutes at 37°C. After washing four times for 30 seconds (300 µl), add 100 µL of conjugate to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film. The plate is washed again 4 times, as described above. Finally, the substrate is distributed, 100 µL/well. After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µL of Stop Solution.

The adsorbance (O.D.) is read at 450 nm or 450/620 nm within 30 min.

9. SCHEME OF TEST PROCEDURE

Manual Technique

STEP 1 Place 100 µL of diluted sample/controls to the wells of the strips.

-

Incubate for 45 min. at 37°C

-

Wash 4 times (300 µl)

STEP 2 Add 100 µL of conjugate to each well

-

Incubate for 45 min. at 37°C

-

Wash 4 times (300 µl)

STEP 3 Add 100 µL of Substrate to each well

-

Incubate for 15 min. at R.T.

STEP 4 Add 100 µL of Stop Solution

-

Read absorbance at 450 nm within 30 min

10. TEST VALIDATION

Subtract the value of the blank (≤ 0.150) from all the other readings. The O.D. values of the control Cut-off serum must be within 25% of the mean value if tested in triplicate. Disregard any abnormal value and recalculate the mean. The Positive control must have an O.D. at least 1.5 times that of the cut-off serum. The ratio between Negative Control and Cut-off must be ≤ 0.6 . The OD of the cut off must be > 0.2 at 450 nm and ≥ 0.16 at 450/620 nm.

11. INTERPRETATION OF THE RESULTS

Qualitative results

If the adsorbance of the sample is higher than that of the cut-off, the sample is positive for the presence of specific IgG. Calculate the ratio between the O.D. value of the sample and that of the Cut-off (INDEX).

The sample is considered:

Positive: if the ratio is > 1.2 .

Doubtful: $\pm 20\%$ of the cut-off.

Negative: if the ratio is < 0.8 .

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

12. LIMITATIONS OF THE PROCEDURE

A serum sample obtained during the acute phase of infection, when only IgM antibodies are present, may be negative by this procedure.

The Herpes Simplex Virus IgM level should be determined using Enzywell Herpes IgM kit. Alternatively, a second serum sample obtained 8-14 days later, should be tested in parallel to determine an increase in the IgG antibody level.

The test result should be used in conjunction with information available from the evaluation of other diagnostic procedures.

13. ANALYTICAL SPECIFICITY

17 samples which were negative for Herpes Simplex, but which contained IgG antibodies against other viruses such as Rubella, Epstein Barr, Cytomegalovirus or Mumps, were tested. None of these antibodies influenced the test.

14. SENSITIVITY AND SPECIFICITY

In a clinical trial performed in a hospital laboratory, 348 samples were analysed, 61 of which proved negative and 287 positive. The samples were analysed with another commercial immunoenzymatic method: 100% agreement was found between the two methods.

The Enzywell kit offers 100% sensitivity and specificity.

15. PRECISION

In-run Precision

<i>Cut Off</i> $n=15$	<i>Lot N. 084</i>	<i>Lot N. 085</i>	<i>Lot N. 086</i>
O.D.	0.446	0.388	0.598
CV%	6	10	6

Between-run Precision

	INDEX						
Sample	I run	II run	III run	IV Run	V Run	Average	CV%
HSG 1	0.4	0.5	0.5	0.4	0.4	0.4	12
HSG 2	1.9	2.4	2.2	2.0	1.7	2.0	13
HSG 3	3.7	4.9	5.5	4.7	3.4	4.4	20
Pos. Control	6.3	6.7	7.9	8.0	5.8	6.9	14

16. TROUBLE SHOOTING GUIDE

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert point 4 for correct code).
		Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow). Repeat test
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.
	Inadequate washing	Ensure that wash apparatus works well
Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well

	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate Color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

17. REFERENCES

- 1.G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
- 2S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
- 3.B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
- 4.C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
- 5.M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
- 6.D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgM antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).
- 7.R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).
- 8.J.E. Kuhn et al.: Analysis of the IgM and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).

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17. LITERATUR

9. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
10. S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
11. B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
12. C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
13. M. Morgan and T. : Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
14. D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgM antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).
15. R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).
16. J.E. Kuhn et al.: Analysis of the IgM and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).

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- 1.G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
- 2.S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
- 3.B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
- 4.C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
- 5.M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
- 6.D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgM antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).
- 7.R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).

8.J.E. Kuhn et al.: Analysis of the IgM and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).

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- 1.G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
- 2.S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
- 3.B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
- 4.C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
- 5.M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
- 6.D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgM antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).
- 7.R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).
- 8.J.E. Kuhn et al.: Analysis of the IgM and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).

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17. BIBLIOGRAPHIE

17. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
18. S. Land et él.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
19. B. Gonik et él.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
20. C. Gleaves et él.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
21. M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
22. D. Ho et él.: Indirect ELISA for the detection of HSV-2 specific IgG and IgM antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).
23. R. Eberle et él.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).
24. J.E. Kuhn et él.: Analysis of the IgM and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).

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