CHORUS

snRNP-C

REF 86092

REF 86092/12

Prodotto da/
Manufactured by/
Fabricado por:

DIESSE

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INSTRUCTIONS FOR USE

CHORUS snRNP-C

For the semiquantitative determination of anti-snRNP-C antibodies

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the semiquantitative determination of anti- snRNP complex antibodies in human serum, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

The U1-snRNP complex (for short snRNP) is a small nuclear ribonucleoprotein particle composed of uridine rich (thus U) small nuclear RNA and a set of proteins, the 70 kDa U1-specific protein plus protein A and C (all formerly summarized as RNPs) and the Sm (Smith) antigen. Because of its components Sm and RNPs the complex has been often called RNP/Sm complex. U1-snRNP is a part of the splicosomal complex, facilitating the processing of pre-mRNA to mature mRNA in the nucleus.

Antibodies against snRNP belong to the heterogeneous group of anti-nuclear antibodies (ANA), which are associated with various autoimmune diseases. They are directed against proteins of the nucleus.

Indirect immunofluorescence test (IFT) on eukaryotic cells has been the established method for the detection of ANAs. Single antibody specificities are distinguished by fluorescence patterns but more specific testing by ELISA employing the target antigen has been established for simple and reliable differentiation of ANAs.

Antibodies against the snRNP complex are directed against Sm as well as the 70 kDa U1-specific protein plus protein A and C. They typically occur in SLE and mixed connective tissue diseases (MCTD), but are also found in Sjögren's syndrome, scleroderma and polymyositis.

The differentiation of the autoantibodies by ELISA employing the specific antigen is an aid for the diagnosis of SLE and MCTD. Antibodies against the 70 kDa U1 RNP protein are found in 95% of MCTD and in 40% of SLE. The isolated occurrence of anti-U1 RNP 70 kDa antibodies is typical for the Sharp syndrome.

In contrast antibodies against Sm are highly specific for SLE and thus are included in diagnostic and classification criteria for SLE. 20-30% of patients with SLE display antibodies against Sm.

3. PRINCIPLE OF THE METHOD

The snRNP-C devices are ready to use for the assay of antibodies against U1-sn-RNP in the CHORUS instruments. The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay).

The antigen, highly purified, is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with diluted human serum.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of monoclonal antibodies, highly specific for anti-human immunoglobulins conjugated to horse radish peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added.

The blue colour which develops is proportional to the concentration of specific antibodies present in the serum sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus instruments.

The Chorus snRNP-C is calibrated against reference sera from the CDC (Atlanta). The results are expressed in AU/ml.

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HBsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the CHORUS instrument.
4. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The conjugate contains phenol
 - b) The substrate is acid

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.

- Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

- Discard devices which show the substrate (well 4) blue colored.**
- Adding the sample into the well verify that it is perfectly distributed on the bottom.
- Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
- The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
- Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
- Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
- Avoid using self-defrosting freezers for the storage of the samples.
- Defective barcodes can be inserted manually in the instrument.
- Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
- The use of strongly hemolyzed samples or samples presenting microbial contamination may all constitute sources of error.
- Before inserting the devices in the instrument, check that the reaction well does not contain foreign bodies.
- Pipette the test serum (50 µl) in well 1 of the device (see figure).
- Do not use the device after the expiry date.
- Make sure that the instrument is connected to the Washing Buffer Autoimmunity REF 86004.**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 86092).


The kit is sufficient for 12 tests (REF 86092/12).

DD DEVICES

6 packages each containing 6 devices (REF 86092).

2 packages each containing 6 devices (REF 86092/12).

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with highly purified snRNP-C

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE 0.35 ml.

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT 0.35 ml

Contents: Proteic solution containing Tween-20 0.2% and Proclin 0.1%.

Position 2: CONJUGATE 0.35 ml

Contents: anti-human IgG monoclonal antibodies labeled with horse radish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

In which the operator must place the undiluted serum

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.175 ml

Contents: Diluted human serum containing IgG antibodies anti-snRNP-C and preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.425 ml

Contents: Diluted human serum containing IgG antibodies anti-snRNP-C and preservative. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY (REF 86004)
- CLEANING SOLUTION 2000 (REF 83609)
- SANITIZING SOLUTION (REF 83604 - 83608)
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT (REF 83607)
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the control serum (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C

POSITIVE CONTROL 8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use. Heat-inactivation can rise to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results. Strongly lipemic, icteric, hemolyzed or contaminated samples cannot be used.

The test cannot be applied to plasma.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 50 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the Chorus instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the control serum to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the control serum has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected. If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
Fax: 0039 0577 366605
email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus instrument expresses the result in AU/ml calculated on the basis of a batch-dependent curve stored in the instrument.

The test serum can be interpreted as follows:

POSITIVE: when the result is > 18.0 AU/ml

NEGATIVE: when the result is < 12.0 AU/ml

DOUBTFUL/EQUIVOCAL: for all values between 12.0 and 18.0 AU/ml

In the case of a doubtful result, repeat the test. If the test remains doubtful, collect a new sample.

11. LIMITATIONS

All the positive results require a careful interpretation.

The test cannot be used as only method for a clinical diagnosis.

Negative results may not exclude an eventual infection.

The test results should be interpreted in conjunction with information available from the clinical evaluation and other diagnostic procedures.

12. CALIBRATION RANGE

Calibration range 3.0-100.0 AU/ml.

For samples > 100.0 AU/ml repeat the test pre-diluting the sample with Negative Control/Sample Diluent (PF83607 – **not supplied with the kit**).

13. REFERENCE RANGE

Among the normal population the expected values, which have been determined by examining 120 sera from healthy donors, were between 3.0 and 11.8 AU/ml.

14. ANALYTICAL SPECIFICITY

5 samples (2 negative, 1 Cut-Off and 2 Positive) containing the following interfering substances were tested:

Rheumatoid Factor (44 – 220 UI/ml)
Bilirubin (4.5 mg/dl – 45 mg/dl)
Triglycerides (10 mg/dl – 250 mg/dl)
Hemoglobin (5 mg/ml – 30 mg/ml)

The presence in the tested serum of the interfering substances listed above did not alter the results of the test.

15. CROSS-REACTIONS

Samples positive to PR-3, MPO, Cardiolipin, Gliadin, AMA-M2, Transglutaminase, SS-A, SS-B, Sm, Scl-70, Cenp-B, Jo-1, dsDNA IgG and ASCA were tested.

No significant cross-reactions were found.

16. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

In an experimentation 114 samples were tested with the Diesse kit and with another commercial kit.

The results are summarized in the following table:

		Reference		
		+	-	Total
Diesse	+	21	1	22
	-	2	90	92
	Total	23	91	114

Percent Positive Agreement (~Diagnostic Sensitivity):

91.3% CI_{95%}: 73.2- 97.4

Percent Negative Agreement: (~Diagnostic Specificity):

98.9% CI_{95%}: 94.0- 99.8

17. PRECISION

Sample	Within-run Precision		Between-run precision	
	Mean (AU/ml)	CV%	Mean (AU/ml)	CV%
1	5.5	9.5	7.5	9.3
2	7.1	8.3	7.4	9.2
3	13.4	7.1	12.9	8.1
4	13.9	8.3	14.3	9.4
5	21.0	8.5	19.0	4.8
6	31.2	9.6	42.2	8.1
7	54.3	5.5	48.1	8.4
8	94.4	7.5	93.3	8.6

Sample	Precision between batches		Precision between instruments	
	Mean (AU/ml)	CV%	Mean (AU/ml)	CV%
1	5.1	4.9	5.1	2.0
2	7.7	7.7	7.6	5.0
3	13.1	5.8	13.1	5.1
4	12.9	7.1	12.9	5.4
5	23.8	6.1	23.7	6.3
6	40.2	6.0	40.2	4.7
7	50.4	6.6	50.4	4.1
8	85.1	4.6	85.1	3.5

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1. Peter JB. et al. (1996) Sciences B.V., Amsterdam.
DIESSE Diagnostica Senese S.p.A.