CHORUS

**Legionella Urinary
Antigen**

DIESSE Diagnostica Senese
S.p.A.
Via delle Rose, 10
53035 Monteriggioni (Siena)
Italy

REF 81301

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CHORUS Legionella Urinary Antigen

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CHORUS Legionella Urinary Antigen

For the qualitative determination of Legionella pneumophila urinary antigen

For *In Vitro* Diagnostic Use Only

1. INTENDED USE

Immunoenzymatic method for the qualitative determination of *Legionella pneumophila* urinary antigen in human urine, using a disposable device applied on the Chorus and Chorus TRIO instruments.

2. INTRODUCTION

Legionnaire's disease is a serious infectious disease with a high mortality rate.

The species most frequently involved in human pathologies is *Legionella pneumophila*, although other species have been isolated. The most serious form of *Legionella* infection is the so-called Legionnaire's disease, characterized by acute pneumonia which becomes evident 2-10 days after contact with the bacteria. Another characteristic clinical picture is Pontiac fever, an acute extra-pulmonary infection.

The infection can affect the whole population. The subjects considered to be most at risk are elderly males who are smokers, consumers of alcohol, affected by chronic diseases (obstructive bronchopulmonary infections, cardiovascular and renal disorders, diabetes, etc.), and those with acquired immunodeficiency following therapy (organ transplant, treatment with steroids and anti-tumoral drugs, etc.), or affected by HIV.

80% of all the cases of legionellosis reported can be attributed to *Legionella pneumophila* serotype 1.

In 1979 Berdal demonstrated the presence of a highly specific soluble antigen in the urine of patients affected by Legionnaire's disease. Therefore the urine of these patients represents the ideal sample for the detection of the antigen in the precocious or late stages of the disease.

3. PRINCIPLE OF THE METHOD

The Chorus Legionella Urinary Antigen device is ready to use for the detection of *Legionella pneumophila* urinary antigen, in the Chorus/ Chorus TRIO instruments.

The test is based on the ELISA principle (Enzyme linked ImmunoSorbent Assay). Monoclonal antibodies, highly specific for the antigen, are bound to the solid phase.

The antigen is bound to the well through incubation with the human urine.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed

of anti-*Legionella pneumophila* monoclonal antibodies conjugated to horse radish peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added. The blue colour which develops is proportional to the concentration of specific antibodies present in the sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus/Chorus TRIO instruments.

The results are expressed in Index (OD sample/OD cut-off).

4. WARNINGS AND PRECAUTIONS

FOR *IN VITRO* DIAGNOSTIC USE ONLY

As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth.
2. Wear disposable gloves and eye protection while handling specimens.
3. Wash hands thoroughly after placing the devices in the Chorus/Chorus TRIO instrument.
4. Consult the relative Material Safety Data Sheet (available on request) for all the information on safety concerning the reagents contained in the kit.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical Precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. Adding the sample into the well verify that it is perfectly distributed on the bottom.


3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent and/or present foreign bodies in the reaction well when visually inspected.
4. The devices are for use with the Chorus/Chorus TRIO instrument; the Instructions for Use must be carefully followed and the Instrument Operating Manual must be consulted.
5. Check that the Chorus/Chorus TRIO instrument is set up correctly (see Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument.
7. Avoid using self-defrosting freezers for the storage of the samples.
8. Defective barcodes can be inserted manually in the instrument (see Operating Manual).
9. Do not expose the devices to strong light or to hypochlorite vapors during storage and use.
10. Do not use the device after the expiry date.
11. **Make sure that the instrument is connected to the Washing Buffer (Ref. 83606).**

5. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 12 tests.

DD DEVICES 2 packages each containing 6 devices

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL

Coated with highly purified anti-*Legionella pneumophila* monoclonal antibodies

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE

Contents: Tetramethylbenzidine and H₂O₂ stabilized in buffer.

Position 3: MONOCLONAL ANTIBODIES

Contents: anti-*Legionella pneumophila* antibodies labelled with Biotin in TRIS buffer.

Position 2: CONJUGATE

Contents: Streptavidin solution labelled with horseradish peroxidase, in phosphate buffer containing phenol 0.05% and Bronidox 0.02%.

Position 1: EMPTY WELL

In which the operator must place the undiluted sample

Use: equilibrate a package at room temperature, open the package and remove the required devices; replace the others in the bag with the silica gel, expel the air and seal by pressing the closure. Store at 2-8°C.

CALIBRATOR CALIBRATOR 1 x 0.350 ml

Contents: Inactivated *Legionella pneumophila* cells in phosphate buffer and preservative. Liquid, ready for use.

CONTROL + POSITIVE CONTROL 1 x 0.850 ml

Contents: Inactivated *Legionella pneumophila* cells in phosphate buffer and preservative. Liquid, ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER **REF** 83606
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- Chorus/Chorus TRIO Instrument
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

6. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the positive control (see section 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening:

DEVICES	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C

7. SPECIMEN COLLECTION AND STORAGE

The sample is composed of human urine to be handled with all precautions dictated by good laboratory practice.

The sample must be collected in a common sterile container and can be maintained for 24 hours at 2/8°C, or frozen for longer periods at -20°C.

The sample be thawed a maximum of 3 times.

Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use.

The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

8. ASSAY PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in chapter 4, Analytical Precautions.
3. Dispense 125 µl of undiluted sample in well no. 1 of each device; at each change of batch, use a device for the calibrator.

4. Place the devices in the Chorus/Chorus TRIO instrument. Perform the calibration (if necessary) and the test as reported in the Instrument Operating Manual.

9. TEST VALIDATION

Use the positive control to check the validity of the results obtained. It should be used as reported in the Instrument Operating Manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the positive control continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus/Chorus TRIO instrument expresses the result in Index (OD sample/OD cut-off).

The test on the examined sample can be interpreted as follows:

POSITIVE: when the result is > 1.0

NEGATIVE: when the result is < 1.0

DOUBTFUL/EQUIVOCAL: when the result is $= 1.0$

If the result is doubtful/equivocal, repeat the test. If it remains doubtful/equivocal, collect a new sample.

11. LIMITATIONS

All the values obtained require a careful interpretation that must consider other indicators relative to the patient.

The test, indeed, can not be used alone for a clinical diagnosis and the test result should be evaluated together with the patient history and other clinical diagnostic evaluation.

12. ANALYTICAL SPECIFICITY

5 samples (2 Negative, 1 Cut-Off and 2 Positive) were spiked with the following potentially interfering factors and then tested:

Bilirubin (0.15 mg/dl – 0.6 mg/dl)
 Hemoglobin (0.01 mg/ml – 0.04 mg/ml)
 Glucose (30 mg/dl – 120 mg/dl)
 Ascorbic Acid (5 mg/dl – 20 mg/dl)
 Calcium Oxalate (2.5 mg/dl – 10 mg/dl)
 Albumin (40 mg/dl – 160 mg/dl)

The presence in the sample of the interfering substances described above does not affect the test result.

13. CROSS-REACTIONS

20 samples, positive to *Escherichia coli*, *Proteus*, *Streptococcus faecalis*, *Staphylococcus* spp, *Streptococcus B*, *Streptococcus D*, *Lactobacillus*, KES bacteria, *Pseudomonas*, *Streptococcus pneumoniae* were tested.

No significant cross-reactions were found.

14. METHOD COMPARISON

In an experimentation 104 samples have been tested with Diesse kit and with a competitor kit.

Data are summarized in the following table :

		Reference		
		+	-	Total
Diesse	+	44	0	44
	-	8	52	60
	Total	52	52	104

Percent Positive Agreement (-Diagnostic Sensitivity):

84.6% CI_{95%}: 72.4 - 91.9

Percent Negative Agreement: (-Diagnostic Specificity): 100%

CI_{95%}: 93.1 – 99.9

The agreement between the two methods is excellent with a Cohen's Kappa of 0.84.

15. PRECISION AND REPEATABILITY

Sample	Within-run Precision		Between-run precision	
	Mean (Index)	CV%	Mean (Index)	CV%
1	0.4	20.0	0.3	16.7
2	0.3	20.0	0.3	16.7
3	0.6	6.7	0.6	13.3
4	0.8	10.0	0.6	11.7
5	1.2	9.2	1.1	10.0
6	2.1	2.4	2.1	5.2
7	8.3	8.8	7.7	6.4
8	4.0	3.3	3.8	5.0
9	4.4	4.3	4.3	3.7

Sample	Precision between batches		Precision between instruments	
	Mean (Index)	CV%	Mean (Index)	CV%
1	0.4	15.0	0.4	15.0
2	0.3	-	0.3	-
3	0.5	-	0.5	-
4	0.6	10.0	0.6	10.0
5	1.2	5.0	1.2	5.0
6	2.4	9.6	2.3	6.5
7	9.3	8.2	9.4	3.4
8	3.4	14.4	3.4	1.8
9	4.4	10.2	4.4	7.0


16. REFERENCES

1. Fraser DW, Tsai TR, Orenstein W, Parkin WE, Beecham HJ, Sharrar RG, Harris J, Mallison GF, Martin SM, McDade JE, Shepard CC, Brachman PS. 1977. Legionnaires' disease: description of an epidemic of pneumonia. *N. Engl. J. Med.* 297:1189-1197.
2. Reingold AL, Thomason BM, Brake BJ, Thacker L, Wilkinson HW, Kuritsky JN. 1984. Legionella pneumonia in the United States : the distribution of serogroups and species causing human illness. *J Infect Dis.* 149(5):819.
3. Berdal BP, Farshy CE, Feeley JC. 1979. Detection of Legionella pneumophila antigen in urine by enzyme-linked immunospecific assay. *J Clin Microbiol.* May; 9(5) :575-8.

DIESSE Diagnostica Senese S.p.A.
Via delle Rose 10
53035 Monteriggioni (Siena)
Italy

	EN Date of manufacture ES Fecha de fabricación IT Data di fabbricazione	FR Date de fabrication GR Ημερομηνία Παραγωγής PT Data de fabrico
	EN Use By ES Fecha de caducidad IT Utilizzare entro	FR Utiliser jusque GR Ημερομηνία λήξης PT Prazo de validade
	EN Do not reuse ES No reutilizar IT Non riutilizzare	FR Ne pas réutiliser GR Μην κάνετε επαναληπτική χρήση PT Não reutilizar
	EN Caution, consult accompanying documents ES Atención, ver instrucciones de uso IT Attenzione, vedere le istruzioni per l'uso	FR Attention voir notice d'instructions GR Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα PT Atenção, consulte a documentação incluída
	EN Manufacturer ES Fabricante IT Fabricante	FR Fabricant GR Κατασκευαστής PT Fabricante
	EN Contains sufficient for <n> tests ES Contenido suficiente para <n> ensayos IT Contenuto sufficiente per "n" saggi	FR Contenu suffisant pour "n" tests GR Περιεχόμενο επαρκές για «n» εξετάσεις PT Conteúdo suficiente para "n" ensaios
	EN Temperature limitation ES Límite de temperatura IT Limiti di temperatura	FR Limites de température GR Περιορισμοί θερμοκρασίας PT Limites de temperatura
	EN Consult Instructions for Use ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso	FR Consulter les instructions d'utilisation GR Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização
	EN Biological risks ES Riesgo biológico IT Rischio biologico	FR Risques biologiques GR Βιολογικοί κίνδυνοι PT Risco biológico
	EN Catalogue number ES Número de catálogo IT Numero di catalogo	FR Référence du catalogue GR Αριθμός καταλόγου PT Referência de catálogo
	EN In Vitro Diagnostic Medical Device ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro	FR Dispositif médical de diagnostic in vitro GR In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro
	EN Batch code ES Código de lote IT Codice del lotto	FR Code du lot GR Αριθμός Παρτίδας PT Código do lote