ENZYWELL

VARICELLA IgM

REF 91079 (48 tests)

REF 91089 (96 tests)

Prodotto da/Manufactured by/ Hersteller: Κατασκευάζεται
από/ Fabricado por/ Fabriqué par:
DIESSE Diagnostica Senese
Strada dei Laghi 39
53035 Monteriggioni (SIENA) - Italy

INDICE / INDEX / INHALT/ ΠΕΡΙΕΧΟΜΕΝΑ / INDICE / SOMMAIRE /CONTEÚDO

1. UTILIZZAZIONE / INTENDED USE / VERWENDUNGSZWECK/ ΥΠΟΔΕΙΚΝΥΟΜΕΝΗ ΧΡΗΣΗ / INDICACIONES DE USO/ UTILISATION/UTILIZAÇÃO PRETENDIDA
2. INTRODUZIONE / SUMMARY AND EXPLANATION OF TEST / PHYSIOLOGIE/ ΕΙΣΑΓΩΓΗ / RESUMEN Y EXPLICACIÓN DEL TEST / INTERET CLINIQUE / SUMÁRIO E EXPLICAÇÃO DO TESTE
3. PRINCIPIO DEL METODO / PRINCIPLE OF THE TEST / TESTPRINZIP/ ΑΡΧΗ ΤΗΣ ΜΕΘΟΔΟΥ / PRINCIPIO DEL MÉTODO / PRINCIPE/ PRINCÍPIOS DO TESTE
4. COMPOSIZIONE DEL KIT E PREPARAZIONE DEI REAGENTI / KIT CONTENTS AND REAGENT PREPARATION / INHALT DES KITS UND VORBEREITUNG DER REAGENZIEN / ΣΥΝΘΕΣΗ ΤΟΥ ΚΙΤ ΚΑΙ ΠΡΟΕΤΟΙΜΑΣΙΑ ΤΩΝ ΑΝΤΙΔΡΑΣΤΗΡΙΩΝ / COMPONENTES DEL KIT Y PREPARACIÓN DEL REACTIVO / CONTENU DU KIT ET PREPARATION DES REACTIFS / CONTEÚDA DO KIT E PREPARAÇÃO DO REAGENTE /
5. MODALITA' DI CONSERVAZIONE E STABILITA' DEI REAGENTI / STORAGE AND STABILITY OF REAGENTS / LAGERUNG UND STABILITÄT DER REAGENZIEN / ΤΡΟΠΟΣ ΑΠΟΘΗΚΕΥΣΗΣ ΚΑΙ ΣΤΑΘΕΡΟΤΗΤΑ ΤΩΝ ΑΝΤΙΔΡΑΣΤΗΡΙΩΝ/ CONSERVACIÓN Y ESTABILIDAD DE LOS REACTIVOS / MODALITES DE CONSERVATION ET STABILITE DES REACTIFS / ARMAZENAMENTO E ESTABILIDADE
6. PRECAUZIONI / PRECAUTIONS / HINWEISE UND VORSICHTSMASSNAHMEN/ ΠΡΟΦΥΛΑΞΕΙΣ / PRECAUCIONES DE USO / PRECAUTIONS D'UTILISATION / PRECAUÇÕES
7. TIPO DI CAMPIONE E CONSERVAZIONE / TYPE AND STORAGE OF SAMPLE / ART UND LAGERUNG DER PROBEN / ΕΙΔΟΣ ΔΕΙΓΜΑΤΟΣ ΚΑΙ ΑΠΟΘΗΚΕΥΣΗ / TIPO DE MUESTRA Y CONSERVACION / TYPE D'ECHANTILLON ET CONSERVATION / TIPO E ARMAZENAMENTO DAS AMOSTRAS
8. PROCEDIMENTO / TEST PROCEDURE / TESTDURCHFÜHRUNG / ΔΙΑΔΙΚΑΣΙΑ / PROCEDIMIENTO DEL TEST / PROCEDURE / PROCEDIMENTO DE TESTE /
9. SCHEMA DEL SAGGIO / SCHEME OF TEST PROCEDURE / TESTSCHEMA /ΔΙΑΓΡΑΜΜΑ ΤΗΣ ΑΝΑΛΥΣΗΣ/ ESQUEMA DEL PROCEDIMIENTO DEL TEST / SCHEMA D'ANALYSE/ ESQUEMA DE PROCEDIMENTO DE TESTE /
10. VALIDAZIONE DEL TEST / TEST VALIDATION / TESTVALIDIERUNG / ΕΓΚΥΠΟΤΗΤΑ ΤΟΥ TEST VALIDACIÓN DEL TEST / VALIDATION DU TEST / VALIDAÇÃO DO TESTE /
11. INTERPRETAZIONE DEL TEST / INTERPRETATION OF RESULTS / INTERPRETATION DER ERGEBNISSE / ΕΡΜΗΝΕΙΑ ΤΩΝ ΑΠΟΤΕΛΕΣΜΑΤΩΝ / INTERPRETACIÓN DE RESULTADOS / INTERPRETATION DU TEST / INTERPRETAÇÃO DOS RESULTADOS
12. LIMITAZIONI DELLA PROCEDURA / LIMITATIONS OF THE PROCEDURE / GRENZEN DES VERFAHRENS / ΠΕΡΙΟΡΙΣΜΟΙ ΣΤΗΝ ΜΕΘΟΔΟ / LIMITACIONES DEL TEST / LIMITES DU TEST/ LIMITAÇÕES DO PROCEDIMENTO /
13. SPECIFICITA' ANALITICA / ANALYTICAL SPECIFICITY / ANALYTISCHE SPEZIFITÄT /ΑΝΑΛΥΤΙΚΗ ΕΙΔΙΚΟΤΗΤΑ/ ESPECIFICIDAD ANALÍTICA / SPECIFICITE ANALYTIQUE / ESPECIFICIDADE ANALÍTICA /
14. SENSIBILITA' E SPECIFICITA' DIAGNOSTICA / DIAGNOSTIC SENSITIVITY AND SPECIFICITY / DIAGNOSTISCHE SENSITIVITÄT UND SPEZIFITÄT /ΕΥΑΙΣΘΗΣΙΑ ΚΑΙ ΔΙΑΓΝΩΣΤΙΚΗ ΕΙΔΙΚΟΤΗΤΑ / SENSIBILIDAD Y ESPECIFICIDAD DIAGNOSTICAS / SENSIBLITE ET SPECIFICITE DIAGNOSTIQUES / ESPECIFICIDADE E SENSIBILIDADE DO DIAGNÓSTICO
15. PRECISIONE / PRECISION / PRÄZISION / ΑΚΡΙΒΕΙΑ / PRECISIÓN / PRECISION / PRECISÃO
16. GUIDA DEI PROBLEMI DI UTILIZZO / “TROUBLE SHOOTING” / “TROUBLE SHOOTING” (PROBLEMLÖSUNGEN) /ΟΔΗΓΟΣ ΕΠΙΛΥΣΗΣ ΠΡΟΒΛΗΜΑΤΩΝ /GUIA DE RESOLUCION DE PROBLEMAS / RESOLUTION DES PROBLEMES / RESOLUÇÃO DE PROBLEMAS
17. BIBLIOGRAFIA / REFERENCES / LITERATUR / ΒΙΒΛΙΟΓΡΑΦΙΑ / BIBLIOGRAFÍA / BIBLIOGRAPHIE / REFERÊNCIAS

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ENZYWELL VARICELLA ZOSTER IgM

REF 91079 (48 tests) **REF** 91089 (96 tests)

(Italiano)

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REF 91079**(48 tests)****REF 91089****(96 tests)****(English)**

1. INTENDED USE

IMMUNOENZYMATIC CAPTURE METHOD FOR THE QUALITATIVE DETERMINATION OF IgM-CLASS ANTIBODIES TO VARICELLA ZOSTER IN HUMAN SERUM, AS AN AID IN THE DIAGNOSIS OF VARICELLA VIRUS INFECTION.

2. SUMMARY AND EXPLANATION OF THE TEST

Varicella and Herpes Zoster are two clinical manifestations resulting from infection by Varicella Zoster Virus (VZV).

Varicella, or chicken pox, is a highly contagious disease which is generally a consequence of primary infection by VZV and normally affects children. Infection caused by VZV during pregnancy can cause disease or malformation of the fetus; if it occurs at the end of pregnancy it may lead to the death of the neonate.

Herpes Zoster is a disease which chiefly affects adults and appears to be caused by a reactivation of the virus, which can remain latent for long periods in the spinal sensorial ganglia following primary infection. The infection causes painful cutaneous eruptions along the path of the nerves affected.

Serological methods are usually adopted to determine the immune state of subjects at risk (chiefly immunodepressed patients) and in pre- and post-natal diagnosis of infected subjects.

3. PRINCIPLE OF THE TEST

The test for the assay of anti-Varicella Zoster IgM is based on the principle of the capture of these immunoglobulins by anti-human IgM monoclonal antibodies found on the solid phase. A subsequent incubation with varicella antigen in a complex with monoclonal antibodies conjugated to horse radish peroxidase selects the IgM antibodies specific for the antigen and is revealed by the addition of the peroxidase substrate. When the enzymatic reaction is stopped by the addition of a sulphuric acid solution, a yellow colouring forms. The colour, which is proportional to the amount of specific antibodies present in the sample, can be read in an ELISA microplate reader.

4. KIT CONTENTS AND REAGENT PREPARATION

- Reagents are sufficient for 48 determinations (REF 91079)
- Reagents are sufficient for 96 determinations (REF 91089)
- **Bring reagents to room temperature before use.**

MT PLATE MICROPLATE 3 x 2 strips (REF 91079); 6 x 2 strips (REF 91089) coated with anti-human IgM monoclonal antibodies.

Use: open the package at the opposite end from the code (M followed by lot number) which is useful for identification purposes, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

CONTROL + POSITIVE CONTROL 1 x 1.6 mL

Contents: Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour of the control sera is proportional to the relative antibody titer.

CONTROL CUT-OFF CUT-OFF CONTROL 1 x 2.5 mL

Contents: Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour of the control sera is proportional to the relative antibody titer.

CONJ CONJUGATE 1 x 10 mL (REF 91079); 1 x 18 mL (REF 91089)

Contents: monoclonal antibodies labelled with peroxidase, in phosphate buffer with mouse ascitic fluid, phenol 0.05% and Bronidox 0.02%.

Preparation: Ready for use without further dilution. The immunocomplex must be prepared 45 minutes before use.

Ag ANTIGEN Freeze-dried x 3 vials (REF 91079), x 6 vials (REF 91089)

Contents: Partially purified Varicella Zoster virus inactivated by treatment with Beta-propiolactone, in phosphate buffer 0.04 mol/L and lactose, pH 7.2.

Preparation: reconstitute with the conjugate volume shown on the label, mixing by inversion.

CONTROL IgM - IgM NEGATIVE CONTROL (PF93900) 1 x 1.6 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

WASH BUF 10x WASH BUFFER 10X (PF93603) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5% .

Preparation: dilute the required volume 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SAMP DIL 2 DILUENT 2 (PF93611) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

To be used to dilute samples.

Contents: Proteic solution in phosphate buffer with sodium azide 0,09% containing methyl orange as dye.

SUBS TMB SUBSTRATE (PF93619) 1 x 12 mL Ready for use **INTERCHANGEABLE BETWEEN LOTS**

Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilised in citrate buffer 0.05 mol/L (pH 3.8).

H₂SO₄ 0.3 M STOP SOLUTION (PF93602) 1x16 mL **INTERCHANGEABLE BETWEEN LOTS**

H₂SO₄ 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

MATERIALS REQUIRED BUT NOT PROVIDED.

- Incubator at 37°C
- Microplate reader (wave length 450 nm or 450/620 nm, with linearity up to OD >= 2000)
- Microplate washer (preferable) able to dispense volumes in the range 225-375 µL
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µL solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
Microplate	5 weeks at 2/8°C, polythene bag
Control sera	5 weeks at 2/8°C
Conjugate	5 weeks at 2/8°C
Reconstituted antigen	5 days at 2/8°C if reconstituted with conjugate; (-20°C if reconstituted with Wash Buffer. Avoid repeated freezing/thawing. See “Analytical Precautions” no. 1).
Substrate	up to the expiry date at 2/8°C, 1 week at 15-30°C; store in the dark
Sample Diluent	up to the expiry date at 2/8°C
Wash Buffer	ready for use 2 weeks at 2/8°C, 5 days at 15/30°C.
Stop Solution	up to the expiry date at 2/8°C

6. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

Caution:

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can

offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The Wash Buffer contains detergents
 - b) The conjugate contains phenol
 - c) The substrate is acid
 - d) The controls contain 0.9% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. **The antigen reconstituted with conjugate is not stable after freezing. In the case of a limited consumption of antigen, proceed as follows: Reconstitute the antigen in 1/10 of the volume reported on the label with Wash Buffer ready for use (eg. volume reported on the label 3 ml: reconstitute with 0.3 ml of Wash Buffer). Take the amount of antigen necessary for immediate use and mix with 10 parts of conjugate. Aliquot and freeze the remaining antigen. At the time of use, thaw and mix with 10 parts of conjugate.**
2. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
3. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
4. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
5. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
6. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
7. Do not allow wells to become dry during the assay procedure.
8. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
9. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
10. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO₂ incubators must not be used.
11. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
12. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.
13. Use of the kit with automatic instruments must be validated by the user.
14. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:

- installation and particular requisites
- operating principles, instructions, precautions and risks
- manufacturer's specifications and instrument performance
- servicing and maintenance

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be carefully shaken before performing the test. Heat inactivation does not lead to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric or contaminated samples should be avoided.

The test cannot be performed on human plasma.

8. TEST PROCEDURE

Manual Technique

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Wash Buffer 10x (100 mL + 900 mL H₂O).
- Prepare the antigen by reconstituting the freeze-dried product directly with the conjugate (volume shown on label). In the case of reduced consumption of the Ag, reconstitute with Wash Buffer ready for use (1/10 of the volume shown on the label) and then 1/11 in the conjugate.

Dilute samples 1:101 distributing 10 µL of serum into 1 mL of diluent; dispense 100 µL of each diluted sample per well (duplicate testing is recommended). Place UNDILUTED controls 1 in a strip (100 µL in each well). The minimum requisite is 1 negative control, 2 cut-off and 1 positive control. Leave one well for the blank, performed using 100 µL of the substrate mixture.

Wells are covered with protective film and incubated for 45 minutes at 37°C. After washing four times for 30 seconds (300 µL), add 100 µL of immunocomplex (antigen/labelled monoclonal antibodies) to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film. The plate is washed again 4 times, as described above. Finally, the substrate is distributed, 100 µL/well.

After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µL of Stop Solution.

The adsorbance (O.D.) is read at 450 nm or 450/620 nm within 30 min.

9. SCHEME OF TEST PROCEDURE

Manual Technique

STEP 1 Place 100 µL of diluted sample/controls in the wells of the strips.

—
Incubate for 45 min. at 37°C

—
Wash 4 times (30" soak time, 300 µL)

STEP 2 Add 100 µL of immunocomplex to each well

—
Incubate for 45 min. at 37°C

—
Wash 4 times (30" soak time, 300 µL)

STEP 3 Add 100 µL of Substrate to each well

—
Incubate for 15 min. at R.T.

STEP 4 Add 100 µL of Stop Solution

—
Read absorbance at 450 nm within 30 min

10. VALIDATION OF THE TEST

Subtract the value of the blank (≤ 0.150) from all the other readings. The OD value of the Cut-off Control must be within 25% of the average value if tested in triplicate. Discard any anomalous values and recalculate the average. The Positive Control must have an OD of at least 1.5 times the Cut-off value. The ratio between Negative Control and Cut-off must be ≤ 0.6 . The O.D. Cut-off must be ≥ 0.2 at 450 nm, and ≥ 0.16 at 450/620 nm.

11. INTERPRETATION OF THE RESULTS

Qualitative results

If the OD of the sample is higher than the Cut-off the sample is positive for the presence of specific IgM.

Calculate the ratio between the OD of the sample and that of the Cut-off. The sample will be considered:

Positive: when the ratio is > 1.2

Doubtful: $\pm 20\%$ of the Cut-off

Negative: when the ratio is < 0.8 .

If the result is doubtful, repeat the test. If it remains doubtful, take a new blood sample.

12. LIMITATIONS OF THE PROCEDURE

The results must always be interpreted together with other clinical and diagnostic data.

13. ANALYTICAL SPECIFICITY

63 serum samples containing potentially interfering substances were tested:

- Serum from pregnant women (n=12)
- Parvovirus IgM (n=3)
- CMV IgM (n=12)
- HSV IgM (n=5)
- VCA IgM (heterophyl Ab) (n=5)
- Rubella IgM (n=5)
- Measles IgM (n=5)
- Mumps IgM (n=5)
- Rheumatoid Factor (up to 1080 UI/dl) (n=5)
- Bilirubin (up to 11 mg/dl)(n=5)
- Triglycerides (up to 1281 mg/dl) (n=5)
- Strongly hemolyzed samples (n=3).

No interference was observed in any case.

14. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

In an external experimentation, 154 samples were tested with the present kit in parallel with the routine method in use in the laboratory.

The samples were divided into 5 different panels:

- Panel 1: 56 samples from patients with a recent varicella infection
- Panel 2: 10 samples with reinfection as demonstrated by serum-conversion or by an increase in the antibody titer in CFT
- Panel 3: 15 samples from patients with a recent infection caused by Cytomegalovirus, containing IgM anti-CMV.
- Panel 4: 8 samples which were positive for infectious mononucleosis, characterized by the presence of IgG and IgM anti-VCA and by the absence of anti-EBNA antibodies.
- Panel 5: 65 samples from the general population, 61 of which contained anti-Varicella IgG.

A comparison of the performance of the two kits shows 99.4% agreement between the two (153/154), with a sensitivity of 100% (63/63) and a specificity of 98.8% (90/91).

15. PRECISION

Precision “in run”

Sample	VZC 1 (Negative< Cut Off)	VZC 2 (Positive>Cut Off)	VZC 3 (Positive)	Cut Off	Positive Control
n (replicates)	24	24	24	12	12
O.D.	0.115	0.511	1.769	0.292	1.755
CV%	10	5	5	9	6

Precision “between runs”

Sample	Index	
	Average	CV%
Positive Control	6.1	2
VZC1	0.4	6
VZC2	1.8	16
VZC3	6.9	4

Precision between batches

Sample	Index			Average	CV%
	Lot n. 034	Lot n. 035	Lot n. 036		
Positive Control	5.7	5.8	6.1	5.9	4
VZC1	0.3	0.4	0.4	0.4	16
VZC2	1.3	2.0	2.0	1.8	23
VZC3	4.9	6.2	6.5	5.9	15

16. TROUBLE SHOOTING GUIDE

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert point 4 for correct code).
		Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow).Repeat test
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.
	Inadequate washing	Ensure that wash apparatus works well
Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate Color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

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