CHORUS

TG

REF 86076

REF 86076/12

Prodotto da/Manufactured by/Fabricado por:
DIESSE Diagnostica Senese
 Strada dei Laghi 39
 53035 Monteriggioni (SIENA) – Italy

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(English)

1. INTENDED USE

QUANTITATIVE DETERMINATION OF THYROGLOBIN (TG) IN HUMAN SERUM, USING A DISPOSABLE DEVICE APPLIED ON THE CHORUS INSTRUMENTS.

2. INTRODUCTION

Thyroglobin (Tg) is a glycoprotein of high molecular weight (660 kDa) localized within the colloid of the thyroid follicle where it is synthesized under the influence of thyrotropin. It plays an essential role in the storage of iodine and acts as precursor for the synthesis of iodinated thyroid hormones: thyroxine (T4) and 3,5,3'-triiodothyronine (T3).

Elevated thyroglobulin serum concentrations have been reported in various thyroid diseases, such as hyperthyroidism, non-toxic goitre, thyroiditis and differentiated thyroid carcinoma.

The main clinical application for the Tg determination, however, is represented by the post-operative monitoring of differentiated thyroid carcinoma. The assay of Tg is used for the early detection or exclusion of metastases, tumour relapses and the follow-up of radioiodine treatments. Serum Tg is non detectable in patients who underwent total thyroidectomy and are free of metastases and tumour. These patients in true complete remission will not display Tg levels, even by endogenous TSH stimulation.

Consequently detectable Tg values in this group of patients are an important indication for still existing or newly developed neoplasia. Particularly if these detectable Tg values are increasing under a TSH-suppressive thyroid hormone treatment (Tg profiles).

In contrast, thyroglobulin levels in patients with medullary carcinoma or undifferentiated tumors remain within the normal range. Since thyroglobulin levels could also be elevated in other benign thyroid diseases, this test is not a criteria for the diagnosis of malignant thyroid tumours.

Determination of thyroglobulin is of prognostic value in Graves' disease patients undergoing therapy. Significantly elevated Tg levels at the end of a thyrostatic therapy are indication for a higher risk of relapse, whereas patients with continuous low thyroglobulin concentrations tend to continual recovery.

In case of patients that need an extended period monitoring, it is recommended to keep using the same analytical test, inserting in the same run the samples previously tested as well.

The presence of anti-Tg antibodies in the serum of the patient can influence the Tg assay.

3. PRINCIPLE OF THE METHOD

The Tg devices are ready to use for the assay of thyroglobulin (Tg), in automation in the CHORUS instruments. The test is based on the ELISA principle (Enzyme linked Immunosorbent Assay).

Monoclonal antibodies highly specific for Tg are bound to the solid phase. The Tg is bound to the well through incubation with the human serum examined.

After washings to eliminate the proteins which have not reacted, incubation is performed with the enzyme tracer conjugate, composed of anti-Tg monoclonal antibodies conjugated to horse radish peroxidase.

The unbound tracer conjugate is eliminated, and the peroxidase substrate added.

The colour which develops is proportional to the concentration of Tg present in the serum sample.

The disposable devices contain all the reagents to perform the test when applied on the Chorus instruments.

4. KIT COMPOSITION AND REAGENT PREPARATION

The kit is sufficient for 36 tests (REF 86076).


The kit is sufficient for 12 tests (REF 86076/12).

DD DEVICES

6 packages each containing 6 devices (REF 86076).

2 packages each containing 6 devices (REF 86076/12).

Use : **equilibrate a package at room temperature**, open the package and remove the devices required; replace the others in the bag with the silica gel, expel the air and **seal** by pressing the closure. Store at 2-8°C.

Description:

Position 8: Space for application of bar code label

Position 7: Empty

Position 6: MICROPLATE WELL coated with anti-Tg monoclonal antibodies

Position 5: Uncoated MICROPLATE WELL

Position 4: TMB SUBSTRATE 0.4 ml.

Contents: Tetramethylbenzidine 0.26 mg/mL and H₂O₂ 0.01% stabilized in 0.05 mol/L citrate buffer (pH 3.8)

Position 3: SAMPLE DILUENT 0.2 ml

Contents: Proteic solution in phosphate buffer, with sodium azide.

Position 2: CONJUGATE 0.35 ml

Contents: anti-Tg monoclonal antibodies labeled with peroxidase.

Position 1: EMPTY WELL in which the operator must place the **undiluted serum (140 µl)**

CALIBRATOR CALIBRATOR

Contents: Human thyroglobulin purified in proteinic stabilizing buffer. Ready for use.

CONTROL + POSITIVE CONTROL

Contents: Human thyroglobulin purified in proteinic stabilizing buffer. Ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- WASHING BUFFER AUTOIMMUNITY **REF** 86004
- CLEANING SOLUTION 2000 **REF** 83609
- SANITIZING SOLUTION **REF** 83604 - 83608
- CHORUS NEGATIVE CONTROL/SAMPLE DILUENT **REF** 83607
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 50-200 µl solution
- Disposable gloves
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C. In the case of storage at an incorrect temperature the calibration must be repeated and the run validated using the positive control (see point 9, Test validation).

The expiry date is printed on each component and on the kit label.

Reagents have a limited stability after opening and/or preparation

DEVICE	8 weeks at 2/8°C
POSITIVE CONTROL	8 weeks at 2/8°C
CALIBRATOR	8 weeks at 2/8°C

6. PRECAUTIONS

ONLY FOR DIAGNOSTIC USE IN VITRO.

Attention:

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples, calibrators and strips once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens. Wash hands thoroughly after placing the devices in the CHORUS instrument.
2. If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
4. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

Bring the devices to room temperature (18-30°C) before use; use within 60 min.

1. **Discard devices which show the substrate (well 4) blue colored.**
2. When adding the sample to the well, verify that it is perfectly distributed on the bottom.
3. Check for the presence of the reagents in the device and that the device is not damaged; do not use devices which are lacking a reagent.
4. The devices are for use with the Chorus instrument; the instructions for use must be carefully followed and the instrument operating manual must be consulted.
5. Check that the Chorus instrument is set up correctly (see Chorus Operating Manual).
6. Do not alter the bar code placed on the handle of the device in order to allow correct reading by the instrument
7. Defective barcodes can be inserted manually in the instrument.
8. Do not expose the devices to strong light or to hypochlorite vapours during storage and use.
9. The use of strongly hemolyzed samples, or samples presenting microbial contamination may all constitute sources of error.
10. Before inserting the devices in the instrument, check that the reaction well does not contain foreign bodies.
11. Pipette the test serum (**140 µl**) in well 1 of the device (see figure).
12. Do not use the device after the expiry date.
13. **Make sure that the instrument is connected to the Washing Buffer Autoimmunity REF 86004.**

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric, hemolyzed or contaminated samples should be avoided.

8. PROCEDURE

1. Open the package (on the side containing the pressure-closure), remove the number of devices required and seal the rest in the bag after expelling the air.
2. Check the state of the device according to the indications reported in the Analytical Warnings points 1 and 8.
3. Dispense 140 µl of undiluted test serum in well no. 1 of each device; at each change of batch, use a device for the calibrator.
4. Place the devices in the instrument. Perform the calibration (if necessary) and the test as reported in the Chorus Operating Manual.

9. TEST VALIDATION

Use the POSITIVE CONTROL SERUM to check the validity of the results obtained. It should be used as reported in the operating manual. If the instrument signals that the positive control has a value outside the acceptable range, the calibration must be repeated. The previous results will be automatically corrected.

If the result of the control serum continues to be outside the acceptable range, contact the Scientific Support.

Tel: 0039 0577 319554
 Fax: 0039 0577 366605
 email: scientificsupport@diesse.it

10. INTERPRETATION OF THE RESULTS

The Chorus instrument gives a result in ng/ml (BCR, Brussels), calculated on the basis of a lot-dependent graph stored in the instrument. **Attention: do not consider the qualitative report N; D; P on the instrument.**

In euthyroidic patients with normal TSH: 50 ng/ml (expected average superior limit of normality)

After total thyroidectomy with suppressed TSH: TG <1 ng/ml TG (1).

11. LIMITATIONS OF THE PROCEDURE

The test result should be used in conjunction with information available from the evaluation of the case history or other diagnostic procedures. The presence of anti-Tg antibodies can influence the test results; for this reason the parallel assay of such auto-antibodies is recommended.

12. ANALYTICAL SENSITIVITY : 0.2 ng/ml

13. FUNCTIONAL SENSITIVITY : 0.4 ng/ml

14. CALIBRATION RANGE : 0.4-100 ng/ml

Standardized on the basis of the European Standard for Human Thyroglobulin CRM-457 according to the Community Bureau of Reference of the European Commission.

In case of out-range samples (>100 ng/ml) repeat the test after diluting with Sample Diluent **REF** 83607.

15. ANALYTICAL SPECIFICITY

5 samples (3 negative and 2 positive) containing the following interfering substances were tested:

Rheumatoid Factor (1.25 – 20 IU/ml)

Bilirubin (5 mg/dl – 100 mg/dl)

Triglycerides (50 mg/dl – 3000 mg/dl)

Hemoglobin (0.6 mg/dl – 10 mg/dl)

The presence in the tested serum of the interfering substances listed above did not alter the results of the test.

16. LINEARITY

3 chosen sera have been tested with this kit and found to dilute linearly. However, due to the heterogeneous nature of human autoantibodies there might be samples that do not follow this rule.

Sample A (European Standard CRM-457)		
Dilution factor	Expected concentration (ng/ml)	Result (ng/ml)
1 :1	127.4	127.4
1 :2	63.7	60.6
1 :10	12.7	12.2
1 :40	3.2	3.0
1 :100	1.3	1.3
1 :200	0.6	0.5

R²= 0.9995

Sample B		
Dilution factor	Expected concentration (ng/ml)	Result (ng/ml)
1 :1	75.5	75.5
1 :2	37.8	44.9
1 :4	18.9	22.4
1 :8	9.4	11.9
1 :16	4.7	5.6
1 :32	2.4	2.4

R²=0.9907

Sample C		
Dilution factor	Expected concentration (ng/ml)	Result (ng/ml)
1 :1	75.3	75.3
1 :2	37.7	42.1
1 :4	18.8	20.5
1 :8	9.4	10.7
1 :16	4.7	5.0
1 :32	2.4	2.1

R²=0.9963

17. METHOD COMPARISON

Cut-off 1 ng/ml:

In an experimentation 67 samples were tested with the Diesse kit and with another commercial kit. The results are summarized in the following table:

		Reference		Total
		+	-	
Diesse	+	28	2	30
	-	1	36	37
	Total	29	38	67

Percent Positive Agreement (~Diagnostic Sensitivity) = 96.6 CI_{95%}= 82.8 99.2
 Percent Negative Agreement (~Diagnostic Specificity)= 94.7 CI_{95%}= 82.7 98.4

Cut-off a 50 ng/ml:

In an experimentation 151 samples were tested with the Diesse kit and with another commercial kit. The results are summarized in the following table:

		Reference		Total
		+	-	
Diesse	+	63	5	68
	-	2	81	83
	Total	65	86	151

Percent Positive Agreement (~Diagnostic Sensitivity) = 96.9 CI_{95%}= 89.4 99.1
 Percent Negative Agreement (~Diagnostic Specificity)= 94.2 CI_{95%}= 87.1 97.4

18. PRECISION

The precision has been evaluated on the basis of 3 internal standards (low, medium and high assay) in replicates: within the same run (Intra-Assay), between different runs performed during different days (Inter-Assay) and between different instruments.

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