DIESE Diagnostica Senese SpA

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Tel. ++39/0577/ 587111 - Fax ++39/0577/318690

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When the germ count is performed on urine samples, the presence of antibacterial substances in the sample can determine false negative results in 20-27% of cases (1,2).

The test can also be performed using other liquids such as drinking water, milk or meat homogenate.

3. PRINCIPLE OF THE TEST

Bacillus Subtilis ATCC 6633 spores, chosen for their sensitivity to all antibiotics, are dried in wells containing medium and growth indicator.

The test urine dissolves the medium. After incubation at 37°C the negative samples produce a red pigmentation, while in the positive ones the medium remains completely yellow.

4. KIT CONTENTS AND REAGENT PREPARATION

- Reagents are sufficient for 96 determinations.
- Bring reagents to room temperature before use.

MT PLATE PLATES CONTAINING SPORES 2 plates with 6 x 8 wells**Contents:** Each plate contains B. subtilis spores, eugon medium and growth indicator, in freeze-dried form.**OTHER MATERIALS**

- 2 adhesive films
- 1 sealable envelope for use as incubation chamber
- 1 diagram of the wells.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
PLATE	1 month

6. PRECAUTIONS**FOR IN VITRO DIAGNOSTIC USE ONLY. STORE AT 2-8°C**

All materials of biological origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling reagents and samples.

Waste disposal: urine samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. If a reagent comes in contact with the skin or eyes, wash thoroughly with water. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should

not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. **Allow all reagents and samples to come to room temperature (18-30°C) before use.** Immediately after use return reagents to the recommended storage temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable.
4. Do not use undiluted urine samples which have been stabilized by the presence of preservatives such as boric acid or other substances.
5. Do not prolong the incubation over 24 hours.

7. TYPE AND STORAGE OF SAMPLE

Urine collected in containers without preservatives. The sample must be tested within 8 hours of collection. For longer storage, freeze at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers.

8. PROCEDURE

Use scissors to cut the number of wells required for the test. Dispense one drop (50 µl) of sample in each well and close with the adhesive film pressing lightly on the sides of the device. Incubate for at least 18 hours at 37°C in the sealed envelope or in a humidity chamber. Read the results and observe whether pigmentation has developed.

9. TEST VALIDATION

Distribute 50 µl of saline or sterile water in at least 1 well. Bacterial development must be found, with a colour change to red.

10. INTERPRETATION OF THE TEST

POSITIVE PAR

The wells contains a clear yellow solution.

NEGATIVE PAR

Diffuse dark red colouring or development of several coloured colonies adhered to the plastic.

11. TEST LIMITATIONS

The test may give a false negative result in the presence of samples of patients treated with antibiotics but infected by resistant strains.

The test gives false positive results in case of urine containing bacteriostatics as borate.

12. ANALYTICAL SENSITIVITY

The test sensitivity depends on the amount of sample. Use samples without dilution to obtain the MAXIMUM TEST SENSITIVITY. For LOWER SENSITIVITY, 25 µl of sterile water and 25 µl of sample may be added to each well.

When samples are collected in the presence of stabilizing agents (e.g. using the Becton Dickinson tubes cod. 4948) add to the wells 90 µl of sterile water and 10 µl of test sample.

13. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

In an external experimentation in which 537 undiluted urine samples collected without the use of preservatives were tested in comparison with a test on plates of agar germs (E. coli, S. aureus), a correlation of 90.1% was found, with 23 false positive and 26 false negative results.

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